



For more information  
on this product please  
scan QR code

---

## Instructions For Use

---

Návod k použití • Gebrauchsanweisung • Instrucciones de uso  
Mode d'emploi • Οδηγίες χρήσης • Istruzioni per l'uso  
Gebruiksaanwijzing • Instruksjoner for bruk • Instrukcja użycia  
Instruções de uso • Användningsinstruktioner

## Disposable Insufflator Filters

STERILE EO



Single Use



CE  
0805



Keep away  
from sunlight



Do not use  
if packaging  
is opened or  
damaged



Keep Dry

---

**Caution - Read all instructions prior to use**

The following instructions relates to the Fairmont's Disposable Insufflator Filters

**Indications / Intended Use – DIF Range**

The Disposable Insufflator Filter (DIF) range is indicated for filtering CO<sub>2</sub> Gas from insufflator to the patient during laparoscopic cases.

1. Peel DIF pouch open using standard aseptic techniques and pass on inner bag to sterile nurse.
2. Sterile nurse to remove from inner bag, screw to entry port device Luer lock.
3. Pass insufflator filter and connector to circulating nurse using standard aseptic technique.
4. Circulating nurse to connect insufflator filter to insufflator.

DIF7002 / 7300 – Push silicone tubing firmly over outlet barb.

DIF7004 – Push plastic 22mm female filter fitting directly over 22mm outlet.

DIF7005 – Screw female Luer connector directly into male Luer outlet.

DIF7007 – (Sopro insufflators only) Connect Luer lock and universal adaptor / tubing to outlet.

DIF7009 - Screw male Luer connector directly into female Luer outlet.

DIF7201 - Push universal connector onto outlet.

DIF7205 - Push plastic 15mm male filter fitting directly into 15mm outlet.

DIF7206 - Push male connector into female outlet.

*Operate insufflator only after both the insufflator connector and patient entry port connector are fully connected.*

1. Peel DIF pouch open using standard aseptic techniques and pass on filter to sterile nurse.
2. Sterile nurse to push silicone tubing firmly over filter outlet barb.
3. Pass insufflator filter and connector to circulating nurse using standard aseptic technique.
4. Circulating nurse to connect insufflator filter to insufflator.

DIF7001 – Push silicone tubing firmly over outlet barb.

DIF7003 – Push plastic 22mm female filter fitting directly over 22mm outlet.

**Precautions**

Do not use if package is compromised/damaged.

This device may only be used by suitably trained healthcare personnel.

Not to be used for any other purpose.

All devices within the DIF range are single use only, do not re-sterilize.

This device is sterile (ETO). Please note, this IFU pertains to Sterile devices sold by Fairmont Medical. If the device has code "XXXX NS", please refer to 3rd party manufacturer for sterility specifications.

**Materials**

All filters within the DIF range are sterile and the filter membrane is hydrophobic PTFE.

All components and packaging are Latex free.

All PVC tubing is Phthalate-free PVC.





Head Office  
11-13 Scoresby Road  
Bayswater Victoria  
3153 Australia  
Ph +61 3 9720 8840  
Fx +61 3 9720 8860  
[www.fairmontmedical.com.au](http://www.fairmontmedical.com.au)



Fairmont Medical Products Ltd  
Elizabeth House  
28 Baddow Road Chelmsford  
Ph +44 (0) 1245 206800  
[www.fairmontmedical.com](http://www.fairmontmedical.com)

Fairmont Medical Products Pty Ltd is certified to ISO13485:2003

PM026 Version 2.1

Issued by E. NITEROS 27/03/2019