

Advazorb® Sacral

Hydrophilic foam sacral dressing with film backing


Fact sheet

Product description

Advazorb® Sacral is a hydrophilic foam dressing for wounds with low to medium levels of exudate around the sacral area, with a soft silicone wound contact layer. Advazorb® Sacral will only adhere to dry skin surrounding the wound and can be easily lifted for adjustment without losing its adherent properties.

Use

Apply to wound (white silicone face down) allowing a minimum of a 2cm border overlap around the wound area. Do not cut Advazorb® Sacral. Can be left in place for up to 7 days but should be changed immediately if strike through is apparent or when dressing reaches capacity. If exudate is visible around the edges of the dressing this is a clear indication that a dressing change is required. Clinical observation is necessary to determine required frequency of change.



Contra-indications

Arterial bleeds, heavily bleeding wounds and vascular fungating tumours.

Storage

Store at room temperature, out of direct sunlight and keep dry.

Advazorb® Sacral ordering information

Size	Stock code	Pack size
17cm x 19cm	CR4401	10
22cm x 23cm	CR4402	10

TOP TIP

Ideal for use on low-moderately exuding wounds, with the soft silicone contact layer extending across the entire surface area of the dressing, not just the border, providing secure but gentle adhesion.



Features

- ✓ Sterile Hydrophilic PU foam
- ✓ Low friction thin film backing
- ✓ Bacteria proof backing
- ✓ Soft and conformable
- ✓ High fluid retention
- ✓ Anatomically shaped

Indications

Suitable for exuding on the sacral area, for example pressure wounds including:

- ✓ For all exuding wounds on the sacral area, for example pressure ulcers.

Care and use symbols

- Read instructions for use
- Single use only
- Do not use if packaging is broken or damaged
- Sterilised using ethylene oxide
- CE marked medical device
- Keep dry
- Store out of direct sunlight
- Do not re-sterilize

Distributed throughout Australia by:

Defries Industries
Hospitals Trust Defries

T: 03 9706 3600 | Freecall: 1300 550 278
Email: sales@defries.com.au www.defries.com.au

