SPIGGLE & THEIS Medizintechnik GmbH
 Burghof 14
 51491 Overath / Germany
 Tel.: + 49 (0) 2206 9081-0
 Fax: + 49 (0) 2206 9081-13
 www.spiggle-theis.com
 info@spiggle-theis.com

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ENGLISH

1. Product description and intended use

ENT-single-use suction tubes:

ENT-single-use suction tubes are designed to aspirate fluids/tissue, e.g. from the surgical site, for the ENT area, and for use in head and neck surgery. The product is connected to a suction pump via suction tubing and is vacuum operated. Different versions of the ENT-single-use suction tubes are available to suit the individual surgical site and anatomy of the patient.

The ENT-single-use suction tubes for otology are used together with the Single use suction handle and are available with two different types of connection caps that allow it to be attached to a Luer connection and a Luer lock.

The connectors that attach to a Luer connection are color coded, making it easier to differentiate between the various suction cannula diameters (Ø 0.4 mm: grey, Ø 0.6 mm: blue, Ø 0.7 mm: black, Ø 0.8 mm: green, Ø 1.0 mm: off-white, Ø 1.2 mm: pink, Ø 1.4 mm: red, Ø 1.6 mm: white and Ø 2.0 mm: green).

ENT disposable suction cannulas for otology in bulk are designed to aspirate fluids/tissue in the ENT area. They are not suitable for use in sterile environments!

ENT disposable suction cannulas for otology with different suction cannula diameters (Ø 1.0 mm: off-white, Ø 1.2 mm, Ø 1.4 mm, Ø 1.6 mm, Ø 2.0 mm, Ø 2.5 mm, Ø 3.0 mm and Ø 4.0 mm) are available in bulk in a sterile pack of 25.

ENT-single-use suction tubes for use in rhinology are available in different lengths and diameters and have a handle with a conical Luer connection that attaches to the suction tubing. The suction power can be adjusted using the aspiration control hole in the handle of the product. Product variants for use on paranasal sinuses are available with an olive at the distal end to protect mucous membranes and highly sensitive structures.

ENT-single-use suction tubes for use in laryngology are available in different lengths and diameters and have a handle with a conical Luer connection that attaches to the suction tubing. The suction power can be adjusted using the aspiration control hole in the handle of the product. The product also has an opening at the distal end for making additional adjustments to the suction power in order to protect tissues.

ENT-single-use suction tubes are designed for single use.

Single use suction handle:

Single use suction handle with aspiration control hole for adjusting the suction power as a connecting piece between the suction device and suction tubing. The single use suction handle is used to regulate the suction power of the suction device for use in the ENT area and for head and neck surgery.

The single use suction handle has a conical Luer connector at the distal end that attaches to ENT-single-use suction tubes of different sizes and a connector for the suction tubing at the proximal end. The aspiration control hole for adjusting the suction power is in embedded in a finger recess.

The product is designed for single use.

Material:

ENT-single-use suction tubes:

Cannula: rustproof stainless steel (1.4301)

Connector: Polypropylene (PP) or styrene-butadiene copolymer (ABS)


Protective tubing: LDPE and PVC

Single use suction handle: ABS

ENT-single-use suction tube, otology.

REF:

301004, 302006, 302014, 303014, 304010, 302010U, 303014U,
301007, 302007, 302016, 303020, 304014, 302012U, 303020U, 304014U,
301010, 302008, 302020, 303025, 304020, 302014U, 303025U, 304020U,
301014, 302010, 303004, 303030, 304025, 302016U, 303030U, 304025U,
302004, 302012, 303007, 303040, 304030, 302020U, 303040U, 304030U

ENT-single-use suction tube, rhinology.

REF:

3012020, 3012035, 3014020, 3014035, 3016020, 3016035,
3012025, 3012040, 3014025, 3014040, 3016025, 3016040,
3012030, 3012050, 3014030, 3014050, 3016030, 3016050

ENT-single-use suction tube, rhinology (with olive).

REF:

3010025, 3010030, 3010040

ENT-single-use suction tube, Laryngology.

REF:

3025020, 3025030

Single use suction handle.

REF:

301000

2. Indication**Indications in rhinology:**

- Septum deviations requiring septoplasty
- Septorhinoplasty
- Plastic surgery on nasal soft tissues
- Surgery on the orbital cavity
- Fractures of the facial skeleton and the nasal bone
- Trauma to the facial skeleton
- Procedures on the paranasal sinuses
- Tumours
- Turbinate surgery
- Conservative procedures such as the removal of foreign bodies

Indications in laryngology:

- Laryngoscopic procedures such as
 - Hoarseness
 - Suspected acute or chronic laryngitis
 - Changes to the vocal cords such as vocal fold polyps (benign regeneration)
 - Malformations of the laryngeal region
 - Tumours
 - Injuries to the larynx
 - Suspected vocal fold paralysis
- Laryngeal surgery (e.g. vocal cord operation and tumour surgery)
- Removal of foreign bodies

Indications in otology:

- Tympanoplasty
- Stapedectomy

- Mastoidectomy
- Facial nerve decompression
- Procedures on the inner and outer ear canal
- Cochlear implant
- Acoustic neuroma
- Vestibular neurectomy
- Traumas
- Tumours
- Removal of cerumen and foreign bodies

Indications for use in otology for bulk products:

- Removal of cerumen and foreign bodies
- External ear canal interventions in non-sterile environment

3. Contraindication

- Allergic reactions to the manufacturing materials.
- There are no other known contraindications when the product is correctly used.

4. Side effects

- Injuries to the nerves, vessels and tissues
- Bleeding
- Iatrogenic injuries
- Temporary hearing problems

5. Directions for use**5.1. General directions for use**

This product should only be used by experienced medical professionals (e.g. a qualified medical specialist).

The accompanying information, such as the instructions for use, is not a substitute for basic medical and technical skills. Such skills must be acquired by the user through specialised training, if necessary. The acquisition of medical skills and their diagnostic and therapeutic consequences are the sole responsibility of the user of the product.

5.2. Special directions for use

- Remove the ENT-single-use suction tube from its protective sleeve and connect the proximal end to a suction pump via suction tubing. Verify that all parts are securely connected before use.
- The ENT-single-use suction tube can be connected to a commercially available syringe and used to gently remove secretions.
- The ENT-single-use suction tube and the single use suction handle are designed for aspiration and should only be used to for this purpose, not for irrigation.
- Hygiene regulations should always be observed in order to prevent infections.
- Do not force the ENT-single-use suction tube by twisting or levering, which can lead to damage or breakage of the product and/or injury to the patient.
- If there is a blockage in the ENT-single-use suction tube or the single use suction handle, the blockage can be removed by connecting it to a commercially available syringe and flushing it out with sterile NaCl solution. If this does not work, the obstructed product must be discarded.
- Spare products should be readily available in order to prevent delays in use in the event of a blockage.
- ENT disposable suction cannulas in bulk are not suitable for use in a sterile environment.

Caution: Please exercise extreme caution when using an ENT-single-use suction tube with a drill. There is a risk of damage to the ENT-single-use suction tube.

Should damage occur to the product as a result of use of the ENT-single-

use suction tube with drills (e.g. abrasion of the material), SPIGGLE & THEIS Medizintechnik GmbH recommends that you remove any remaining dust from the abrasion in the OP site through irrigation.

Furthermore, the damage to the product may lead to the formation of sharp edges. Should this occur, the ENT-single-use suction tube must be replaced.

6. General handling

It is very important to avoid product contamination. The product must be stored in its sealed protective packaging. The packaging must not show any signs of damage. The protective packaging should only be opened just before removing the product. Before opening the package, it should be checked for damage and whether it is properly sealed, since its sterility may be compromised. In addition, the product must be inspected for any visible or functional damage insofar as possible.

The product is intended for single use only. If the sealed seam of the sterile packaging is open, the contents cannot be returned to the manufacturer.

ENT disposable suction cannulas for otology in bulk are no longer sterile after the sterile packaging is opened. Contamination of the products that are not used immediately should be avoided to the extent possible. The product must be inspected before use for any visible or functional damage insofar as possible. The product is intended for single use only.

7. Packaging and sterility

The product packaging consists of:

- Primary packaging/sterile barrier system
- Secondary packaging/protective packaging

The packaging complies with the standards in force. When intact, the packaging protects the product from external influences and ensures product sterility during storage.

7.1. Handling of sterile packaging

The appropriate aseptic technique should be used when removing the product from the packaging.

7.2. Sterilisation

The product is EO (ethylene oxide) sterilised in its protective packaging (see package illustration).

7.3. Resterilisation

The reprocessing and/or reuse of a previously used or non-sterile product is not permitted.

The processing of single-use products is inconsistent with patient welfare. It is impossible to determine whether cleaning and sterilisation will result in defects such as kinks, tears and material changes, microorganisms, endotoxins and/or chemical residues, and therefore changes in functionality.

The user bears full responsibility for processing and reusing the product. When the user processes the product, he becomes its manufacturer, according to the MDD 93/42/EEC and/or MDR 2017/745. He therefore assumes liability for any damages resulting from the use of a single-use product that has been damaged through processing.

8. Storage

The product should be stored in a dust-free area protected from moisture and contamination and away from direct sunlight. The product may no longer be used after the expiry date.

9. Disposal

When disposing of the product, care should be taken to prevent the risk of injury and infection. Contaminated products should be placed in a hazardous

waste container and handled in such a way as to eliminate the risk of contaminating third parties.

10. Obligation to report

Any serious incidents that occur in conjunction with the product must be reported to the manufacturer and the competent local authority of the EU Member State where the user and/or patient is resident.

11. Declaration of conformity

We declare under our sole responsibility that the products as supplied comply with the requirements of the MDD 93/42/EEC and/or MDR 2017/745. The products are in accordance with the relevant provisions of the MDD 93/42/EEC and/or MDR 2017/745 and bear the CE marking.

Revisions and typographical errors reserved.

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© SPIGGLE & THEIS Medizintechnik GmbH, 02_210727.
© SPIGGLE & THEIS Medizintechnik GmbH · 02_210727.