	DO NOT REUSE
	MEDICAL DEVICE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE C USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION
 Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INDICATION AND USAGE
The DeRoyal® angiographic control syringes are intended to be used for the intra-arterial or intravenous administration and aspiration of radiographic contrast media, saline and bodily fluids.

- ⚠️ WARNINGS**
- To eliminate the potential for air aspiration, severe back and forth bending motions should not be applied to this device during use.
 - All connections should be hand tightened, without over-tightening, to ensure their security and prevent damage to the device.
 - Examine the products carefully for entrapped air and fully de-bubble prior to injection.
 - DO NOT** store fluid in the syringe. Inject fluid immediately after filling.
 - DO NOT** use this syringe to inflate a balloon catheter.
 - DO NOT** use product for biopsy procedure.
 - DO NOT** use an expired device.
 - Contents supplied non-sterile should be sterilized prior to use.
 - The tip of the syringe should be kept pointing down and the handle of the manifold should be slightly raised so any air bubbles can float up to avoid injecting air into the system.
 - Connections should be made only with devices compliant with ISO 80369-7:2016 for luer connections.

- CAUTIONS**
-  Product should not be used if sterile packaging is damaged or opened. If damage is found, remove the syringe from the operative area to prevent unintended use, and return to the manufacturer.
 - This device is for single patient use only. It is not intended for disinfection and subsequent re-use, which may result in device failure or create the risk of contamination.
 - This device has not been evaluated for reprocessing or re-sterilization. Reprocessing or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.

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- INSTRUCTION FOR USE**
- Remove from packaging using aseptic technique.
 - Inspect device for damage.
 - Prep syringe by moving the plunger back and forth several times to ensure it is moving appropriately within the barrel.
 - To fill the syringe, pull the plunger back.
 - Connect syringe to other appropriate ISO 80369-7 compliant devices. **DO NOT** overtighten.
 - To inject, use plunger and graduation markings to control injection speed.
 - Dispose of according to facility protocol.

CLEANING AND/OR MAINTENANCE
Please refer to hospital approved guidelines for cleaning and disinfecting of the device.

STERILIZATION INFORMATION
DeRoyal intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization. These products are to be sterilized by ethylene oxide, utilizing the following paramaters: Temperature: 120°F - 135°F (48.9°C – 57.2°C), Ethylene Oxide Gas Concentration: 625 – 800 mg/l, Exposure Time: 4 hrs minimum (Preconditioning: Temperature 100°F -120°F (37.8°C – 48.9°C), Humidity 55-75%, Time 18hrs minimum), (Aeration: Temperature 100°F -135°F (37.8°C – 42.2°C), Time 18hrs minimum). Other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

STORAGE AND TRANSPORT CONDITIONS
Do not store product at extreme temperatures or in a moist/damp environment.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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PART#971923G I REVISED 10/23

©2023 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

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