DEROYAL® MANIFOLDS AND STOPCOCKS

|                                                                                   |                                                                                                                              |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
|  | DO NOT REUSE                                                                                                                 |
| STERILE/EO                                                                        | STERILIZED WITH ETHYLENE OXIDE                                                                                               |
|                                                                                   | NOT MADE WITH NATURAL RUBBER LATEX                                                                                           |
| RX ONLY                                                                           | FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER. |

**INTENDED USE**  
DeRoyal® manifolds and stopcocks are intended to be used during cardiac catheterization procedures for intra-arterial and intravenous administration of contrast, saline or radiographic contrast media. The product is intended for use by a properly licensed practitioner.

| PRODUCT                                                                                         | MAXIMUM PRESSURE RATING* (PSI) | MAXIMUM PRESSURE RATING (KPA) |
|-------------------------------------------------------------------------------------------------|--------------------------------|-------------------------------|
| Low Pressure Manifolds                                                                          | 250 PSI                        | 1723 kPa                      |
| Medium Pressure Manifolds                                                                       | 550 PSI                        | 3792 kPa                      |
| High Pressure Manifolds (designed for procedures requiring small catheters and power injection) | 800 PSI                        | 5515 kPa                      |

-  **WARNINGS**
- All connections should be hand tightened, without over-tightening, to ensure their security and prevent damage to the device.
  - Check all connections to ensure they are tight and secure to prevent air from entering the system.
  - Examine the products carefully for entrapped air and fully de-bubble prior to injection.
  - Contents supplied non-sterile should be sterilized prior to use. Sterilization should be by Ethylene Oxide according to the manufacturer's validated sterilization parameters.
  - Connections should be made only with devices compliant with ISO 80369-7 for luer connections.
  - Low Pressure Manifolds are not intended for power injection as this may result in leaking and/or allow air bubbles into the system.
  - The PSI ratings listed are for injections made through the end port. The side port(s) are not rated and should not be used for high pressure injections.
  - Inspect packaging for any damage and product for any damage, contaminate or missing parts prior to use.
- PRECAUTIONS**
- DO NOT** reuse. This device is not intended for disinfection and subsequent re-use, which may result in device failure or create the risk of contamination.
  - The device has not been evaluated for re-processing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure which could result in patient illness, injury or death.

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**STERILIZATION INFORMATION**  
DeRoyal intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization. These products are to be sterilized by ethylene oxide, utilizing the following paramaters: Temperature: 120°F - 135°F (48.9°C – 57.2°C), Ethylene Oxide Gas Concentration: 625 – 800 mg/l, Exposure Time: 4 hrs minimum (Preconditioning: Temperature 100°F -120°F (37.8°C – 48.9°C), Humidity 55-75%, Time 18hrs minimum), (Aeration: Temperature 100°F -135°F (37.8°C – 42.2°C), Time 18hrs minimum). Other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

**\*PLEASE NOTE:** All pressure ratings are based on DeRoyal's sterilization cycle. Other sterilization may affect product performance.

|                                                                                   |                         |
|-----------------------------------------------------------------------------------|-------------------------|
|  | KEEP DRY                |
|  | KEEP AWAY FROM SUNLIGHT |

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

**WARRANTY**  
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.