



ATTENTION:

See instructions for use.



Caution.

Storage: Exposure to elevated temperature and ultraviolet light should be avoided during storage.



Sterilized using ethylene oxide.



Not made with DEHP.



Not made with natural rubber latex.



Medical device.



Authorised representative in the European Community.



Catalogue number.



Batch code.



Date of Manufacture.



Do not resterilize.



Do not use if package is damaged.



Do not reuse.



Use by.



Manufacturer.



Keep away from sunlight.



Keep dry.



Upper limit of temperature: 49°C.



Pyrogen-free.

SUCTION CATHETER INSTRUCTIONS FOR USE

These instructions apply to the use of the following suction catheters:

DSC7008	Y-Type Light Blue 8FG 100mm
DSC7010	Y-Type Black 10FG 100mm
DSC7012	Y-Type White 12FG 100mm
DSC7410	Y-Type Black 10FG 430mm
DSC7412	Y-Type White 12FG 430mm
DSC7510	Y-Type Black 10FG 560mm
DSC7512	Y-Type White 12FG 560mm

Description

Manufactured from clear polyvinyl chloride (PVC) incorporating the following features:

- Polished eyelets.
- Available in two different tip configurations (elephant foot tip and standard tip). The elephant foot tip is designed to enlarge the area touching patients and reduce puncture wound on patients.
- Colour coded funnel identifies.

Indications for use

The suction catheters, in combination with aspirator, are intended for suction of the phlegm or secretory juice in the respiratory tract.

Contraindications

The product shall be forbidden to insert from nose if the patient has skull base fracture.

Warnings

To help reduce the potential risk of infection and/or other complications, do not re-use. Dispose of appropriately after procedure.

Precautions

1. For single patient use only.
2. Use the product during the validity period.
3. The product is forbidden for use if the individual package is damaged.
4. Not intended for reprocessing.
5. Avoid freezing and excessive heat.
6. The inflammation patients shall cautiously use.

Instructions

1. Open the individual package, take away Suction catheters.
2. Insert the patient end into the cavity, or trachea.
3. Connect the connector with the suction machine/device (aspirator).
4. Adjust the pressure and take suction.
5. Control the suction pressure by press the control button (if any).



Fairmont Medical Products Pty Ltd
11-13 Scoresby Road
Bayswater Vic 3153
Australia



Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands



APPROVED ARTWORK/批准版面	
Controlled Doc. /受控文件	SC-005IFU 001 Rev: A 22/08 吸痰管说明书客户005版面
Dimension(mm) /尺寸(毫米)	148.5 × 210 (正反两面)
Color/颜色	Black (黑色)
Validation /有效期	四年
Approval Ref. /批准索引	Mail dated on Aug. 4, 2022 2022年8月4日确认