LIST NO.:	011-OP560-01
COMMODITY NO.:	DF-4565
DATE PREPARED:	8-24-2012
ARTIST:	KENT

MATERIAL
Paper: 60 pound White if printed in-house or
Paper: 80 grams Bond (50 pound White) if printed outside vendor
Ink: Black on uncoated surface
SPECIFICATIONS:
Panel Length: 10 1/2 inches
Sheet Width: 7 1/2 inches
Folded Size: 3 3/4 inches x 10 1/2 inches
Center Folded Fold Style: 2 Panel
MADE FROM: EN-1455, DF-4191

VENDOR TO VERIFY PAPER STOCK:

VENDOR:
NOSCO, INC.
651 SOUTH UTICA STREET
WAUKEGAN, IL. 60085

ALTERNATE:
ICU MEDICAL, INC.

CARDIAC CATHETERISATION KIT

icumedical

LIST NO.:	011-OP560-01
COMMODITY NO.:	DF-4565
DATE PREPARED:	8-24-2012
ARTIST:	KENT

MATERIAL
 Paper: 60 pound White if printed in-house or
 Paper: 80 grams Bond (50 pound White) if printed outside vendor
 Ink: Black on uncoated surface
 SPECIFICATIONS:
 Panel Length: 10 1/2 inches
 Sheet Width: 7 1/2 inches
 Folded Size: 3 3/4 inches x 10 1/2 inches
 Center Folded Fold Style: 2 Panel
 MADE FROM: EN-1455, DF-4191

VENDOR TO VERIFY PAPER STOCK:

VENDOR:
 NOSCO, INC.
 651 SOUTH UTICA STREET
 WAUKEGAN, IL. 60085

ALTERNATE:
 ICU MEDICAL, INC.

CARDIAC CATHETERISATION KIT

Instructions for Use	3-4
Does Not Contain DEHP	
Mode d'emploi (FR).....	5-6
Ne contient pas de DEHP	
Instrucciones de uso (ES)	7-8
No Contiene DEHP	
Instruções de uso (PT)	9-10
Não Contém DEHP	
Bruksanvisning (SV).....	11-12
Innehåller Inte DEHP	
Käyttöohjeet (FI).....	13-14
Ei Sisältää DEHP	
Brugsvejledning (DA)	15-16
Indeholder ikke DEHP	
Gebruiksaanwijzing (NL)	17-18
Bevat geen DEHP	
Gebrauchsanweisung (DE)	19-20
Enthält kein DEHP	
Istruzioni per l'uso (IT).....	21-22
Non Contiene DEHP	
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ (EL).....	23-24
Δεν περιέχει DEHP	

LIST NO.:	011-OP560-01
COMMODITY NO.:	DF-4565
DATE PREPARED:	8-24-2012
ARTIST:	KENT

MATERIAL

Paper: 60 pound White if printed in-house or
 Paper: 80 grams Bond (50 pound White) if printed outside vendor
 Ink: Black on uncoated surface
SPECIFICATIONS:
 Panel Length: 10 1/2 inches
 Sheet Width: 7 1/2 inches
 Folded Size: 3 3/4 inches x 10 1/2 inches
 Center Folded Fold Style: 2 Panel
 MADE FROM: EN-1455, DF-4191

VENDOR TO VERIFY PAPER STOCK:

VENDOR:
 NOSCO, INC.
 651 SOUTH UTICA STREET
 WAUKEGAN, IL. 60085

ALTERNATE:
 ICU MEDICAL, INC.

CARDIAC CATHETERISATION KIT

Instructions for Use

Critical Care Systems

EN

I. INSTALLATION OF TRANSDUCER CABLES

Connect the transducer cables to the monitor.

II. MANIFOLD KIT SET UP

A. Set up the manifold kit using aseptic technique

1. Open package and remove kit. Check all fittings to ensure tight connections.
2. Attach pressure line to desired transducer system (transducer/dome, disposable transducer, stopcock etc).
CAUTION: All electrical connectors must be kept dry or erratic readings may result.
3. Prepare a collapsible flush solution bag. Heparinise as desired and remove air from the bag by using a sterile needle and syringe or by inverting bag during initial filling of administration set.
CAUTION: If an air-free solution source is not used (i.e., air is not extracted from the bag) air may be forced into the line when solution is exhausted.
4. Close the CAIR® clamp on the administration set. Remove the protective cover from the administration set spike and insert the spike into the solution bag.
CAUTION: To prevent inadvertent puncture of the solution bag, invert and insert the spike carefully, using a downward twisting motion.
5. With the bag inverted, open CAIR® clamp halfway, open 'flush' stopcock and fill drip chamber to halfway mark. Close CAIR® clamp.
6. Insert the solution bag into the pressure cuff then hang the cuff from the IV pole. Close 'flush' stopcock to the administration line and open CAIR® clamp.
7. Remove the protective cover from the contrast spike and insert spike into contrast bottle. Hang bottle from IV pole.

B. Purging air from the system.

- CAUTION: The manifold system must be thoroughly flushed prior to patient usage.**
1. Remove any caps on manifold or transducer system.
 2. Turn 'contrast' or 'proximal' stopcock "off" to the 'distal' end of the manifold (toward the rotator). Flush the contrast line, contrast port and injection port of manifold with contrast fluid using gravity flow. Clear all air and bubbles from the system (tap as needed to dislodge bubbles).
 3. Cap 'proximal' (injection) port or attach syringe. Turn 'contrast' stopcock "off" to proximal port and flush with contrast fluid toward 'distal' (rotator) end of manifold until fluid column passes the 'flush' (center) stopcock. Close contrast stopcock "off" to contrast administration set.
 4. Close 'flush' stopcock "off" to proximal end of manifold and flush administration line and manifold through the 'distal' (rotator) port, clearing all air and bubbles from the system.
 5. Close 'pressure' stopcock to the distal (rotator) port and flush pressure line, transducer dome, and stopcocks clearing all air and bubbles from the system. Close 'flush' stopcock to administration line and close 'pressure' stopcock to pressure line. Close and cap transducer stopcock(s).
NOTE: Take special care to ensure no air is trapped in the fluid pathway. The monitoring system must be totally airfree for maximum performance, i.e., optimal dynamic response.
 6. Pressurise cuff to 300 mmHg. Pressurise contrast bottle if desired.

III. BALANCING AND CALIBRATING

A. After the system has been flushed and mounted, zero balance the transducer.

LIST NO.:	011-OP560-01
COMMODITY NO.:	DF-4565
DATE PREPARED:	8-24-2012
ARTIST:	KENT

MATERIAL
 Paper: 60 pound White if printed in-house or
 Paper: 80 grams Bond (50 pound White) if printed outside vendor
 Ink: Black on uncoated surface
 SPECIFICATIONS:
 Panel Length: 10 1/2 inches
 Sheet Width: 7 1/2 inches
 Folded Size: 3 3/4 inches x 10 1/2 inches
 Center Folded Fold Style: 2 Panel
 MADE FROM: EN-1455, DF-4191

VENDOR TO VERIFY PAPER STOCK:

VENDOR:
 NOSCO, INC.
 651 SOUTH UTICA STREET
 WAUKEGAN, IL. 60085

ALTERNATE:
 ICU MEDICAL, INC.

EN

1. Turn 'pressure' stopcock of the manifold "off" to the transducer.
2. Remove cap of stopcock on the transducer and open to air.
NOTE: The air-fluid interface of the zero reference or 'balancing' stopcock should be at or near the right atrial (mid-axillary) level.
3. Zero balance transducer according to monitor manufacturer's instructions.
4. Close zero reference stopcock and replace cap.

---of---

1. Attach desired catheter to manifold and flush, purging all air bubbles from catheter.
2. Open pressure stopcock on the manifold to the catheter. (The catheter tip is now the system air-fluid interface.)
3. Place manifold and transducer in the position (horizontal plane) they will maintain during pressure measurement.
4. Place the catheter tip at the right atrial (mid-axillary) level.
5. Zero balance transducer according to monitor manufacturer's instructions.

NOTE: As long as the relative level between catheter tip and transducer are maintained, zero balance will be appropriate. Any significant change in relative level necessitates re-zeroing.

- B. Calibrate each transducer to the monitor. Follow monitor manufacturer's printed calibration procedures.

NOTE: An accurate calibration can only be achieved by using a known pressure source, e.g., mercury manometer.

NOTE: The above procedure should be used for periodic balancing and calibration checks.

IV. CHECKING FOR LEAKS

The system should be checked periodically for proper fluid source pressure and leaks.

V. CONNECTION TO PATIENT

- A. Before connecting manifold rotator to catheter, make certain system is free of air and that all air is back-flushed from catheter.
- B. Connect rotator to catheter and proceed as needed.

VI. INDICATIONS AND CONTRAINDICATIONS

A. Indications

- Cardiac catheterisation.
- Radiology special procedure.

B. Contraindications

- Monitoring where continuous flush is indicated.

THE ICU MEDICAL CARDIAC CATHETERISATION KIT IS A SINGLE USE DEVICE. DO NOT REUSE. REUSE NEGATIVELY IMPACTS PERFORMANCE / STERILITY POTENTIALLY RESULTING IN PRODUCT FAILURE / CONTAMINATION.

LIST NO.:	011-OP560-01
COMMODITY NO.:	DF-4565
DATE PREPARED:	8-24-2012
ARTIST:	KENT

MATERIAL

Paper: 60 pound White if printed in-house or
 Paper: 80 grams Bond (50 pound White) if printed outside vendor
 Ink: Black on uncoated surface
SPECIFICATIONS:
 Panel Length: 10 1/2 inches
 Sheet Width: 7 1/2 inches
 Folded Size: 3 3/4 inches x 10 1/2 inches
 Center Folded Fold Style: 2 Panel
 MADE FROM: EN-1455, DF-4191

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ICU Medical, Inc.