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Copyright © 2016, Medtronic Inc.
All rights reserved. Printed in the U.S.A.
M967754A001 Rev.1A

Manufacturer: Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
Tel.: 763/514-4000
800/328-2518
Fax: 763/514-4879

EC Authorized Representative:
Medtronic Ireland
Parkmore Business Park West
Galway, Ireland
Tel.: +353-91-708000
Fax: +353-91-757524

Australian Sponsor:
Medtronic Australasia Pty Ltd
5 Alma Road
Macquarie Park, NSW 2113 Australia
1800 668 670

M967754A001

Rev.1A

ENGLISH

EVEREST 20cc DISPOSABLE INFLATION DEVICE & SURVIVAL KIT: INFORMATION FOR USE
READ INSTRUCTIONS PRIOR TO USE

Refer to the end of this document for symbol definitions.

Contents of Package:

Contains one or more of the following: Everest20 or Everest30 20cc Disposable Inflation Device, Stopcock, Y/Tri-Adapter with Hemostasis Valve, Guide Wire Insertion Tool, and Steering Handle.

Warnings:

- For single patient use only. Do not reuse, reprocess or resterilize this product. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device which could result in patient injury, illness or death. Cleaning, disinfection and resterilization may compromise the essential material and design characteristics of the device leading to device failure.
- These instruments are designed for use by physicians engaged in the practice of a specialized branch of medicine. Use of these devices should be restricted to those specialists trained to perform this procedure.
- Always follow the manufacturer's directions accompanying the balloon catheter for instructions for use, maximum inflation pressure, precautions, and warnings for use.
- Improper tightening of the connection between the Inflation Device and Y/Tri-Adapter and accompanying equipment may result in introduction of air into the vascular system. Do not inject any fluid if air bubbles are visible within the fluid path of the Inflation Device or the Y/Tri-Adapter.
- The pressure displayed on the inflation device must not exceed maximum recommended balloon catheter inflation pressure.
- Prolonged exposure to contrast medium may cause components of the device to malfunction and thereby create a safety hazard.

Precautions:

- The components of the package have been delivered STERILE. Thoroughly inspect the components prior to use to verify that no damage has occurred during handling. If the package and/or its contents have been damaged, or contamination or damage is suspected for any reason, contact Medtronic (see address page).

Description:

The Medtronic Everest Disposable Inflation Device is a sterile 20cc inflation device with a locking mechanism that is operated via the Trigger. Normally, the locking mechanism is engaged. Once the Trigger is pulled back, the locking mechanism is released and the piston can be manually manipulated.

The Everest20 Device is outfitted with a manometer with measuring pressures ranging from vacuum to 20 bars in 0.5 bar increments. The Everest30 Device is outfitted with a manometer with measuring pressures ranging from vacuum to 30 bars in 1 bar increments. The manometers are accurate to within +/- 3% of the gauges full scale. A high pressure connecting tube with a male rotating adapter and a disposable 3-way stopcock are also included to aid in preparation of the Device. When purchased as a "Survival Kit", the package will include a Y/Tri-Adapter with hemostasis valve, a Guide Wire Insertion Tool and a Steering Handle.

Indications:

The Everest 20cc Inflation Device/Survival Kit is to be used to facilitate the use of catheters and guide wires during interventional procedures. The Everest 20cc Inflation Device is designed to be used to inflate/deflate balloon catheters as well as to monitor pressure within the balloon. The Y/Tri-Adapter with Hemostasis Valve is designed to be used on a guiding catheter or dilatation catheter to control backbleeding and to provide a port for introduction of fluids into the interventional system. The Guide Wire Insertion Tool is designed to facilitate placement of a guide wire tip through the Y/Tri-Adapter and into the wire lumen of an interventional catheter. The Guide Wire Steering Handle is designed to hold a small diameter guide wire and provide a handle for manipulating the wire.

Preparation:

1. Prepare solution of contrast medium and normal saline. Check balloon catheter instructions for recommendations, if any, for specific mixture requirements.
2. Pull the Trigger to release the lock mechanism and manually advance the piston via the knob forward to the 0 (zero) cc position.
3. Submerge the male rotating adapter into the prepared contrast solution. With the Trigger pulled back, aspirate an appropriate volume of contrast solution within the syringe. Release the Trigger.
4. Install the Stopcock onto the male rotating adapter. Adjust the valve so that an open fluid path to the syringe is established.
5. While holding the device upright, purge the air from the syringe and connecting tube by advancing the piston. If necessary, tap the syringe lightly to remove all entrapped air within the system.
6. Once certain that all entrapped air has been removed, adjust the valve of the Stopcock so that the fluid path to the syringe is closed.

Connecting the Inflation Device to the Balloon Dilatation Catheter:

Note: Prior to connecting the Inflation Device, the balloon dilatation catheter should be prepared and tested according to the manufacturers instructions.

1. Remove the Stopcock from the male rotating adapter of the Inflation Device.

Caution: Failure to remove the stopcock could compromise the integrity of the fluid path allowing air to enter the system and/or contrast solution to leak out. Inflation pressures may not be maintained.

2. Create a fluid to fluid connection between the male rotating adapter of the Inflation Device and the Balloon inflation port and tighten the hubs securely to one another.

Operating Inflation Device:

1. Release the piston by squeezing the trigger below the syringe. In this position the piston can be pushed in freely for inflation or pulled back easily for deflation.
2. To inflate the balloon, first squeeze and hold the trigger while slowly pushing the piston knob forward. Next, release the trigger and gradually turn the piston knob clockwise to achieve desired pressure. To decrease pressure, turn the piston knob counter-clockwise.
3. To rapidly deflate balloon, squeeze and hold the trigger to release the piston and then pull back on the piston knob. Release the trigger to hold the piston in negative pressure.

Instructions for Use of Accessories:

If the Y/Tri-Adapter is to be used with a guiding catheter:

1. Thread the dilatation catheter through the hole in the cap of the Y/Tri-Adapter. Flush thoroughly to remove air.
2. Advance the dilatation catheter into the guiding catheter.
3. To join the rotating connector of the Y/Tri-Adapter (with the dilatation catheter passing through it) to the guiding catheter, vigorously flush forward through the adapter to remove air while allowing backbleeding through the guiding catheter. This will permit a fluid-fluid interface while the adapter is attached to the guiding catheter.

If the Y/Tri-Adapter is to be used with a dilatation catheter:

1. Attach the rotating connector to the hub of the dilatation catheter.
2. Flush thoroughly to remove any air that may be trapped in the device.
3. Insert the distal tip of the appropriate guide wire into the hub of the guide wire insertion tool. Insert the guide wire insertion tool through the Y/Tri-Adapter and advance the guide wire into the desired device.
4. Remove the guide wire insertion tool, and tighten down the Y/Tri-Adapter hemostasis valve until no back bleeding is evident.

Guide Wire Steering Handle:

1. Advance the appropriate guide wire into the lumen of the dilatation catheter.
2. Slide the steering handle over the proximal end of the guide wire and tighten the handle on the wire.

Storage:

Store accessories in a clean, cool and dark area.

Warranty Disclaimer – United States

Note: Although the Disposable Inflation Device and Survival Kit, hereafter referred to as “product”, has been manufactured under carefully controlled conditions, Medtronic, Inc., and affiliates (collectively, “Medtronic”) have no control over conditions under which this product is used. Medtronic therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty or merchantability or fitness for a particular purpose, shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Warranty Disclaimer – Outside United States

Although the Disposable Inflation Device and Survival Kit, hereafter referred to as “Product” has been carefully designed, manufactured, and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of this Disclaimer of Warranty. Medtronic therefore disclaims

all warranties, both expressed and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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