
Electrosurgical Handpiece Electrosurgical Pencil Single Use

Instructions for Use

MD-IFU-PTD.17A

(EN)

Description/Intended Use:
These devices are designed to be used as accessories in conjunction with the electrosurgical units and electrodes with which they are known to be compatible. Their use enables the operator to remotely conduct an electrosurgical current from the output connector of an electrosurgical unit to the operative site for the desired surgical effect. Maximum voltage is not to exceed 6 KV peak. These devices are compatible with Medstar disposable blade/needle/ball/leep active electrodes.

Contraindications:
These devices should never be used when:
 • There is visible evidence of damage to the exterior of the device such as cracked or damaged plastic or connector damage.
 • These devices fail the inspection described herein.
 • In the presence of flammable gases, flammable prep solutions or drapes, oxidizing gases such as Nitrous Oxide(N2O), or in oxygen-enriched environments.
 • Do not use monopolar electrosurgery on small appendages, as in circum-cision or finger surgery.
 Electrosurgery is potentially hazardous for patients with pacemakers, AICDs, neuro-stimulators or other active implants. High frequency energy may damage or adversely affect active implants. Ventricular fibrillation may occur.

Warning:
 • Keep active accessory away from the patient when not in use. An Accessory Holster may be used to hold active accessories safely.
 • Do not use in the presence of combustible or explosive materials!
 • The ability to effectively clean and re-sterilize this single patient use device has not been established and subsequent re-use may adversely affect the performance, safety and/or sterility of the device.

Caution:
 • The use-by-date is indicated on the packaging. Do not use this product if the use-by-date has been exceeded!
 • Make sure to remove tip protector before using.
 • Use lowest possible power setting on the associated electrosurgical unit capable of achieving desired surgical effect. Activation time should be as short as possible.
 • Never allow the associated accessories connected to these devices to be in contact with skin of the patient or operator except during intended activations.
 • Do not permit the cables connected to the associated electrosurgical accessories to be parallel and in close proximity to the leads of other electrical devices.
 • Always place unused associated electrosurgical accessories in a safe insulated location such as a holster when not in use. Inspect and test each device before each use.
 • After use, this product may be a potential biohazard. Handle and dispose of as required by hospital policy and applicable law.

Inspection:
These devices should be inspected before each use. Visually examine the devices for obvious physical damage including:
 • Cracked, broken or otherwise distorted plastic parts.
 • Broken or significantly bent connector contacts.
 • Damage including cuts, punctures , nicks, abrasions, unusual lumps, significant discoloration.
 • Verify that the electrode is fully and securely seated in the handpiece before use.

Instructions for Use:
 • Use aseptic techniques when removing the product from the packaging!
 • Pull the protective sleeve off the electrode. Ensure that the electrode is firmly seated in the handle. If you change the electrode, make sure that it is inserted fully into the handle.
 • Pull the securing collar off the cable.
 • When using Fingerswitch Handpiece :
 Plug the 3-prong plug into the appropriate receptacle on the electrosurgical unit. Press the yellow button on the handpiece to activate cutting mode. Press the blue button on the handpiece to activate coagulation mode. Make sure that the electrosurgical unit responds correctly.
 • When using Footswitch Handpiece :
 Plug the 1-pin-plug into the appropriate receptacle on the electrosurgical unit. Step on the yellow pedal connects with electrosurgical unit to activate cutting mode. Step on the blue pedal connects with electrosurgical unit to activate coagulation mode. Make sure that the electrosurgical unit responds correctly.

Symbols used on labeling:

REFERENCE NUMBER	BATCH CODE	USE BY	DATE OF MANUFACTURE
DO NOT REUSE	LATEX FREE	BF	KEEP AWAY FROM SUNSHINE
KEEP DRY	CONSULT INSTRUCTIONS FOR USE	CAUTION	DO NOT USE IF PACKAGE IS DAMAGED
PRODUCT COMPLIES WITH REQUIREMENTS OF DIRECTIVE 93/42/EEC MEDICAL DEVICES	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY	MANUFACTURER	RELATIVE HUMIDITY LIMITS
		TEMPERATURE LIMITS	

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MEDSTAR TECHNOLOGY
 Zhejiang Medstar Technology Co., Ltd.
 1/F, Bldg.1, No.239 Yatai Rd.,
 Nanhu District, Jiaxing, 314050, China

OBELIS S.A.
 BOULEVARD GENERAL WAHIS 53
 1030 BRUSSELS, BELGIUM

EC REP

V2.01