



UK Smith & Nephew Healthcare Limited, Healthcare House, 101 Hessele Road, Hull HU3 2BN. Tel: 01482 222200.

AUSTRALIA Smith & Nephew Pty Ltd., 315 Ferntree Gully Road, Mount Waverley, PO Box 242, Victoria 3149. Tel: (03) 8540 6777.

AUSTRIA Smith & Nephew GmbH, Concorde Business Park C3, A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

BELGIUM Smith & Nephew NV-S.A., Hector Henneslaan 366 - 1930 Zaventem. Tel: (02) 702 29 11.

CANADA Smith & Nephew Inc., 2280 Arpenfia Road, Mississauga, ON L5N 6H8. Tel: 1 800 463 7439.

CHINA Smith & Nephew Ltd, 12F Headquarters Building, 168 Xizang Road (M), Huangpu District, Shanghai 200001. Tel: +86 21 2330 3000.

DENMARK Smith & Nephew A/S, Slotsmarken 14, 2970 Hørsholm. Tel: +45 45 80 61 00.

EASTERN EUROPE Smith & Nephew GmbH, Concorde Business Park C3, A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

FRANCE Smith & Nephew SAS, 25 bd Alexandre Oyon 72019 LE MANS Cedex 2. Tel: 02 43 83 23 23.

FINLAND Smith & Nephew Oy, Ayrille 12 C, 01510 Vantaa. +358 207 866 300.

GERMANY Smith & Nephew GmbH, Friesenweg 4/ Haus 21, 22763 Hamburg. Tel: 040 879744-0.

GREECE Costas A. Papaellinas (Hellas) SA, 26th km, Paeanias - Markopoulou Ave, PO Box 182, 19400 Koropi. Tel: 662 6201-7.

HONG KONG Smith & Nephew Ltd., Unit 813-818, 8/F, Delta House, 3 On Yu Street, Shatin, N.T. Tel: (852) 2648 7700.

INDIA Smith and Nephew Healthcare Pvt. Ltd. B-501-509, Dynasty Business Park, Andheri-Kurla Road, Andheri (E), Mumbai-400059. Tel: 022 4005 5090.

IRELAND Smith & Nephew Limited, Carriga Court, George's Avenue, Blackrock, Co. Dublin. Tel: 01-27 0444.

ITALY Smith & Nephew S.r.l. - Via De Capitani 2A, 20864 Agrate Brianza (MB). Tel: +39 039 60941.

JAPAN Smith & Nephew Wound Management KK, 2-4-1 Shiba Park, Minato-ku, Tokyo 105-0011 Japan. Tel: 03-5403-8930.

KOREA Smith & Nephew Ltd, 13F, ASEM Tower 517, Yeongdong-daero, Gangnam-gu, Seoul, Korea (135-798). Tel: 2 6007 7575.

MALAYSIA Smith & Nephew Healthcare Sdn. Bhd. Suite 11.01B & 11.02, Level 11, Menara AmFIRST, No 1 Jalan 19/3, 46300 Petaling Jaya, Selangor. Tel: +603-7958 7103.

MALTA Cherubino Ltd, Delf Building, Sliema Road, Gzira. Tel: 21 343270.

MEXICO Smith & Nephew, SA de CV, Patricio Sanz No. 442, Colonia del Valle, C.P. 03100, México D.F., Tel. 53 40 22 60.

MIDDLE EAST Smith & Nephew FZE, 4th Floor National Bank of Abu Dhabi Building, Al Jazari Building 47, Dubai Healthcare City - District 8, PO Box 9715, Dubai, U.A.E. Tel: 971 4 4299111.

NETHERLANDS Smith & Nephew B.V., PO Box 525, 2130 AM Hoofddorp, Kruisweg 637, 2132 NB Hoofddorp. Tel: (+31) (0) 20 6543999.

NEW ZEALAND Smith & Nephew Limited, PO Box 442 Auckland 1140, 621 Rosebank Road Avondale. Tel: 9 828 4059.

NORWAY Smith & Nephew A/S, PO Box 224, Nye Vakas vei 64, 1379 Nesbru. Tel: 66842020.

PAKISTAN BSN Medical (PVT) Ltd, A/69 SITE, Manghopir Road, PO Box 3659, Karachi 75700. Tel: +92 21 32561176-77.

POLAND Smith & Nephew Sp z o.o., Ul. Marynarska 15, 02-674 Warszawa. Tel: +48 22 360 41 20.

PORTUGAL Smith & Nephew Lda., Estrada Nacional nº 10, Km 131, Parque Tejo - Bloco C, 2625-445 Forte da Casa. Tel: +351 214 460 650.

SINGAPORE Smith & Nephew Pte Ltd., 1A International Business Park, #09-03 Tolaram 609933. Tel: (+65) 6270 0552.

SOUTH AFRICA Smith & Nephew (Pty) Limited, 30 The Boulevard, West End Office Park, Westville, 3629. Tel: 031 2428111.

SPAIN Smith & Nephew SAU, Fructuós Gelabert, 2-4, 08970 Sant Joan Desp, Barcelona. Tel: 93 37373 01.

SWEDEN Smith & Nephew AB, Box 143, Kråkelörspgatan 20, 43122 Mölndal. Tel: 031-746 58 00.

SWITZERLAND Smith & Nephew AG, Glutz Blozheim-Str.1, CH-4502 Solothurn. Tel: 32 624-5660.

TAIWAN (R.O.C.) Smith & Nephew (Overseas) Ltd. Taiwan Branch, 9F-2, No. 50, Sec.1, Xinsheng South Road Taipei City. Tel: 886 2 3322 5020.

THAILAND Smith & Nephew Ltd. 16th Floor Building A, The 9th Towers Grand Rama 9, 33/4 Rama 9 Road, Huaykwang Sub-District, Huaykwang District, Bangkok 10310. Tel: 66(0)2-116-5581.

TURKEY Smith ve Nephew Medikal Cizlar Tic. LTD. Şti. Bahçeköy Mah. Orkide Sokak No:8/E 84473 Sarıyer/ İSTANBUL. Tel: +90 212 226 22 55.

UNITED STATES Smith & Nephew, Inc., Advanced Wound Management, 3909 Hulen Street, Fort Worth, TX 76107. Customer Care Center: 1 800 876-1261.

VENEZUELA Ace Sales Corp, 8085 NW 68th St, Miami, Florida 33166. Tel: 305.835.0310



EN Do not use if pack is open or damaged

FR Ne pas utiliser si l'emballage est ouvert ou endommagé

DE Nicht verwenden wenn Verpackung bereits geöffnet oder beschädigt ist

ES No utilizar si el envase está abierto o deteriorado

IT Non utilizzare se la confezione è aperta o danneggiata

NL Niet gebruiken indien verpakking geopend of beschadigd is

GR Να μην χρησιμοποιείται εάν η συσκευασία είναι ανοικτή ή κατεστραμμένη

DK Má ikke anvendes hvis pakningen er åben eller beskadiget

NO Má ikke brukes hvis emballasjen er åpen eller skadet

FI Ei saa käyttöä, jos pakkaus on avattu tai vahingoittunut

SE Använd inte om förpackningen är öppen eller skadad

PT Não utilizar se a embalagem estiver aberta ou danificada

TR Paket açılmış ya da hasar görmüş ise kullanmayınız

EN

smith&nephew

IV3000^o STANDARD

Moisture responsive catheter dressing

Description

IV3000 Standard is a dressing based on Smith & Nephew's REACTICTM moisture responsive transparent film, specifically designed to meet the needs of catheter fixation.

Mode of Action

The REACTIC film in IV3000 Standard has a unique molecular structure which is significantly more permeable to water vapour than ordinary film. The film is not porous and it therefore helps prevent bacterial contamination to the IV site and is waterproof. These properties result in drier catheter sites, reduced skin maceration and lower skin colonisation, thus reducing the risk of catheter related infection.

Further, it can be applied easily and aseptically, and due to the pattern spread adhesive, pain on removal is reduced. IV3000 Standard is spread with a non-irritating non-sensitising adhesive. Due to these properties IV3000 Standard offers overall advantages over ordinary films, tapes or gauze, and is believed to be a 'technological advance for IV site care'.

Instructions for use

Clean and thoroughly dry application area.

1. Remove approximately half of the protector paper and centre the dressing over the insertion site.
2. Smooth the dressing down onto the patient's skin working from the catheter outwards and remove the remaining protector paper. Be careful not to apply the dressing under tension.
3. The protector paper may be removed from the orange tape handles and the entire handle adhered to the patient's skin.
4. Alternatively, a portion of the tape handles can be torn off with the protector paper still in place. This

can then be removed and the tape handles used to help secure the catheter hub or tubing. They may also be used for documentation.

5. If using more than one dressing, do not allow the dressings to overlap.
6. Remove dressing by loosening the orange handles and stretching the dressing parallel to the patient's skin while stabilising the patient's skin and catheter with the other hand. Alternatively, use an alcohol wipe to break the adhesive seal. If the skin is particularly fragile then warm soap and water may be used to break the adhesive bond and to ease the dressing off.

Indications

The above mentioned products are suitable for dressing peripheral, central line and epidural catheters directly onto the skin.

Contraindications

IV3000 Standard should be used on IV sites only. If a film dressing is required for open wounds such as pressure sores, burns or donor sites, OPSITE^o FLEXIGRID^o should be used.

Precautions

Use on prepared dry skin only. Do not stack dressings, or allow dressings to overlap. Periodically monitor the dressing and catheter site to confirm secure attachment and continued proper infusion, especially after bathing, showering, or if the dressing and catheter site becomes wet. If the dressing comes off, evaluate to ensure proper catheter placement, and then apply a new dressing. As with all adhesive products, apply and remove carefully from sensitive or fragile skin.

Product Availability

Code	Size	Indication	Carton
4924	6cm x 8.5cm/2 3/8" x 3 3/8"	Non ported peripheral catheters	100 (not available in the USA)
4925	10cm x 14cm/4" x 5 1/2"	Central line catheters	10
4973	10cm x 14cm/4" x 5 1/2"	Central line catheters	50
4649	10cm x 20cm/4" x 8"	Epidural catheters	50

Smith & Nephew Medical Limited, 101 Hessele Rd, Hull HU3 2BN England. *Trade Marks of Smith & Nephew ©Smith & Nephew Date of issue 10/2015 www.smith-nephew.com



0086

STERILE EO



10601264 65174