



Push-Click Y Connector Kit Instructions For Use

Please read instructions carefully before use.

Ref No.	Push-Click Y Connector	Guidewire Torque Device	Guidewire Insertion Tools
DPY01	Push-Click Y Connector	N/A	N/A
DPY01B		Guidewire Torque	Guidewire Introducer
DPY02	Push-Click Y Connector with Fixed Extension Tube and 3-Way Stopcock	N/A	N/A
DPY02 B		Guidewire Torque	Guidewire Introducer
DPY03	Push-Click Y Connector with Removable Extension Tube and 3-Way Stopcock	N/A	N/A
DPY03B		Guidewire Torque	Guidewire Introducer
DPY04B	N/A	Guidewire Torque	N/A
DPY05	N/A	N/A	Guidewire Introducer



(a) Push-Click Y Connector with Fixed Extension Tube and 3-Way Stopcock



(b) Guidewire Introducer



(c) Guidewire Torque


Graph 1

Indications:

Y connector intends to provide two tunnels into the blood vessel, this can introduce the contrast agent and guidewire or catheters during the interventional surgery. Y connector valve intends to prevent blood loss. Guidewire Introducer inserts the valve of Y connector, making a tunnel for guidewire to get through and protecting guidewire tip. Guidewire Torquer holds guidewire and makes torque guidewire move easier.

Instruction for use (Push-Click Y Connector Kit):

1. Check the package and products prior to use to verify that no damage or contamination.
2. Flush the devices with normal saline before use.
3. Connect the removable 3-way stopcock with bypass if necessary.




Graph 2


1

2

4. Securely connect the distal luer lock with catheter or other devices.
5. Press the knob to the end to keep the valve in "Open" position. ("click" sound should be heard) (see Graph 2)
6. When the valve is in fully "Open" position, operator can put through the operating device such as guidewire, balloon catheter. Do not try to put any devices through the Y connector, when valve is in "Close" position.
7. Press the release key (small blue key) will reset the valve to "Close" position. (see Graph 3)



Graph 3



Graph 4

8. Press knob by fingers between "Open" and "Close" can adjust the valve from full open to full close which may help the device through with minimum blood loss. (see Graph 4)
9. If try to decrease bleeding, operator can try not push knob to fully open position and make the Y connector's valve half open. But in this way, operator must be strongly careful to make sure that the device is not damaged.
10. After put the device in the catheter through the Y connector valve, operator can close the Y connector's valve, but the device still can be pushed and rotated.
11. Guidewire Introducer can help the insertion of guidewire.
12. Using Guidewire Torquer is recommended to help controlling guidewire.

Warnings:

1. Read the instructions before use.
2. The device is designed for single use only. Do not de-sterilize and /or reuse.
3. Do not use it after use before date (expiry date).
4. Do not use it if the package is damaged or not well sealed.
5. Luer connector must be well connected, and the valve must be closed. Incorrect connection, lose connection or the valve not being closed properly may cause leaking or air getting into blood vessel.
6. Do not inject air into the blood vessel.
7. Suggested working pressure is under 58psi (400kPa).
8. For qualified physicians use only.
9. Please dispose the device after use, according to the local legislation.
10. The products with extension tubes contain DEHP.
11. The intactness of the content should be inspected after opening.
12. Strict aseptic operation should be carried out during use to avoid bacterial contamination.
13. It is important to make that the valve of the push-pull Y-connector is open before using it to connect any device.
14. It is recommended to use the full set of PTCA accessories manufactured by our company, as mixing with products of other brands may cause poor mismatch. The company bears no responsibility for any problems arising therefrom.



Caution: Please dispose the device after use according to local laws and regulations.

3

4

Signs and Legends:



5

Manufacture Information:

Beijing Demax Medical Technology Co., Ltd
Registration Address: A13-7, Jingshengnansi Street, Tongzhou District, 101102 Beijing, PEOPLE'S REPUBLIC OF CHINA
Manufacturing Address: A13-7, Jingshengnansi Street, Tongzhou District, 101102 Beijing, PEOPLE'S REPUBLIC OF CHINA
Tel.: +86-10-59771799
Fax: +86-10-59771883
Email: marketing@demaxmedical.com

European Representative:

Shanghai International Holding Corp. GmbH(Europe)
Address: Eiffestrasse 80, 20537 Hamburg, Germany
Tel.: +49-40-2513175
Fax: +49-40-255726

Storage

Store in a dry, dark, cool place.

The use before date (expiry date) see label.

Transportation condition:

keep away from heavy and sunshine, keep dry

Temperature limit:-18 C~55 C

Humidity limitation: 0~85%

Customer Service:

Australia and New Zealand:

Demax Medical PTY LTD

Address: Suite 241,117 Old Pittwater Road, Brookvale, NSW 2100 (Sydney), Australia

Tel.: +61 2 9939 7093

Fax: +61 2 9939 0397

Email: csau@demaxmedical.com

China Mainland:

Beijing Demax Medical Technology Co., Ltd

Address: A13-7, Jingshengnansi Street, Tongzhou District, 101102 Beijing, PEOPLE'S REPUBLIC OF CHINA

Tel.: +86-10-59771799

Fax: +86-10-59771883

Email: cscn@demaxmedical.com

6

7