

ECO Basic Podiatry Kit

Sterile and Single Use



Podiatry instruments can be hard to clean between patients and time consuming to collect, sterilise and set up. Importantly, minimising the risk of cross infection between clients represents a major challenge in healthcare clinics.

Our single-use Podiatry Kits from the UK are sterile, contain quality medical grade instruments and are supplied conveniently in an easy to open peel down pouch.

- + Saves setup time and post procedure clean-up
- + Ensures all components are readily available
- + Protects patients against cross-contamination infection risks
- + Cost efficient alternative to reprocessing reusable instruments



Contact us today on 1300 550 278 for further information or order direct from our online shop at defries.com.au

CODE	DESCRIPTION	BOX QTY.
------	-------------	----------

ECO BASIC PODIATRY KIT

PP16015	Each sterile Podiatry Kit includes: 1 x Nail Nipper Roll Spring Concaved Blade 14.5cm 1 x BP Scalpel Handle No. 3 1 x Emery Board 1 x Blacks File Medium Single End Sealed in a peel down pouch	100
PP16015-5		5

