

PRODUCT CATALOGUE

Effective October 2018





Defries Industries is a 100% Australian family owned medical business. As a world class manufacturer, our vision has remained clear over the years:

“WE ARE COMMITTED TO IMPROVING PATIENT CARE BY OFFERING SINGLE-USE, STERILE MEDICAL PRODUCTS AND CUSTOMISED PACKS TO HEALTHCARE PROFESSIONALS.”

CONTENTS

Introduction	2
Clinical drapes and equipment covers	4
Clinical protective apparel	13
Consumables	20
Custom procedure and standard packs	43
Instruments and instrument sets	53
Maternity products	68
Medical imaging and cath lab products	72
Medical labels	76
Plasticware	84
Wound care	90
Index	97

INTRODUCTION

The Defries Industries team share a vision of high standards and a commitment to superior patient care. We continually strive to create outstanding products that excel in design, ease of use, safety and price. Defries has built a strong national reputation since 1981 and continues to expand and develop our product range. We are proud to be a 100% Australian family owned medical business and aim for ongoing growth and investment in Australian and New Zealand healthcare into the future.



FREECALL 1300 550 278

WE ARE COMMITTED TO IMPROVING PATIENT CARE WITH OUR SAFE AND EFFECTIVE SINGLE-USE, STERILE MEDICAL PRODUCTS AND CUSTOMISED PROCEDURE PACKS

Our range of over 3000 different medical and surgical products are a combination of single sterile items and procedure packs. This catalogue features all our single sterile items and standard general packs. For customised procedure packs, please contact our Customer Service team for your local Representative's details.

AT DEFRIES WE COLLABORATE WITH INDUSTRY LEADING SUPPLIERS

We regularly meet with our primary suppliers to explore new technologies and products available, so that we can offer our customers quality products at competitive prices. We visit our suppliers' facilities to check their Quality Management Systems and ensure the ongoing supply of components to Defries.

We are proud to distribute innovative products from Advancis Medical, DeRoyal Industries, OBP Medical, Degania Medical and Tava Surgical. These medical companies are leaders in healthcare from around the world.

REGULATORY COMPLIANCE IS VITALLY IMPORTANT AT DEFRIES INDUSTRIES

Our quality management system, which is based on ISO13485:2003, ensures we meet stringent medical device regulatory and quality requirements. Our conformity assessment certificate number is AU Q00187/01.

The key role of our Quality Assurance is to deliver products that are superior in safety, quality and efficacy by achieving ongoing compliance with approved specifications and relevant industry standards.



SUSTAINABILITY AND ENVIRONMENTAL DECISIONS

We aim to apply and monitor the highest environmental standards in all areas of our operations at Defries. We are committed to identifying and improving environmental performance and constantly review how the design, manufacture and distribution of Defries products affects the environment. Our corporate sustainable strategy includes environmental values, impacts and commitments which include the following:

- Energy use and Carbon emissions
- Ethical company suppliers
- Waste management
- Water resources
- Transport and fuel emissions
- General and end user waste management





CLINICAL DRAPES AND EQUIPMENT COVERS

FREECALL 1300 550 278

CLINICAL DRAPES AND EQUIPMENT COVERS

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BAGS

CLIP SEAL BAG



DEF640

DEF640	20cm x 25cm	100
DEF3286	20cm x 25cm, Packet of 20	50
DEF1703	28cm x 38cm	25

FREEZER BAG



DEF639

DEF639	30cm x 45cm	100
DEF66	30cm x 45cm with 2 x Rubber Bands #34	100
DEF1336	30cm x 45cm with 2 x Rubber Bands #34, Double Packed	100

DRAWSTRING BAG



DEF846

DEF846	25cm x 50cm	25
DEF864	44.5cm x 46.5cm	25
DEF861	44.5cm x 46.5cm	100
DEF1687	44.5cm x 46.5cm, Packet of 2	50
DEF848	50cm x 50cm	25
DEF847	50cm x 85cm	25

PLASTIC BAG



DEF1878


DEF1878	14cm x 30cm	100
DEF1276	14cm x 30cm with a Blank Adhesive Label	100
DEF638	28cm x 32.5cm	100
DEF650	Burns 30.5cm x 39cm	100
DEF642	45.5cm x 51cm with Cotton Tape 60cm Long	100
DEF808	Heavy Duty 57.5cm x 62cm	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

EQUIPMENT COVERS

BACK TABLE COVER		
DEF5015	112cm x 228cm (44" x 90")	40
DEF5016	152cm x 228cm (60" x 90")	40
DEF5016AQ	152cm x 228cm (60" x 90"), Aqua	40
DEF1920A	152cm x 228cm (60" x 90"), Fan Folded	40
DEF3544	150cm x 190cm, Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3545	150cm x 230cm, Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF5017	152cm x 228cm (60" x 90"), Heavy Duty	40
DEF5018	178cm x 280cm (70" x 110"), Heavy Duty	40

CAMERA & CORD COVER		
		
DEF1164		
DEF1164	Camera Cord Cover 15cm x 230cm, Sealed End	50
DEF1201	Camera Cord Cover 15cm x 230cm with Cuff, Sealed End & Tapes, Double Packed	50
DEF3591	Camera Light Handle Cover 66mm	25

CHAIR COVER		
		
DEF772		
DEF772	2 x 57.5cm x 62cm with Sealed Ends, with 2 x Rubber Bands #33	100

DRILL SLEEVE COVER		
		
DEF887		
DEF1224	5cm x 80cm with 3 Labels	100
DEF1368	5cm x 110cm with Tape	100
DEF300	7.5cm x 75cm with Tape	100
DEF1516	7.5cm x 100cm with Adhesive both Ends	100
DEF1154	7.5cm x 150cm with 2 x Rubber Bands #31	100
DEF318	10cm x 60cm, Packet of 3	100
DEF887	10cm x 100cm	100
DEF889	10cm x 150cm	100
DEF618	10cm x 152cm with Tape, Double Packed	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HANDLE SLEEVE COVER		
		
DEF1053		
DEF1053	14cm x 30cm with Adhesive and Sealed End	100

MAYO STAND COVER		
DEF5019AQ	61cm x 135cm with Sealed End	30
DEF5019	58.5cm x 142cm, Reinforced with Absorbent pad 50cm x 65cm, Cuff and Sealed End	50
DEF761	62cm x 91cm, Clear Plastic with Sealed End	100
DEF775	70cm x 137cm, Blue Plastic with Cuff and Sealed End	50
DEF3488	81cm x 165cm, Reinforced, Heavy Duty with Sealed End	25

PLASTIC SLEEVE COVER		
		
DEF454		
DEF454	15cm x 100cm, with 2 x Rubber Bands #31, Double Packed	100

PROTECTIVE COVER		
DEF613	15cm x 100cm, Plastic with Sealed End	50
DEF1744	15cm x 2m Sealed End, with 1 x Rubber Band #31, Double Packed	100
DEF763	104cm x 150cm, Protective	40
DEF1573	150cm x 208cm, Protective	50
DEF2582	27cm x 92cm, IV Pole Cover with Sealed End, Double Packed	100

TROLLEY COVER		
		
DEF641		
DEF641	Trolley Bag Cover 61cm x 102cm with Sealed End	50
DEF1239	Shower Trolley Drape 2m x 3m	25



CUSTOMISED PROCEDURE DRAPES

THE FLEXIBILITY OF CHOICE

At Defries we have a comprehensive range of patient drapes which cover all procedures in the Operating Room and we can also design drapes in conjunction with Surgeons, Nurses and Theatre Managers. We provide our customers with superior products that match their clinical needs.

Our drapes are made to specific requirements such as size, fabric, fenestrations and extra details like velcro taping and clear plastic anaesthetic screens.

CHOICE OF FABRIC IS IMPORTANT

At Defries we use a variety of drape fabrics:

Azure drapes use a viscose impervious fabric with a polyethylene laminate, suitable for general surgery with minimal fluid or short procedures. This range includes the Azure minor procedure drapes.

Sapphire drapes use a bi-laminate or tri-laminate fabric with extra reinforcement. This range includes back table covers, mayo covers, some procedure drapes and bespoke designed drapes to meet specific customer requirements.

Aqua drapes are our newest drape range available and use an innovative drape fabric that is a tri-laminate – a three layer drape which offers an absorbent top level, an impervious middle coating and a soft and comfortable patient layer.

This high performance fabric is soft and drapeable with additional features such as reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. We have an excellent standard range of designs covering all procedures from general surgery and universal packs, to orthopaedics and specialised surgery such as neuro.

Supreme drapes are made from a tri-laminate fabric and can be tailored and customised to your surgeons' requirements. They are specifically designed for high level performance, coupled with great absorbency from fabric reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. The Supreme range includes drape sheets and specialty procedure drapes.

Our drape sheets, universal packs and specialty drapes include the instructions for correct aseptic application with markings and directions on the drape to assist ease of use.

All drape fabrics are conformable to patients' shapes. They are low linting and have excellent barrier properties to prevent strikethrough of airborne and blood borne bacteria.

Please contact your Defries representative for further information.

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GENERAL DRAPES AND PACKS

AZURE MINOR PROCEDURE DRAPE		
Soft viscose fabric combined with a polyethylene laminate to create a drape ideal for minor procedures.		
DEF5094	38cm x 45cm, Azure Circumcision Drape with a 3cm Circular Fenestration	50
DEF5090	50cm x 60cm, Azure Minor Procedure Drape with a 6.5cm Circular Fenestration with Full Adhesive	25
DEF5091	91cm x 91cm, Azure Minor Procedure Drape with a 10cm x 15cm Oval Fenestration with Full Adhesive	40
DEF5092	200cm x 300cm, Azure Minor Procedure Drape with a 15cm x 25cm Fenestration with Full Adhesive	12

AQUA DRAPE SHEETS		
Innovative tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.		
DEF5024AQ	Half Drape, 100cm x 150cm	20
DEF5025AQ	Medium Drape, 100cm x 180cm	20
DEF5026AQ	3/4 Drape, 150cm x 200cm	20
DEF5028AQ	Large Drape, 183cm x 215cm	20
DEF5029AQ	XLarge Drape, 178cm x 254cm	12
DEF5030AQ	Top/Head Drape, 135cm x 250cm with Reinforcement and Adhesive	15
DEF5031AQ	Side Drape, 100cm x 130cm with Reinforcement and Adhesive	20
DEF5032AQ	Bottom/Foot Drape, 190cm x 185cm with Reinforcement and Adhesive	20

SUPREME DRAPE SHEETS

Supreme tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.		
DEF5024	Half Drape, 112cm x 145cm	20
DEF5025	Medium Drape, 112cm x 180cm	20
DEF5027	3/4 Drape, 142cm x 193cm with Reinforcement	20
DEF5028	Large Drape, 183cm x 215cm	20
DEF5030	Top/Head Drape, 135cm x 250cm with Reinforcement and Adhesive	30
DEF5031	Side Drape, 100cm x 110cm with Reinforcement and Adhesive	20
DEF5032	Bottom/Foot Drape, 185cm x 190cm with Reinforcement and Adhesive	20

CLEAR PLASTIC DRAPE

		
DEF742		
DEF742	45cm x 61cm	200
DEF741	61cm x 91cm	25
DEF079	61cm x 91cm	200
DEF749	70cm x 120cm	100
DEF745	76cm x 106cm	100
DEF151	91cm x 115cm	100
DEF2430	Procedure Adhesive Drape, 91cm x 115cm with 2.5cm x 115cm adhesive down length of drape	50
DEF748	120cm x 140cm	100
DEF752	120cm x 140cm, Double Packed	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
UNIVERSAL DRAPE PACKS					
UNIVERSAL PACK SHORT					
DEF5212	1 x Top Drape, 135cm x 250cm with Reinforcement 40cm x 70cm and Adhesive 2 x Side Drapes, 100cm x 100cm with Reinforcement 40cm x 60cm and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement 40cm x 65cm and Adhesive 1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement 50cm x 65cm Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement 80cm x 230cm	10			
			DEF1481		
DEF1481	76cm x 106cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100			
DEF748FA-S	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	10			
DEF3394	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100			
UNIVERSAL PACK LONG					
DEF5213	1 x Top Drape, 135cm x 250cm with Reinforcement 40cm x 70cm and Adhesive 2 x Side Drapes, 100cm x 195cm with Reinforcement 40cm x 75cm and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement 40cm x 65cm and Adhesive 1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement 50cm x 65cm Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement 75cm x 230cm	5			
POLYBACKED DRAPE					
An absorbent drape with a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor procedures.					
					
			DEF1146		
DEF1146	Polybacked Towel, 33cm x 46cm	25			
DEF833	Polybacked Dressing Towel, 33cm x 46cm Large	750			
DEF1287	Polybacked Towel, 40cm x 60cm	100			
DEF1288	Polybacked Towel, 40cm x 60cm with a 10cm Circular Fenestration	100			
DEF1289	Polybacked Towel, 40cm x 60cm with a 10cm x 15cm Fenestration and Adhesive on each corner of Drape 2.5cm x 2.5cm	100			
DEF2072	Polybacked Procedural Drape, 75cm x 105cm with a 10cm x 15cm Fenestration and a 5cm Adhesive length along the 10cm sides of	25			

EPIDURAL DRAPE

			DEF1481		
DEF1481	76cm x 106cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100			
DEF748FA-S	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	10			
DEF3394	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100			

POLYBACKED DRAPE

An absorbent drape with a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor procedures.

			DEF1146		
DEF1146	Polybacked Towel, 33cm x 46cm	25			
DEF833	Polybacked Dressing Towel, 33cm x 46cm Large	750			
DEF1287	Polybacked Towel, 40cm x 60cm	100			
DEF1288	Polybacked Towel, 40cm x 60cm with a 10cm Circular Fenestration	100			
DEF1289	Polybacked Towel, 40cm x 60cm with a 10cm x 15cm Fenestration and Adhesive on each corner of Drape 2.5cm x 2.5cm	100			
DEF2072	Polybacked Procedural Drape, 75cm x 105cm with a 10cm x 15cm Fenestration and a 5cm Adhesive length along the 10cm sides of Fenestration	25			

CODE

DESCRIPTION

BOX QTY.

CODE

DESCRIPTION

BOX QTY.

FABRIC SMS DRAPE

SMS fabric is a 3 layer polypropylene material bonded together to form a drape that is fluid repellent. Ideal for minor procedures.



DEF1594

DEF1122	30.5cm x 61cm	100
DEF1594	61cm x 61cm	50

STOCKINETTE



DEF5060

DEF5060	Medium 20cm x 122cm	25
DEF5061	Large 23cm x 122cm	25
DEF5062	XLarge 30cm x 122cm	25

LEGGINGS

DEF2481AQ	Aqua Lithotomy Leggings, 80cm x 140cm, Packet of 2	26
-----------	--	----

SUPPLEMENTARY

HAND DRYING AND STERILE FIELDS

BLUE HUCK TOWEL

High absorbency cotton towels.



DEF1300

DEF1976	40cm x 68cm	30
DEF1300	40cm x 68cm	100
DEF2018	40cm x 68cm Blue, Packet of 5	40
DEF2219	40cm x 68cm with X-Ray, Blue	25

SUPPLEMENTARY

HAND DRYING AND STERILE FIELDS

BLUE HUCK TOWEL


High absorbency cotton towels.

			DEF1300		
DEF1976	40cm x 68cm	30			
DEF1300	40cm x 68cm	100			
DEF2018	40cm x 68cm Blue, Packet of 5	40			
DEF2219	40cm x 68cm with X-Ray, Blue	25			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAND DRYING AND STERILE FIELDS (CONTINUED)

GREEN HUCK TOWEL




DEF2619

DEF2619	40cm x 68cm, Packet of 2	100
DEF2501	40cm x 68cm, Packet of 4	50
DEF2140	40cm x 68cm, Packet of 5	40

WHITE HUCK TOWEL

DEF2753	40cm x 68cm	30
---------	-------------	----

DRESSING T TOWEL



DEF865

DEF865	32cm x 32.5cm	100
DEF243	32cm x 32.5cm	1000
DEF862	32cm x 32.5cm, Packet of 2	50
DEF879	32.5cm x 60cm Large	100
DEF857	32.5cm x 60cm Large, Packet of 2	50

DRESSING D TOWEL



DEF656


DEF656	32cm x 34cm, Packet of 2	50
DEF831	32cm x 34cm, Packet of 2	300
DEF832	32cm x 34cm, Packet of 10	100
DEF894	34cm x 60cm Large	50
DEF471	34cm x 60cm	500
DEF912	34cm x 60cm Large, Packet of 2	50
DEF910	34cm x 60cm Large, Packet of 5	25

DRESSING TOWEL

DEF470	40cm x 42cm, White	50
--------	--------------------	----



DEF470



CLINICAL
PROTECTIVE
APPAREL

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GOWNS

PATIENT GOWNS

Patient Gowns are lightweight, soft and comfortable and are available in a range of sizes to accommodate all body shapes. They are ideal for patients in hospitals, day procedures, GP clinics, radiology or any setting where a patient is required to wear a disposable gown.



PATIENT-Gown-SS-MED	Medium Short Sleeves	50
PATIENT-Gown-SS-LRG	Large Short Sleeves	50
PATIENT-Gown-SS-XL	X-Large Short Sleeves	50

SCRUB SUITS

Scrub Suits are made from a comfortable and soft fabric, and offer good protection. Ideal for nurses in the operating room, ER Departments, GP Clinics, Dentists, Vets and Students.



SCRUB-SUIT-SML-TOPS

SCRUB-SUIT-SML-TOPS	Top – Size Small	100
SCRUB-SUIT-MED-TOPS	Top – Size Medium	100
SCRUB-SUIT-LRG-TOPS	Top – Size Large	100
SCRUB-SUIT-XL-TOPS	Top – Size XL	100
SCRUB-SUIT-XXL-TOPS	Top – Size XXL	100



SCRUB-SUIT-SML-PANTS

SCRUB-SUIT-SML-PANTS	Pants – Size Small	100
SCRUB-SUIT-MED-PANTS	Pants – Size Medium	100
SCRUB-SUIT-LRG-PANTS	Pants – Size Large	100
SCRUB-SUIT-XL-PANTS	Pants – Size XL	100
SCRUB-SUIT-XXL-PANTS	Pants – Size XXL	100

SURGICAL GOWNS

Our single-use surgical gowns have specifically been developed for the Australian and New Zealand environment. Our range encompasses gowns for short and dry procedures through to high fluid and lengthy procedures.

OUR COMFORT RANGE – OUR SOFTEST AND MOST COMFORTABLE GOWNS

- Comfort Surgical Gown – Level 2.
- Comfort Unreinforced with Glove Assist – Level 3 gown with a reinforced sleeve to allow extra strength during gloving, but still cool, soft and ideal for short and dry procedures.
- Comfort Plus Fabric Reinforced Surgical Gown – Level 3 gown includes fluid resistant fabric reinforced panels yet remains lightweight and durable.
- Comfort Plus Poly Reinforced Surgical Gown – Level 3 gown which includes poly reinforced impervious panels for added protection.
- Comfort Ultimate Protection – Fully impervious high tech fabric, but still soft, comfortable and cool which meets AAMI Level 4 fabric standards.

OUR CLASSIC RANGE:

- Classic Standard Surgical Gown – suitable for most surgical procedures.
- Classic Plus Reinforced Surgical Gown – Includes a reinforced panel on the front and sleeves, for surgical cases requiring extra protection.
- Coolmax Surgical Gown – breathable and cool, these gowns are notable for their superb comfort and being fully impervious. With AAMI Level 4 fabric.

Designed in close co-operation with hospital staff, our Surgical Gowns are cool and breathable with an easy fit, while offering excellent protection against strikethrough from microbial and liquid penetration.




CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COMFORT SURGICAL GOWNS

LEVEL 2 COMFORT 35 UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 35gsm T-cut Surgical Gown ideal for dry and short procedures.




Level 2 Comfort 35 Unreinforced

DEF5316	Medium 119cm	24
DEF5317	Large 127cm	24
DEF5318	X-Large 137cm	24
DEF5319	XL-XLong 149cm	24
DEF5320	XX-Large 155cm	24

LEVEL 3 COMFORT UNREINFORCED WITH GLOVE ASSIST

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41gsm T-cut Surgical Gown ideal for dry and short procedures. Has innovative Glove Assist fabric reinforcement to sleeve area to ensure extra strength. Comfortable and cool to wear.




Level 3 Comfort Unreinforced With Glove Assist

DEF5301	Medium 119cm	24
DEF5302	Large 127cm	24
DEF5303	X-Large 137cm	24
DEF5304	XL-XLong 149cm	24
DEF5305	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS FABRIC REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41gsm T-cut Surgical Gown with fabric reinforcement on front and sleeves. Fluid resistant on panels. Comfortable to wear.




Comfort Plus Fabric Reinforced

DEF5306	Medium 119cm	24
DEF5307	Large 127cm	24
DEF5308	X-Large 137cm	24
DEF5309	XL-XLong 149cm	24
DEF5310	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41 gsm T-cut Surgical Gown which includes poly reinforced panels. Durable and impervious to fluid on panels. Comfortable to wear.




Comfort Plus Poly Reinforced

DEF5311	Medium 119cm	24
DEF5312	Large 127cm	24
DEF5313	X-Large 137cm	24
DEF5314	XL-XLong 149cm	24
DEF5315	XX-Large 155cm	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 4 COMFORT ULTIMATE PROTECTION

Surgical Gown T-cut wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile Level 4, 67gsm trilaminate material. High tech, fully impervious and breathable fabric which is extremely comfortable and cool for surgeons and or nurses.




Comfort Ultimate Protection made with Level 4 fabric

DEF5321	Medium 119cm	24
DEF5322	Large 127cm	24
DEF5323	X-Large 137cm	24
DEF5324	XL-XLong 149cm	24
DEF5325	XX-Large 155cm	24

CLASSIC SURGICAL GOWNS

LEVEL 3 CLASSIC UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown ideal for short and dry procedures.




Level 3 Classic Unreinforced

DEF5000A	Medium 110cm	25
DEF5001A	Large 130cm	25
DEF5002A	X-Large 140cm	25
DEF5004A	XL-XLong 157cm	25

LEVEL 3 CLASSIC PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown with impervious poly reinforcement. Colour coded neck band – purple for level 3 Gown.



Level 3 Classic Plus Poly Reinforced

DEF5005A	Medium 110cm	25
DEF5006A	Large 130cm	25
DEF5007A	X-Large 140cm	25
DEF5009A	XL-XLong 157cm	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------


WEARABLE CONSUMABLES

CLASSIC SURGICAL GOWNS (CONTINUED)

LEVEL 4 COOLMAX

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile Level 4, trilaminate material.

High quality T-cut Surgical Gown ideal for high fluid procedures or lengthy procedures. Colour Coded neck band – blue for level 4 Gown.



Coolmax Surgical Gown Level 4

DEF5010A	Medium 110cm	25
DEF5011A	Large 130cm	25
DEF5012A	X-Large 140cm	25
DEF5013A	XL-XLong 160cm	25
DEF3198A	Coolmax Surgical Gown Pack – 1 Large and 1 X-Large Gown wrapped in a 200cm x 200cm SMS field with two huck towels	10

SPEyes EYEWEAR NON STERILE




Our popular DeRoyal SPEyes eyewear is designed to minimise exposure to pathogens from waste fluid. They combine innovation, quality, value, and comfort for clinicians where personal protective equipment is a vital part of self-care.

- Wrap-around disposable lens for complete protection
- UltraClear™ technology prevents glare, fog, and static
- Frame fits completely against forehead for splash protection from above
- Sealed lens mounting for added fluid protection
- Clip-On/Off design for quick and easy changing
- Coloured top strip allows finger-print-free handling
- Latex free



23-505

23-505	SPEyes Eye Shields pre-assembled assorted colours	100
--------	---	-----



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPEyes EYEWEAR NON STERILE



23-506

23-506	Dispenser box with pre-assembled SPEyes Eye Shields (30 per box 4 boxes per case)	120 total
23-501	SPEyes Eye Shield Replacement Lens	100
23-502	SPEyes Eye Shield Frame	25

FACESHIELDZ™



23-002

23-002	Face Mask Full with Elastic Strap	50
--------	-----------------------------------	----

GLASSES



DEF1755

DEF1755	Disposable Glasses Clear Polycarbonate Non Sterile	1
---------	--	---

BERET CAPS



DEF1142

DEF1142	Beret Cap 53cm (21") Blue Non Sterile	250
---------	---------------------------------------	-----

OVERSHOES NON STERILE



DEF1145

DEF1145	Non Skid Universal Overshoes Polypropylene	100
---------	--	-----

CONSUMABLES

FREECALL 1300 550 278

CONSUMABLES

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

THEATRE CONSUMABLES

All surgical consumables are pre-packed and sterile, ready to be used in the operating room, unless indicated as Non Sterile.

BANDS

RUBBER BANDS (CONTAINS LATEX)

ID indicates Inner Diameter. BW indicates Band Width.



DEF1405






DEF1265	#28 – (ID) 40mm x 1.5mm (BW), Packet of 1	25
DEF1590	#16 – (ID) 60mm x 1.5mm (BW), Packet of 2, Double Packed	100
DEF624	#64 – (ID) 55mm x 6mm (BW), Packet of 1	100
DEF627	#64 – (ID) 55mm x 6mm (BW), Packet of 3, Double Packed	100
DEF626	#64 – (ID) 55mm x 6mm (BW), Packet of 5, Double Packed	100
DEF578	#31 – (ID) 63mm x 3mm (BW), Packet of 1	100
DEF1405	#31 – (ID) 63mm x 3mm (BW), Packet of 2	100
DEF339	#31 – (ID) 63mm x 3mm (BW), Packet of 3	100
DEF625	#31 – (ID) 63mm x 3mm (BW), Packet of 10, Double Packed	100
DEF599	#34 – (ID) 64mm x 3mm (BW), Packet of 2, Double Packed	25
DEF619	#34 – (ID) 64mm x 3mm (BW), Packet of 5	25
DEF1793	#33 – (ID) 87mm x 3mm (BW) Packet of 2	25
DEF669	#33 – (ID) 87mm x 3mm (BW), Packet of 5, Double Packed	25

LATEX FREE BANDS

DEF3542	#33 – (ID) 87mm x 3mm (BW), Packet of 2	100
DEF3560	#33 – (ID) 87mm x 3mm (BW), Packet of 3	100

CLAMPS/CLIPS/HOLDERS

TUBE CLAMPS

DEF1575	Tube Clamp with Lock Blue Plastic 12cm	50		DEF1575
DEF2075	Towel Clamp U-Shaped Cream	50		DEF2075
DEF2190	Towel Clip Green Plastic 11.5cm	25		DEF2190
DEF2488	Towel Clip Stainless Steel 11cm Sharp/Sharp	10		DEF2488
DEF2298	Tube Clamps 13cm Smooth Jaw	10		DEF2298

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CLAMPS/CLIPS/HOLDERS (CONTINUED)

DRAPE CLIP

General purpose, reusable Drape Clip, used during surgery for clamping purposes to ensure items remain where placed. A wide easy to use grip surface ensures heavy duty holding power. Available in 8 different colours. Non sterile.




DRAPE-CLIP

DRAPE-CLIP	Non Sterile Mixed Colours	1
------------	---------------------------	---

TUBE HOLDER

Ideal for securing and holding tubes during a procedure. Adhesive backing for securement onto a drape or towel. Adjustable velcro to allow gentle and easy placement.



DEF1889

DEF1889	Tube Holder with Velcro Adhesive 7.5cm x 7.5cm	25
---------	--	----

COMBINE DRESSING

Highly absorbent layer of cellulose blend enclosed in a soft and conformable non woven fabric. Each non woven Combine is pre-cut and each end is sealed into a variety of sizes.



DEF711

DEF711	9cm x 20cm	25
DEF719	9cm x 20cm, Packet of 5	100
DEF720	9cm x 20cm, Packet of 10	50
DEF1322	9cm x 30cm, Packet of 2	50
DEF697	9cm x 90cm	25
DEF715	20cm x 20cm	25
DEF80	20cmx 60cm	100
DEF1160	20cm x 80cm	50

COMBINE ROLL

DEF714	Non Woven Combine Roll 9cm x 10m, Double Packed	16
DEF1004	Non Woven Combine Roll 20cm x 2.5m	25

COMBINE SPECIALISED



DEF804

DEF2289	Non Woven Amputation Dressing T-Shaped Melolin 40cm x 20cm	35
DEF804	Non Woven Amputation Dressing 40cm x 70cm	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DENTAL

BITE BLOCK



DEF1828

DEF1828	Bite Block 7.5cm X-ray with Tail	100
DEF1246	Bite Block 7.5cm X-ray with Tail Non Sterile	100

DENTAL PELLET

DEF1180	Dental Pellet Size #0, Packet of 5, Double Packed	100
---------	---	-----

TOOTHBRUSH



DEF690

DEF690	Toothbrush	25
--------	------------	----

DEFOGGER

Specially formulated solution that safely prevents lens fogging. Comes with a non-abrasive, x-ray detectable foam pad. Sterile.



28-0101	DeFogger Oval Bottle	20
---------	----------------------	----

DENTAL PRESSURE PACK



DEF799

DEF799	X-ray with Cotton Tail	100
DEF1063	X-ray with 35cm Tape, Packet of 2	100
DEF144	3cm x 4cm X-ray Small with 10cm tail, Packet of 2	50



DEF144

DEF145	4cm x 5cm X-ray Large with 9cm tail, Packet of 2	50
--------	--	----



DEF145

DEF1209	Dental Roll 4mm x 35mm, Packet of 5	100
---------	-------------------------------------	-----



DEF1209

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
ELECTROSURGICAL					
DIATHERMY PENCILS					
Hand activated Electrosurgical Pencils with an ergonomic design.					
88-E0031	Rocker Switch Pencil and Stainless Steel Blade	100			
			88-E0031		
88-E0041	Rocker Switch Pencil and Teflon Blade	100			
			88-E0041		
DIATHERMY PACK WITH STAINLESS STEEL BLADE					
DEF2864	1 x Diathermy Pencil with Rocker Switch and Stainless Steel Blade 1 x Quiver 5cm x 25cm with Metal Ring 1 x Scratch Pad Tip Cleaner	10			
DIATHERMY PACK WITH TEFLON BLADE & QUIVER					
DEF2863A	1 x Diathermy Pencil with Rocker Switch and Teflon Blade 1 x Quiver 5cm x 25cm with Metal Ring	10			

DIATHERMY TIP CLEANER		
This coated scratch pad cleaner quickly and easily removes debris from instrument tips and diathermys. The adhesive backing will secure the Tip Cleaner onto draped surfaces and contains a radiopaque strip to ensure X-ray detectability.		
DEF1886	Tip Cleaner 5cm x 5cm with Adhesive and Radiopaque Strip	50
		
DEF1886		
QUIVER		
A single use disposable polypropylene Quiver used to hold diathermy during surgical procedures. Quiver is recyclable (metal ring needs to be removed before being placed in recycling container).		
DEF2772	5cm x 25cm with Metal Ring	90
		
DEF2772		
DEF3068	5cm x 25cm with Cable Tie	90
		
DEF3068		
DEF3291	Large 9cm x 35cm with Metal Ring	40
DEF3296	Large 9cm x 35cm with Cable Tie	40

DIATHERMY TIP CLEANER

This coated scratch pad cleaner quickly and easily removes debris from instrument tips and diathermys. The adhesive backing will secure the Tip Cleaner onto draped surfaces and contains a radiopaque strip to ensure X-ray detectability.

DEF1886	Tip Cleaner 5cm x 5cm with Adhesive and Radiopaque Strip	50
		
DEF1886		

QUIVER

A single use disposable polypropylene Quiver used to hold diathermy during surgical procedures. Quiver is recyclable (metal ring needs to be removed before being placed in recycling container).

DEF2772	5cm x 25cm with Metal Ring	90
		
DEF2772		

DEF3068	5cm x 25cm with Cable Tie	90
		
DEF3068		

DEF3291	Large 9cm x 35cm with Metal Ring	40
DEF3296	Large 9cm x 35cm with Cable Tie	40

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.			
EYE								
EYE PAD								
								
DEF069								
DEF069	Eye Pad	100						
EYE SPEAR HIGHLY ABSORBENT								
								
DEF484								
DEF484	Eye Spear, Packet of 2	100						
DEF485	Eye Spear, Packet of 5	20						
DEF486	Eye Spear, Packet of 6	20						
DEF487	Eye Spear, Packet of 10	20						
EYE STROLL HIGHLY ABSORBENT								
								
DEF112								
DEF112	Eye Stroll, Packet of 5	50						
DEF081	Eye Stroll, Packet of 10	50						
DEF704	Eye Strolls Bulk Non Sterile	500						

FOAM PAD								
Versatile polyurethane foam pads used in the operating room for positioning of patients or relieving pressure areas for patients.								
DEF792								
DEF792	100mm x 100mm x 13mm	100						
DEF769	100mm x 100mm x 25mm	50						
DEF783	100mm x 180mm x 13mm	50						
DEF784	100mm x 180mm x 25mm	40						
DEF738	150mm x 150mm x 13mm	50						
FOAM PAD								
DEF738								
DEF1738	150mm x 315mm x 25mm	50						
DEF914	Plastics Foam 150mm x 200mm x 25mm	25						
DEF1549	200mm x 25mm x 25mm	100						
FOAM PAD								
DEF1549								
DEF739	250mm x 250mm x 13mm	50						
DEF1128	250mm x 250mm x 13mm with 4 x #31 Rubber Bands	50						
DEF786	300mm x 600mm x 13mm	20						

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
IDENTILOOPS					
These radio-opaque Identiloops are used in surgical procedures for retraction, identification or occlusion of arteries, veins, nerves or tendons. A soft and flexible, inert medical grade silicone which is oval in cross-section, allows for even pressure distribution to minimise damage to tissues. Available in 4 colours – Blue, White, Red and Yellow in mini and maxi sizes. All Double Packed.					
2200220-2	MINI-RED 40cm, Packet of 2	10	2202220-2	MAXI-RED 40cm, Packet of 2	10
2200220-4	MINI-RED 40cm, Packet of 4	10	2202220-4	MAXI-RED 40cm, Packet of 4	10
		2200220			2202220
2200420-2	MINI-YELLOW 40cm, Packet of 2	10	2202420-2	MAXI-YELLOW 40cm, Packet of 2	10
2200420-4	MINI-YELLOW 40cm, Packet of 4	10	2202420-4	MAXI-YELLOW 40cm, Packet of 4	10
		2200420			2202420
2200620-2	MINI-BLUE 40cm, Packet of 2	10	2202620-2	MAXI-BLUE 40cm, Packet of 2	10
2200620-4	MINI-BLUE 40cm, Packet of 4	10	2202620-4	MAXI-BLUE 40cm, Packet of 4	10
		2200620			2202620
2200920-2	MINI-WHITE 40cm, Packet of 2	10	2202920-2	MAXI-WHITE 40cm, Packet of 2	10
2200920-4	MINI-WHITE 40cm, Packet of 4	10	2202920-4	MAXI-WHITE 40cm, Packet of 4	10
		2200920			2202920
		2200920			DEF2182

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LARYNGOSCOPE

SURESCOPE

SURESCOPE is a single-use one piece laryngoscope with a built-in LED light source. This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting by obsoleting reusable lights sources, handles and blades.



C090120	MAC 2	20
C090130	MAC 3	14
C090140	MAC 4	14
C090205	MILLER 00	20
C090200	MILLER 0	20
C090210	MILLER 1	20
C090220	MILLER 2	14
C090230	MILLER 3	14

LIGHT HANDLE COVERS – CAMERA

Designed to be used in the operating room to create a sterile barrier between the user and the camera. Contain a lens inserted into the cover that is ultra clear and durable.



DEF3673

DEF3591	Camera Light Handle Cover 66mm with Lens Box	25
DEF3673	Camera Light Handle Cover 75mm with Lens Box	25

LIGHT HANDLE COVERS

Our range of flexible light handle covers are single sterile and come either coloured as blue or green. These adaptable covers are twisted or clicked onto the light handle adapter and are firmly held for duration required. We can supply the compatible light handle adapters for each operating room light. Easy and versatile to use.



DEF2082

DEF2082	Light Handle Cover Green – friction fit	50
DEF2919	Light Handle Cover Green – friction fit, packet of 2	40
DEF2675	Light Handle Cover Blue – click on	50



DEF2675

DEF3170	Light Handle Cover Blue – click on, packet of 2	40
DEF2557	Light Handle Cover Foam Blue	50




DEF2557

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
NEEDLE COUNTERS DEROYAL					
Excellent selection of Needle Counters from DeRoyal with easy to read numbers and description, facilitating correct counting of blades and needles used in a procedure. Full surface magnetic face provides extra space for sharps, pins, blade or needles. 'Y' codes indicate Needle Counter case is Yellow.					
DELUXE NEEDLE COUNTERS					
MAGNET					
25-0103Y	Needle Counter Single Magnetic/Clear Top, 10 Count, 20 Capacity	50			
			25-0103Y		
MAGNET/FOAM COMBINATION					
25-0401	Needle Counter Magnetic/Foam, 20 Count, 40 Capacity	20			
			25-0401		
25-0701	Needle Counter Foam Strips/Magnetic 15 Count, 30 Capacity	50			
			25-0701		
25-0702	Needle Counter Magnetic/Foam Strips, 30 Count, 60 Capacity	27			
			25-0702		
FOAM					
25-0504	Needle Counter Foam/Adhesive, 40 Count, 80 Capacity	48			
			25-0504		
25-0803	Needle Counter Foam Strips, 60 Count, 120 Capacity	27			
			25-0803		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
DEBLADE NEEDLE COUNTERS					
This range of Needle Counters all feature a blade removal system which provides safe and simple removal of contaminated blades.					
MAGNET					
25-1104Y	Needle Counter Double Magnetic, 30 Count, 60 Capacity with DeBlade	50	25-1402	Needle Counter Magnetic/ Foam block 10 Count, 20 Capacity with DeBlade	25
		25-1104Y			25-1701
25-1108	Needle Counter Double Magnetic 60 Count, 120 Capacity with DeBlade	27	25-1702	Needle Counter Magnetic/ Foam Strips 30 Count, 60 Capacity with DeBlade	27
		25-1108			25-1702
MAGNET/FOAM COMBINATION					
25-1401	Needle Counter Magnet/ Foam Block 20 Count, 40 Capacity with DeBlade	20	FOAM		
		25-1401	25-1801Y	Needle Counter Foam Strips 30 Count, 60 Capacity with DeBlade	25
					25-1801Y


CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

NEURO SKIN HOOKS



DEF405

DEF405	5 on a Tray Simmons Needle	20
DEF336	5 on a Tray Small Hook, Double Packed	20
DEF1032	5 on a Tray Size: 3/0 Large Hook	20



DEF1032

NEURO YARD STRIPS

DEF858	Neuro 1/2 (53cm) Yard Strips X-ray, Packet of 5	50
--------	---	----



DEF858

PREP APPLICATORS

DEF921	Prep Foam 50mm x 60mm x 25mm Cream No Handle	100
--------	--	-----



DEF921

DEF1846	Prep Foam on a Handle 15cm Blue	100
---------	---------------------------------	-----



DEF1846

DEF1847	Prep Foam on a Handle 15cm Blue, Packet of 2	100
DEF2035	Prep Foam on a Handle 20cm Orange	100



DEF2035

SAFETY PINS

DEF071SP	#2 (38mm) Single Packed	50
DEF072SP	#3 (51mm) Single packed	50
DEF073SP	#4 (57mm) Single packed	50


NON-STERILE

DEF635	#3 (51mm) Non Sterile	144
--------	-----------------------	-----



DEF635

DEF636	#4 (57mm) Non Sterile	144
--------	-----------------------	-----



DEF636

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SKIN MARKING PENS

Skin Markers are used during surgical procedures and contain medical grade, non irritating gentian violet ink. They mark a smooth, clear line on a patient's skin, and the ink is non toxic, non tattooing and quick drying.

26-001	Skin Markers Standard Tip with Ruler 15cm long	50
--------	--	----



26-001

26-002	Skin Marker Fine Point with Ruler 15cm	50
--------	--	----



26-002

26-019	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip	50
--------	--	----



26-019

26-020	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip with Ruler & Labels	50
DEF464	Prepacked Marking Pencil – Omnichrom 108-9	24



DEF464


NON-STERILE

DEF2669	Mini pre-op Skin Markers Non Sterile	100
DEF2183	Skin Markers Standard Tip	50
26-018PRNS	Prep-Resistant Ink Mini Skin Marker, Pre-Op	100
26-005PRNS	Prep-Resistant Ink Skin Marker, Standard Tip	100

SPONGES & SWABS


SPONGES ABDOMINAL X-RAY

Soft highly absorbent cotton sponges. Each complete with an X-ray radiopaque element to ensure detection.



DEF1350

DEF1350	23cm x 23cm, Packet of 5	100
DEF162	23cm x 23cm with Tape, Packet of 1	100
DEF160	23cm x 23cm with Tape, Packet of 3	100
DEF168	23cm x 23cm with Tape, Packet of 5	100
DEF1351	30cm x 30cm, Packet of 5, Double Packed	75
DEF1352	30cm x 48cm, Packet of 5, Double Packed	60



DEF162

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPONGES & SWABS (CONTINUED)

SPONGES PURCOTTON ABDOMINAL X-RAY

PurCotton X-ray Abdominal sponges are made from a unique 100% cotton sourced from the US. They are soft and comfortable, low linting, have superior absorption and wet tensile strength, and are 100% biodegradable. The PurCotton is the next evolution of cotton gauze.



DEF3741

DEF3741	30cm x 30cm X-ray 50g, Packet of 5	80
DEF3745	30cm x 45cm X-ray 50g, Packet of 5	60
DEF3743	45cm x 45cm X-ray 50g, Packet of 5	30

GAUZE SWABS X-RAY

Excellent absorbency and quality swabs with 8 ply cotton blend. Includes an X-ray radiopaque element for easy detection.



DEF044

DEF202	7.5cm x 7.5cm X-ray, Packet of 5, Double Packed	100
DEF044	10cm x 10cm X-ray, Packet of 5, Double Packed	100
DEF164	10cm x 10cm X-ray, Packet of 10, Double Packed	50

GAUZE SWABS PURCOTTON X-RAY



DEF3737

DEF3737	10cm x 10cm X-ray 40g, Packet of 5, Double Packed	500
DEF3738	10cm x 10cm X-ray 50g, Packet of 5, Double Packed	200

GAUZE SWABS GREEN



DEF190

DEF190	10cm x 10cm Green, Packet of 2	100
DEF192	10cm x 10cm Green, Packet of 3	100
DEF207	10cm x 10cm Green, Packet of 5, Double Packed	100
DEF196	10cm x 10cm Green, Non Sterile	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GAUZE SWABS FOLDED

DEF347

Gauze Swab 7cm x 21cm X-ray, Folded, Packet of 5

50



DEF347

DEF185

Gauze Swab Neuro 10cm x 10cm 32ply X-ray, Packet of 10, Double Packed

100



DEF185

DEF1074

Gauze Swab Neuro 10cm x 20cm X-ray, Folded, Packet of 5

25



DEF1074

GAUZE HEAD DRESSINGS

DEF611

Head Dressing Kit includes:
8 x Gauze Swabs 10cm x 20cm
20 x Gauze Swabs 10cm x 10cm
2 x 90cm x 90cm Gauze, Folded

35



DEF611

GAUZE PREP SPONGES X-RAY

Round sponges with X-ray radiopaque element. Excellent absorbency. Used as a prep sponge prior to a surgical procedure.



DEF251

DEF255	Sponges Prep X-ray Large, Packet of 1	100
DEF251	Sponges Prep X-ray Large, Packet of 2	100

GAUZE PEANUT SPONGES

Round sponges ideal for delicate surgical sponging or soft tissue dissection. X-ray radiopaque element for easy detection.

DEF315

X-ray, 5 x Small in an Instapack

50



DEF315

DEF046

X-ray, 5 x Small on a Safety Pin

50



DEF046

DEF313

X-ray, 5 x Small in a Bag

50



DEF313

DEF316

X-ray, 5 x Large in an Instapack

50



DEF316

SPECIALISED GAUZE



DEF169

DEF986	Donor Site Gauze Dressing 5cm x 10cm, Double Packed	100
DEF169	Recovery Gauze 5cm x 5cm Plain 32 Ply	100

32 | Hospitals Trust Defries

33

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL RETRACTOR

ONETRAC

ONETRAC is a single-use cordless surgical retractor with a built-in LED light source and smoke evacuation channel. It provides illumination in the depths of a surgical tissue pocket or cavity. No additional parts or components or fibre optic cables are required. Ideal for breast or reconstructive procedures, as well as other surgical procedures requiring access, illumination and retraction of soft tissue. Reduce cost, increase efficiency and eliminate unnecessary risk of cross-contamination with ONETRAC.




ONETRAC

C100110-5	135mm x 30mm	5
C100140-5	90mm x 22mm	5
C100150-5	40mm x 20mm	5

SNUGGERS PVC TUBING

Medical grade, clear and transparent Tubing Snuggers used during cardiac surgery. A radiopaque strip along length of tubing ensures X-ray detectability.

DEF3051	PVC Tubing Snuggers 16Fr x 17cm X-ray, 6 per Packet, Double Packed	25
---------	--	----



DEF3051

DEF3052	PVC Tubing Snuggers 16Fr x 22cm X-ray, 6 per Packet, Double Packed	25
---------	--	----

SCOPE WARMER

Long, double sided Scope Warmer from DeRoyal that pre-warms endoscopes to body temperate. Easy to use design provides 1 hour of warming time on each side.

TM500D	Scope Warmer	10
--------	--------------	----



TM500D

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ESMARCH BANDAGE

Offers excellent compression and elastic properties to achieve exsanguination, whilst delivering protection for latex sensitive patients and health professionals.

DEF1006	6.5cm x 2m Green Silicone	10
---------	---------------------------	----



DEF1006

DEF2764	10cm x 2.7m Blue Neoprene	20
---------	---------------------------	----



DEF2764

DEF1081	10cm x 3m Green Silicone	5
---------	--------------------------	---



DEF1081

DEF2668	15cm x 2.7m Blue Neoprene	20
---------	---------------------------	----



DEF2668

SILICONE TUBING

Our 100% medical grade, clear silicone tubing range is available pre-cut and sterile or as non-sterile rolls. Please contact our customer service to discuss sizing.

Sizing below is represented as:
ID (Internal Diameter) x OD (Outer Diameter) x Length of tube in cm or m.


DEF2307	0.30mm x 0.64mm x 20cm	75
DEF2522	1mm x 2mm x 50cm	50
DEF477	1.57mm X 2.41mm X 30cm, 2 per Packet, Double Packed	25
DEF1564A	1.57mm x 3.18mm x 50cm	25
DEF1545	1.57mm x 3.18mm x 1m	15
DEF1428	1.57 mm x 2.41mm x 15cm, Double Packed	100
DEF473	1mm x 3mm x 50cm, 4 per Packet	20
DEF1944A	3mm x 5mm x 2.3m	15
DEF2314	3.18mm x 6.35mm x 25cm	60
DEF2038	4mm x 6mm x 45cm	55
DEF1665	5mm x 7mm x 20cm, 2 per Packet, Double Packed	25
DEF1268	5mm x 7mm x 42cm	25
DEF481	5mm x 8mm x 45cm	50
DEF1253	5mm x 8mm x 1.6m, Double Packed	50
DEF2060	5mm x 8mm x 3.1m	8
DEF1261	5mm x 9mm x 1.5m	50
DEF2786	6.4mm x 9.5mm x 2cm, 2 per Packet	375
DEF2032	7mm x 10mm x 3m	8
DEF482	7mm x 12mm x 85cm, 2 per Packet, Double Packed	100
DEF1204	8mm x 14mm x 3m	5
DEF2468	10mm x 16mm x 1.5m	16

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE TUBING (CONTINUED)

SILICONE TUBING WITH CONNECTOR

Flexible pre-cut and sterile silicone tubing with connector attached.




DEF1186

DEF735	5mm x 8mm x 1m, Double Packed	50
DEF1186	5mm x 9mm x 30cm, Double Packed	50
DEF306	5mm x 9mm x 50cm Neuro	50
DEF1712	5mm x 9mm x 60cm	40
DEF1337	5mm x 9mm x 1m Neuro	50
DEF1231A	5mm x 9mm x 1.5m	30
DEF305	7mm x 11mm x 1m, Double Packed	50

SUTURE BOOTS WITH ADHESIVE TRAY

Suture boots offer a secure way to hold delicate sutures without slippage or damage during a surgical procedure. The tray comes with an adhesive strip for securing and easy removal. Radio-opaque element ensures X-ray detection. Available in a Yellow colour.




30-702

DEF1093	2 per Packet	25
DEF1094	4 per Packet	25
30-702	6 per Packet	50
DEF1096	8 per Packet	25
30-701	10 per Packet	50
DEF2284	12 per Packet, Double Packed	25


TUBING CONNECTOR

DEF465	Connector Small 4mm – 7mm	100
--------	---------------------------	-----



DEF465

DEF466	Connector Medium 7mm – 10mm	100
--------	-----------------------------	-----



DEF466

WOUND DRAINAGE

Wounds drains and evacuators provide drainage for various surgical procedures, allowing for the removal of blood and excess fluid to improve patient outcome. Our wound drainage systems are made up of evacuators, flat drains, fluted drains, round drains, capillary action drains and easy flow drains. They are all made from 100% Medical Grade Silicone, which is inert, soft and flexible and designed to reduce the potential of post – operative infections.

WOUND EVACUATOR BULB

Our wound evacuator is a clear bulb designed to slowly draw fluid and blood from the wound site. It creates an ideal vacuum to a maximum of 125mmHg pressure and is compatible with our flat drains, fluted drains and round drains.




31300380DS

31300380DS	Evacuator Bulb 100ml/cc	10
------------	-------------------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

WOUND FLAT DRAIN

Soft and comfortable silicone flat drains with multiple drainage holes (36 x 2) to enhance the fluid removal in a wound. They have an ergonomic flat shape which enables an easier placement. X-ray detectable.



11100200780DS

11100200780DS	Flat Drain White 7mm x 190mm	10
11100201080DS	Flat Drain White 10mm x 190mm	10

WOUND FLUTED DRAIN

Fluted drains are made of soft, flexible silicone and available in a round and flat shape, with or without attached trocar. Designed with a multi channel to allow for great drainage and assisting in minimal tissue in-growth. The inner core of the H design offers drain strength and flexibility.

FLAT FLUTED DRAIN


The flat fluted drains are available as a white, radio-opaque X-ray detectable with a clear tube. Available in 7mm and 10mm.

11101310780DS	Flat white drain 3/4 fluted 7mm with Trocar and Hub	10
11101410780DS	Flat white drain fully fluted 7mm with Trocar and Hub	10
11101311080DS	Flat white drain 3/4 fluted 10mm with Trocar and Hub	10
11101411080DS	Flat white drain fully fluted 10mm with Trocar and Hub	10


ROUND FLUTED DRAIN

The round fluted drains come in 2 versions, the white radio-opaque with a clear tube and the hubless drain which is clear with a blue X-ray detectable line. Available in 10,15, 19 & 24 Fr.

11003131080DS	Round clear full fluted drain hubless 10Fr with Trocar	10
---------------	--	----



11003131080DS



11003831080DS

11003831080DS	Round white full fluted drain with Hub 10Fr	10
11003131580DS	Round clear full fluted drain hubless 15Fr with Trocar	10
11003131980DS	Round clear full fluted hubless 19Fr with Trocar	10
11003811980DS	Round white full fluted drain hubless 19Fr	10
11003122480DS	Round clear full fluted hubless 24Fr	10

ROUND CLEAR DRAIN

Clear soft silicone drain 500:5 which consists of 5 side holes with a total drain length of 500mm. X-ray detectable.

11000301280DS	500:5 12Fr 2mm x 4mm	10
11000301880DS	500:5 18Fr 4mm x 6mm	10
11000302480DS	500:5 24Fr 5mm x 8mm	10
11000302780DS	500:5 27Fr 6mm x 9mm	10
11000303080DS	500:5 30Fr 7mm x 10mm	10
11000303680DS	500:5 36Fr 8mm x 12mm	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL CONSUMABLES

PASSIVE WOUND DRAINS

MULTILUMEN CAPILLARY ACTION DRAINS

These clear radio-opaque drains have varying lumens according to drainage requirements. The lumens can be separated insitu and positioned as necessary. The drainage takes place both through, under and over the lumens.



DEF1750

DEF1110	Multilumen drain 2 Lumens x 12.5cm	10
DEF1076	Multilumen drain 5 Lumens x 25cm	10
DEF1024	Multilumen drain 9 Lumens x 25cm	10
DEF1750	Multilumen drain 16 Lumens x 20cm	10
DEF1751	Multilumen drain 16 Lumens x 40cm	10

EASY FLOW DRAINS

These are a 30cm radio-opaque clear drain, which feature inner teeth to assist in maintaining an open drain.




21740632

1120010680DS	6mm X-ray, Double Packed	10
21740632	40mm X-ray, Double Packed	10

COTTON PRODUCTS

COTTON BALLS


Packaged and sterilised for ease of entry into a sterile field. Ideal for prepping, wound cleaning, perineal care and fluid application.



DEF074

DEF725	Large 3cm, Packet of 4	100
DEF074	Large 3cm, Packet of 5	100

COTTON WOOL




DEF859

DEF859	10cm x 10cm	25
DEF615	10cm x 15cm Double Packed	100

COTTON WOOL ROLL

Soft and absorbent cotton wool roll used for padding in patient care. Sterilised and ready to be used.



DEF874


DEF874	15cm x 55cm, Double Packed	50
DEF875	30cm x 90cm, Double Packed	25
DEF880	38cm x 165cm, Double Packed	15

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPONGES

SPONGES NASAL BOLSTER

Used after nasal or sinus surgery to absorb any discharge or ooze from the nose. High quality non woven net filled with absorbent cotton padding, and stretch elastic sides to secure around the ears. Extremely soft and comfortable for post surgery patients. Latex free and available as sterile or non sterile.



DEF1382


DEF1381	Nasal Bolster	20
DEF1382	Nasal Bolster	100

SPONGES NASAL BOLSTER – NON STERILE

DEF1380	Nasal Bolster Bulk	20
---------	--------------------	----


GAUZE SWABS PLAIN

Sterile and ready to be used for a wide variety of applications, including as a dressing, for prepping procedures and wound cleaning. Available as a sterile Plain 8 ply gauze.



DEF181

DEF183	7.5cm x 7.5cm, Packet of 3	100
DEF181	7.5cm x 7.5cm, Packet of 5	100
DEF371	7.5cm x 7.5cm, Packet of 20	50
DEF189	10cm x 10cm, Packet of 2	100
DEF186	10cm x 10cm, Packet of 3	100



DEF045

DEF045	10cm x 10cm, Packet of 5	100
DEF208	10cm x 10cm, Packet of 10	50
DEF213	10cm x 10cm, Packet of 10	300
DEF178	10cm x 10cm, Packet of 20	25
DEF177	10cm x 10cm, Packet of 25	25
DEF214	10cm x 10cm, Packet of 25	200
DEF1696	10cm x 10cm, Packet of 50	50
DEF1033	10cm x 20cm, Packet of 4	100
DEF198	10cm x 20cm, Packet of 5	250
DEF195	10cm x 20cm, Packet of 10	50
DEF853	10cm x 80cm, Packet of 4, 16 ply Folded	50

PURCOTTON GAUZE SWABS PLAIN




DEF3559	10cm x 10cm, Packet of 25	100
---------	---------------------------	-----




GAUZE SWABS PLAIN – NON STERILE

DEF136	10cm x 10cm Plain Bulk	100
--------	------------------------	-----

NONWOVEN GAUZE SWABS


DEF1131	5cm x 5cm, Packet of 5	100
DEF197	7.5 x 7.5cm, Plain Bulk Non Sterile	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
SWAB STICKS/COTTON TIPPED APPLICATORS					
Packaged and sterilised, ready to use for taking a swab or specimen from a patient. Paper, plastic or wooden shaft in different length sizes with a single or double ended soft cotton head.					
SWAB STICKS WITH WOODEN SHAFT					
					
DEF057					
DEF057	7.5cm Single Ended, Packet of 3	100			
DEF442	7.5cm Single Ended, Packet of 5	100			
DEF439	15cm Single Ended, Packet of 2	100			
DEF059	15cm Single Ended, Packet of 3	100			
DEF443	15cm Single Ended, Packet of 5	100			
					
DEF439					
SWAB STICKS WITH WOODEN SHAFT – NON STERILE					
DEF062	15cm Single Ended Bulk	1000			
SWAB STICKS WITH PAPER SHAFT					
					
DEF421					
DEF424	7.5cm Double Ended, Packet of 5	100			
DEF421	7.5cm Double Ended, Packet of 10	100			
DEF104	7.5cm Double Ended, Packet of 50	80			

SWAB STICKS WITH PAPER SHAFT – NON STERILE					
DEF061	7.5cm Double Ended Bulk	1250			
SWAB STICKS WITH PLASTIC SHAFT					
					
DEF058					
DEF058	7.5cm Double Ended, Packet of 3	100			
SWAB STICK MOUTH & THROAT JUMBO					
					
DEF401					
DEF400	15cm, Packet of 2	150			
DEF401	15cm, Packet of 3	150			
DEF402	15cm, Packet of 5	150			
DEF455	20cm, Packet of 2	150			
DEF407	20cm, Packet of 3	150			
DEF1402	20cm, Packet of 4	100			
DEF1593	20cm, Packet of 1	200			
					
DEF407					
SWAB STICK MOUTH & THROAT JUMBO – NON STERILE					
DEF055	15cm Bulk	200			
DEF054	20cm Bulk	200			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
BURNS					
BURNS DRESSING					
Specialised products for patients with burns. Non adherent, highly absorbent wound dressings.					
					
DEF987					
DEF199	10cm x 20cm Plain, Packet of 20	50			
DEF968	10cm x 100cm Plain, Packet of 5	100			
DEF987	40cm x 50cm Melolin, Double Packed	50			
DEF679	50cm x 70cm, Double Packed	35			
DEF3570	50cm x 175cm Melolin, Double Packed	25			
DEF974	90cm x 2m Plain	50			
					
DEF679					
BURNS GAUZE					
Soft roll suitable for wounds and burns.					
DEF965	10cm x 80cm, Packet of 16	50			
DEF2436	25cm x 90cm, Packet of 2	350			
DEF1039	60cm x 90cm, Packet of 6	100			
DEF342	90cm x 2m	50			
DEF344	90cm x 2m, Packet of 4	25			
DEF1662	90cm x 5m	10			
BURNS PLASTIC BAG					
DEF650	305mm x 390mm	100			

SANITARY DISPOSAL BAG		
DEF709	Packet of 100, Non Sterile	20
		
DEF709		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.																																																																																																																						
<div>SILICONE TUBING ROLL – NON STERILE</div> <div>Excellent range of medical grade silicone tubing available in a variety of different metric sizes. Available as non sterile rolls of 15 or 25 metres.</div> <div></div> <div>2110030239</div> <table><tr><td>2110030239</td><td>0.30mm (ID) x 0.64mm (OD) 15m Roll</td></tr><tr><td>2110050134</td><td>0.5mm (ID) x 1mm (OD) 25m Roll</td></tr><tr><td>203005009</td><td>0.51mm (ID) x 0.94mm (OD) 15m Roll</td></tr><tr><td>203006011</td><td>0.64mm (ID) x 1.19mm (OD) 15m Roll</td></tr><tr><td>2110070439</td><td>0.7mm (ID) x 1.6mm (OD) 15m Roll</td></tr><tr><td>201010020</td><td>1mm (ID) x 2mm (OD) 25m Roll</td></tr><tr><td>2110100234</td><td>1mm (ID) x 3mm (OD) 25m Roll</td></tr><tr><td>2110100439</td><td>1.016mm (ID) x 2.159mm (OD) 15m Roll</td></tr><tr><td>203010021</td><td>1.02mm (ID) x 2.16mm (OD) 15m Roll</td></tr><tr><td>203014019</td><td>1.47mm (ID) x 1.96mm (OD) 15m Roll</td></tr><tr><td>2110150434</td><td>1.5mm (ID) x 3mm (OD) 25m Roll</td></tr><tr><td>2110150939</td><td>1.57mm (ID) x 2.41mm (OD) 15m Roll</td></tr><tr><td>2110151039</td><td>1.57mm (ID) x 3.18mm (OD) 15m Roll</td></tr><tr><td>2110190239</td><td>1.98mm (ID) x 3.18mm (OD) 15m Roll</td></tr><tr><td>201020030</td><td>2mm (ID) x 3mm (OD) 25m Roll</td></tr><tr><td>2110200334</td><td>2mm (ID) x 4mm (OD) 25m Roll</td></tr><tr><td>2110200534</td><td>2mm (ID) x 6mm (OD) 25m Roll</td></tr><tr><td>2110260239</td><td>2.65mm (ID) x 4.88mm (OD) 15m Roll</td></tr><tr><td>2110300234</td><td>3mm (ID) x 5mm (OD) 25m Roll</td></tr><tr><td>2110300334</td><td>3mm (ID) x 6mm (OD) 25m Roll</td></tr><tr><td>2110300434</td><td>3mm (ID) x 7mm (OD) 25m Roll</td></tr><tr><td>2110310239</td><td>3.18mm (ID) x 6.35mm (OD) 15m Roll</td></tr><tr><td>2110330239</td><td>3.35mm (ID) x 4.65mm (OD) 15m Roll</td></tr><tr><td>2110400434</td><td>4mm (ID) x 6mm (OD) 25m Roll</td></tr></table> <table><tr><td>2110400634</td><td>4mm (ID) x 7mm (OD) 25m Roll</td></tr><tr><td>2110400734</td><td>4mm (ID) x 8mm (OD) 25m Roll</td></tr><tr><td>2110470139</td><td>4.76mm (ID) x 7.93mm (OD) 15m Roll</td></tr><tr><td>20304707901</td><td>4.77mm (ID) x 7.9mm (OD) 15m Roll</td></tr><tr><td>20304709501</td><td>4.77mm (ID) x 9.50 (OD) 15m Roll</td></tr><tr><td>2110500234</td><td>5mm (ID) x 7mm (OD) 25m Roll</td></tr><tr><td>2110500334</td><td>5mm (ID) x 8mm (OD) 25m Roll</td></tr><tr><td>2110500534</td><td>5mm (ID) x 9mm (OD) 25m Roll</td></tr><tr><td>2110500734</td><td>5mm (ID) x 10mm (OD) 25m Roll</td></tr><tr><td>2110600234</td><td>6mm (ID) x 9mm (OD) 25m Roll</td></tr><tr><td>2110600334</td><td>6mm (ID) x 10mm (OD) 25m Roll</td></tr><tr><td>2110600434</td><td>6mm (ID) x 11mm (OD) 25m Roll</td></tr><tr><td>2110600534</td><td>6mm (ID) x 12mm (OD) 25m Roll</td></tr><tr><td>2110630139</td><td>6.4mm (ID) x 9.5mm (OD) 15m Roll</td></tr><tr><td>2110630239</td><td>6.4mm (ID) x 11mm (OD) 15m Roll</td></tr><tr><td>2110630339</td><td>6.4mm (ID) x 12.7mm (OD) 15m Roll</td></tr><tr><td>2110700334</td><td>7mm (ID) x 10mm (OD) 25m Roll</td></tr><tr><td>2110700534</td><td>7mm (ID) x 11mm (OD) 25m Roll</td></tr><tr><td>2110700634</td><td>7mm (ID) x 12mm (OD) 25m Roll</td></tr><tr><td>201070130</td><td>7mm (ID) x 13mm (OD) 25m Roll</td></tr><tr><td>2110790239</td><td>7.9mm (ID) x 12.7mm (OD) 15m Roll</td></tr><tr><td>2110800334</td><td>8mm (ID) x 12mm (OD) 25m Roll</td></tr><tr><td>2110800534</td><td>8mm (ID) x 14mm (OD) 25m Roll</td></tr><tr><td>2110800734</td><td>8mm (ID) x 16mm (OD) 25m Roll</td></tr><tr><td>2110900334</td><td>9mm (ID) x 13mm (OD) 25m Roll</td></tr><tr><td>2110950139</td><td>9.5mm (ID) x 12.7mm (OD) 15m Roll</td></tr><tr><td>2110950239</td><td>9.5mm (ID) x 14.3mm (OD) 15m Roll</td></tr><tr><td>2110950339</td><td>9.5mm (ID) x 15.9mm (OD) 15m Roll</td></tr><tr><td>2111000334</td><td>10mm (ID) x 14mm (OD) 25m Roll</td></tr><tr><td>2111000434</td><td>10mm (ID) x 15mm (OD) 25m Roll</td></tr><tr><td>201100160</td><td>10mm (ID) x 16mm (OD) 25m Roll</td></tr><tr><td>201100180</td><td>10mm (ID) x 18mm (OD) 25m Roll</td></tr><tr><td>2111200334</td><td>12mm (ID) x 17mm (OD) 25m Roll</td></tr><tr><td>2111250134</td><td>12.5mm (ID) x 19mm (OD) 25m Roll</td></tr><tr><td>2111270239</td><td>12.7mm (ID) x 19.05mm (OD) 15m Roll</td></tr></table>						2110030239	0.30mm (ID) x 0.64mm (OD) 15m Roll	2110050134	0.5mm (ID) x 1mm (OD) 25m Roll	203005009	0.51mm (ID) x 0.94mm (OD) 15m Roll	203006011	0.64mm (ID) x 1.19mm (OD) 15m Roll	2110070439	0.7mm (ID) x 1.6mm (OD) 15m Roll	201010020	1mm (ID) x 2mm (OD) 25m Roll	2110100234	1mm (ID) x 3mm (OD) 25m Roll	2110100439	1.016mm (ID) x 2.159mm (OD) 15m Roll	203010021	1.02mm (ID) x 2.16mm (OD) 15m Roll	203014019	1.47mm (ID) x 1.96mm (OD) 15m Roll	2110150434	1.5mm (ID) x 3mm (OD) 25m Roll	2110150939	1.57mm (ID) x 2.41mm (OD) 15m Roll	2110151039	1.57mm (ID) x 3.18mm (OD) 15m Roll	2110190239	1.98mm (ID) x 3.18mm (OD) 15m Roll	201020030	2mm (ID) x 3mm (OD) 25m Roll	2110200334	2mm (ID) x 4mm (OD) 25m Roll	2110200534	2mm (ID) x 6mm (OD) 25m Roll	2110260239	2.65mm (ID) x 4.88mm (OD) 15m Roll	2110300234	3mm (ID) x 5mm (OD) 25m Roll	2110300334	3mm (ID) x 6mm (OD) 25m Roll	2110300434	3mm (ID) x 7mm (OD) 25m Roll	2110310239	3.18mm (ID) x 6.35mm (OD) 15m Roll	2110330239	3.35mm (ID) x 4.65mm (OD) 15m Roll	2110400434	4mm (ID) x 6mm (OD) 25m Roll	2110400634	4mm (ID) x 7mm (OD) 25m Roll	2110400734	4mm (ID) x 8mm (OD) 25m Roll	2110470139	4.76mm (ID) x 7.93mm (OD) 15m Roll	20304707901	4.77mm (ID) x 7.9mm (OD) 15m Roll	20304709501	4.77mm (ID) x 9.50 (OD) 15m Roll	2110500234	5mm (ID) x 7mm (OD) 25m Roll	2110500334	5mm (ID) x 8mm (OD) 25m Roll	2110500534	5mm (ID) x 9mm (OD) 25m Roll	2110500734	5mm (ID) x 10mm (OD) 25m Roll	2110600234	6mm (ID) x 9mm (OD) 25m Roll	2110600334	6mm (ID) x 10mm (OD) 25m Roll	2110600434	6mm (ID) x 11mm (OD) 25m Roll	2110600534	6mm (ID) x 12mm (OD) 25m Roll	2110630139	6.4mm (ID) x 9.5mm (OD) 15m Roll	2110630239	6.4mm (ID) x 11mm (OD) 15m Roll	2110630339	6.4mm (ID) x 12.7mm (OD) 15m Roll	2110700334	7mm (ID) x 10mm (OD) 25m Roll	2110700534	7mm (ID) x 11mm (OD) 25m Roll	2110700634	7mm (ID) x 12mm (OD) 25m Roll	201070130	7mm (ID) x 13mm (OD) 25m Roll	2110790239	7.9mm (ID) x 12.7mm (OD) 15m Roll	2110800334	8mm (ID) x 12mm (OD) 25m Roll	2110800534	8mm (ID) x 14mm (OD) 25m Roll	2110800734	8mm (ID) x 16mm (OD) 25m Roll	2110900334	9mm (ID) x 13mm (OD) 25m Roll	2110950139	9.5mm (ID) x 12.7mm (OD) 15m Roll	2110950239	9.5mm (ID) x 14.3mm (OD) 15m Roll	2110950339	9.5mm (ID) x 15.9mm (OD) 15m Roll	2111000334	10mm (ID) x 14mm (OD) 25m Roll	2111000434	10mm (ID) x 15mm (OD) 25m Roll	201100160	10mm (ID) x 16mm (OD) 25m Roll	201100180	10mm (ID) x 18mm (OD) 25m Roll	2111200334	12mm (ID) x 17mm (OD) 25m Roll	2111250134	12.5mm (ID) x 19mm (OD) 25m Roll	2111270239	12.7mm (ID) x 19.05mm (OD) 15m Roll
2110030239	0.30mm (ID) x 0.64mm (OD) 15m Roll																																																																																																																										
2110050134	0.5mm (ID) x 1mm (OD) 25m Roll																																																																																																																										
203005009	0.51mm (ID) x 0.94mm (OD) 15m Roll																																																																																																																										
203006011	0.64mm (ID) x 1.19mm (OD) 15m Roll																																																																																																																										
2110070439	0.7mm (ID) x 1.6mm (OD) 15m Roll																																																																																																																										
201010020	1mm (ID) x 2mm (OD) 25m Roll																																																																																																																										
2110100234	1mm (ID) x 3mm (OD) 25m Roll																																																																																																																										
2110100439	1.016mm (ID) x 2.159mm (OD) 15m Roll																																																																																																																										
203010021	1.02mm (ID) x 2.16mm (OD) 15m Roll																																																																																																																										
203014019	1.47mm (ID) x 1.96mm (OD) 15m Roll																																																																																																																										
2110150434	1.5mm (ID) x 3mm (OD) 25m Roll																																																																																																																										
2110150939	1.57mm (ID) x 2.41mm (OD) 15m Roll																																																																																																																										
2110151039	1.57mm (ID) x 3.18mm (OD) 15m Roll																																																																																																																										
2110190239	1.98mm (ID) x 3.18mm (OD) 15m Roll																																																																																																																										
201020030	2mm (ID) x 3mm (OD) 25m Roll																																																																																																																										
2110200334	2mm (ID) x 4mm (OD) 25m Roll																																																																																																																										
2110200534	2mm (ID) x 6mm (OD) 25m Roll																																																																																																																										
2110260239	2.65mm (ID) x 4.88mm (OD) 15m Roll																																																																																																																										
2110300234	3mm (ID) x 5mm (OD) 25m Roll																																																																																																																										
2110300334	3mm (ID) x 6mm (OD) 25m Roll																																																																																																																										
2110300434	3mm (ID) x 7mm (OD) 25m Roll																																																																																																																										
2110310239	3.18mm (ID) x 6.35mm (OD) 15m Roll																																																																																																																										
2110330239	3.35mm (ID) x 4.65mm (OD) 15m Roll																																																																																																																										
2110400434	4mm (ID) x 6mm (OD) 25m Roll																																																																																																																										
2110400634	4mm (ID) x 7mm (OD) 25m Roll																																																																																																																										
2110400734	4mm (ID) x 8mm (OD) 25m Roll																																																																																																																										
2110470139	4.76mm (ID) x 7.93mm (OD) 15m Roll																																																																																																																										
20304707901	4.77mm (ID) x 7.9mm (OD) 15m Roll																																																																																																																										
20304709501	4.77mm (ID) x 9.50 (OD) 15m Roll																																																																																																																										
2110500234	5mm (ID) x 7mm (OD) 25m Roll																																																																																																																										
2110500334	5mm (ID) x 8mm (OD) 25m Roll																																																																																																																										
2110500534	5mm (ID) x 9mm (OD) 25m Roll																																																																																																																										
2110500734	5mm (ID) x 10mm (OD) 25m Roll																																																																																																																										
2110600234	6mm (ID) x 9mm (OD) 25m Roll																																																																																																																										
2110600334	6mm (ID) x 10mm (OD) 25m Roll																																																																																																																										
2110600434	6mm (ID) x 11mm (OD) 25m Roll																																																																																																																										
2110600534	6mm (ID) x 12mm (OD) 25m Roll																																																																																																																										
2110630139	6.4mm (ID) x 9.5mm (OD) 15m Roll																																																																																																																										
2110630239	6.4mm (ID) x 11mm (OD) 15m Roll																																																																																																																										
2110630339	6.4mm (ID) x 12.7mm (OD) 15m Roll																																																																																																																										
2110700334	7mm (ID) x 10mm (OD) 25m Roll																																																																																																																										
2110700534	7mm (ID) x 11mm (OD) 25m Roll																																																																																																																										
2110700634	7mm (ID) x 12mm (OD) 25m Roll																																																																																																																										
201070130	7mm (ID) x 13mm (OD) 25m Roll																																																																																																																										
2110790239	7.9mm (ID) x 12.7mm (OD) 15m Roll																																																																																																																										
2110800334	8mm (ID) x 12mm (OD) 25m Roll																																																																																																																										
2110800534	8mm (ID) x 14mm (OD) 25m Roll																																																																																																																										
2110800734	8mm (ID) x 16mm (OD) 25m Roll																																																																																																																										
2110900334	9mm (ID) x 13mm (OD) 25m Roll																																																																																																																										
2110950139	9.5mm (ID) x 12.7mm (OD) 15m Roll																																																																																																																										
2110950239	9.5mm (ID) x 14.3mm (OD) 15m Roll																																																																																																																										
2110950339	9.5mm (ID) x 15.9mm (OD) 15m Roll																																																																																																																										
2111000334	10mm (ID) x 14mm (OD) 25m Roll																																																																																																																										
2111000434	10mm (ID) x 15mm (OD) 25m Roll																																																																																																																										
201100160	10mm (ID) x 16mm (OD) 25m Roll																																																																																																																										
201100180	10mm (ID) x 18mm (OD) 25m Roll																																																																																																																										
2111200334	12mm (ID) x 17mm (OD) 25m Roll																																																																																																																										
2111250134	12.5mm (ID) x 19mm (OD) 25m Roll																																																																																																																										
2111270239	12.7mm (ID) x 19.05mm (OD) 15m Roll																																																																																																																										

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------



CUSTOM
PROCEDURE AND
STANDARD PACKS

CUSTOM PROCEDURE PACKS

COST EFFECTIVE AND TIME EFFICIENT:
A PACK FOR EVERY PROCEDURE

Our Custom Procedure Pack range covers all areas of the hospital from Cardiac, Orthopaedic, Vascular, ENT and more. They contain all the individual consumables for a surgical procedure, wrapped in a 'field' and sterilised. All products are ready to use when opened, unlike individually wrapped consumables that are touched many times prior to use and require opening separately. Plus there is reduced wastage. Using a Procedure pack saves on setup time for procedures and ensures that all components are readily available for every case. Ordering and inventory control is easier.

CUSTOM TAILOR A PROCEDURE PACK TO MEET YOUR SPECIFIC NEEDS

Select individual consumables such as patient drapes and surgical gowns, to diathermy, suction tubing, abdominal packs, gauze and instruments. Your pack will only include items specifically chosen by you, the surgeon and theatre nurses. Defries assembles these customised packs in cleanrooms according to hospital specifications.

Custom procedure packs can be varied to include cost-effective components as they become available, and allow easy accountability of costs per procedure. As a local Melbourne manufacturer we are able to respond quickly to any changes requested by a hospital to the configuration of your custom packs.

Designing and implementing Defries custom procedure packs is an easy process:



We'll provide you with ongoing education, listen to feedback and are always available to discuss your specific needs.

READY-FOR-USE STANDARD PROCEDURE PACKS

We have a range of standard procedure packs such as the anaesthetic packs, epidural packs, regional trays, catheter packs, Bronny IV insertion pack, lumbar puncture packs, ICU packs including central line packs, arterial line packs, all designed, made and ready to be ordered.



THEATRE & WARD ANAESTHETIC PACKS

All our anaesthetic packs come sterile, prepacked and ready to use. Indications for use include:

- General and Local Anaesthetic procedures
- Spinal/Epidural/Regional/Lumbar Puncture procedures
- IV Access.

Or please contact us to customise your own anaesthetic pack.

ANAESTHETIC PACKS		
DEF1409	1 x Dressing Tray 24cm x 15cm x 4cm 2 x Gauze Swabs 10cm x 10cm Green	36
		
DEF1409		
DEF1162	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 3 x Gauze Swabs 10cm x 10cm Green 1 x Towel 2 Ply 40cm x 42cm	60
		
DEF1162		

ANAESTHETIC TRAY		
DEF1524	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 2 x Gauze Swabs 10cm x 10cm Green 1 x Dressing Towel 40cm x 42cm	80

ANAESTHETIC PACKS

DEF1618	1 x Drape Blue Plastic 75cm x 75cm with 10cm x 15cm Fenestration with Adhesive Tabs on Corners of Fenestration and Corners of Drape 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Forcep Plastic Green 1 x Syringe Luer Slip 5ml 1 x Syringe Luer Slip 3ml 1 x Dressing Towel Quick Dry 33cm x 33cm 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt Drawing up 1 x Needle 25g x 1" Wrapped in a Blue Plastic Field 75cm x 75cm	60
		
DEF1618		
DEF4117B	1 x Dressing Tray 20.5cm x 26.5cm x 3.3cm with 4 compartments 1 x Syringe 10ml Luer Slip Eccentric Tip 1 x Syringe 30ml Luer Slip Eccentric Tip 1 x Syringe 2ml Luer Slip 1 x Syringe 5ml Luer Slip 2 x Gauze Swabs 10cm x 10cm Green Wrapped in a Blue Plastic Field 75cm x 75cm	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
ANAESTHETIC PROCEDURE PACK			EPIDURAL PACKS		
DEF2255A	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 91cm x 115cm 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 3ml Luer Lock 1 x Syringe 5ml Luer Lock 2 x Forceps Pearman Plastic 1 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 25g x 1.5" Hypodermic 2 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 1 x Gallipot 60ml 1 x Particle Filter 5 Micron Luer Lock Wrapped in a Blue Plastic Field 91cm x 115cm	60	DEF1796	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 2 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swab 10cm x 10cm Green 1 x Dressing Towel T Large 32cmx 60cm 2 x Forceps Plastic Green Wrapped in a Blue Plastic Field 91cm x 115cm	50
REGIONAL TRAY WITH LARGE DRAPE			DEF1988B	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 2 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Syringe Luer Slip 5ml 1 x Syringe Luer Slip 3ml 2 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 25g x 1" Hypodermic 1 x Needle 25g x 1.5" Hypodermic 1 x Dressing Towel T Small 32cm x 32cm 3 x Absorbent Foam Applicator 20cm Orange Wrapped in a Blue Plastic Field 75cm x 75cm	50
DEF1479	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 4 x Gauze Swabs 10cm x 10cm Green 2 x Forceps Plastic Green 1 x Towel White 2 Ply 40cm x 42cm Wrapped in a Blue Plastic Field 91cm x 115cm	60	DEF3085	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Syringe Luer Slip 5ml 1 x Syringe Luer Slip 3ml 2 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 25g x 1" Hypodermic 1 x Needle 25g x 1.5" Hypodermic 1 x Dressing Towel T Small 30cm x 32cm Wrapped in a Blue Plastic Field 75cm x 75cm	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
COMPLETE EPIDURAL PACKS			EPIDURAL/SPINAL TRAY		
DEF1480	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swab 10cm x 10cm Green 1 x Forceps Plastic Green 1 x Syringe Luer Slip 5ml 1 x Syringe Luer Slip 3ml 1 x Needle 19g x 1.5" Hypodermic 1 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 25g x 1" Hypodermic 1 x Dressing Towel T Small 30cm x 32cm Wrapped in a Blue Plastic Field 75cm x 75cm	60	DEF2345A	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 20cm x 25cm x 2.5cm with 3 compartments 1 x Syringe Luer Lock 10ml 1 x Syringe Luer Lock 5ml 1 x Syringe Luer Slip 3ml 1 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 25g x 1.5" Hypodermic 1 x Particle Filter 5 Micron Luer Lock 2 x Dressing Towel T Small 32cm x 32cm 5 x Gauze Swabs 10cm x 10cm Green 1 x Gallipot 60ml 1 x Foam Prep Stick Blue 15cm	24
DEF2464	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Syringe Luer Slip 5ml 1 x Syringe Luer Slip 3ml 1 x Needle 19g x 1.5" Hypodermic 1 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 23g x 1.25" Hypodermic 1 x Dressing Towel T Small 32cm x 32cm 2 x Forceps Plastic Green 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	60			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

THEATRE & WARD CATHETER PACKS

All Catheter Procedure Packs have been designed in consultation with the healthcare professionals who use them. They include all the consumables used to insert or change a catheter aseptically. Alternatively please contact us to customise a catheter pack to your own specific requirements.

CATHETER PACKS



DEF1601	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 2 x Dressing Tray 23cmx 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 1 x Blue Huck Towel 40cm x 68cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF2074	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Mediwipe Towel 35cm x 60cm 2 x Forceps Plastic Pearman 1 x Syringe 10ml Luer Slip 5 x Gauze Swabs 10cm x 10cm X-ray 2 x Kidney Dishes 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24

DEF2061A	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Mediwipe Towel 35cm x 60cm 1 x Syringe 10ml Luer Slip Concentric Tip 2 x Gallipot 150ml 1 x Forcep Plastic Pearman 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Green 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Sponge Holder 19cm Plastic Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF2649A	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 2 x Kidney Dishes 700ml Clear 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Syringe 10ml Luer Slip 2 x Gallipots 150ml 1 x Sponge Holder 24cm Plastic Wrapped in a SMS Field 130cm x 150cm	18
DEF1262	1 x Drape SMS 61cm x 61cm with a 10cm x 10cm Fenestration 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Gallipot 60ml 1 x Kidney Dish 700ml Clear 1 x Forcep Plastic Pearman 1 x Syringe Luer Slip 10ml Concentric Tip 1 x Dressing Towel T Small 32cm x 32cm 1 x Tube Clamp with Lock 5 x Gauze Swabs 10cm x 10cm X-ray Wrapped in a Blue Plastic Field 75cm x 75cm	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DEF1333T	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe Luer Slip 10ml Concentric Tip 2 x Forceps Plastic Green 1 x Dressing Towel T Small 32cm x 32cm 1 x Water Sachet 10ml Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1650A	1 x Drape Azure Minor Procedure 50cm x 60cm with 6.5cm Circular Fenestration and Adhesive around Fenestration 2 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 2 x Forceps Plastic Green 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Dressing Towel White 40cm x 42cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1335	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Dressing Tray 23cmx 9cm x 3cm with 3 Compartments 5 x Cotton Balls Large 2 x Forceps Plastic Green 1 x Dressing Towel T Small 32cm x 32cm Wrapped in a Blue Plastic Field 91cm x 115cm	50
DEF1522	3 x Drape Polybacked 33cm x 46cm 2 x Dressing Tray 23cmx 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Plain 2 x Forceps Green Plastic 1 x Syringe 10ml Luer Slip Wrapped in a Blue Plastic Field 75cm x 75cm	60

DEF2253C	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Dressing Towel T Small 32cm x 32cm 2 x Forceps Plastic Green 2 x Gloves Nitrile Large 2 x Syringe 10ml Luer Slip 5 x Cotton Balls Large Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1036	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Kidney Dishes 700ml Clear 2 x Gallipots 150ml 2 x Forceps Plastic Green 1 x Dressing Towel T Large 32cm x 60cm 6 x Cotton Balls Large Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1038A	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 1 x Kidney Dish 700ml Clear 1 x Gallipot 150ml 1 x Tube Clamp with Lock Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1881	3 x Blue Huck Towels 40cm x 68cm 2 x Kidney Dishes 700ml Clear 1 x Gallipot 150ml 7 x Non Woven Balls 2 x Non Woven Swabs 7.5cm x 7.5cm 2 x Forceps Plastic Pearman Wrapped in a Blue Plastic Field 75cm x 75cm	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL PACKS

THEATRE CATHETER PACK		
DEF1932	1 x Drape PolyBacked 50cm x 60cm with 10cm Circular Fenestration 2 x Gallipot 150ml 1 x Kidney Dish 700ml Clear 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Sponge Holder 24cm Plastic 1 x Syringe 30ml Eccentric Tip 1 x Towel D 32cm x 34cm Wrapped in a Blue Plastic Field 91cm x 115cm	30
DEF1999A	1 x SMS Drape 91cm x 91cm with 10cm x 15cm Fenestration 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Kidney Dish 700ml Clear 1 x Blue Huck Towel 40cm x 68cm 1 x Sponge Holder 24cm Plastic 2 x Gallipot 150ml 1 x Syringe Luer Slip 10ml Concentric Tip Wrapped in a Blue Plastic Field 75cm x 75cm	24

IV INSERTION PACKS		
For the insertion of an IV line.		
BRONNY IV PACK		
DEF920	2 x Gauze Swabs 5cm x 5cm Plain 1 x Plastic Drape 45cm x 60cm	200
BRONNY IV PACK WITH 4 SWABS		
DEF1676	4 x Gauze Swabs 5cm x 5cm Plain 1 x Plastic Drape 45cm x 60cm	200

MSU (MID STREAM URINE) PACK		
DEF1866	1 x Kidney Dish 400ml Green 1 x Gallipot 60ml 3 x Non Woven Swabs 7.5cm x 7.5cm 7 x Non Woven Balls 2 x Vinyl Gloves Large Wrapped in Blue Plastic Field 50cm x 50cm	50

MOUTH TOILET PACK		
DEF399	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 10 x M & T Swab Sticks 15cm Paper 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	50

FINE PROCEDURE PACK		
DEF2453A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Kidney Dish 700ml Clear 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Scalpel Handle # 3 4 x Gauze Swabs 7.5cm x 7.5cm Plain 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

THYROID PACK		
DEF1834A	1 x Towel Blue Cotton 40cm x 68cm 2 x Forcep Spencer Wells 13cm Straight 1 x Scissors Sharp/Blunt 13cm 1 x Staple Remover Scissor Action with Plastic Handle 10 x Gauze Swabs 10cm x 10cm Plain Wrapped in a SMS Blue Field 61cm x 61cm	20

PLASTICS PACK		
DEF219	1 x Dressing Tray Long 23cm x 9cm x 3cm 1 x Gallipot 60ml 12 x Gauze Swabs 10cm x 10cm Plain 2 x Towel White 2 Ply 40cm x 42cm 2 x Swab Sticks 15cm with Wooden Shaft Wrapped in a Blue Plastic Field 50cm x 50cm	60

LARGE DRESSING PACK		
DEF2234A	1 x Dressing Tray 23.5cm x 29.5cm x 5cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml 2 x Dressing Towel T Small 32cm x 32cm 10 x Non Woven Balls 3 x Forceps Plastic Green 10 x Non Woven Swabs 7.5cm x 7.5cm Wrapped in a Blue Plastic Field 91cm x 115cm	10

LUMBAR PUNCTURE PACK		
DEF2746B	1 x Drape SMS 61cm x 61cm with 14cm Circular Fenestration and with Corner Tabs 1 x Dressing Tray 25.5cm x 20.5cm x 2.5cm with 3 Compartments 5 x Gauze Swabs 5cm x 5cm Plain 2 x Forcep Plastic Green 1 x Spinal Manometer 0-34 1 x Hypodermic Needle 25g x 5/8" 1 x Syringe 5ml Luer Lock 1 x IV 3000 Dressing 6cm x 8.5cm 1 x Specimen Jar 5ml 4 x Labelled Specimen Tube 10ml with 'V' Bottom 1 x Dressing Towel T Small 32cm x 32cm Wrapped in a Blue Plastic Field 75cm x 75cm	24

NEONATAL LUMBAR PUNCTURE PACK		
DEF4973	1 x Drape Clear Plastic 50cm x 60cm with 5cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray 25.5cm x 20.5cm x 2.5cm with 3 Compartments 5 x Gauze Swabs 7.5cm x 7.5cm Plain 10 x Cotton Balls Large 1 x Set of 3 Labelled Specimen Tubes 10ml 'V' Bottom enclosed in a Plastic Bag 2 x Gallipot 60ml 2 x Forcep Plastic Green 1 x Spinal Manometer 0-34 1 x Syringe 3ml Luer Lock 1 x Stopcock 3 way Low Pressure-Off 1 x Hypodermic Needle 23g x 1.25" 1 x Hypodermic Needle 23g x 5/8" 1 x Dressing Tegaderm 4.4cm x 4.4cm 1 x PurCotton Blue Towel 40cm x 68cm 1 x Needle Spinal Quinkce 22g x 38mm (added after sterilisation) in a clipseal Bag 4"x 7" Wrapped in a SMS Field 91cm x 91cm	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INTENSIVE CARE PACKS

TRACHY PROCEDURE PACK		
DEF2445	1 x Tracheal Dilator Trousseau 2 Prong 1 x Retractor Berman 16.5cm	1

ARTERIAL TRAY		
DEF2161A	1 x Drape Polybacked 33cm x 46cm with 8cm Square Fenestration 1 x Dressing Tray 25.5cm x 20.5cm x 2.5cm with 3 Compartments 6 x Gauze Swabs 7.5cm x 7.5cm Plain 2 x Blue Huck Towel 48cm x 60cm 1 x Scissors 13cm Sharp/Blunt Straight 1 x Scissors Iris 11.5cm Sharp/Sharp 1 x Needle Holder 13cm 2 x Forceps Artery Mosquito 12.5cm Straight 2 x Forceps Artery Mosquito 12.5cm Curved 1 x Forcep Iris Plain 10cm 1 x Forcep Iris Toothed 10cm 1 x Probe Lachrymal Size 1-2 1 x Probe Lachrymal Size 3-4 Wrapped in a Clear Plastic Field 75cm x 75cm	5

CVC INSERTION PACK		
DEF2393A	1 x Kidney Dish 700ml Clear 1 x Scissors 13cm Sharp/Blunt Straight 5 x Gauze Swabs 10cm x 10cm Plain 2 x Blue Huck Towel 48cm x 60cm 2 x Gallipot 150ml 5 x Non Woven Balls 1 x Needle Holder 13cm 1 x Forcep Adson Micro Toothed 2 x Forcep Spencer Wells 13cm Straight Wrapped in a Blue Plastic Field 91cm x 115cm	25

INTERCOSTAL PACK		
DEF1609	2 x Kidney Dish 700ml Clear 2 x Gallipot 150ml 6 x Cotton Balls Large 4 x Gauze Swabs 10cm x 10cm Plain 4 x Dressing Towel T Small 32cm x 32cm 1 x Needle Holder 16cm 2 x Forcep Spencer Wells 13cm Straight 2 x Forcep Spencer Wells 13cm Curved 1 x Scissors 12.5cm Sharp/Sharp Straight 1 x Scalpel Handle # 4 1 x Forcep Gillies Toothed 15.5cm Wrapped in a SMS Blue Field 91cm x 91cm	24

INSTRUMENTS AND
INSTRUMENT SETS

Our single use sterile instruments will equip you with quality tools for your surgical and medical procedures. The technical design of each instrument is modelled on reusable instruments, whilst providing you the benefit of no repeat sterilisation and a reduced risk of cross infection.

Each instrument has a green dipped handle to visually identify it as single use. Medical grade stainless steel instruments are made with an exacting mechanical finish. They can be recycled into meltdown metal or incinerated to become non-toxic file ash.

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FORCEPS

CLAMPS



DEF1781

DEF2298	13cm Smooth Jaw	10
DEF2304	13cm Smooth Jaw	100
DEF1781	18cm Smooth Jaw	1
DEF2275	18cm Smooth Jaw, Packet of 2	5
DEF2775	18cm Smooth Jaw	10
DEF1575	Tube Clamp 12cm Blue Plastic with Lock	50



DEF1575

DENTAL



DEF1450

DEF1450	Mirror	100
DEF3474	Mirror with Handle #4	50
DEF1457	1 x Mirror, 1 x Probe, 1 x Forcep	100



DEF1457

ADSON FORCEP



DEF2142

DEF2147	Plain 12.5cm	10
DEF2148	Plain 12.5cm	100
DEF2142	Toothed 12.5cm	10
DEF2143	Toothed 12.5cm	100

MICRO ADSON FORCEP



DEF2346A

DEF2135-10	Plain 12.5cm	10
DEF2135	Plain 12.5cm	100
DEF2346A	Toothed Super Fine 12.5cm	10
DEF2134-10	Toothed 12.5cm	10
DEF2134	Toothed 12.5cm	100



DEF2134-10

ROCHESTER PEAN FORCEP



DEF1782A

DEF1782A	22cm Straight	10
DEF3233	22cm Straight, Packet of 2	25
DEF2920	20cm Curved	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MOSQUITO ARTERY FORCEP



DEF2413-10

DEF2413-10	Straight 12.5cm	10
DEF2413	Straight 12.5cm	100
DEF2414-10	Curved 12.5cm	10



DEF2414-10

DEF2414	Curved 12.5cm	100
DEF2108	Curved 12.5cm Micro	10



DEF2108

SPENCER WELLS ARTERY FORCEP



DEF2277-10

DEF2277-10	Straight 13cm	10
DEF2277	Straight 13cm	100
DEF2415-10	Curved 13cm	10



DEF2415-10

DEF2415	Curved 13cm	100
DEF2417-10	Straight 16cm	10



DEF2417-10

DEF2417	Straight 16cm	100
DEF2416-10	Curved 16cm	10



DEF2416-10

DEF2416	Curved 16cm	100
---------	-------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CRILE ARTERY FORCEP




DEF2441

DEF2441	Straight 14cm	10
DEF1073	Straight 14cm	100
DEF2442	Curved 14cm	10



DEF2442

DRESSING FORCEP




DEF2152

DEF2152	Blunt 13cm	10
DEF2153	Blunt 13cm	100
DEF1072	Blunt 16cm	100

BASIC POINTED FORCEP

DEF1304	Dressing Forcep 13cm Pointed with Tip Protector	10
DEF1264A	Basic Pointed 13cm with Tip Protector	100



DEF1304


GILLIES FORCEP



DEF2418-10

DEF2418-10	Gillies Toothed 15.5cm	10
DEF2418	Gillies Toothed 15.5cm	100

MCINDOE FORCEP



DEF2677

DEF2677	McIndoe 15.5cm Non Toothed	10
DEF2948	McIndoe 15.5cm Non Toothed	100

MAGILL FORCEP



DEF1807

DEF1807	Magill 20cm Child	10
DEF2083	Magill 25cm Adult	10



DEF2083

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GREEN PLASTIC FORCEP



DEF686

DEF686	Green Plastic 11cm	100
--------	--------------------	-----

SPECIALISED FORCEPS

DEF1578	Forcep Wilde 15cm Angled	10
---------	--------------------------	----



DEF1578

DEF2578	Kelman-McPherson Forcep 11cm	10
---------	------------------------------	----



DEF2578

DEF2443	Lister Sinus Forcep 13cm	10
---------	--------------------------	----




DEF2443

DEF1579	Bozman Forcep Straight 25cm	5
---------	-----------------------------	---




DEF1579

DEF2174	Splinter Forcep Feichenfeldt 11cm	10
---------	-----------------------------------	----




DEF2174

DEF2750	Hartmann Croc-Action Forcep 14cm	10
---------	----------------------------------	----




DEF2750

DEF3141	Fine Croc-Action Forcep 15.5cm	10
---------	--------------------------------	----




DEF3141

DEF2695	Iris Forcep 10cm Plain	10
---------	------------------------	----




DEF2695

DEF3205	Iris Forcep Fine 10cm Curved	10
---------	------------------------------	----




DEF3205

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------




DEF2683

DEF2683	Tilley Ear Forcep 15cm	10
DEF2238	Packing Forcep 15cm	5
DEF1304	Dressing Forcep 13cm, Pointed	10
DEF2875	Tenaculum Forcep 25cm	10




DEF2875

DEF3550	Forcep Allis 4/5 Teeth 15.5cm Basic	20
---------	-------------------------------------	----




DEF3550

DEF2353	Moorfield Forcep 10.5cm	25
---------	-------------------------	----




DEF2353

NEEDLE HOLDERS




DEF2136

DEF2149	Needle Holder 13cm	10
DEF2136	Needle Holder 13cm	100
DEF2150	Needle Holder 16cm	10




DEF2150

DEF2151	Needle Holder 16cm	100
DEF2707	Nievert Needle Holder 13cm	10



DEF2707

RETRACTORS



DEF2402

CODE

DESCRIPTION

BOX QTY.


CODE

DESCRIPTION

BOX QTY.

SCALPEL


SCALPEL HANDLE



DEF2145


DEF2145	#3	10
DEF2137	#3	50
DEF2146	#3	100

SAFETY SCALPEL AND HANDLE



DEF2700

DEF2700	#10	100
DEF2701	#11	100
DEF2702	#15	100




DEF2702

DEF3935	#23	100
---------	-----	-----

SCALPEL PENBLADE

This pen-style scalpel is ergonomically designed with a high quality Swann Morton blade. The blade can be retracted without looking or repositioning the handle, whilst the safe – trimming groove can be used to trim sutures and PICC lines.




PB-M-10

PB-M-10	#10	10
PB-M-11P	#11P	10
PB-M-15	#15	10

SCISSORS

SHARP/BLUNT SCISSORS




DEF1871

DEF1871	Sharp/Blunt 13cm	10
DEF1870	Sharp/Blunt 13cm	100

NON STERILE

DEF1280	Sharp/Blunt 13cm Non Sterile	12
---------	------------------------------	----


SHARP/SHARP SCISSORS



DEF1311

DEF1311	Sharp/Sharp 12.5cm	10
DEF2276	Sharp/Sharp 12.5cm	100

BLUNT/BLUNT SCISSORS



DEF1637

DEF1637	Blunt/Blunt 12.5cm	100
---------	--------------------	-----



DEF1871

DEF1871	Sharp/Blunt 13cm	10
DEF1870	Sharp/Blunt 13cm	100

NON STERILE

DEF1280	Sharp/Blunt 13cm Non Sterile	12
---------	------------------------------	----



DEF1311

DEF1311	Sharp/Sharp 12.5cm	10
DEF2276	Sharp/Sharp 12.5cm	100



DEF1637

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

IRIS SCISSORS



DEF2141

DEF2144	Sharp/Sharp 7.5cm	10
DEF2141	Sharp/Sharp 11.5cm	10
DEF2133	Sharp/Sharp 11.5cm	100
DEF2460	Sharp/Sharp Curved 11.5cm	10

MAYO SCISSORS



DEF2408

DEF2408	Straight 17cm	10
DEF3361	Straight 17cm	50
DEF2117	Curved 17cm	10
DEF3370	Curved 17cm	50



DEF2117


METZENBAUM SCISSORS



DEF2511

DEF2511	Curved 16cm	10
---------	-------------	----

VANNAS SCISSORS



DEF2972

DEF2972	Straight Fine Tip	10
---------	-------------------	----

UNIVERSAL SCISSORS




DEF2851

DEF2851	Black Handle 18cm	25
---------	-------------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPECULUMS

EAR SPECULUM



DEF2240

DEF2240	Ear Speculum	5
---------	--------------	---

NOSE SPECULUM



DEF1356, DEF1357, DEF1358

DEF1356-10	Thudichum #1 20mm x 70mm	10
DEF1356	Thudichum #1 20mm x 70mm	25
DEF1357-10	Thudichum #2 22mm x 70mm	10
DEF1357	Thudichum #2 22mm x 70mm	25
DEF1358-10	Thudichum #3 25mm x 70mm	10
DEF1358	Thudichum #3 25mm x 70mm	25
DEF2292	Mod Wein Small	10
DEF2293	Medium	10
DEF2294	Large	10



DEF2292, DEF2293, DEF2294

VAGINAL SPECULUM



DEF2398

DEF2398	SIMS Small	20
DEF2399	SIMS Medium	20
DEF2400	SIMS Large	20

VAGINAL SPECULUM WITH LED LIGHT SOURCE

Single use disposable vaginal speculum with built-in LED light source from OBP Medical. Each comes ready for use – no assembly required and no extra parts needed. High powered, battery operated LED light source that produces peak light for 30 plus minutes with no heat generated. Batteries used are landfill friendly and require no special disposal instructions.

STERILE SPECULUMS WITH LIGHT SOURCE

C020300	Small White	12
C020310	Medium Green	12
C020320	Large Blue	12

NON STERILE SPECULUMS WITH LIGHT SOURCE

ER-SPEC-SMALL	Small White	18
ER-SPEC-MEDIUM	Medium Green	18
ER-SPEC-LARGE	Large Blue	16

STERILE SPECULUMS WITHOUT LIGHT SOURCE

DEF3157	Small	20
DEF3158	Medium	20
DEF3159	Large	20

NON STERILE SPECULUMS WITHOUT LIGHT SOURCE

DEF3157NS	Small	20
DEF3158NS	Medium	20
DEF3159NS	Large	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SUCTION INSTRUMENTS

ANAL SPECULUM NON STERILE

ANOSPEC WITH LED LIGHT SOURCE

Single use disposable Anoscope with built-in LED light source from OBP Medical. Designed to aid in the examination of the anus to screen, diagnose and evaluate problems of the anal canal. This ready-to-use device eliminates the time and expense of reprocessing. Each Anospec comes complete with its own obturator.

ANO-SPEC-SMALL	90mm x 18mm – Slotted	20
ANO-SPEC-MEDIUM	96mm x 23mm – Slotted	16
ANO-SPEC-LARGE	103mm x 18mm – Bevelled	20

SPONGE HOLDER

DEF2609	Rampley 24cm Stainless Steel	10
---------	------------------------------	----



DEF2609

DEF1923	Plastic Blue 19cm	10
---------	-------------------	----



DEF1923

DEF2071	Plastic Blue 24cm	10
---------	-------------------	----



DEF2071

STAPLE REMOVER



DEF1826

DEF1826	Staple Remover Plastic Handle	50
---------	-------------------------------	----

NON STERILE

DEF1825	Staple Remover Plastic Handle	20
---------	-------------------------------	----

YANKAUER SUCTION INSTRUMENT

Popular suction instrument with a large lumen providing rapid aspiration whilst resisting clogging. Removes fluid and debris from a surgical site or body orifice.



DEF2329, DEF2925, DEF3030

DEF2329	With Tip 30cm	10
DEF2925	No Tip 30cm	10
DEF3030	Straight No Tip 30cm	10

PAEDIATRIC SUCTION INSTRUMENT



DEF2498, DEF2497

DEF2497	Small 27cm	10
DEF2498	Medium 29cm	10

POOLE SUCTION INSTRUMENT

Designed to evacuate pooled blood, fluid and debris during deep abdominal surgeries and caesarean sections.

DEF2554	25cm	10
---------	------	----



DEF2554

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FRAZIER SUCTION INSTRUMENT

Used in open surgical procedures where fine, delicate suctioning is necessary to remove debris from the surgical site. Comes with an obturator to clear the instrument of obstructions during use.

FRAZIER SUCTION INSTRUMENT STANDARD



DEF2325, DEF2326, DEF2327, DEF2328, DEF3477

DEF2325	6 Fr Gauge 17cm long	10
DEF2326	8 Fr Gauge 17cm long	10
DEF2327	10 Fr Gauge 17cm long	10
DEF2328	12 Fr Gauge 17cm long	10
DEF3477	14 Fr Gauge 17cm long	10

FRAZIER SUCTION INSTRUMENT LONG



DEF3482, DEF3483, DEF3484

DEF3562	6 Fr Gauge 20.5cm long	50
DEF3482	8 Fr Gauge 20.5cm long	50
DEF3483	10 Fr Gauge 20.5cm long	50
DEF3484	12 Fr Gauge 20.5cm long	50
DEF3485	14 Fr Gauge 20.5cm long	50

FRAZIER SUCTION INSTRUMENT EXTRA LONG

The suction tube connector opening allows for a secure fit of a 20ml Luer Slip syringe for flushing of the suction tip.



DEF3080A, DEF3081A, DEF3082A, DEF3486

DEF3080A	8 Fr Gauge 23cm long	50
DEF3081A	10 Fr Gauge 23cm long	50
DEF3082A	12 Fr Gauge 23cm long	50
DEF3486	14 Fr Gauge 23cm long	50

FUJITA SUCTION INSTRUMENT SHORT LENGTH



DEF3500, DEF3502, DFE3501

DEF3500	4 Fr Gauge 18cm long	50
DEF3501	5 Fr Gauge 18cm long	50
DEF3502	7 Fr Gauge 18cm long	50

RHOTON SUCTION INSTRUMENT



DEF3787, DEF3786, DEF3785, DEF2807


DEF3787	5 Fr Gauge 10cm long	10
DEF3786	7 Fr Gauge 10cm long	10
DEF3785	10 Fr Gauge 10cm long	10
DEF2807	10 Fr Gauge 13cm long	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INSTRUMENT PACKS


TRACHEAL DILATOR

DEF2373	3 Prong	1
---------	---------	---



DEF2373

DEF1467	Trousseau 2 Prong	1
---------	-------------------	---



DEF1467

DEF2585	Trousseau 2 Prong	5
---------	-------------------	---

WIRE CUTTER



DEF1800

DEF1800	Universal 12cm	50
DEF2003	Universal 12cm	10

PLASTIC SET

DEF1328A	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep McIndoe 15.5cm	10
----------	--	----

SUTURE KIT


DEF1811	1 x Safety Scalpel # 11 with Handle 1 x Needle Holder 16cm	25
---------	---	----

FINE SUTURE SET

DEF1656	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain	25
---------	--	----

SUTURE SET SMALL

DEF1560	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 13cm	10
---------	--	----



DEF1560

SUTURE SET LARGE

DEF1561	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 16cm	10
---------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INSTRUMENT PACK

DEF1913B	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Forcep Adson 12.5cm Plain 1 x Needle Holder 13cm	10
----------	---	----

FINE FACIAL MICRO SUTURE SET

DEF1636-10	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm	10
DEF1636-25	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm	25
DEF1636	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm	100

SUTURE KIT

DEF2270C	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel White 2 Ply 40cm x 42cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Green Plastic 10 x Gauze Swabs X-ray 10cm x 10cm Wrapped in a Blue Plastic Field 50cm x 50cm	10
----------	--	----

FINE PLASTIC SUTURE TRAY

DEF2051	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Scalpel Handle # 3 2 x Non Woven Swabs 7.5cm x 7.5cm 5 x Gauze Sponges Plain 10cm x 10cm Wrapped in a Blue Plastic Field 50cm x 50cm	10
---------	---	----

SUTURE KIT

DEF2179	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 2 x Dressing Towel White 2 Ply 40cm x 42cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Scissor Sharp/Sharp Straight 12.5cm 1 x Needle Holder 16cm 1 x Forcep Gillies Toothed 15.5cm 1 x Forcep Green Plastic 5 x Gauze Sponges Plain 10cm x 10cm Wrapped in a Blue Plastic Field 75cm x 75cm	50
---------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
FINE PLASTIC SUTURE TRAY					
DEF2507B	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm 1 x Scalpel Handle # 3 5 x Gauze Sponges Plain 10cm x 10cm 2 x Non Woven Swabs 7.5cm x 7.5cm Wrapped in a Blue Plastic Field 75cm x 75cm	10			
SUTURE TRAY					
DEF2050A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Scissor 12.5cm Sharp/Blunt 1 x Forcep Adson 12.5cm Plain 1 x Needle Holder 13cm 5 x Gauze Sponges Plain 10cm x 10cm 2 x Non Woven Swabs 7.5cm x 7.5cm Wrapped in a Blue Plastic Field 50cm x 50cm	10			
FINE PLASTIC SUTURE TRAY WITH TOWEL					
DEF3468	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm 5 x Gauze Sponges Plain 10cm x 10cm Placed into a plastic bag 20cm x 20cm	40			
FINE SUTURE SET					
DEF2226	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel T Small 32cm x 32cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Forcep Green Plastic 10 x Gauze Sponges Plain 10cm x 10cm 1 x Syringe 5ml Luer Slip 1 x Hypodermic Needle 25g x 1" Wrapped in a Blue Plastic Field 75cm x 75cm	10			
SUTURE KIT					
DEF2835A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Scissor 12.5cm Sharp/Blunt 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Green Plastic 1 x Gauze Swab Plain 10cm x 10cm Wrapped in a Blue Plastic Field 75cm x 75cm	10			
SUTURE SET					
DEF2422	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel T Small 32cm x 32cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 16cm 1 x Forcep Green Plastic 5 x Non Woven Swabs 7.5cm x 7.5cm 5 x Non Woven Balls Wrapped in a Blue Plastic Field 50cm x 50cm	60			







CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
SUTURE KIT WITH COTTON BALLS					
DEF1671A	1 x Dressing Tray with 3 compartments 23cm x 9cm x 3cm deep 1 x Dressing Towel White 2 Ply 40cm x 42cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 16cm 1 x Forcep Blunt 13cm 1 x Forcep Green Plastic 1 x Scalpel Handle # 3 3 x Gauze Sponges Plain 10cm x 10cm 5 x Cotton Balls Large Wrapped in a Blue Plastic Field 50cm x 50cm	10			
INSTRUMENT DRESSING PACK					
DEF1677A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Tray Long 23cm x 9cm x 3cm 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissors Sharp/Blunt 12.5cm 1 x Forcep Blunt 13cm 1 x Forcep Pointed 13cm 1 x Forcep Green Plastic 5 x Gauze Sponges Plain 7.5cm x 7.5cm 1 x Dressing Towel T Small 32cm x 32cm Wrapped in a Blue Plastic Field 75cm x 75cm	60			
PAEDIATRIC SUTURE SET					
DEF2955	1 x Scissors Sharp/Blunt 12.5cm 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm	25			
MEDIUM DRESSING TRAY					
DEF2502A	1 x Dressing Tray Long 23cm x 9cm x 3cm 1 x Towel Blue Cotton 40cm x 68cm 14 x Gauze Sponges Plain 7.5cm x 7.5cm 1 x Kidney Dish 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	30			
PHARMACY INSTRUMENT PACK					
DEF2486	1 x Scissors Sharp/Blunt 12.5cm 1 x Forcep Blunt 13cm 1 x Forcep Crile 14cm 1 x Kidney Dish 700ml Clear Wrapped in a Blue Plastic Field 75cm x 75cm	10			

MATERNITY PRODUCTS



FREECALL 1300 550 278

MATERNITY PRODUCTS

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
AMNIOTIC HOOK					
An amniotic hook is a thin, long piece of plastic with a hook on one end, used to break the amniotic sac to help initiate labour.					
					
AMNIOTIC-HOOK					
AMNIOTIC-HOOK	Amniotic Hook 26cm	100			
BABY TEATS					
					
DEF1059					
DEF1059	Teat Cherry Slow Flow	100			
BABY CUP WITH LID					
DEF2129	Baby Cup with Lid	100			
					
DEF2129					
BABY RECEIVING BLANKET					
Soft, cotton flannel blanket used for wrapping a baby after birth.					
					
DEF3227					
DEF3227	Flannel 75cm x 100cm	40			
DEF3621	Flannel 75cm x 100cm, Packet of 2	40			
MATERNITY PAD					
					
DEF721					
DEF721	10cm x 30cm	50			
NON STERILE					
DEF700	Bulk 10cm x 30cm	300			
DEF238	10cm x 30cm, Packet of 10 in a Clipseal Bag	25			
NIPPLE SHIELD					
					
DEF410					
DEF410	Nipple Shield Gentle Feed, Packet of 2, Non Sterile	1			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MATERNITY PACKS

PLACENTA BASIN

DEF2822	26cm x 18.5cm x 7cm	12
---------	---------------------	----



DEF2822

UMBILICAL CORD CUTTING DEVICE

DEF2535	Umbilical Cord Cutting Scissor	10
---------	--------------------------------	----



DEF2535

DEF2168	Umbilical Cord Clamp Cutter Non Sterile, with Plastic Handle	10
---------	---	----



DEF2168

MATERNITY SUTURE SET

DEF2131A	1 x Dressing Tray Long 23cm x 9cm x 3cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Needle Holder 18cm 1 x Forcep Gillies Toothed 15.5cm 1 x Forcep Mosquito Artery 12.5cm Straight 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	10
----------	---	----

DISPOSABLE DELIVERY SUTURE SET

DEF2118	1 x Scissors Mayo 17cm Straight 1 x Needle Holder 16cm 1 x Forcep Bonney 18cm ½ tooth 2 x Forcep Crile Curved 14cm	25
---------	---	----

EPISIOTOMY SUTURE SET

DEF2671A	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Needle Holder 16cm 1 x Forcep Gillies Toothed 15.5cm Wrapped in a Blue SMS field 91cm x 91cm	10
----------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DELIVERY SET

DEF2752	1 x Blue SMS Drape 76cm x 106cm 1 x Maternity Pad 10cm x 30cm 1 x Baby Receiver Flannel Blanket 70cm x 100cm 20 x Non Woven Swabs 10cm x 10cm 3 x Towels White 2 Ply 40cm x 42cm 1 x Kidney Dish Clear 700ml 1 x Bowl Blue 500ml 2 x Forcep Rochester Pean 22cm Straight 1 x Mayo Scissor 17cm Curved 2 x Umbilical Cord Clamps Wrapped in a Blue SMS Field 91cm x 91cm	15
---------	---	----



DEF2752

LABOUR PACK

DEF2526	2 x Forcep Rochester Pean 22cm Straight 1 x Mayo Scissor 23cm Straight 1 x Bowl 600ml Green 1 x Bowl 2.5L 1 x Gallipot 150ml Wrapped in a Blue SMS Field 91cm x 91cm	10
---------	---	----

BIRTHROOM BUNDLE

DEF2670A	1 x Placental Basin 26cm x 18.5cm x 7cm deep 1 x Blue bowl 500ml 1 x Scissor Sharp/Blunt 12.5cm 1 x Forcep Rochester Pean 22cm Straight Wrapped in a Blue SMS Field 91cm x 91cm	10
----------	--	----

IUD INSERTION KIT

This procedure Pack is for an IUD insertion procedure.

The kit is sterile and disposable, ensuring all components are available for every insertion.

Reduced Packaging waste, plus the benefit of no repeat sterilisation costs.

DEF3256	2 x Swab Sticks 20cm Paper Shaft 2 x Gallipots 150ml 1 x Kidney Dish 1000ml Blue 1 x Vaginal Speculum Large 1 x Scissor Mayo 23cm Straight 1 x Forcep Tenaculum 25cm 1 x Sponge Holder 19cm Plastic 1 x Uterine Sound (Hysterometer CH14) 4.67mm with Nozzle Wrapped in a Blue Plastic Field 75cm x 75cm	10
---------	--	----



DEF3256

MEDICAL IMAGING AND CATH LAB PRODUCTS

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANGIO PROCEDURAL DRAPE

ANGIOGRAPHY PROCEDURE DRAPE WITH BENDABLE TABS		
DEF2847	Size 205cm x 346cm, with Supersorb Fabric 82cm x 306cm Dual 15cm Round Fenestrations with Cut-Out Incise 5 x Bendable Formers Clear Side and Bottom Panels Wrapped in a SMS 91cm x 91cm Field	10

FEMORAL RADIAL ANGIOGRAPHY DRAPE		
DEF5187	Size 224cm x 342cm, with Sapphire Fabric and Clear Side Panels Dual 20cm Round Fenestration with Built-in Cut-Out Incise Dual 18cm Round Fenestration with Built-in Cut-Out Incise Wrapped in a SMS 91cm x 91cm Field	10

EQUIPMENT COVERS

IMAGE INTENSIFIER COVER		
		
DEF764		
DEF764	Cap Small 61cm x 91cm, Double Packed	40
DEF770	Cover Small 76cm x 106cm Double Packed	50
DEF767	Cover Medium 91cm x 115cm	50
DEF1393	Cover 104cm x 150cm, Non Sterile	50
DEF774	Cover Large 120cm x 140cm	40
DEF1045	Cover 105cm x 150cm with Sealed End, 1 Securing Tape, Cuff and Labels, Double Packed	15
DEF1348	Cover 105cm x 110cm with Sealed End, 1 Securing Tape, Cuff and Labels, Double Packed	15

PROBE COVER		
		
DEF301		
DEF301	Probe Cover 10cm x 100cm with a Sealed End	100
DEF648	Drawstring Probe Cover 23cm x 100cm	25
DEF1482	Probe Cover 23cm x 100cm with a Sealed End	100
DEF2630	Drawstring Probe Cover 23cm x 150cm with a Sealed End	25
DEF1088	Ultrasound Probe Cover 23cm x 80cm with a Sealed End	100
DEF643	Ultrasound Probe Cover 23cm x 150cm with a Sealed End	25
DEF2409	Ultrasound Probe Cover 25cm x 150cm with a Sealed End and a Rubber Band #64	100

LEAD SCREEN COVER		
		
DEF1245		
DEF1245	80cm x 80cm with Elastic Opening	50
DEF3281	75cm x 90cm with a Sealed End	50
DEF2073	Fluoroscope/Lead Screen cover 76cm x 91cm with a Sealed End	50
DEF1762	104cm x 105cm with a Sealed End	50
DEF805	90cm x 105cm with a Sealed End	50
DEF1161	80cm x 110cm with a Sealed End	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GUIDEWIRE BOWL

Patented design from DeRoyal with built up lid and recessed bottom. Will hold guidewire in place and allows for easy removal.

DEF2579	Guidewire Bowl	10
---------	----------------	----



DEF2579

SPLASH STOP

An innovative DeRoyal safe alternative to messy, open basins.

A closed bowl with a secure, non-splatter lid which minimizes staff and patient exposure to blood-borne pathogens.

Contains a safe, one-step biohazard warning label.



77-301371

77-301370	Waste Collection Bowl	20
77-301371	Waste Collection Bowl with Solidifier	20

MEDICAL IMAGING & CATH LAB PACKS

Please contact Defries sales representative to discuss your custom pack requirements.

X-RAY INSTRUMENT PACK

DEF2485	1 x Scissors Sharp/Blunt 12.5cm 1 x Forcep Gillies Toothed 15.5cm 1 x Forcep Crile Straight 14cm 1 x Needle Holder 16cm 1 x Forcep Artery Mosquito 12.5cm Curved 1 x Forcep McIndoe 15.5cm Plain	10
---------	---	----


ANGIO INSTRUMENT SET

DEF2768	1 x Scissors Sharp/Sharp 12.5cm 1 x Forcep Artery Mosquito 12.5cm Straight 1 x Forcep Artery Mosquito 12.5cm Curved 1 x Sponge Holder 24cm Plastic Blue	10
---------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

PACEMAKER PACK

DEF3586	1 x Dressing Tray 3 Cavity 255mm x 205mm x 25mm 1 x Angio Drape 346cm x 205cm with Femoral Fenestration 1 x Clipseal Bag 8" x 10" (With small components enclosed) 2 x Absorbent Foam Applicator 15cm Blue 1 x Forcep Adson Plain with Green Handle 1 x Dressing Forcep 16cm Green Handle 2 x Forcep Crile Straight 14cm with Green Handle 1 x Gillies Forcep 15.5cm Toothed Green Handle 4 x Forcep Artery Mosquito 12.5cm Curved Green Handle 3 x Forcep Spencer Wells Curved Green Handle 16cm 1 x Forcep Adson Toothed with Green Handle 1 x Senn Muller Retractor 16cm Blunt 1 x Weitlaner Retractor 18cm Blunt 4/3 Teeth 1 x Iris Scissor 11.5cm Sharp/Sharp with Green Handles 1 x Mayo Scissor 17cm Curved Green Handle 1 x Needle Holder 16cm with Green Handle 5 x Gauze Swabs 10cm x 10cm X-ray Double Packed 1 x Bowl 1 Litre Blue 2 x Gallipot 150ml 1 x Kidney Dish 700ml Clear Wrapped in a Back Table Cover 60" x 90"	4
---------	---	---



DEF3586

MEDICAL LABELS

Please contact Customer Service on 1300 550 278 to discuss your label requirements. Our durable pre-printed Sterile and Non Sterile labels are used throughout the hospital in different wards such as the Emergency Department, General Wards, Sterilising Dept, Operating Theatre and Intensive Care.

Our labels are designed to clearly and easily identify items in a hospital on a professional, easy to read label.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GENERAL USE NON STERILE LABELS

WHITE BACKGROUND

These labels are available with a white background and a coloured print in either blue, black, red, magenta, yellow and green. Rolls of either 100 labels or 1000 labels.

ANAESTHETIC SIZE

- 15mm x 300mm
- 21mm x 61mm
- 32mm x 50mm

STANDARD SIZE

- 32mm x 70mm
- 50mm x 70mm
- 40mm x 102mm
- 40mm x 170mm
- 80 x 102mm – provided as 2 rolls of 500

FULL COLOUR BACKGROUND

Any colour combination can be made either for the background or print for this type of label. Rolls of either 100 labels or 1000 labels.

ANAESTHETIC SIZE

- 15mm x 300mm
- 21mm x 61mm
- 32mm x 50mm

STANDARD SIZE

- 32mm x 70mm
- 50mm x 70mm
- 40mm x 102mm
- 40mm x 170mm
- 80 x 102mm – provided as 2 rolls of 500

FLUORESCENT LABELS

These labels are fluorescent in colour with black print. Fluorescent colours available are Orange, Yellow, Red, Green, Blue or Pink.

- 32mm x 500mm – ONLY Fluorescent Red
- 30mm x 70mm
- 50mm x 70mm – ONLY Fluorescent Green
- 80mm x 102mm – PLUS Fluorescent Purple
- ME9 Tray Incomplete Label – Orange colour 400

PIGGY BACK LABELS

These labels are made up of two labels where one labels sits on top of (or Piggy backs) the other label. So the top label can be removed and adhered. These versatile labels can be used through out the hospital.

PBLABELS – B (BLANK) IN ROLLS OF 100 OR 1000

- 32mm x 50mm
- 32mm x 50mm
- 36mm x 105mm

PBLABELS – P (PRINTED) IN ROLLS OF 100 OR 1000

- 32mm x 50mm
- 32mm x 50mm
- 36mm x 105mm

PBLABELS – S (STERRAD) IN ROLLS OF 100 OR 1000

- 32mm x 50mm
- 32mm x 50mm
- 36mm x 105mm

STERILE ANAESTHESIA DRUG LABELS

Sterile anaesthetic drug labels are used in the operating theatre to label syringes during procedures to help prevent medication mix ups. These easy to use labels come pre-printed and pre-packed in a peel down sterile packet. The sizes of these labels are 15mm x 30mm Anaesthetic size with a white background. The colour of print is either Blue, Black, Red, Magenta, Yellow or Green. Available in boxes of 100.

- Prepacked in 1's
- Prepacked in 2's
- Prepacked in 3's
- Prepacked in 6's
- Prepacked in 12's
- Prepacked in 18's
- Prepacked in 24's

The sizes of these labels are 15mm x 30mm Anaesthetic size with a full colour background. These labels can have any colour combination for either the background or print. Available in boxes of 100.

- Prepacked in 1's
- Prepacked in 2's
- Prepacked in 3's
- Prepacked in 6's
- Prepacked in 12's
- Prepacked in 18's
- Prepacked in 24's

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

STERILE PERI-OPERATIVE LABELS

These labels sheets are customisable according to the National Recommendations for user-applied labelling. Please contact our Customer Service on 1300 550 278 for further information.

Labels sheets are printed, folded and packaged in a peel down pouch and then sterilised.

Labels sheets are waterproof and can be written on by surgical marking pens such as the DEF3219 Surgical Marker Permanent Fine Tip.

LAB9009 and LAB9010 are currently used in hospitals and are available to order from Defries Industries. Alternatively, you can design your own label sheets for your department requirements.

Economical price and easy to use and read label sheets, ready for every procedure.

Sterile and latex free.



LAB9009A

LAB9009	Sterile Medical Label Sheet – specialties include General, Ortho, Neuro, Cardiac, Vascular, Urology, Obstetrics & Gynaecology, Breast and Endocrine	100
LAB9010	Sterile Medical Label Sheet – specialties include ENT, Max/Facial, Plastic, Ophthalmology and Neuro Surgery	100



PERSPEX DISPENSER

Our Anaesthesia roll labels are used in operating rooms and day surgery units. Easy to use and apply drug labels that conform and are compliant to the National Recommendations for user-applied labelling. These labels have a large font and with an easy to read print. Can be used with our perspex dispenser boxes.

PD10	Perspex dispenser to hold 10 rolls	1
PD15	Perspex dispenser to hold 15 rolls	1
PD20	Perspex dispenser to hold 20 rolls	1
PD15-2T	Perspex dispenser to hold 30 rolls	1

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHESIA LABEL ROLL – NON STERILE

ANAESTHESIA LABEL ROLL

AL/AA-1	AN-LABELS-ATROPINE MG/ML (GREEN)	500
AL/AA-2	AN-LABELS-GLYCOPYROLATE (GREEN)	500
AL/AA-3	AN-LABELS-BUSCOPAN MG/ML (GREEN)	500
AL/AAB	AN-LABELS – BLANK (NO TEXT) (GREEN)	500
AL/AE	AN-LABELS-ANTI-EMETIC (SALMON)	500
AL/AE-2	AN-LABELS-METOCLOPRAMIDE MG/ML (SALMON)	500
AL/AE-3	AN-LABELS-STEMETIL MG/ML (SALMON)	500
AL/AE-5	AN-LABELS-ONDANSETRON HCL (SALMON)	500
AL/AE-7	AN-LABELS-TROPISETRON (SALMON)	500
AL/AE-8	AN-LABELS-GRANISETRON	500
AL/AI-4	AN-LABELS-MAXALON (SALMON)	500
AL/BW	AN-LABELS – BLANK (NO TEXT) (WHITE)	500
AL/HA-1	AN-LABELS-NITROGLYCERINE MG/ML (VIOLET/WHITE)	500
AL/HA-10	AN-LABELS-NITROPRUSSIDE	500
AL/HA-2	AN-LABELS-NITROPRUSSIDE MG/ML (VIOLET/WHITE)	500
AL/HA-3	AN-LABELS-TRIMEHAPHAN MG/ML (VIOLET/WHITE)	500
AL/HA-4	AN-LABELS-PENTOLAMINE MG/ML (VIOLET/WHITE)	500
AL/HA-5	AN-LABELS-HYDRALLAZINE (VIOLET/WHITE)	500
AL/HA-6	AN-LABELS-CLONIDINE (VIOLET/WHITE)	500
AL/HA-7	AN-LABELS-GTN (VIOLET/WHITE)	500
AL/HA-8	AN-LABELS-VERAPAMIL (VIOLET/WHITE)	500
AL/HA-9	AN-LABELS-PHENOXYBENZAMINE HCL (VIOLET/WHITE)	500
AL/HAB	AN-LABELS-BLANK (NO TEXT) (VIOLET&WHITE)	500

AL/IA-1	AN-LABELS-THIOPENTONE MG/ML (YELLOW)	500
AL/IA-2	AN-LABELS-PROPOFOL MG/ML (YELLOW)	500
AL/IA-3	AN-LABELS-METHOHEXITONE MG/ML (YELLOW)	500
AL/IA-4	AN-LABELS-KETAMINE MG/ML (YELLOW)	500
AL/IA-6	AN-LABELS-TRAMADOL (YELLOW)	500
AL/IAB	AN-LABELS – BLANK (NO TEXT) (YELLOW)	500
AL/LA	AN-LABELS-LOCAL (GREY)	500
AL/LA-09	AN-LABELS-XYLOCAINE MG/ML	500
AL/LA-1	AN-LABELS-LIGNOCAINE MG/ML (GREY)	500
AL/LA-1B	AN-LABELS-LIGNOCAINE 1% MG/ML (GREY)	500
AL/LA-1C	AN-LABELS-LIGNOCAINE 2% MG/ML (GREY)	500
AL/LA-2	AN-LABELS-PROCAINE (GREY)	500
AL/LA-3	AN-LABELS-BUPIVACAINE MG/ML (GREY)	500
AL/LA-3B	AN-LABELS-BUPIVACAINE 0.25% MG/ML (GREY)	500
AL/LA-3C	AN-LABELS-BUPIVACAINE 0.5% MG/ML (GREY)	500
AL/LA-4	AN-LABELS-ROPIVACAINE MG/ML (GREY)	500
AL/LA-5	AN-LABELS-LOCAL MEPIVACAINE (GREY)	500
AL/LAB	AN-LABELS – BLANK (NO TEXT) (GREY)	500
AL/MES-1	AN-LABELS-RESTRICTED DRUG (RED)	500
AL/MES-11	AN-LABELS-ALERT (BLUE)	500
AL/MES-18	AN-LABELS-CHANGE ON (BLACK ON WHITE)	500
AL/MES-19	AN-LABELS-MAINTENANCE (BLACK ON GREEN)	500
AL/MES-2	AN-LABELS-EPIDURAL (YELLOW)	500
AL/MES-21	AN-LABELS-ALLERGY (RED)	500
AL/MES-3	AN-LABELS-PERIPHERAL LINE/DATE (ORANGE)	500
AL/MES-30	AN-LABELS-VOLUME	500

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
ANAESTHESIA LABEL ROLL (CONTINUED)					
AL/MES-33	AN-LABELS-PCA (LIGHT BLUE)	500	AL/MIS-27	AN-LABELS-GENTAMYCIN (BLACK ON WHITE)	500
AL/MES-34	AN-LABELS-CVP (BLUE)	500	AL/MIS-28	AN-LABELS-WATER FOR INJECT. (BLACK/WHITE)	500
AL/MES-35	AN-LABELS-REGIONAL INFUSION (VIOLET)	500	AL/MIS-3	AN-LABELS-0.9% SALINE (BLACK ON WHITE)	500
AL/MES-4	AN-LABELS-ATERIAL (RED)	500	AL/MIS-3/GW	AN-LABELS-0.9% SALINE (GREEN ON WHITE)	500
AL/MES-5	AN-LABELS-C.V.C DATE INIT (BLUE)	500	AL/MIS-30	AN-LABELS-THROMBIN (BLACK ON WHITE)	500
AL/MES-54	AN-LABELS-CONTRAST (BLACK ON WHITE)	500	AL/MIS-31	AN-LABELS-FRUSEMIDE (BLACK ON WHITE)	500
AL/MES-57	AN-LABELS-MEDICATION ADDED	500	AL/MIS-31B	AN-LABELS-LASIX 50MG/45ML NORMAL SALINE	500
AL/MES-58	AN-LABELS-IV FLUSH (BLACK ON WHITE)	500	AL/MIS-32	AN-LABELS-POT.CHLORID (BLACK ON WHITE)	500
AL/MES-7	AN-LABELS-ORAL MEDICATION (RED & WHITE STRIPE)	500	AL/MIS-33	AN-LABELS-ERGOMETRINE (BLACK ON WHITE)	500
AL/MES-83	CYTOTOXIC LABEL (BLACK ON VIOLET)	500	AL/MIS-34	AN-LABELS-SYNTOMERINE MG/ML (BLACK/WHITE)	500
AL/MES-9	AN-LABELS-PULMONARY ARTERY (YELLOW)	500	AL/MIS-36	AN-LABELS-MILRINONE (BLACK ON WHITE)	500
AL/MES36	AN-LABELS-B.S.S (BLACK ON WHITE)	500	AL/MIS-37	AN-LABELS-SODIUM BIOG. (BLACK ON WHITE)	500
AL/MIS-1	AN-LABELS-HEPARIN (BLACK ON WHITE)	500	AL/MIS-39	AN-LABELS-METHADONE (BLACK/WHITE)	500
AL/MIS-10	AN-LABELS-DEXAMETHASONE MG/ML (BLACK/WH)	500	AL/MIS-4	AN-LABELS-ANTIBIOTIC (BLACK ON WHITE)	500
AL/MIS-106	AN-LABELS-PARACETAMOL (RED ON WHITE)	500	AL/MIS-43	AN-LABELS-PROMETHAZINE (BLACK ON WHITE)	500
AL/MIS-11	AN-LABELS-INSULIN (BLACK ON WHITE)	500	AL/MIS-47	AN-LABELS-ABCIXIMAB (BLACK ON WHITE)	500
AL/MIS-12	AN-LABELS-AMIODARONE (BLACK ON WHITE)	500	AL/MIS-5	AN-LABELS-STEROID (BLACK ON WHITE)	500
AL/MIS-14	AN-LABELS-MIOCHOL (BLACK ON WHITE)	500	AL/MIS-50	AN-LABELS-KETOROLAC (BLACK ON WHITE)	500
AL/MIS-19	AN-LABELS-CALCIUM CHLORIDE (BLACK/WHITE)	500	AL/MIS-59	AN-LABELS-PARECOXIB SODIUM (BLACK/WHITE)	(DYNASTAT) 500
AL/MIS-1A	AN-LABELS-HEPARIN (RED ON WHITE)	500	AL/MIS-6	AN-LABELS-OCXYTOCIN UNITS/ML (BLACK/WHITE)	500
AL/MIS-2	AN-LABELS-PROTAMINE (BLACK ON WHITE)	500	AL/MIS-60	AN-LABELS-DIGOXIN (BLACK ON WHITE)	500
AL/MIS-20	AN-LABELS-ADENOSINE (BLACK ON WHITE)	500	AL/MIS-62	AN-LABELS-TRAMADOL (BLACK ON WHITE)	500
AL/MIS-21	AN-LABELS-MORPHINE/ MIDAZOLAM (BLACK/WHITE)	500	AL/MIS-67	AN-LABELS-HYDROCORTISONE (WHITE)	500
AL/MIS-22	AN-LABELS-AMINOPHYLLINE (BLACK ON WHITE)	500	AL/MIS-7	AN-LABELS-CALCIUM GLUCONATE (BLACK/WHITE)	500

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
AL/MIS-70	AN-LABELS-POTASSIUM (BLACK ON WHITE)	500	AL/N-3	AN-LABELS-PETHIDINE MG/ML (BLUE)	500
AL/MIS-70/RW	AN-LABELS-POTASSIUM (RED ON WHITE)	500	AL/N-4	AN-LABELS-ALFENTANIL MICROG/ML (BLUE)	500
AL/MIS-72	AN-LABELS-MAGNESIUM SULFATE (BLCK/WHITE)	500	AL/N-6	AN-LABELS-REMIFENTANYL MICROG/ML (BLUE)	500
AL/MIS-73	AN-LABELS-METOPROLOL (BLACK ON WHITE)	500	AL/NA-01	AN-LABELS-NALOXONE (BLUE/WH STRIPE)	500
AL/MIS-86	AN-LABELS-TPN (WHITE LABEL)	500	AL/NAB	AN-LABELS – BLANK (NO TEXT) (BLUE&WHITE)	500
AL/MIS-9	AN-LABELS-HEPARINIZED SALINE (BLACK/WHITE)	500	AL/NB	AN-LABELS – BLANK (NO TEXT) (BLUE)	500
AL/MIS-90	AN-LABELS-SALINE 3%	500	AL/PATH01	PATHOLOGY LABELS-TERMINAL ILEUM (BLACK/WHITE)	500
AL/MIS-91	AN-LABELS-DEXTROSE	500	AL/PATH03	PATHOLOGY LABELS-CAECUM (BLACK/WHITE)	500
AL/MIS-93	AN-LABELS-PROSTIN	500	AL/PATH04	PATHOLOGY LABELS- TRANSVERSE (BLACK/WHITE)	500
AL/MR	AN-LABELS-MUSCLE RELAX (RED)	500	AL/PATH07	PATHOLOGY LABELS-L. COLON (BLACK/WHITE)	500
AL/MR-1	AN-LABELS-SUXAMETHONIUM MG/ML (RED)	500	AL/PATH08	PATHOLOGY LABELS-R. COLON (BLACK/WHITE)	500
AL/MR-2	AN-LABELS-PANCURONIUM MG/ML (RED)	500	AL/PATH09	PATHOLOGY LABELS- ASCENDING COLON (BLACK/WHITE)	500
AL/MR-3	AN-LABELS-ATRACURIUM MG/ML (RED)	500	AL/PATH10	PATHOLOGY LABELS- DESCENDING COLON (BLACK/WHITE)	500
AL/MR-4	AN-LABELS-VECURONIUM MG/ML (RED)	500	AL/PATH11	PATHOLOGY LABELS-T.I	500
AL/MR-5	AN-LABELS-ROCURONIUM MG/ML (RED)	500	AL/PATH12	PATHOLOGY LABELS-SIGMOID (BLACK/WHITE)	500
AL/MR-6	AN-LABELS-MIVACURIUM MG/ML (RED)	500	AL/PATH13	PATHOLOGY LABELS-RECTUM (BLACK/WHITE)	500
AL/MR-7	AN-LABELS-CISATRACURIUM (RED)	500	AL/PATH15	PATHOLOGY LABELS-GASTRIC (RED/WHITE)	500
AL/MR-8	AN-LABELS-PAPAVERINE MG/ML (RED)	500	AL/PATH16	PATHOLOGY LABELS- OESOPHAGUS (RED/WHITE)	500
AL/MRB	AN-LABELS – BLANK (NO TEXT) (RED)	500	AL/PATH18	PATHOLOGY LABELS-ANTRUM (BLACK/WHITE)	500
AL/MT-01	AN-LABELS-HALOPERIDOL (SALMON)	500	AL/PATH21	PATHOLOGY LABELS-S.B.B. (RED/WHITE)	500
AL/MT-02	AN-LABELS-CHLORPROMAZINE (SALMON)	500	AL/PATH22	PATHOLOGY LABELS- DUODENUM (RED/WHITE)	500
AL/MTB	AN-LABELS – BLANK (NO TEXT) (SALMON)	500	AL/PATH24	PATHOLOGY LABELS-RANDOM COLONIC (BLACK/WHITE)	500
AL/N	AN-LABELS-NARCOTIC (BLUE)	500	AL/PATH30	PATHOLOGY LABELS-GASTRIC BODY	500
AL/N-1	AN-LABELS-MORPHINE MG/ML (BLUE)	500			
AL/N-2	AN-LABELS-FENTANYL MCGM/ML (BLUE)	500			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
ANAESTHESIA LABEL ROLL (CONTINUED)					
AL/PATHBR01	PATHOLOGY LABELS-UPPER LOBE (BLACK/BLUE)	500	AL/V-3	AN-LABELS-NORADRENALINE (VIOLET)	500
AL/PATHBR02	PATHOLOGY LABELS-MIDDLE LOBE (BLACK/BLUE)	500	AL/V-4	AN-LABELS-EPHEDRINE MG/ML (VIOLET)	500
AL/PATHBR03	PATHOLOGY LABELS-LOWER LOBE (BLACK/BLUE)	500	AL/V-5	AN-LABELS-ISOPRENALINE (VIOLET)	500
AL/PATHC01	PATHOLOGY LABELS-MAIN CARINA (BLACK/ORANGE)	500	AL/V-6	AN-LABELS-DOPAMINE MG/ML (VIOLET)	500
AL/PATHC02	PATHOLOGY LABELS-UPPER LOBE (BLACK/ORANGE)	500	AL/V-7	AN-LABELS-DOBUTAMINE MG/ML (VIOLET)	500
AL/PATHC03	PATHOLOGY LABELS-MIDDLE LOBE(BLACK/ORANGE)	500	AL/V-8	AN-LABELS-PHENYLEPHRINE (VIOLET)	500
AL/PATHC04	PATHOLOGY LABELS-LOWER LOBE (BLACK/ORANGE)	500	AL/V-9	AN-LABELS-VASOPRESSIN (VIOLET)	500
AL/PATHTB01	PATHOLOGY LABELS-TRACHEA (BLACK/SALMON)	500	AL/VB	AN-LABELS – BLANK (NO TEXT) (VIOLET)	500
AL/PATHTB02	PATHOLOGY LABELS-UPPER LOBE (BLACK/SALMON)	500	AL/VB-1	AN-LABELS-INOTROPE	500
AL/PATHTB03	PATHOLOGY LABELS-INTERMEDIATE (BLACK/SALMON)	500			
AL/RRA	AN-LABELS-NMB – REVERSAL	500			
AL/RRA-1	AN-LABELS-NEOSTIG/ATROP (RED & WHITE)	500			
AL/RRA-2	AN-LABELS-EDROPHONIUM MG/ML (RED/WHITE)	500			
AL/RRA-3	AN-LABELS-NEOSTIGMINE (RED & WHITE)	500			
AL/RRA-4	AN-LABELS-PHYSOSTIGMINE MG/ML(RED & WHITE)	500			
AL/RRAB	AN-LABELS – BLANK (NO TEXT) (RED & WHITE)	500			
AL/T-1	AN-LABELS-MIDAZOLAM MG/ML (ORANGE)	500			
AL/T-2	AN-LABELS-DIAZEPAM MG/ML (ORANGE)	500			
AL/TA-01	AN-LABELS-FLUMAZENIL (ORANGE/WHITE STRIPE)	500			
AL/TB	AN-LABELS – BLANK (NO TEXT) (ORANGE)	500			
AL/V-1	AN-LABELS-ADRENALINE (VIOLET)	500			
AL/V-2	AN-LABELS-METARAMINOL (VIOLET)	500			
AL/V-2B	AN-LABELS-ARAMINE 10MG/10ML NORMALSALINE	500			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

POLY INSTRUMENT REPAIR TAG – NON STERILE

POLY INSTRUMENT REPAIR TAG		
PT001	POLY TAG – FAST TRACK (RED ON WHITE)	5cm x 9cm pack of 100
PT013/BLANK	POLY TAG – BLANK (NO TEXT)	50mm x 90mm with nylon tie pack of 100
PT013/BLUE	POLY TAG – BLUE – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/GOLD	POLY TAG – GOLD – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/GREEN	POLY TAG – GREEN – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/ORANGE	POLY TAG – ORANGE – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/PINK	POLY TAG – PINK – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/RED	POLY TAG – RED – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/VIOLET	POLY TAG – VIOLET – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT013/YELLOW	POLY TAG – YELLOW – BLANK (NO TEXT) N/S	50mm x 90mm with nylon tie pack of 100
PT015	INSTRUMENT REPAIR TAG – POLY	Pack of 100, tray or specialty
PT019/VIOLET	POLY TAG – LOAN SET VIOLET	50mm x 90mm with ties pack of 100
PT076-D	PRIORITY RED ON WHITE	7cm x 9cm roll 100
PT080/BLUE	POLY TAG – LOAN SET BLUE	50mm x 90mm with ties pack of 100
PT080/GREEN	POLY TAG – LOAN SET GREEN	50mm x 90mm with ties pack of 100
PT080/ORANGE	POLY TAG – LOAN SET ORANGE	50mm x 90mm with ties pack of 100
PT080/PINK	POLY TAG – LOAN SET PINK	50mm x 90mm with ties pack of 100
PT080/RED	POLY TAG – LOAN SET RED	50mm x 90mm with ties pack of 100
PT080/VIOLET	POLY TAG – LOAN SET VIOLET	50mm x 90mm with ties pack of 100
PT080/YELLOW	POLY TAG – LOAN SET YELLOW	50mm x 90mm with ties pack of 100
PT107	POLY TAG – FAST TRACK	5cm x 7cm pack of 100
PT108	INSTRUMENT REPAIR TAG	5cm x 7cm pack of 100
PT153/01G	THEATRE NUMBER 1 GREEN	Polytag with nylon tie pack of 100
PT153/02O	THEATRE NUMBER 2 ORANGE	Polytag with nylon tie pack of 100
PT153/03B	THEATRE NUMBER 3 BLUE	Polytag with nylon tie pack of 100
PT153/04Y	THEATRE NUMBER 4 YELLOW	Polytag with nylon tie pack of 100
PT153/05P	THEATRE NUMBER 5 PURPLE	Polytag with nylon tie pack of 100
PT153/06R	THEATRE NUMBER 6 RED	Polytag with nylon tie pack of 100
PT153/07S	THEATRE NUMBER 7 GOLD	Polytag with nylon tie pack of 100
PT153/08PUR	THEATRE NUMBER 8 PURPLE	Polytag with nylon tie pack of 100
PT153/09W	POLY TAG – THEATRE NUMBER 9 (WHITE)	With nylon tie pack of 100
PT166-5045MM	INSTRUMENT LABELS – NO TIES	50mm x 45mm with a hole
PT224	POLY TAG – LOAN SET BLACK ON WHITE	
PT237	POLY TAG – LIMITED REUSE	50mm x 90mm pack of 100

PLASTICWARE



FREECALL 1300 550 278

PLASTICWARE

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEASURING CUP



DEF459

DEF459	60ml, Height 5cm, Diameter 5cm	50
DEF460	60ml, Height 5cm, Diameter 5cm	250

NON STERILE

DEF1526	60ml, Height 5cm, Diameter 5cm	100
---------	--------------------------------	-----

GALLIPOT



DEF867

DEF867	60ml, Height 2cm, Diameter 7.5cm	50
DEF860	60ml, Height 2cm, Diameter 7.5cm	250
DEF1147	150ml, Height 5cm, Diameter 7.5cm	25
DEF888	150ml, Height 5cm, Diameter 7.5cm	150
DEF2792	150ml, Height 5cm, Diameter 7.5cm, Packet of 3	25

NON STERILE

DEF1114	150ml	500
DEF2362	150ml	200

DISPOSABLE JUG

These jugs are a clear recyclable plastic for easy viewing of contents. Comfortable handles and are graduated on both sides in ml and cc for easy measurement.

DEF2773	1L Tall Jug Increments in 100ml to 1L	12
---------	---------------------------------------	----



DEF2773

DEF2092	1.2L Jug Increments in 100ml to 1L	12
---------	------------------------------------	----



DEF2092

SPECIMEN JAR














DEF842 (Left), DEF758 (Middle), DEF840 (Right)


DEF758	70ml Yellow Lid	50
DEF843	70ml Yellow Lid	150
DEF278	70ml Yellow Lid, Double Packed, Packet of 4	40
DEF842	120ml Yellow Lid, Height 10.8cm, Diameter 4.2cm	150
DEF840	250ml Yellow Lid, Height 10cm, Diameter 6.2cm	80
DEF1109	250ml Yellow Lid	25
DEF1139	500ml White Lid	36




DEF1139

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
SCALPEL HOLDER					
DEF2820	Scalpel Holder	50			
			DEF2820		
SAFETY SHARPS TRANSFER DEVICE					
Device ensures the safe transfer and passing of sharps between scrub nurses and surgeons by shielding the sharp section of the blade.					
DEF2631	Bladesafe Safety Scalpel Transfer Device	24			
			DEF2631		
DISPOSABLE TRAY					
DEF2546	Square Fully Deep Tray 27cm x 25cm x 4cm Deep	20			
			DEF2546		
DEF321	Drug Tray 3 compartments 23cmx 9cm x 3cm Deep	80			
			DEF321		
DEF2690A	Dressing Tray 3 Compartment 20cm x 25cm x 2.5cm Deep	30			
			DEF2690A		
MIXING BOWL AND SPATULA					
					
			DEF1530		
DEF1530	Mixing Bowl 600ml and Spatula Double Packed	10			
DEF1532	Plastic Spatula 3cm x 17cm x 2.5cm	10			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
DISPOSABLE BOWLS AND RECEIVER SETS					
A wide selection of bowls, kidney dishes and containers that are recyclable and sturdy to use through out the hospital setting. Sterile and latex free.					
DEF2901	500ml Blue	20			
			DEF2901		
DEF2902	1L Blue	20			
			DEF2902		
DEF2489	1.5L Blue Double Packed	20			
			DEF2489		
DEF2550	1.5L Blue	20			
			DEF2550		
DEF2553	2.5L Blue Double Packed	10			
			DEF2553		

SPLASH BOWL		
		
DEF2650A		
DEF2650A	Splash Bowl 6L wrapped in 1 x SMS Field 130cm x 150cm	6
DEF2434A	Splash Bowl 6L wrapped in 2 x SMS Field 130cm x 150cm	6

BOWL DRAPE		
This flexible bowl drape is made from a tough plastic material that fits over a re-usable Splash Bowl. Size is 137cm x 137cm, bowl size reinforcement is 36cm x 36cm and wrapped in a laminated field 50cm x 50cm.		
80-13506-S	Bowl Drape	30
		
80-13506-S		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BOWL AND KIDNEY DISH SETS

KIDNEY DISH

Our disposable and recyclable polypropylene Kidney Dishes are flexible and easy to hold, with graduated markings up to either 700ml or 1L. Available sterile or non sterile for any area of the hospital.

DEF2791	Blue 700ml	35
---------	------------	----



DEF2791

DEF1174	Clear 700ml	35
---------	-------------	----



DEF1174

DEF2109	Yellow 700ml	35
---------	--------------	----



DEF2109

DEF2551	Blue 1L	20
---------	---------	----



DEF2551

NON STERILE

DEF1661	Green 400ml	640
---------	-------------	-----

PERI-OPERATIVE KIDNEY DISH SET

DEF2541A	2 x Kidney Dish Green 700ml 1 x Kidney Dish Clear 700ml 3 x Gallipot 150ml Wrapped in a blue SMS Field 75cm x 75cm	20
----------	---	----

KIDNEY DISH SET

DEF2957	1 x Kidney Dish Yellow 700ml 1 x Kidney Dish Blue 700ml 2 x Gallipot 150ml Wrapped in a blue plastic field 75cm x 75cm	20
---------	---	----

KIDNEY DISH AND GALLIPOT PACK

DEF409A	1 x Kidney Dish Clear 700ml 1 x Gallipot 150ml Wrapped in a blue plastic field 75cm x 75cm	35
---------	--	----

KIDNEY DISH SET

DEF2837	1 x Kidney Dish Clear 700ml 1 x Kidney Dish Blue 700ml 2 x Gallipot 150ml Wrapped in a blue plastic field 75cm x 75cm	20
---------	--	----

MAJOR BOWL SET

DEF2549	1 x Bowl 6L 1 x Bowl 2.5L 2 x Kidney Dish Clear 700ml 2 x Gallipot 150ml 1 x Tray 2 Compartments Wrapped in a blue SMS Field 130cm x 150cm	5
---------	---	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DISPOSABLE BOWL SET

DEF2433E	1 x Bowl 6L 1 x Bowl 1.5L 1 x Bowl 500ml 3 x Kidney Dish Green 700ml 4 x Gallipot 150ml 1 x Sharps Tray Rectangular Yellow 500ml 1 x Label "Has Time Out Been Done?" 1 x Skin Marker Standard Tip Wrapped in 2 blue SMS Field 130cm x 150cm	5
----------	---	---

STANDARD BOWL SET

DEF3180A	1 x Bowl 6L 2 x Bowl 250ml 1 x Kidney Dish Clear 700ml 1 x Kidney Dish Yellow 700ml 1 x Kidney Dish Blue 1000ml Wrapped in a blue SMS Field 150cm x 150cm	5
----------	--	---

BOWL SET

DEF2993	1 x Bowl 1.5L 1 x Bowl 1L 2 x Kidney Dish Blue 700ml 1 x Kidney Dish Yellow 700ml 4 x Gallipot 150ml Wrapped in a blue SMS Field 91cm x 91cm	12
---------	---	----

NEURO RECEIVER SET

DEF2490A	6 x Gallipots 150ml 1 x Jug 1.2L 1 x Label "Has Time Out Been Done?" 1 x Plastic Bag	8
----------	---	---



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAEMOSTATIC BANDAGES

First aid and Haemostatic products developed in Japan to simplify haemostasis treatment. These adhesive bandages protect and secure punctured skin after an IV injection or catheter removal. Used for various situations such as covering a punctured site following an injection, to stop the flow of blood from an intravenous line using pressure immobilisation, and to treat bleeding of a peripheral vein.

CHUSHAVAN MEDIUM SIZE – FIRST AID ADHESIVE BANDAGE

These adhesive bandages are ideal as a protective covering after an injection. They are 27mm x 27mm adhesive size with a pad size of 13mm x 13mm with 1mm thickness.

DEF1250	Chushavan Pressure Pad Medium	Box of 60 packets – each packet contains 8 pressure pads
---------	-------------------------------	--



DEF1250

IV PRESSURE PAD

These IV Pressure pads are suitable to use following the removal of an IV, after giving an intravenous injection or after drawing blood or immunisation injections. Comprises of a 36mm round adhesive with a 16mm haemostatic pad. The 3mm thick round pad is made of highly absorbent non-woven cellulose and creates a compression effect, resulting in bleeding stopping quickly. After bleeding has stopped, remove the pressure pad slowly and carefully.

DEF929	IV Pressure Pad White	100
--------	-----------------------	-----



DEF929

STEPTY HAEMOSTATIC BANDAGE

These pressure pads comprise of a rectangular 4 way stretch adhesive with a thick compression haemostatic pad. Designed as a protective covering of peripheral veins after blood sampling, a blood transfusion or an infusion. Two products available with different pad thicknesses.

DEF1251	Stepty Pressure Pad 6mm thick, length 39mm x 80mm, white	50
---------	--	----



DEF1251

DEF1252	Stepty P Pressure Pad 9mm thick, length 40mm x 120mm, white	30
---------	---	----



DEF1252

ISLAND AND PAD DRESSING

Quality island and pad dressing by Nichiban, Japan. These dressings are gentle, secure and have a long lasting urethane adhesive. The pad is a high quality absorbent cotton, whilst the clear adhesive is breathable, vapour permeable and skin friendly. Ideal for use as a post surgical dressing or for cuts and abrasions. Sterile and latex free.

CPSP0508	Cathereepus Film & Pad Island Dressing 5x8cm (Pad 2.5x5cm)	20
CPSP0915	Cathereepus Film & Pad Island Dressing 9x15cm (Pad 5x10cm)	20
CPSP0921	Cathereepus Film & Pad Island Dressing 9x21cm (Pad 5x16cm)	20

CATHEREEPLUS ROLL

A new urethane gel adhesive which offers high permeability and breathability. Comfortable and gentle on the skin, and easy to remove. Pull out necessary length and cut using scissors.

CPSR1010	100mm x 10m	20
----------	-------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BANDAGES

Excellent range of sterile bandages that can be used in the operating room.

MEDIUM WEIGHT

Medium weight elastic crepe bandage provides support for soft tissue injuries and will secure a dressing in place. Non fraying edges and a soft, comfortable cotton.



DEF268

DEF266	5cm x 1.6m Single Packed	48
DEF267	7.5cm x 1.6m Single Packed	48
DEF268	10cm x 1.6m Single Packed	48
DEF269	15cm x 1.6m Single Packed	48

HOSPITAL WEIGHT

Made from hospital quality, porous cotton with non-fraying edges. Ideally suited for light compression and holding a dressing in place.



DEF822

DEF819	2.5cm x 1.5m Double Packed	48
DEF820	5cm x 1.5m Double Packed	48
DEF821	7.5cm x 1.5m Double Packed	48
DEF822	10cm x 1.5m Double Packed	48
DEF823	15cm x 1.5m Double Packed	36

COTTON/WOOL CREPE BANDAGE (VIC)

Blended wool & cotton elastic crepe bandage with light support properties. Ideal for providing warmth and absorbency.



DEF517

DEF516	7.5cm x 2.2m Double Packed	48
DEF517	10cm x 2.2m Double Packed	48
DEF518	15cm x 2.2m Double Packed	24

HEAVY DUTY

Sterile, heavy weight supporting crepe bandage that provides a firm support to joint injuries. Has a centre guide for assisting application and is comfortable to wear.



DEF272

DEF270	5cm x 2.3m Double Packed	48
DEF271	7.5cm x 2.3m Double Packed	48
DEF272	10cm x 2.3m Double Packed	48
DEF273	15cm x 2.3m Double Packed	48

COHESIVE BANDAGE

Sterile, lightweight, highly conformable bandage that is non slipping. It will grip against the skin, without any adhesive, so ideal for sensitive skin.

DEF1056	6cm x 1.5m Double Packed	25
---------	--------------------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ORTHOPAEDIC PADDING

NATURAL UNDERCAST PADDING

Soft, undercast orthopaedic padding made from rayon fleece and used under casts to provide comfort and cushioning for a fracture. Pulls apart easily for a quick application.



DEF825N

DEF825N	15cm x 2.7m Double Packed	24
---------	---------------------------	----

WATER RESISTANT UNDERCAST PADDING

Comfortable and smooth orthopaedic padding featuring high cushioning and protection for a fracture. Easily conformable and is resistant to moisture absorbency.



DEF824

DEF829	5cm x 2.7m Double Packed	48
DEF830	7.5cm x 2.7m Double Packed	48
DEF824	10cm x 2.7m Double Packed	48
DEF825	15cm x 2.7m Double Packed	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ABSORBENT UNDERCAST PADDING

Absorbent padding which is multilayered to keep the skin dry and comfortable. Provides excellent comfort and cushioning.



DEF826

DEF828	7.5cm x 2.75m Double Packed	48
DEF826	10cm x 2.75m Double Packed	48
DEF827	15cm x 2.75m Double Packed	24

COTTON UNDERCAST PADDING

100% cotton undercast padding to provide soft cushioning and excellent absorbency. Has a unique crinkle weave pattern for loft and bulk. Can be used with either synthetic or plaster casting material.



DEF675

DEF673	5cm x 3.6m Double Packed	48
DEF674	7.5cm x 3.6m Double Packed	48
DEF675	10cm x 3.6m Double Packed	48
DEF676	15cm x 3.6m Double Packed	24

SURGICAL TAPE

Surgical or medical tapes are used to secure a dressing onto a wound. They have a low allergenic adhesive which is designed to hold firmly onto the skin or dressing materials but removes easily without damaging the skin.

ADHESIVE SURGICAL FIXATION TAPE – F

Sterile surgical tape that is used for easy fixation in the operating theatre. Coated with a hypoallergenic adhesive whilst the backing is a stretch non woven polyester. Suitable for fixation of wound dressings or fixation of tubes, catheters or instruments.



DEF1504

DEF1504	10cm x 12cm	100
DEF1505	10cm x 20cm	100
DEF1506	10cm x 45cm	100
DEF1614	15cm x 70cm	50
DEF1507	20cm x 20cm	100
DEF3943	30cm x 45cm	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ADHESIVE SURGICAL FIXATION TAPE – H

Excellent fixation tape which offers stretch and conformability and a low allergenic adhesive. Non woven fabric dressing sheet designed for rapid and quick fixation. Pre-cut and sterilised to a variety of sizes for immediate use in the operating theatre.



DEF1044

DEF1630	2.5cm x 15cm	100
DEF1623	2.5cm x 30cm	50
DEF3174	2.5cm x 100cm	100
DEF1647	5cm x 8cm	50
DEF1995	5cm x 15cm	66
DEF1461	5cm x 24cm	50
DEF1599	5cm x 25cm	50
DEF1624	5cm x 30cm	50
DEF1797	5cm x 50cm	60
DEF1044	10cm x 10cm	100
DEF943	10cm x 20cm	50
DEF944	10cm x 25cm	50
DEF1057	10cm x 36.5cm	100
DEF1798	10cm x 50cm	60
DEF2195	15cm x 7cm	140
DEF1600	15cm x 30cm	50
DEF946	15cm x 40cm	25
DEF1026	15cm x 60cm	100
DEF1025	15cm x 100cm	100
DEF2160	20cm x 8cm	125
DEF937	20cm x 25cm	25
DEF678	20cm x 30cm	100
DEF475	20cm x 75cm Double Packed	50
DEF3942	30cm x 45cm	100

SURGICAL TAPE NICHIMESH

Meshpore Tape from Japan, it is made of an elastic non-woven mesh with a gentle adhesive. Properties include excellent flexibility, softness, breathability and moisture permeability. Ideal for protective covering and securing.



DEF947

DEF947	5cm x 30cm, packet of 10, Double Packed	50
--------	---	----

COTTON TAPE FOR RETENTION

100% pure cotton tape pre-cut and sterilised ready for use in the operating room. Ideal for retention of tubes or retractors during a procedure.



DEF167

DEF173	1.5cm x 60cm Double Packed	100
DEF2106	1.25cm x 1.5m	100
DEF167	2.5cm x 2m Rondo	100

COTTON TAPE FOR RETENTION ROLL NON STERILE



DEF919

DEF918	1.2cm x 500m	1
DEF919	2.5cm x 250m	1
DEF1591	6mm x 1000m	1

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE TAPE – SILTAPE ROLL

Soft silicone perforated tape roll by Advancis Medical.
Gentle on the skin and easy removal and inspection of wound or dressings. Ideal use on fragile or sensitive skin.

CR3938	Siltape 2cmx 3m	1
--------	-----------------	---



CR3938

CR3939	Siltape 4cm x 1.5m	1
--------	--------------------	---



CR3939

TRACHEOSTOMY SPONGE

Soft, absorbent gauze ideal to use around any drainage area to prevent chaffing.

DEF053	10cm x 10cm Plain	100
--------	-------------------	-----



DEF053

DEF2165	10cm x 10cm Y Cut non woven, packet of 2	100
---------	--	-----



DEF2165

CR3900	Ultra Thin Hydrophilic Foam Tracheostomy Dressing 8cm x 8cm	20
CR4212	Standard Hydrophilic Foam Tracheostomy Dressing 8cm x 8cm	20

INDEX



A

ABDOMINAL SPONGES 31

AMNIOTIC HOOK 69

ANAESTHETIC PACKS

 Epidural pack 46

 Procedure pack 46

 Regional tray 46

 Spinal tray 47

ANAESTHETIC TRAY 45

ANAL SPECULUM 62

ANGIO INSTRUMENT SET 74

ANGIO PROCEDURAL DRAPES 73

ARTERIAL LINE TRAY 52

B

BACK TABLE COVER 6

BAGS

 Clip seal bag 5

 Burns plastic bag 5

 Drawstring bag 5

 Freezer bag 5

 Plastic bag 5

 Sanitary disposal bag 41

BANDAGES

 Crepe 92

 Cohesive 92

 Conforming 93

 Esmarch 35

 H conforming 93

 Haemostatic 91

 Heavy duty 92

 Hospital weight 92

 Medium weight 92

 Undercast 93-94

BANDS

 Latex free bands 21

 Rubber latex bands 21

BERET CAP 19

BIRTHROOM BUNDLE 71

BITE BLOCK 23

BLADESAFE SCALPEL TRANSFER DEVICE 86

BLANKET – BABY RECEIVING 69

BOWL

 Blue plastic 87

 Drape 87

 Guidewire 74

 Kidney dish 88

 Mixing bowl and spatula 86

 Set 88

 Splash 87

BOWL DRAPE 87

BOWL SETS

 Disposable bowl set 88

 Receiver set 89

 Standard bowl set 89

BRONNY IV PACKS 50

BURNS

 Dressing 41

 Gauze 41

 Plastic bag 5, 41

C

CAMERA & CORD COVERS 6, 27

CAPS BERET 19

CATH LAB PACKS 74

CATHETER PACKS

 Plain 48-50

 X-ray 48-50

CATS PAW RETRACTOR 58

CHAIR COVER 6

CHUSHAVEN PRESSURE PAD 91

CLAMPS

 Towel clamp U-shaped 21

 Tube clamp blue plastic 21, 54

 Tube clamp smooth jaw 21, 54

CLINICAL PROTECTIVE APPAREL 13

CLIPS

 Drape clip 22

 Towel clip green plastic 21

 Towel clip stainless steel 21

CLIP SEAL BAG 5

COMBINE

 Dressing 22

 Roll 22

 Specialised 22

COMBINE DRESSING 22

CORD CLAMP CUTTING DEVICE 70

COTTON TIPPED APPLICATORS/SWAB STICKS

 Mouth and throat jumbo 40

 Paper shaft 40

 Plastic shaft 40

 Wooden shaft 40

COTTON

 Balls 38

 Roll 38

 Wool 38

COVERS

 Back table 6

 Camera 6, 27

 Chair 6

 Cord 6

 Drill sleeve 6

 Equipment 6

 Handle sleeve 7

 Image intensifier 73

 Lead screen 73

 Light Handle 27

 Mayo stand 7

 Plastic sleeve 7

 Probe 73

 Protective 7

 Trolley 7

 Ultrasound 73

CUPS

 Baby cup with lid 69

 Measuring cups 85

CUSTOM PROCEDURE PACKS 44

CVC INSERTION PACK 52

D

DELIVERY SUTURE SET 70

DELIVERY BUNDLE 71

DEFOGGER 23

DELIVERY SET 70

DENTAL

 Bite block 23

 Dental instruments plastic 54

 Mirror 54

 Pellets 23

 Pressure pack 23

 Roll 23

 Toothbrush 23

DIATHERMY

 Pencil 24

 Quiver 24

 Tip cleaner 24

DRAINS

 Easy flow drain 38

 Evacuator bulb 36

 Flat drain 37

 Fluted drain 37

 Multilumen drain 38

 Round clear drain 37

DRAPES

 Angio procedural 73

 Aqua drape sheets 9

 Azure minor procedure 9

 Bowl 87

 Clear plastic 9

 Epidural 10

Fabric SMS	11	DRUG LABELS – AUSTRALIAN STANDARDS	77
Femoral radial angiography	73	DRUG TRAY PLASTIC	86
Leggings	11	E	
Minor procedure	9	EAR SPECULUM	61
Polybacked	10	EPIDURAL DRAPES	10
Stockinette		EPIDURAL/REGIONAL/ANAESTHETIC PACKS	46
Supreme drape sheets	9	EPISIOTOMY SUTURE SET	70
DRAPE CLIP	22	ESMARCH BANDAGE	
DRAPES/EQUIPMENT COVERS		Blue neoprene	35
Back table cover	6	Green silicone	35
Clear plastic	9	EVACUATOR BULB	36
Camera & cord cover	6	EYE	
Chair cover	6	Eye pad	25
Drill sleeve cover	6	Eye spear	25
Handle sleeve cover	7	Eye stroll	25
Mayo stand cover	7	Glasses disposable	19
Minor procedure	9	Face mask full	19
Plastic sleeve cover	7	SPEyes protective eyewear	18
Protective cover	7	F	
Polybacked	10	FACE SHIELD FULL	19
Tray liner	11	FEMORAL RADIAL ANGIOGRAPHY DRAPE	73
Trolley cover	7	FINE PROCEDURE PACK	50
Universal packs	10	FOAM PADS	25
DRAPE PACKS		FORCEPS	
Aqua universal drape pack long sides	10	Allis 4/5 teeth	58
Universal pack long	10	Adson	54
Universal pack short	10	Bozman	57
DRAWSTRING BAG	5	Crile artery	56
DRESSING PACK LARGE	51	Dressing blunt	56
DRESSINGS		Dressing pointed	56
Amputation	22	Fine croc action	57
Burns	41	Gillies	56
Combine	22	Hartmann croc action	57
Combine roll	22	Iris forcep	57
Island and pad	91		
DRESSING TRAY	86		
DRILL SLEEVE COVER	6		

Kelman-Mcpherson	57	GOWNS	
Lister sinus	57	Classic surgical gowns	17
Magill	56	Comfort surgical gowns	16
McIndoe	56	Coolmax surgical gowns	18
Micro Adson	54	Patient	14
Moorfield	58	Scrub suit pants	14
Mosquito	55	Scrub suit tops	14
Packing	58	Surgical	15
Plastic green	57	GUIDEWIRE BOWL	74
Rochester Pean	54	H	
Spencer Wells	55	HAND DRYING	
Splinter Feichenfeldt	57	Blue huck towel	11
Tenaculum	58	Dressing D towel	12
Tilley ear	58	Dressing T towel	12
Wilde	57	Dressing towel	12
FRAZIER SUCTION INSTRUMENT		Green huck towel	11
Standard	63	White huck towel	12
Long	63	HOOKS	
Extra long	63	Amniotic	69
FREEZER BAG	5	Neuro skin	30
FUJITA SUCTION INSTRUMENT	63	I	
G		IDENTILOOPS	26
GALLIPOTS	85	IMAGE INTENSIFIER COVER	73
GAUZE SWABS/SPONGES		INSTRUMENTS	54
Gauze burns	41	INSTRUMENT DRESSING PACK	67
Gauze folded	33	INSTRUMENT PACKS	64
Gauze green	32	INTERCOSTAL PACK	52
Gauze head dressing	33	IUD INSERTION KIT	71
Gauze non woven	39	IV INSERTION PACKS	50
Gauze plain	39	IV PRESSURE PAD	91
Gauze PurCotton plain swabs	39	J	
Gauze recovery	33	JUGS PLASTIC	85
Gauze specialised	33		
Gauze X-ray swabs	32		
Gauze X-ray PurCotton swabs	32		
Sponges abdominal	31		
Sponges peanut	33		

K	
KIDNEY DISHES PLASTIC	88
KIDNEY DISH SETS	88
L	
LABELS	
Anaesthesia roll	79
Customised	77
Drug.....	77
Peri-operative label sheets.....	78
Piggybacked	77
LABOUR PACK	71
LARYNGOSCOPE	
Surescope	27
LEAD SCREEN COVERS	73
LIGHT HANDLE COVERS	27
LUMBER PUNCTURE PACK – ADULT/CHILD	51
M	
MARKING PENS	31
MATERNITY PADS	69
MATERNITY SUTURE SET	70
MEDICAL IMAGING PACKS	74
MEDICAL LABELS	77
MID STREAM URINE PACK	50
MOUTH & THROAT JUMBO SWAB STICKS	40
MOUTH TOILET PACK.....	50
N	
NASAL BOLSTER	39
NASAL SPECULUM THUDICHUM	61
NASAL SPECULUM MOD WEIN	61
NEEDLE COUNTERS	28
NEEDLE HOLDERS	58
NEURO	
Gauze	33
Head dressing	33
Receiver set	89
Skin hooks	30
Yard strips	30
NEONATAL LUMBER PUNCTURE PACK	51
NIEVERT NEEDLE HOLDER	58
NIPPLE SHIELDS	69
O	
OVERSHOES	19
P	
PACEMAKER PACK	75
PADDING UNDERCAST	93-94
PHARMACY INSTRUMENT PACK	67
PERSPEX DISPENSER	78
PLACENTA BASIN	70
PLASTIC BAGS	5
PLASTICWARE	85
PLASTICS PACK	51
POLY INSTRUMENT REPAIR TAGS	83
PREP APPLICATORS	
Gauze prep sponges X-ray	33
Prep foam	30
Prep foam on a handle	30
PROBE COVERS	73
PROTECTIVE COVER	7
PRESSURE PAD	
Chushavan	91
I.V.	91
Stepty.....	91
Stepty P	91

Q	
QUIVER	24
R	
RETRACTORS	
Cats paw	58
Onetrac	34
RUBBER BANDS	21
S	
SAFETY PINS	30
SANITARY DISPOSABLE BAG	41
SCALPEL	
Handle	59
Holder	86
Penblade	59
Scalpel blade & handle	59
SCISSORS	
Blunt/blunt	59
Dressing sharp/blunt	59
Iris	60
Mayo	60
Metzenbaum	60
Sharp/sharp	59
Umbilical cord cutting.....	70
Universal	60
Vannas	60
SCOPE WARMER	34
SILICONE TUBING	
Pre-cut sterile	35
Roll non sterile	42
With connector	36
SILTAPE	96
SKIN MARKING PENS	31
SNUGGERS	
PVC tubing.....	34
SPATULA	86
SPECIMEN JARS	85
SPECULUMS	
Anal.....	62
Ear.....	61
Nasal	61
Vaginal.....	61
SPEYES EYEWEAR PROTECTION	18
SPLASH BOWL.....	87
SPLASH STOP	74
SPONGE HOLDERS	62
SPONGES	
Abdominal X-ray.....	31
Abdominal X-ray with tape	31
Abdominal PurCotton X-ray.....	31
Nasal bolster	39
Peanut	33
Tracheostomy	96
STAPLE REMOVER	62
STEPTY & STEPTY P PRESSURE PADS	91
STOCKINETTE.....	11
SUCTION INSTRUMENTS	
Frazier.....	63
Fujita	63
Paediatric	62
Poole	62
Rhoton	63
Yankauer	62
SURESCOPE	27
SURGICAL GOWNS	
Classic surgical gowns	17
Comfort surgical gowns	16
Coolmax surgical gowns.....	18
SURGICAL RETRACTOR	
Onetrac.....	34

SURGICAL TAPE		TOWEL	
Cotton	95	Blue huck towel	11
Nichimesh.....	95	Dressing D towel	12
Retraction	95	Dressing T towel	12
Siltape silicone	96	Dressing towel.....	12
Surgical tape F	94	Green huck towel.....	12
Surgical tape H	95	White huck towel.....	12
SUTURE BOOTS.....	36	TRACHEAL DILATOR.....	64
SUTURE SETS	64	TRACHY PROCEDURE PACK	52
SWABS/GAUZE		TRACHEOSTOMY SPONGE	96
Gauze head dressing	33	TRAY DISPOSABLE	86
Gauze non-woven.....	39	TUBE HOLDER	22
Gauze prep	33	TUBE CLAMP.....	21,54
Peanut	33	TUBING CONNECTOR	36
Swabs green	32		
Swabs folded	33	U	
Swabs plain	39	ULTRASOUND PROBE COVER	73
Swabs plain PurCotton	39	UMBILICAL CORD CLAMP CUTTING DEVICE.....	70
Swabs X-ray	32	UMBILICAL CORD CUTTING SCISSORS	70
Swabs X-ray PurCotton	32	UNIVERSAL PACKS	10
SWAB STICKS/COTTON TIPPED APPLICATOR		V	
Mouth and throat jumbo	40	VAGINAL SPECULUM	61
Paper shaft.....	40		
Plastic shaft	40	W	
Wooden shaft	40	WOUND CARE	91
T		WOUND DRAINS	
TAPE		Easy flow drain.....	38
Cotton	95	Evacuator bulb	36
Nichimesh.....	95	Flat drain	37
Retraction	95	Fluted drain	37
Siltape silicone	96	Multilumen drain.....	38
Surgical tape F	94	Round clear drain	37
Surgical tape H	95	WIRE CUTTER.....	64
TEATS BABY	69	X	
THYROID PACK	51	X-RAY INSTRUMENT PACK	74
TIP CLEANER	24		
TOOTHBRUSH	23		

Australia

Head Office
1-7 Vision Street
Dandenong South VIC 3175
t: 61 3 9706 3600
f: 61 3 9706 3619
freecall: 1300 550 278
e: sales@defries.com.au
www.defries.com.au

New Zealand

Level 4, 21 Queens Street
Auckland CBD 1010
New Zealand
toll free ph: 0800 511 740
toll free fax: 0800 511 840
e: sales@defries.co.nz
www.defries.co.nz

Manufacturing Site

18-24 Villas Road
Dandenong South VIC 3175