

CLINICAL PROTECTIVE APPAREL

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GOWNS

PATIENT GOWNS

Patient Gowns are lightweight, soft and comfortable and are available in a range of sizes to accommodate all body shapes. They are ideal for patients in hospitals, day procedures, GP clinics, radiology or any setting where a patient is required to wear a disposable gown.



PATIENT-Gown-SS-MED	Medium Short Sleeves	50
PATIENT-Gown-SS-LRG	Large Short Sleeves	50
PATIENT-Gown-SS-XL	X-Large Short Sleeves	50

SCRUB SUITS

Scrub Suits are made from a comfortable and soft fabric, and offer good protection. Ideal for nurses in the operating room, ER Departments, GP Clinics, Dentists, Vets and Students.



SCRUB-SUIT-SML-TOPS

SCRUB-SUIT-SML-TOPS	Top – Size Small	100
SCRUB-SUIT-MED-TOPS	Top – Size Medium	100
SCRUB-SUIT-LRG-TOPS	Top – Size Large	100
SCRUB-SUIT-XL-TOPS	Top – Size XL	100
SCRUB-SUIT-XXL-TOPS	Top – Size XXL	100



SCRUB-SUIT-SML-PANTS

SCRUB-SUIT-SML-PANTS	Pants – Size Small	100
SCRUB-SUIT-MED-PANTS	Pants – Size Medium	100
SCRUB-SUIT-LRG-PANTS	Pants – Size Large	100
SCRUB-SUIT-XL-PANTS	Pants – Size XL	100
SCRUB-SUIT-XXL-PANTS	Pants – Size XXL	100

SURGICAL GOWNS

Our single-use surgical gowns have specifically been developed for the Australian and New Zealand environment. Our range encompasses gowns for short and dry procedures through to high fluid and lengthy procedures.

OUR COMFORT RANGE – OUR SOFTEST AND MOST COMFORTABLE GOWNS

- Comfort Surgical Gown – Level 2.
- Comfort Unreinforced with Glove Assist – Level 3 gown with a reinforced sleeve to allow extra strength during gloving, but still cool, soft and ideal for short and dry procedures.
- Comfort Plus Fabric Reinforced Surgical Gown – Level 3 gown includes fluid resistant fabric reinforced panels yet remains lightweight and durable.
- Comfort Plus Poly Reinforced Surgical Gown – Level 3 gown which includes poly reinforced impervious panels for added protection.
- Comfort Ultimate Protection – Fully impervious high tech fabric, but still soft, comfortable and cool which meets AAMI Level 4 fabric standards.

OUR CLASSIC RANGE:

- Classic Standard Surgical Gown – suitable for most surgical procedures.
- Classic Plus Reinforced Surgical Gown – Includes a reinforced panel on the front and sleeves, for surgical cases requiring extra protection.
- Coolmax Surgical Gown – breathable and cool, these gowns are notable for their superb comfort and being fully impervious. With AAMI Level 4 fabric.

Designed in close co-operation with hospital staff, our Surgical Gowns are cool and breathable with an easy fit, while offering excellent protection against strikethrough from microbial and liquid penetration.




CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COMFORT SURGICAL GOWNS

LEVEL 2 COMFORT 35 UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 35gsm T-cut Surgical Gown ideal for dry and short procedures.




Level 2 Comfort 35 Unreinforced

DEF5316	Medium 119cm	24
DEF5317	Large 127cm	24
DEF5318	X-Large 137cm	24
DEF5319	XL-Long 149cm	24
DEF5320	XX-Large 155cm	24

LEVEL 3 COMFORT UNREINFORCED WITH GLOVE ASSIST

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41gsm T-cut Surgical Gown ideal for dry and short procedures. Has innovative Glove Assist fabric reinforcement to sleeve area to ensure extra strength. Comfortable and cool to wear.




Level 3 Comfort Unreinforced With Glove Assist

DEF5301	Medium 119cm	24
DEF5302	Large 127cm	24
DEF5303	X-Large 137cm	24
DEF5304	XL-Long 149cm	24
DEF5305	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS FABRIC REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41gsm T-cut Surgical Gown with fabric reinforcement on front and sleeves. Fluid resistant on panels. Comfortable to wear.




Comfort Plus Fabric Reinforced

DEF5306	Medium 119cm	24
DEF5307	Large 127cm	24
DEF5308	X-Large 137cm	24
DEF5309	XL-Long 149cm	24
DEF5310	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile 41 gsm T-cut Surgical Gown with impervious poly reinforcement. Durable and impervious to fluid on panels. Comfortable to wear.




Comfort Plus Poly Reinforced

DEF5311	Medium 119cm	24
DEF5312	Large 127cm	24
DEF5313	X-Large 137cm	24
DEF5314	XL-Long 149cm	24
DEF5315	XX-Large 155cm	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 4 COMFORT ULTIMATE PROTECTION

Surgical Gown T-cut wrapped in a 91cm x 91cm SMS field with two hand towels 40cm x 58cm. Sterile Level 4, 67gsm trilaminate material. High tech, fully impervious and breathable fabric which is extremely comfortable and cool for surgeons and or nurses.




Comfort Ultimate Protection made with Level 4 fabric

DEF5321	Medium 119cm	24
DEF5322	Large 127cm	24
DEF5323	X-Large 137cm	24
DEF5324	XL-Long 149cm	24
DEF5325	XX-Large 155cm	24

CLASSIC SURGICAL GOWNS

LEVEL 3 CLASSIC UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown ideal for short and dry procedures.




Level 3 Classic Unreinforced

DEF5000A	Medium 110cm	25
DEF5001A	Large 130cm	25
DEF5002A	X-Large 140cm	25
DEF5004A	XL-Long 157cm	25

LEVEL 3 CLASSIC PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown with impervious poly reinforcement. Colour coded neck band – purple for level 3 Gown.



Level 3 Classic Plus Poly Reinforced

DEF5005A	Medium 110cm	25
DEF5006A	Large 130cm	25
DEF5007A	X-Large 140cm	25
DEF5009A	XL-Long 157cm	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------


WEARABLE CONSUMABLES

CLASSIC SURGICAL GOWNS (CONTINUED)

LEVEL 4 COOLMAX

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile Level 4, trilaminate material.

High quality T-cut Surgical Gown ideal for high fluid procedures or lengthy procedures. Colour Coded neck band – blue for level 4 Gown.



Coolmax Surgical Gown Level 4

DEF5010A	Medium 110cm	25
DEF5011A	Large 130cm	25
DEF5012A	X-Large 140cm	25
DEF5013A	XL-Long 160cm	25
DEF3198A	Coolmax Surgical Gown Pack – 1 Large and 1 X-Large Gown wrapped in a 200cm x 200cm SMS field with two huck towels	10

SPEyes EYEWEAR NON STERILE




Our popular DeRoyal SPEyes eyewear is designed to minimise exposure to pathogens from waste fluid. They combine innovation, quality, value, and comfort for clinicians where personal protective equipment is a vital part of self-care.

- Wrap-around disposable lens for complete protection
- UltraClear™ technology prevents glare, fog, and static
- Frame fits completely against forehead for splash protection from above
- Sealed lens mounting for added fluid protection
- Clip-On/Off design for quick and easy changing
- Coloured top strip allows finger-print-free handling
- Latex free



23-505

23-505	SPEyes Eye Shields pre-assembled assorted colours	100
--------	---	-----



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPEyes EYEWEAR NON STERILE



23-506

23-506	Dispenser box with pre-assembled SPEyes Eye Shields (30 per box 4 boxes per case)	120 total
23-501	SPEyes Eye Shield Replacement Lens	100
23-502	SPEyes Eye Shield Frame	25

FACESHIELDZ™



23-002

23-002	Face Mask Full with Elastic Strap	50
--------	-----------------------------------	----

GLASSES



DEF1755

DEF1755	Disposable Glasses Clear Polycarbonate Non Sterile	1
---------	--	---

BERET CAPS



DEF1142

DEF1142	Beret Cap 53cm (21") Blue Non Sterile	250
---------	---------------------------------------	-----

OVERSHOES NON STERILE



DEF1145

DEF1145	Non Skid Universal Overshoes Polypropylene	100
---------	--	-----