

# CONSUMABLES

FRECALL 1300 550 278

CONSUMABLES

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## THEATRE CONSUMABLES

All surgical consumables are pre-packed and sterile, ready to be used in the operating room, unless indicated as Non Sterile.

### BANDS

#### RUBBER BANDS (CONTAINS LATEX)

ID indicates Inner Diameter. BW indicates Band Width.



DEF1405

DEF1265	#28 – (ID) 40mm x 3.5mm (BW), Packet of 1	25
DEF1590	#16 – (ID) 60mm x 1.5mm (BW), Packet of 2, Double Packed	100
DEF624	#64 – (ID) 90mm x 6mm (BW), Packet of 1	100
DEF627	#64 – (ID) 90mm x 6mm (BW), Packet of 3, Double Packed	100
DEF626	#64 – (ID) 90mm x 6mm (BW), Packet of 5, Double Packed	100
DEF578	#31 – (ID) 60mm x 3mm (BW), Packet of 1	100
DEF1405	#31 – (ID) 60mm x 3mm (BW), Packet of 2	100
DEF339	#31 – (ID) 60mm x 3mm (BW), Packet of 3	100
DEF625	#31 – (ID) 60mm x 3mm (BW), Packet of 10, Double Packed	100
DEF599	#34 – (ID) 100mm x 3mm (BW), Packet of 2, Double Packed	25
DEF619	#34 – (ID) 100 mm x 3mm (BW), Packet of 5	25
DEF1793	#33 – (ID) 90mm x 3mm (BW) Packet of 2	25
DEF669	#33 – (ID) 90mm x 3mm (BW), Packet of 5, Double Packed	25

#### LATEX FREE BANDS

DEF3542	#33 – (ID) 90mm x 3mm (BW), Packet of 2	100
DEF3560	#33 – (ID) 90mm x 3mm (BW), Packet of 3	100

### CLAMPS/CLIPS/HOLDERS

#### TUBE CLAMPS

DEF1575	Tube Clamp with Lock Blue Plastic 12cm	50
---------	--	----



DEF1575

DEF2190	Towel Clip Green Plastic 11.5cm	25
---------	---------------------------------	----



DEF2190

DEF2488	Towel Clip Stainless Steel 11cm Sharp/Sharp	10
---------	---	----



DEF2488

DEF2298	Tube Clamps 13cm Smooth Jaw	10
---------	-----------------------------	----



DEF2298

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**CLAMPS/CLIPS/HOLDERS (CONTINUED)**

**DRAPE CLIP**

General purpose, reusable Drape Clip, used during surgery for clamping purposes to ensure items remain where placed. A wide easy to use grip surface ensures heavy duty holding power. Available in 8 different colours. Non sterile.




DRAPE-CLIP

DRAPE-CLIP	Non Sterile Mixed Colours	1
------------	---------------------------	---

**TUBE HOLDER**

Ideal for securing and holding tubes during a procedure. Adhesive backing for securement onto a drape or towel. Adjustable velcro to allow gentle and easy placement.



DEF1889

DEF1889	Tube Holder with Velcro Adhesive 7.5cm x 7.5cm	25
---------	--	----

**COMBINE DRESSING**

Highly absorbent layer of cellulose blend enclosed in a soft and conformable non woven fabric. Each non woven Combine is pre-cut and each end is sealed into a variety of sizes.



DEF711

DEF711	9cm x 20cm	25
DEF719	9cm x 20cm, Packet of 5	100
DEF720	9cm x 20cm, Packet of 10	50
DEF1322	9cm x 30cm, Packet of 2	50
DEF697	9cm x 90cm	25
DEF715	20cm x 20cm	25
DEF80	20cm x 60cm	100
DEF1160	20cm x 80cm	50

**COMBINE ROLL**

DEF714	Non Woven Combine Roll 9cm x 10m, Double Packed	16
DEF1004	Non Woven Combine Roll 20cm x 2.5m	25

**COMBINE SPECIALISED**



DEF804

DEF2289	Non Woven Amputation Dressing T-Shaped Melolin 40cm x 20cm	35
DEF804	Non Woven Amputation Dressing 40cm x 70cm	50

**DENTAL**

**BITE BLOCK**



DEF1828

DEF1232	Bite Block 7.5cm X-ray without Tail	100
DEF1232C	Bite Block 7.5cm X-ray without Tail Non Sterile	100
DEF1828	Bite Block 7.5cm X-ray with Tail	100
DEF1246	Bite Block 7.5cm X-ray with Tail Non Sterile	100

**DENTAL PELLET**

DEF1180	Dental Pellet Size #0, Packet of 5, Double Packed	100
---------	---	-----

**TOOTHBRUSH**



DEF690

DEF690	Toothbrush	25
--------	------------	----

**DEFOGGER**

Specially formulated solution that safely prevents lens fogging. Comes with a non-abrasive, x-ray detectable foam pad. Sterile.



28-0101	DeFogger Oval Bottle	20
---------	----------------------	----

**DENTAL PRESSURE PACK**



DEF799

DEF799	X-ray with Cotton Tail	100
DEF1063	X-ray with 35cm Tape, Packet of 2	100
DEF144	3cm x 4cm X-ray Small with 10cm tail, Packet of 2	50



DEF144

DEF145	4cm x 5cm X-ray Large with 9cm tail, Packet of 2	50
--------	--	----



DEF145

DEF1209	Dental Roll 4mm x 35mm, Packet of 5	100
---------	-------------------------------------	-----



DEF1209

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**ELECTROSURGICAL**

**DIATHERMY PENCILS**

Hand activated Electrosurgical Pencils with an ergonomic design.

88-E0031	Rocker Switch Pencil and Stainless Steel Blade	100
----------	--	-----



88-E0031

88-E0041	Rocker Switch Pencil and Teflon Blade	100
----------	---------------------------------------	-----



88-E0041

**DIATHERMY PACK WITH STAINLESS STEEL BLADE**

DEF2864	1 x Diathermy Pencil with Rocker Switch and Stainless Steel Blade 1 x Quiver 5cm x 25cm with Metal Ring 1 x Scratch Pad Tip Cleaner	10
---------	---	----


**DIATHERMY PACK WITH TEFLON BLADE & QUIVER**

DEF2863A	1 x Diathermy Pencil with Rocker Switch and Teflon Blade 1 x Quiver 5cm x 25cm with Metal Ring	10
----------	---	----

**DIATHERMY TIP CLEANER**

This coated scratch pad cleaner quickly and easily removes debris from instrument tips and diathermys. The adhesive backing will secure the Tip Cleaner onto draped surfaces and contains a radiopaque strip to ensure X-ray detectability.

DEF1886	Tip Cleaner 5cm x 5cm with Adhesive and Radiopaque Strip	50
---------	--	----



DEF1886

**QUIVER**

A single use disposable polypropylene Quiver used to hold diathermy during surgical procedures. Quiver is recyclable (metal ring needs to be removed before being placed in recycling container).

DEF2772	5cm x 25cm with Metal Ring	90
---------	----------------------------	----



DEF2772

DEF3068	5cm x 25cm with Cable Tie	90
---------	---------------------------	----



DEF3068

DEF3291	Large 9cm x 35cm with Metal Ring	40
---------	----------------------------------	----

DEF3296	Large 9cm x 35cm with Cable Tie	40
---------	---------------------------------	----

**EYE**

**EYE SPEAR HIGHLY ABSORBENT**



DEF484

DEF484	Eye Spear, Packet of 2	100
DEF485	Eye Spear, Packet of 5	20
DEF486	Eye Spear, Packet of 6	20
DEF487	Eye Spear, Packet of 10	20

**EYE STROLL HIGHLY ABSORBENT**



DEF112

DEF112	Eye Stroll, Packet of 5	50
DEF081	Eye Stroll, Packet of 10	50
DEF704	Eye Strolls Bulk Non Sterile	500

**FOAM PAD**

Versatile polyurethane foam pads used in the operating room for positioning of patients or relieving pressure areas for patients.



DEF792

DEF792	100mm x 100mm x 13mm	100
DEF769	100mm x 100mm x 25mm	50
DEF783	100mm x 180mm x 13mm	50
DEF784	100mm x 180mm x 25mm	40
DEF738	150mm x 150mm x 13mm	50



DEF738

DEF1738	150mm x 315mm x 25mm	50
DEF914	Plastics Foam 150mm x 200mm x 25mm	25
DEF1549	200mm x 25mm x 25mm	100



DEF1549

DEF739	250mm x 250mm x 13mm	50
DEF786	300mm x 600mm x 13mm	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**IDENTILOOPS**

These radio-opaque Identiloops are used in surgical procedures for retraction, identification or occlusion of arteries, veins, nerves or tendons. A soft and flexible, inert medical grade silicone which is oval in cross-section, allows for even pressure distribution to minimise damage to tissues. Available in 4 colours – Blue, White, Red and Yellow in mini and maxi sizes. All Double Packed.

2200220-2	MINI-RED 40cm, Packet of 2	10
2200220-4	MINI-RED 40cm, Packet of 4	10



2200220

2200420-2	MINI-YELLOW 40cm, Packet of 2	10
-----------	-------------------------------	----



2200420

2200620-2	MINI-BLUE 40cm, Packet of 2	10
2200620-4	MINI-BLUE 40cm, Packet of 4	10



2200620

2200920-2	MINI-WHITE 40cm, Packet of 2	10
2200920-4	MINI-WHITE 40cm, Packet of 4	10



2200920

2202220-2	MAXI-RED 40cm, Packet of 2	10
2202220-4	MAXI-RED 40cm, Packet of 4	10



2202220

2202420-2	MAXI-YELLOW 40cm, Packet of 2	10
2202420-4	MAXI-YELLOW 40cm, Packet of 4	10



2202420

2202620-2	MAXI-BLUE 40cm, Packet of 2	10
-----------	-----------------------------	----



2202620

2202920-2	MAXI-WHITE 40cm, Packet of 2	10
-----------	------------------------------	----



2202920

DEF1193	SUPER MAXI BLUE 40cm, Packet of 2	10
---------	-----------------------------------	----

DEF1318	MIXED MAXI 40cm, Packet of 5	10
---------	------------------------------	----

DEF2182	SUPER MAXI BLUE 60cm, Packet of 2	10
---------	-----------------------------------	----




DEF2182

**LARYNGOSCOPE**

**SURESCOPE**

SURESCOPE is a single-use one piece laryngoscope with a built-in LED light source. This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting by obsoleting reusable lights sources, handles and blades.



C090120	MAC 2	20
C090130	MAC 3	14
C090140	MAC 4	14
C090205	MILLER 00	20
C090200	MILLER 0	20
C090210	MILLER 1	20
C090220	MILLER 2	14
C090230	MILLER 3	14

**CAMERA COVERS**

Designed to be used in the operating room to create a sterile barrier between the user and the camera. Contain a lens inserted into the cover that is ultra clear and durable.



DEF3673

DEF3591	Camera Light Handle Cover 66mm with Lens Box	25
DEF3673	Camera Light Handle Cover 75mm with Lens Box	25

**LIGHT HANDLE COVERS**

Our range of flexible light handle covers are single sterile and come either coloured as blue or green. These adaptable covers are twisted or clicked onto the light handle adapter and are firmly held for duration required. We can supply the compatible light handle adapters for each operating room light. Easy and versatile to use.



DEF2082

DEF2082	Light Handle Cover Green – friction fit	50
---------	---	----

DEF2919	Light Handle Cover Green – friction fit, packet of 2	40
---------	--	----

DEF2675	Light Handle Cover Blue – click on	50
---------	------------------------------------	----



DEF2675

DEF3170	Light Handle Cover Blue – click on, packet of 2	40
---------	---	----

DEF2557	Light Handle Cover Foam Blue	50
---------	------------------------------	----



DEF2557

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------


**NEEDLE COUNTERS DEROYAL**

Excellent selection of Needle Counters from DeRoyal with easy to read numbers and description, facilitating correct counting of blades and needles used in a procedure. Full surface magnetic face provides extra space for sharps, pins, blade or needles. 'Y' codes indicate Needle Counter case is Yellow.

**DELUXE NEEDLE COUNTERS**

**MAGNET**

25-0103Y	Needle Counter Single Magnetic/Clear Top, 10 Count, 20 Capacity	50
----------	---	----



25-0103Y

**MAGNET/FOAM COMBINATION**

25-0401	Needle Counter Magnetic/Foam, 20 Count, 40 Capacity	20
---------	---	----




25-0401

25-0701	Needle Counter Foam Strips/Magnetic 15 Count, 30 Capacity	50
---------	---	----



25-0701

25-0702	Needle Counter Magnetic/Foam Strips, 30 Count, 60 Capacity	27
---------	--	----



25-0702

**FOAM**

25-0504	Needle Counter Foam/Adhesive, 40 Count, 80 Capacity	12
---------	---	----



25-0504

25-0803	Needle Counter Foam Strips, 60 Count, 120 Capacity	27
---------	--	----



25-0803

**DEBLADE NEEDLE COUNTERS**

This range of Needle Counters all feature a blade removal system which provides safe and simple removal of contaminated blades.

**MAGNET**

25-1104Y	Needle Counter Double Magnetic, 30 Count, 60 Capacity with DeBlade	50
----------	--	----



25-1104Y

25-1108	Needle Counter Double Magnetic 60 Count, 120 Capacity with DeBlade	27
---------	--	----



25-1108

**MAGNET/FOAM COMBINATION**

25-1401	Needle Counter Magnet/Foam Block 20 Count, 40 Capacity with DeBlade	20
---------	---	----



25-1401

25-1402	Needle Counter Magnetic/Foam block 10 Count, 20 Capacity with DeBlade	25
---------	---	----

25-1402Y	Needle Counter Magnetic/Foam block 10 Count, 20 Capacity with DeBlade	25
----------	---	----

25-1701	Needle Counter Magnetic/Foam Strips 15 Count, 30 Capacity with DeBlade	50
---------	--	----



25-1701

25-1702	Needle Counter Magnetic/Foam Strips 30 Count, 60 Capacity with DeBlade	27
---------	--	----



25-1702

**FOAM**


25-1801Y	Needle Counter Foam Strips 30 Count, 60 Capacity with DeBlade	25
----------	---	----



25-1801Y

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**NEURO SKIN HOOKS**



DEF1032

DEF336	5 on a Tray Small Hook, Double Packed	20
DEF1032	5 on a Tray Size: 3/0 Large Hook	20

**NEURO YARD STRIPS**

DEF858	Neuro 1/2 (53cm) Yard Strips X-ray, Packet of 5	50
--------	---	----



DEF858

**PREP APPLICATORS**

DEF921	Prep Foam 50mm x 60mm x 25mm Cream No Handle	100
--------	--	-----



DEF921

DEF1846	Prep Foam on a Handle 15cm Blue	100
---------	---------------------------------	-----



DEF1846

DEF1847	Prep Foam on a Handle 15cm Blue, Packet of 2	100
---------	--	-----

DEF2035	Prep Foam on a Handle 20cm Orange	100
---------	-----------------------------------	-----



DEF2035

**SAFETY PINS**

DEF071SP	#2 (38mm) Single Packed	50
DEF072SP	#3 (51mm) Single packed	50
DEF073SP	#4 (57mm) Single packed	50

**NON-STERILE**

DEF635	#3 (51mm) Non Sterile	144
--------	-----------------------	-----



DEF635

DEF636	#4 (57mm) Non Sterile	144
--------	-----------------------	-----



DEF636

**SKIN MARKING PENS**

Skin Markers are used during surgical procedures and contain medical grade, non irritating gentian violet ink. They mark a smooth, clear line on a patient's skin, and the ink is non toxic, non tattooing and quick drying.

26-001	Skin Markers Standard Tip with Ruler 15cm long	50
--------	--	----



26-001

26-002	Skin Marker Fine Point with Ruler 15cm	50
--------	--	----



26-002

26-019	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip	50
--------	--	----



26-019

26-020	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip with Ruler & Labels	50
--------	--	----

DEF464	Prepacked Marking Pencil – Omnichrom 108-9	24
--------	--	----



DEF464

**NON-STERILE**

DEF2669	Mini pre-op Skin Markers Non Sterile	100
DEF2183	Skin Markers Standard Tip	50
26-018PRNS	Prep-Resistant Ink Mini Skin Marker, Pre-Op	100
26-005PRNS	Prep-Resistant Ink Skin Marker, Standard Tip	100

**SPONGES & SWABS**

**SPONGES ABDOMINAL X-RAY**

Soft highly absorbent cotton sponges. Each complete with an X-ray radiopaque element to ensure detection.



DEF1350

DEF1350	23cm x 23cm, Packet of 5	100
DEF162	23cm x 23cm with Tape, Packet of 1	100
DEF160	23cm x 23cm with Tape, Packet of 3	100
DEF168	23cm x 23cm with Tape, Packet of 5	100
DEF1351	30cm x 30cm, Packet of 5, Double Packed	75
DEF1352	30cm x 48cm, Packet of 5, Double Packed	60



DEF162

**SPONGES & SWABS (CONTINUED)**

**SPONGES PURCOTTON ABDOMINAL X-RAY**

PurCotton X-ray Abdominal sponges are made from a unique 100% cotton sourced from the US. They are soft and comfortable, low linting, have superior absorption and wet tensile strength, and are 100% biodegradable. The PurCotton is the next evolution of cotton gauze.



DEF3741

DEF3741	30cm x 30cm X-ray 50g, Packet of 5	80
DEF3745	30cm x 45cm X-ray 50g, Packet of 5	60
DEF3743	45cm x 45cm X-ray 50g, Packet of 5	30

**GAUZE SWABS X-RAY**

Excellent absorbency and quality swabs with 8 ply cotton blend. Includes an X-ray radiopaque element for easy detection.



DEF044

DEF202	7.5cm x 7.5cm X-ray, Packet of 5, Double Packed	100
DEF044	10cm x 10cm X-ray, Packet of 5, Double Packed	100
DEF164	10cm x 10cm X-ray, Packet of 10, Double Packed	50

**GAUZE SWABS PURCOTTON X-RAY**



DEF3737

DEF3737	10cm x 10cm X-ray 40g, Packet of 5, Double Packed	500
DEF3738	10cm x 10cm X-ray 50g, Packet of 5, Double Packed	200

**GAUZE SWABS GREEN**



DEF190

DEF190	10cm x 10cm Green, Packet of 2	100
DEF192	10cm x 10cm Green, Packet of 3	100
DEF207	10cm x 10cm Green, Packet of 5, Double Packed	100
DEF196	10cm x 10cm Green, Non Sterile	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**GAUZE SWABS FOLDED**

DEF347	Gauze Swab 7cm x 21cm X-ray, Folded, Packet of 5	50			
--------	--	----	--	--	--



DEF347

DEF185	Gauze Swab Neuro 10cm x 10cm 32ply X-ray, Packet of 10, Double Packed	100			
--------	---	-----	--	--	--



DEF185

DEF1074	Gauze Swab Neuro 10cm x 20cm X-ray, Folded, Packet of 5	25			
---------	---	----	--	--	--

**GAUZE HEAD DRESSINGS**

DEF611	Head Dressing Kit includes: 8 x Gauze Swabs 10cm x 20cm 20 x Gauze Swabs 10cm x 10cm 2 x 90cm x 90cm Gauze, Folded	35			
--------	---	----	--	--	--

**GAUZE PREP SPONGES X-RAY**

Round sponges with X-ray radiopaque element. Excellent absorbency. Used as a prep sponge prior to a surgical procedure.



DEF251

DEF255	Sponges Prep X-ray Large, Packet of 1	100			
DEF251	Sponges Prep X-ray Large, Packet of 2	100			

**GAUZE PEANUT SPONGES**

Round sponges ideal for delicate surgical sponging or soft tissue dissection. X-ray radiopaque element for easy detection.

DEF315	X-ray, 5 x Small in an Instapack	50			
--------	----------------------------------	----	--	--	--



DEF315

DEF046	X-ray, 5 x Small on a Safety Pin	50			
--------	----------------------------------	----	--	--	--



DEF046

DEF313	X-ray, 5 x Small in a Bag	50			
--------	---------------------------	----	--	--	--

DEF316	X-ray, 5 x Large in an Instapack	50			
--------	----------------------------------	----	--	--	--



DEF316

**SPECIALISED GAUZE**



DEF169

DEF986	Donor Site Gauze Dressing 5cm x 10cm, Double Packed	100			
--------	---	-----	--	--	--

DEF169	Recovery Gauze 5cm x 5cm Plain 32 Ply	100			
--------	---------------------------------------	-----	--	--	--

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**SURGICAL RETRACTOR**

**ONETRAC**

ONETRAC is a single-use cordless surgical retractor with a built-in LED light source and smoke evacuation channel. It provides illumination in the depths of a surgical tissue pocket or cavity. No additional parts or components or fibre optic cables are required. Ideal for breast or reconstructive procedures, as well as other surgical procedures requiring access, illumination and retraction of soft tissue. Reduce cost, increase efficiency and eliminate unnecessary risk of cross-contamination with ONETRAC.



ONETRAC

C100110-5	135mm x 30mm	5
C100140-5	90mm x 22mm	5
C100150-5	40mm x 20mm	5

**SNUGGERS PVC TUBING**

Medical grade, clear and transparent Tubing Snuggers used during cardiac surgery. A radiopaque strip along length of tubing ensures X-ray detectability.

DEF3051	PVC Tubing Snuggers 16Fr x 17cm X-ray, 6 per Packet, Double Packed	25
---------	--	----



DEF3051

DEF3052	PVC Tubing Snuggers 16Fr x 22cm X-ray, 6 per Packet, Double Packed	25
---------	--	----

**SCOPE WARMER**

Long, double sided Scope Warmer from DeRoyal that pre-warms endoscopes to body temperate. Easy to use design provides 1 hour of warming time on each side.

TM500D	Scope Warmer	10
--------	--------------	----



TM500D

**ESMARCH BANDAGE**

Offers excellent compression and elastic properties to achieve exsanguination, whilst delivering protection for latex sensitive patients and health professionals.

DEF1006	6.5cm x 2m Green Silicone	10
---------	---------------------------	----



DEF1006

DEF2764	10cm x 2.7m Blue Neoprene	20
---------	---------------------------	----



DEF2764

DEF1081	10cm x 3m Green Silicone	5
---------	--------------------------	---



DEF1081

DEF2668	15cm x 2.7m Blue Neoprene	20
---------	---------------------------	----



DEF2668

**SILICONE TUBING**

Our 100% medical grade, clear silicone tubing range is available pre-cut and sterile or as non-sterile rolls. Please contact our customer service to discuss sizing.

Sizing below is represented as:  
ID (Internal Diameter) x OD (Outer Diameter) x Length of tube in cm or m.

DEF2307	0.30mm x 0.64mm x 20cm	75
DEF2522	1mm x 2mm x 50cm	50
DEF477	1.57mm X 2.41mm X 30cm, 2 per Packet, Double Packed	25
DEF1564A	1.57mm x 3.18mm x 50cm	25
DEF1545	1.57mm x 3.18mm x 1m	15
DEF1428	1.57 mm x 2.41mm x 15cm, Double Packed	100
DEF473	1mm x 3mm x 50cm, 4 per Packet	20
DEF1944A	3mm x 5mm x 2.3m	15
DEF2314	3.18mm x 6.35mm x 25cm	60
DEF2038	4mm x 6mm x 45cm	55
DEF1665	5mm x 7mm x 20cm, 2 per Packet, Double Packed	25
DEF1268	5mm x 7mm x 42cm	25
DEF481	5mm x 8mm x 45cm	50
DEF1253	5mm x 8mm x 1.6m, Double Packed	50
DEF2060	5mm x 8mm x 3.1m	8
DEF1261	5mm x 9mm x 1.5m	50
DEF2786	6.4mm x 9.5mm x 2cm, 2 per Packet	375
DEF2032	7mm x 10mm x 3m	8
DEF482	7mm x 12mm x 85cm, 2 per Packet, Double Packed	100
DEF1204	8mm x 14mm x 3m	5
DEF2468	10mm x 16mm x 1.5m	16



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**SILICONE TUBING (CONTINUED)**

**SILICONE TUBING WITH CONNECTOR**

Flexible pre-cut and sterile silicone tubing with connector attached.



DEF1186

DEF735	5mm x 8mm x 1m, Double Packed	50
DEF1186	5mm x 9mm x 30cm, Double Packed	50
DEF306	5mm x 9mm x 50cm Neuro	50
DEF1712	5mm x 9mm x 60cm	40
DEF1337	5mm x 9mm x 1m Neuro	50
DEF1231A	5mm x 9mm x 1.5m	30
DEF305	7mm x 11mm x 1m, Double Packed	50

**SUTURE BOOTS WITH ADHESIVE TRAY**

Suture boots offer a secure way to hold delicate sutures without slippage or damage during a surgical procedure. The tray comes with an adhesive strip for securing and easy removal. Radio-opaque element ensures X-ray detection. Available in a Yellow colour.



30-702

DEF1093	2 per Packet	25
DEF1094	4 per Packet	25
30-702	6 per Packet	50
DEF1096	8 per Packet	25
30-701	10 per Packet	50
DEF2284	12 per Packet, Double Packed	25

**TUBING CONNECTOR**

DEF465	Connector Small 4mm – 7mm	100
--------	---------------------------	-----



DEF465

DEF466	Connector Medium 7mm – 10mm	100
--------	-----------------------------	-----



DEF466

**WOUND DRAINAGE**

Wounds drains and evacuators provide drainage for various surgical procedures, allowing for the removal of blood and excess fluid to improve patient outcome. Our wound drainage systems are made up of evacuators, flat drains, fluted drains, round drains, capillary action drains and easy flow drains. They are all made from 100% Medical Grade Silicone, which is inert, soft and flexible and designed to reduce the potential of post-operative infections.

**WOUND EVACUATOR BULB**

Our wound evacuator is a clear bulb designed to slowly draw fluid and blood from the wound site. It creates an ideal vacuum to a maximum of 125mmHg pressure and is compatible with our flat drains, fluted drains and round drains.



31300380DS

31300380DS	Evacuator Bulb 100ml/cc	10
------------	-------------------------	----

**WOUND FLAT DRAIN**

Soft and comfortable silicone flat drains with multiple drainage holes (36 x 2) to enhance the fluid removal in a wound. They have an ergonomic flat shape which enables an easier placement. X-ray detectable.



11100200780DS

11100200780DS	Flat Drain White 7mm x 190mm	10
11100201080DS	Flat Drain White 10mm x 190mm	10

**WOUND FLUTED DRAIN**

Fluted drains are made of soft, flexible silicone and available in a round and flat shape, with or without attached trocar. Designed with a multi channel to allow for great drainage and assisting in minimal tissue in-growth. The inner core of the H design offers drain strength and flexibility.

**FLAT FLUTED DRAIN**

The flat fluted drains are available as a white, radio-opaque X-ray detectable with a clear tube. Available in 7mm and 10mm.

11101310780DS	Flat white drain 3/4 fluted 7mm with Trocar and Hub	10
11101410780DS	Flat white drain fully fluted 7mm with Trocar and Hub	10
11101311080DS	Flat white drain 3/4 fluted 10mm with Trocar and Hub	10
11101411080DS	Flat white drain fully fluted 10mm with Trocar and Hub	10

**ROUND FLUTED DRAIN**

The round fluted drains come in 2 versions, the white radio-opaque with a clear tube and the hubless drain which is clear with a blue X-ray detectable line. Available in 10, 15, 19 & 24 Fr.

11003131080DS	Round clear full fluted drain hubless 10Fr with Trocar	10
---------------	--	----



11003131080DS



11003831080DS

11003831080DS	Round white full fluted drain with Hub 10Fr	10
11003131580DS	Round clear full fluted drain hubless 15Fr with Trocar	10
11003131980DS	Round clear full fluted hubless 19Fr with Trocar	10
11003811980DS	Round white full fluted drain hubless 19Fr	10
11003122480DS	Round clear full fluted hubless 24Fr	10

**ROUND CLEAR DRAIN**

Clear soft silicone drain 500:5 which consists of 5 side holes with a total drain length of 500mm. X-ray detectable.

11000301280DS	500:5 12Fr 2mm x 4mm	10
11000301880DS	500:5 18Fr 4mm x 6mm	10
11000302480DS	500:5 24Fr 5mm x 8mm	10
11000302780DS	500:5 27Fr 6mm x 9mm	10
11000303080DS	500:5 30Fr 7mm x 10mm	10
11000303680DS	500:5 36Fr 8mm x 12mm	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**MEDICAL CONSUMABLES**

**PASSIVE WOUND DRAINS**

**MULTILUMEN CAPILLARY ACTION DRAINS**

These clear radio-opaque drains have varying lumens according to drainage requirements. The lumens can be separated insitu and positioned as necessary. The drainage takes place both through, under and over the lumens.



DEF1750

DEF1110	Multilumen drain 2 Lumens x 12.5cm	10
DEF1076	Multilumen drain 5 Lumens x 25cm	10
DEF1024	Multilumen drain 9 Lumens x 25cm	10
DEF1750	Multilumen drain 16 Lumens x 20cm	10
DEF1751	Multilumen drain 16 Lumens x 40cm	10

**COTTON PRODUCTS**

**COTTON BALLS**

Packaged and sterilised for ease of entry into a sterile field. Ideal for prepping, wound cleaning, perineal care and fluid application.



DEF074

DEF725	Large 3cm, Packet of 4	100
DEF074	Large 3cm, Packet of 5	100

**COTTON WOOL**



DEF859

DEF859	10cm x 10cm	25
DEF615	10cm x 15cm Double Packed	100

**COTTON WOOL ROLL**

Soft and absorbent cotton wool roll used for padding in patient care. Sterilised and ready to be used.



DEF874

DEF874	15cm x 55cm, Double Packed	50
DEF875	30cm x 90cm, Double Packed	25
DEF880	38cm x 165cm, Double Packed	15

**SPONGES**

**SPONGES NASAL BOLSTER**

Used after nasal or sinus surgery to absorb any discharge or ooze from the nose. High quality non woven net filled with absorbent cotton padding, and stretch elastic sides to secure around the ears. Extremely soft and comfortable for post surgery patients. Latex free and available as sterile or non sterile.



DEF1382

DEF1381	Nasal Bolster	20
DEF1382	Nasal Bolster	100

**SPONGES NASAL BOLSTER – NON STERILE**

DEF1380	Nasal Bolster Bulk	20
---------	--------------------	----

**GAUZE SWABS PLAIN**

Sterile and ready to be used for a wide variety of applications, including as a dressing, for prepping procedures and wound cleaning. Available as a sterile Plain 8 ply gauze.



DEF181

DEF197	7.5 x 7.5cm, Plain Bulk Non Sterile	100
DEF183	7.5cm x 7.5cm, Packet of 3	100
DEF181	7.5cm x 7.5cm, Packet of 5	100
DEF371	7.5cm x 7.5cm, Packet of 20	50
DEF189	10cm x 10cm, Packet of 2	100
DEF186	10cm x 10cm, Packet of 3	100



DEF045

DEF045	10cm x 10cm, Packet of 5	100
DEF208	10cm x 10cm, Packet of 10	50
DEF213	10cm x 10cm, Packet of 10	300
DEF178	10cm x 10cm, Packet of 20	25
DEF177	10cm x 10cm, Packet of 25	25
DEF214	10cm x 10cm, Packet of 25	200
DEF1696	10cm x 10cm, Packet of 50	50
DEF1033	10cm x 20cm, Packet of 4	100
DEF198	10cm x 20cm, Packet of 5	250
DEF195	10cm x 20cm, Packet of 10	50
DEF853	10cm x 80cm, Packet of 4, 16 ply Folded	50

**PURCOTTON GAUZE SWABS PLAIN**

DEF3559	10cm x 10cm, Packet of 25	25
---------	---------------------------	----

**GAUZE SWABS PLAIN – NON STERILE**

DEF136	10cm x 10cm Plain Bulk	100
--------	------------------------	-----

**NONWOVEN GAUZE SWABS**

DEF1131	5cm x 5cm, Packet of 5	100
---------	------------------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**SWAB STICKS/COTTON TIPPED APPLICATORS**

Packaged and sterilised, ready to use for taking a swab or specimen from a patient. Paper, plastic or wooden shaft in different length sizes with a single or double ended soft cotton head.

**SWAB STICKS WITH WOODEN SHAFT**



DEF057

DEF057	7.5cm Single Ended, Packet of 3	100
DEF442	7.5cm Single Ended, Packet of 5	100
DEF439	15cm Single Ended, Packet of 2	100
DEF059	15cm Single Ended, Packet of 3	100
DEF443	15cm Single Ended, Packet of 5	100



DEF439

**SWAB STICKS WITH WOODEN SHAFT – NON STERILE**

DEF062	15cm Single Ended Bulk	1000
--------	------------------------	------

**SWAB STICKS WITH PAPER SHAFT**



DEF421

DEF424	7.5cm Double Ended, Packet of 5	100
DEF421	7.5cm Double Ended, Packet of 10	100
DEF104	7.5cm Double Ended, Packet of 50	80

**SWAB STICKS WITH PAPER SHAFT – NON STERILE**

DEF061	7.5cm Double Ended Bulk	1250
--------	-------------------------	------

**SWAB STICKS WITH PLASTIC SHAFT**



DEF058	7.5cm Double Ended, Packet of 3	100
--------	---------------------------------	-----

**SWAB STICK MOUTH & THROAT JUMBO**



DEF401

DEF400	15cm, Packet of 2	150
DEF401	15cm, Packet of 3	150
DEF402	15cm, Packet of 5	150
DEF455	20cm, Packet of 2	150
DEF407	20cm, Packet of 3	150
DEF1402	20cm, Packet of 4	100
DEF1593	20cm, Packet of 1	200



DEF407

**SWAB STICK MOUTH & THROAT JUMBO – NON STERILE**

DEF055	15cm Bulk	200
DEF054	20cm Bulk	200

**BURNS**

**BURNS DRESSING**

Specialised products for patients with burns. Non adherent, highly absorbent wound dressings.



DEF987

DEF199	10cm x 20cm Plain, Packet of 20	50
DEF968	10cm x 100cm Plain, Packet of 5	100
DEF987	40cm x 50cm Melolin, Double Packed	50
DEF679	50cm x 70cm, Double Packed	35
DEF3570	50cm x 175cm Melolin, Double Packed	25
DEF974	90cm x 2m Plain	50



DEF679

**BURNS GAUZE**

Soft roll suitable for wounds and burns.

DEF965	10cm x 80cm, Packet of 16	50
DEF2436	25cm x 90cm, Packet of 2	350
DEF1039	60cm x 90cm, Packet of 6	100
DEF342	90cm x 2m	50
DEF344	90cm x 2m, Packet of 4	25
DEF1662	90cm x 5m	10

**BURNS PLASTIC BAG**

DEF650	305mm x 390mm	100
--------	---------------	-----

**SANITARY DISPOSAL BAG**

DEF709	Packet of 100, Non Sterile	20
--------	----------------------------	----



DEF709

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

**SILICONE TUBING ROLL – NON STERILE**

Excellent range of medical grade silicone tubing available in a variety of different metric sizes. Available as non sterile rolls of 15 or 25 metres.



2110030239

2110030239	0.30mm (ID) x 0.64mm (OD) 15m Roll
2110050134	0.5mm (ID) x 1mm (OD) 25m Roll
203005009	0.51mm (ID) x 0.94mm (OD) 15m Roll
203006011	0.64mm (ID) x 1.19mm (OD) 15m Roll
2110070439	0.7mm (ID) x 1.6mm (OD) 15m Roll
201010020	1mm (ID) x 2mm (OD) 25m Roll
2110100234	1mm (ID) x 3mm (OD) 25m Roll
2110100439	1.016mm (ID) x 2.159mm (OD) 15m Roll
203010021	1.02mm (ID) x 2.16mm (OD) 15m Roll
203014019	1.47mm (ID) x 1.96mm (OD) 15m Roll
2110150434	1.5mm (ID) x 3mm (OD) 25m Roll
2110150939	1.57mm (ID) x 2.41mm (OD) 15m Roll
2110151039	1.57mm (ID) x 3.18mm (OD) 15m Roll
2110190239	1.98mm (ID) x 3.18mm (OD) 15m Roll
201020030	2mm (ID) x 3mm (OD) 25m Roll
2110200334	2mm (ID) x 4mm (OD) 25m Roll
2110200534	2mm (ID) x 6mm (OD) 25m Roll
2110260239	2.65mm (ID) x 4.88mm (OD) 15m Roll
2110300234	3mm (ID) x 5mm (OD) 25m Roll
2110300334	3mm (ID) x 6mm (OD) 25m Roll
2110300434	3mm (ID) x 7mm (OD) 25m Roll
2110310239	3.18mm (ID) x 6.35mm (OD) 15m Roll
2110330239	3.35mm (ID) x 4.65mm (OD) 15m Roll
2110400434	4mm (ID) x 6mm (OD) 25m Roll

2110400634	4mm (ID) x 7mm (OD) 25m Roll
2110400734	4mm (ID) x 8mm (OD) 25m Roll
2110470139	4.76mm (ID) x 7.93mm (OD) 15m Roll
20304707901	4.77mm (ID) x 7.9mm (OD) 15m Roll
20304709501	4.77mm (ID) x 9.50 (OD) 15m Roll
2110500234	5mm (ID) x 7mm (OD) 25m Roll
2110500334	5mm (ID) x 8mm (OD) 25m Roll
2110500534	5mm (ID) x 9mm (OD) 25m Roll
2110500734	5mm (ID) x 10mm (OD) 25m Roll
2110600234	6mm (ID) x 9mm (OD) 25m Roll
2110600334	6mm (ID) x 10mm (OD) 25m Roll
2110600434	6mm (ID) x 11mm (OD) 25m Roll
2110600534	6mm (ID) x 12mm (OD) 25m Roll
2110630139	6.4mm (ID) x 9.5mm (OD) 15m Roll
2110630239	6.4mm (ID) x 11mm (OD) 15m Roll
2110630339	6.4mm (ID) x 12.7mm (OD) 15m Roll
2110700334	7mm (ID) x 10mm (OD) 25m Roll
2110700534	7mm (ID) x 11mm (OD) 25m Roll
2110700634	7mm (ID) x 12mm (OD) 25m Roll
201070130	7mm (ID) x 13mm (OD) 25m Roll
2110790239	7.9mm (ID) x 12.7mm (OD) 15m Roll
2110800334	8mm (ID) x 12mm (OD) 25m Roll
2110800534	8mm (ID) x 14mm (OD) 25m Roll
2110900334	9mm (ID) x 13mm (OD) 25m Roll
2110950139	9.5mm (ID) x 12.7mm (OD) 15m Roll
2110950239	9.5mm (ID) x 14.3mm (OD) 15m Roll
2110950339	9.5mm (ID) x 15.9mm (OD) 15m Roll
2111000334	10mm (ID) x 14mm (OD) 25m Roll
2111000434	10mm (ID) x 15mm (OD) 25m Roll
201100160	10mm (ID) x 16mm (OD) 25m Roll
201100180	10mm (ID) x 18mm (OD) 25m Roll
2111200334	12mm (ID) x 17mm (OD) 25m Roll
2111250134	12.5mm (ID) x 19mm (OD) 25m Roll
2111270239	12.7mm (ID) x 19.05mm (OD) 15m Roll