



Quality  
Products

Trusted Experts

## Custom Procedure Packs

A Defries Custom Procedure Pack is designed by the healthcare professional and the Defries representative to contain all the consumables used during a surgical procedure



# Defries Custom Procedure Packs

Defries Custom Procedure Packs are easy to use, complete and ready on the shelf for use in the operating theatre. Custom Procedure Packs save set up time and ensure all components are readily available for every procedure. With all components packed together, packaging waste is reduced and infection control maximised.

## Cost Effective and Time Efficient

By using a Custom Procedure Pack in the operating theatre, the amount of time required for set up is reduced by up to 15-20 minutes compared with using individually packed consumables. The sterile Custom Procedure Packs are self-contained with all components ready to use when opened.

This provides efficiencies and saves time, unlike individually wrapped consumables that are touched many times prior to use and require opening separately. Also, the Custom Procedure Packs ensure that every consumable is on hand and allays issues of stock unavailability in the operating theatre



## Infection Control

A Defries Custom Procedure Pack maximises infection control. The inclusion of our single use drapes and gowns in our packs minimise lint that can occur from reusable linen.

# Defries Custom Procedure Packs in 9 Simple Steps



# Trusted Experts

Defries Industries have been in business since 1981 and have been manufacturing Custom Procedure Packs for over 15 years. Our reputation for quality and consistency have made us a number 1 choice for supplying custom packs throughout Australia. We have won numerous tenders and contracts to supply our packs to healthcare providers.

Our Quality Assurance team is integral to Defries Industries. Manufacturing Custom Procedure Packs is overseen by our strict quality management systems. Our manufacturing department, procurement specialist team and sales team also play a vital role in the delivery of high quality products.



## Clean Room

Our new state-of-the-art cleanroom with innovative technology, more than doubles our former manufacturing and processing capacity. In this streamlined facility, our raw materials storage will directly feed into the cleanroom via a purpose built decanting area designed to accommodate both custom procedure packs and consumables. This will provide us with greater efficiencies and further opportunities for growth and development throughout Australia and New Zealand. It will provide a controlled, contaminant free environment and meet certification Class ISO8 in accordance with ISO 14644-1 international cleanroom standards.



## Locally Made

Custom Procedure Packs are packed into Tyvek pouches and sterilised at Steritech with ethylene oxide (EtO) gas. All Custom Procedure Packs undergo post sterilisation checks by the QA team. The cartons are then stored in our warehouse ready for dispatch.

Defries Custom Procedure Packs are manufactured in Australia. Our local company ensures quality made products and supports the Australian manufacturing workforce.



# Delivery and VMI Service

Defries Industries offers a Vendor Managed Impress or VMI service for Custom Procedure Packs to hospitals. VMI involves our employees monitoring usage of packs on a daily basis, delivering cartons directly to the operating theatre, and unpacking and placing Custom Procedure Packs directly onto the theatre shelves ready for picking and use.

A monthly consolidated invoice is generated for the hospital. This service has been implemented and well received at a number of our customers' sites. We are happy to discuss requirements that would best suit your organisation. This service delivers many benefits to hospitals as staff are then alleviated from monitoring stock, placing orders and restocking.



## Quality Products

Defries' Custom Procedure Packs meet all stringent Australian standards and are delivered to hospitals ready for use.

Each pack has a check sheet enclosed for easy component identification:

- Clear labelling. Easy to identify the name of pack, company details and include a component list. A column is included so components can be checked by theatre staff.
- Removable patient labels for easy tracking. These include the date of manufacture, expiry date, reorder code and batch number.
- Bar coding will be included on patient labels.
- Full traceability in case of a product recall.
- Check sheet meets all TGA labelling requirements.
- EO Sterilisation Dot visually confirms pack is sterile.

# Environmental Choices

Environmental sustainability is very important to our company.  
We are implementing the following initiatives:

- Our new head office and cleanroom facilities are built to incorporate the most advanced technologies according to the environmental sustainability guidelines.
- Energy efficiencies, water efficiencies, waste reduction, and operations and maintenance optimisation has been incorporated into these facilities.
- All Defries Industries disposable bowls, kidney dishes and gallipots are recyclable.
- All SMS material is recyclable through Replas.
- Our cartons and plastic inserts are 100% recyclable.
- Tyvek packaging is recyclable.
- We incorporate an environmentally proactive transport strategy.



## Talented Team

Defries Industries has a team of trusted professionals ranging from customer service representatives, account managers, sales managers and a product specialist. Our staff are available to provide an overview of our facility and cleanroom, and discuss our 34 year reputation working in the Australian and New Zealand hospital industry.

Our product knowledge on drapes, gowns and specialised consumables combined with the expertise of our team ensures a smooth, easy development process. As our custom packs are made and sterilised right here in Australia, we have

a quick turnaround time for new packs and any adjustments or changes can be easily arranged. We also support the Australian manufacturing industry by employing local people in a variety of positions.

Defries Industries is a local, community driven company and we thank you for your ongoing support to our 100% owned Australian medical company.

**Please contact your Sales Account Manager to discuss your Custom Procedure Pack requirements.**

### Head Office

1-7 Vision Street, Dandenong South VIC 3175

**freecall:** 1300 550 278 **t:** 61 3 9706 3600

**f:** 61 3 9706 3619

**e:** [sales@defries.com.au](mailto:sales@defries.com.au) **w:** [defries.com.au](http://defries.com.au)

### Manufacturing Site

18-24 Villas Road, Dandenong South VIC 3175