



September 2026



Defries Industries is a 100% Australian owned family medical company. As a world class manufacturer, our vision has remained clear over the years.

"WE ARE COMMITTED TO IMPROVING PATIENT CARE BY OFFERING QUALITY SINGLE-USE STERILE MEDICAL PRODUCTS AND CUSTOMISED PACKS TO HEALTHCARE PROFESSIONALS."



CONTENTS

Introduction	2
Clinical Drapes and Equipment Covers	4
Consumables	13
Custom Procedure and Standard Packs	36
Hollowware	50
Instruments and Instrument Sets	57
Maternity Products	73
Medical Imaging and Cath Lab Products	76
Medical Labels	80
Personal Protective Apparel	82
Surgical Gowns	85
Wound Care	91
Index	101

INTRODUCTION

The Defries Industries team share a vision of high standards and a commitment to superior patient care.

We continually strive to create outstanding products that excel in design, ease of use, effectiveness and price.

Defries has built a strong national reputation since 1980 and continues to expand and develop our product range. We are proud to be a 100% Australian family owned medical business and aim for ongoing growth and investment in Australian and New Zealand healthcare into the future.



WE ARE COMMITTED TO IMPROVING PATIENT CARE WITH OUR QUALITY SINGLE-USE, STERILE MEDICAL PRODUCTS AND CUSTOMISED PROCEDURE PACKS

Our range of over 3000 different medical and surgical products are a combination of single sterile items and procedure packs. This catalogue features all our single sterile items and standard general packs. For customised procedure packs, please contact our Customer Service team for your local Representative's details.

AT DEFRIES WE COLLABORATE WITH INDUSTRY LEADING SUPPLIERS

We regularly meet with our primary suppliers to explore new technologies and products available, so that we can offer our customers quality products at competitive prices. We visit our suppliers' facilities to check their Quality Management Systems and ensure the ongoing supply of components to Defries.

We are proud to distribute innovative products from Advancis Medical, DeRoyal Industries, OBP Medical, Degania Medical and NewGen Surgical. These medical companies are leaders in healthcare from around the world.

REGULATORY COMPLIANCE IS VITALLY IMPORTANT AT DEFRIES INDUSTRIES

Our quality management system, which is based on ISO13485:2016, is regularly assessed by the TGA to ensure ongoing compliance with this ISO standard and the relevant provisions of the Therapeutic Goods (Medical Devices) Regulations 2002.

The key role of our Quality Assurance is to deliver products that are superior in safety, quality and performance by achieving ongoing compliance with approved specifications and relevant industry standards.



SUSTAINABILITY AND ENVIRONMENTAL DECISIONS

We are committed to managing our resources actively and wisely to prevent pollution, minimise our environmental impact and reduce our carbon footprint. We regularly review how the design, manufacture and distribution of Defries products affects the environment. Our corporate sustainable strategy includes environmental values, impacts and commitments which include the following:

- Energy use and Carbon emissions
- Ethical company suppliers
- Waste management
- Water resources
- Transport and fuel emissions
- General and end user waste management

Our company is certified to the ISO14001 Environmental Management System.



CLINICAL DRAPES AND EQUIPMENT COVERS



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

EQUIPMENT COVERS

BACK TABLE COVER

DEF5016AQ	152cm x 228cm (60" x 90") Aqua	40
DEF3544	150cm x 190cm Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3545	150cm x 230cm Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3678	150cm x 300cm Sapphire, Heavy Duty, Impervious and Absorbent	12
DEF5802	152cm x 228cm Aqua, Zoned, Heavy Duty	20

CHAIR COVER

		
DEF772		
DEF772	2 x 57.5cm x 62cm with 2 x Latex-free Bands #33	100

CAMERA CORD COVER

		
DEF1201		
DEF1201	Camera Cord Cover 15cm x 230cm with Cuff Sealed End & Tapes, Double Packed	50

DRILL SLEEVE COVER

		
DEF889		
DEF1224	5cm x 80cm with 3 Labels	100
DEF1368	5cm x 110cm with Tape	100
DEF1516	7.5cm x 100cm with Adhesive both Ends	100
DEF1154	7.5cm x 150cm with 2 x Latex-free Bands #31	100
DEF889	10cm x 150cm	100
DEF618	10cm x 152cm with Tape Double Packed	25

C-ARM DRAPE

Large, one-piece clear poly cover for complete coverage of C-Arm with white poly adhesive strips to secure the drape in place. Packaged with a telescoping fold. Adhesive straps included to secure excess drape material.

DEF5911	104cm x 188cm with Bands	24
---------	--------------------------	----

O.R. TABLE - C-ARM PROTECTOR DRAPE

Made from transparent polyethylene, this drape securely attaches to the patient drape on the operating table to help maintain a sterile field during C-Arm use. The drape is not applied to the C-Arm itself, as the imaging equipment is independently covered with the user's choice of C-Arm cover.

DEF5927	137cm x 100cm	40
---------	---------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HANDLE SLEEVE COVER



DEF1053

DEF1053	14cm x 30cm with Adhesive	100
DEF1466	14cm x 30cm with Adhesive Packet of 2	100

PLASTIC SLEEVE COVER



DEF454

DEF454	15cm x 100cm, with 2 x Latex-free Bands #31, Double Packed	100
DEF1744	15cm x 200cm, with 1 x Latex-free Band #31 and Sealed End, Double Packed	100

MAYO STAND COVER

DEF5019AQ	61cm x 135cm, Aqua	30
DEF761	62cm x 91cm Clear Plastic with Sealed End	100
DEF775	70cm x 137cm, Blue Plastic with Cuff	50
DEF6100	80cm x 145cm, Sapphire, Heavy Duty, Blue Plastic with Reinforcement	25
DEF4232	80cm x 165cm Fully Reinforced Blue PE Film with Absorbent Material	15

PROTECTIVE COVER

DEF1276	14cm x 30cm with a Blank Adhesive Label	100
DEF647	20cm x 20cm, Plastic	100
DEF650	30.5cm x 39cm, Plastic	100
DEF763	104cm x 150cm, Clear Plastic	40
DEF1573	150cm x 208cm, Clear Plastic	50

TROLLEY COVER



DEF641

DEF641	Trolley Bag Cover 61cm x 102cm	50
DEF1239	Shower Trolley Drape 2m x 3m	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------



CUSTOMISED PROCEDURE DRAPES

THE FLEXIBILITY OF CHOICE

At Defries we have a comprehensive range of patient drapes which cover all procedures in the operating room. We provide our customers with superior products that match their clinical needs.

Our drapes are made to specific requirements such as size, fabric, fenestrations and extra details like velcro taping and clear plastic anaesthetic screens.

CHOICE OF FABRIC IS IMPORTANT

At Defries we use a variety of drape fabrics:

Azure drapes use a viscose impervious fabric with a polyethylene laminate, suitable for general surgery with minimal fluid or short procedures. This range includes the Azure minor procedure drapes.

Sapphire drapes use a bi-laminate or tri-laminate fabric with extra reinforcement providing excellent barrier protection. This range includes back table covers, mayo covers, some procedure drapes and bespoke designed drapes to meet specific customer requirements.

Aqua drapes are our newest drape range available and use an innovative drape fabric that is a tri-laminate – a three layer drape which offers an absorbent top level, an impervious middle coating and a soft and comfortable patient layer. Soft and light, this material feels the closest to linen when draping the patient.

This high performance fabric is soft and drapeable with additional features such as reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. We have an excellent standard range of designs covering all procedures from general surgery and universal packs, to orthopaedics and specialised surgery such as neuro.

Supreme drapes are made from a tri-laminate fabric and are specifically designed for high level performance whilst being soft, drapeable and easy to use. These drapes provide great absorbency from fabric reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. The Supreme range includes drape sheets and specialty procedure drapes.

Our drape design includes the instructions for correct aseptic application with markings and directions on the drape to assist ease of use.

All drape fabrics are conformable to patients' shapes. They are low linting and have excellent barrier properties to prevent strike through of airborne and blood borne bacteria, meeting regulatory standards.

Please contact your Defries representative for further information.

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DRAPES FOR EVERY PROCEDURE

Our wide range of drapes include:

- **ANCILLARY**

Absorbent Pad, Adhesive Drape Towel, Drape Tape

- **ANGIOGRAPHY & CATH LAB**

Angio, Femoral Radial Angio, Biliary, Hickman, Pacemaker, Radial

- **CARDIOVASCULAR**

Cardiac & Cardiac Splits, CV Groin, Cardiac Saddle

- **ENT**

Minor Procedure, Oral, Turban, Split, Paediatric ENT Split

- **EQUIPMENT COVERS**

Back Table Cover, Mayo Stand Covers, Utility

- **GENERAL DRAPES**

Drape Sheets, Fenestrated Minor Procedure, Laparoscopy, Laparotomy, PICC/Arterial, Central Line, Universal Drape Packs

- **GYNAECOLOGY & OBSTETRICS**

Caesarean, Endoscopy, Laparoscopy, Lithotomy, Leggings, Under Buttocks

- **NEPHROLOGY**

PERC

- **NEURO SPINAL**

Neuro-craniotomy, Isolation, Transsphenoidal

- **PAEDIATRICS**

Neonatal Split, Neonatal, Paediatric Drape Sheets

- **OPHTHALMIC**

Eye

- **ORTHOPAEDICS**

Extremity Foot, Hand, Hip, Knee, Shoulder Arthroscopy, Shoulder Extremity, Shoulder Beach Chair, Split & U, Stockinette, Vertical Isolation

- **UROLOGY**

TUR

GENERAL DRAPES AND PACKS

CLEAR PLASTIC DRAPE



DEF742

DEF742	45cm x 61cm	200
DEF079	61cm x 91cm	200
DEF749	70cm x 120cm	100
DEF745	76cm x 106cm	100
DEF3910	76cm x 152cm, Double Packed	100
DEF151	91cm x 115cm	100
DEF2430A	Procedure Adhesive Drape 91cm x 115cm with 2.5cm x 115cm adhesive down length of drape	50
DEF748	120cm x 140cm	100
DEF752	120cm x 140cm, Double Packed	100
DEF5713	200cm x 200cm	100

EPIDURAL DRAPE



DEF1481

DEF1481	76cm x 106cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100
DEF748FA-S	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	10
DEF3394	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

UNIVERSAL DRAPE PACKS

AQUA RANGE – LONG SIDES

DEF5400	1 x Top Drape, 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drapes, 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 80cm x 145cm with Reinforcement Wrapped in a Back Table Cover 127cm x 230cm with Reinforcement	5
---------	---	---

AQUA RANGE – HEAVY DUTY - MEDIUM SIDES

DEF5800A	1 x Top Drape, 150cm x 276cm with Reinforcement and Adhesive 2 x Side Drapes, 99cm x 145cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 193cm x 193cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 80cm x 145cm with Reinforcement 1 x Paper Disposal Bag Wrapped in a Back Table Cover 152cm x 248cm with Reinforcement	6
----------	--	---

AQUA RANGE – HEAVY DUTY - LONG SIDES

DEF5801-4	1 x Top Drape, 150cm x 276cm with Reinforcement and Adhesive 2 x Side Drapes, 105cm x 195cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 193cm x 193cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 80cm x 145cm with Reinforcement 1 x Paper Disposal Bag Wrapped in a Back Table Cover 152cm x 248cm with Reinforcement	4
-----------	---	---

AQUA RANGE – LONG SIDES WITHOUT MAYO STAND COVER

DEF5243AQ	1 x Top Drape, 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drapes, 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement and Adhesive Wrapped in a Back Table Cover 127cm x 230cm with Reinforcement	5
-----------	--	---

AZURE MINOR PROCEDURE DRAPE

Soft viscose fabric combined with a polyethylene laminate to create a drape ideal for minor procedures.

DEF5094	38cm x 45cm Circumcision Drape with a 3cm Circular Fenestration	50
DEF5090	50cm x 60cm with a 6.5cm Circular Fenestration and Full Adhesive	25
DEF5201	90cm x 90cm Fenestrated Drape with a 15cm x 20cm Rectangular Fenestration	25
DEF5091	91cm x 91cm with a 10cm x 15cm Oval Fenestration and Full Adhesive	40
DEF2923	112cm x 153cm ICU Procedure Drape with a 8cm x 10cm Rectangular Fenestration and Full Adhesive	10
DEF5203AQ	120cm x 160cm with 10cm x 15cm Fenestration and Reinforcement	50
DEF5092	200cm x 300cm with a 15cm x 25cm Fenestration and Full Adhesive	12


AQUA DRAPE SHEETS

Innovative tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.

DEF5141AQ	Square Drape, 100cm x 100cm	40
DEF808823	100cm x 150cm	27
DEF5025AQ	100cm x 180cm	20
DEF5024AQ	110cm x 150cm	20
DEF5026AQ	150cm x 200cm	20
DEF5028AQ	180cm x 215cm	20
DEF5030AQ	Top/Head Drape 135cm x 250cm with Reinforcement and Adhesive	15
DEF5828	Top/Head Drape 150cm x 276cm with Reinforcement Adhesive and Tube Holders	15
DEF5055A	Top/Head Drape 180cm x 300cm with Reinforcement Adhesive and Cord Tabs	15
DEF5031AQ	Side Drape 100cm x 130cm with Reinforcement and Adhesive	40
DEF5023AQ	Side Drape 100cm x 195cm with Reinforcement and Adhesive	20
DEF5032AQ	Bottom/Foot Drape 190cm x 185cm with Reinforcement and Adhesive	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

UTILITY DRAPE



DEF5130

DEF5130	65cm x 65cm with Adhesive wrapped in a 91cm x 91cm SMS	50
---------	--	----

SMS DRAPE

SMS fabric is a 3 layer polypropylene material bonded together to form a drape that is fluid repellent. Ideal for minor procedures.



DEF1594

DEF1122	30.5cm x 61cm	100
DEF1594	61cm x 61cm	50
DEF3806NS	61cm x 61cm, Non Sterile	100
DEF2300A	91cm x 91cm	100
DEF3807NS	91cm x 91cm, Non Sterile	100
DEF3805NS	150cm x 200cm, Non Sterile	50

POLYBACKED DRAPE

An absorbent drape with a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor procedures.



DEF833

DEF833	Polybacked Drape 33cm x 46cm White	750
DEF1287	Polybacked Drape 40cm x 60cm Blue	100
DEF1288A	Polybacked Drape 40cm x 60cm with a 10cm Circular Fenestration Blue	100
DEF2072	Polybacked Procedure Drape 75cm x 105cm with a 10cm x 15cm Fenestration and a 5cm Adhesive length along the 10cm sides of Fenestration, Blue	25

STOCKINETTE



DEF5060AQ

9503-36	7.6cm x 91cm (3" x 36"), Single Ply	20
DEF5060AQ	20cm x 96cm, Medium, Aqua	25
DEF5061AQ	23cm x 122cm, Large, Aqua	25
DEF5062AQ	30cm x 122cm, XLarge, Aqua	25

LEGGINGS

DEF2481AQ	Aqua Lithotomy Leggings 80cm x 140cm, Packet of 2	26
DEF5681	Clear Lithotomy Leggings 82cm x 125cm, Packet of 2	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SUPPLEMENTARY

HAND DRYING AND STERILE FIELDS

BLUE HUCK TOWEL

High absorbency cotton towels.



DEF1300

DEF1976	40cm x 68cm, Blue	30
DEF1300	40cm x 68cm, Blue	100
DEF2018	40cm x 68cm, Blue, Packet of 5, Double Packed	40

GREEN HUCK TOWEL



DEF2619

DEF2619	40cm x 68cm, Green, Packet of 2	100
DEF2140	40cm x 68cm, Green, Packet of 5	40
DEF5556	60cm x 90cm, Green	100

WHITE HUCK TOWEL

DEF2753	40cm x 68cm, White	30
---------	--------------------	----



DEF2753

PURCOTTON TOWEL

This PurCotton towel is a superior alternative to traditional Operating Room hand towels. They are a non-woven fabric made with a natural cotton, have a soft and comfortable feel, high liquid absorbency and are 100% biodegradable. They are also low linting and have an excellent wet-strength.



DEF3536A

DEF3536A	42cm x 68cm, Blue	100
DEF3674	90cm x 100cm, White	100



DEF5824

DEF5824	90cm x 100cm, White Non Sterile	100
---------	---------------------------------	-----

PURCOTTON DRAPE



DEF2521A

DEF2521A	61cm x 61cm, Blue, Packet of 4	100
DEF2404A	76cm x 106cm, Green	30
DEF2412	76cm x 106cm, Green, Packet of 2	50
DEF2663D	120cm x 140cm with a 7cm circular fenestration	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.																																				
DESSING TOWEL (DETOWEL) Made from airlaid paper with cotton yarn.  DEF5612 <table><tr><td>DEF5612</td><td>30cm x 30cm</td><td>100</td></tr><tr><td>DEF5614</td><td>30cm x 30cm</td><td>1000</td></tr><tr><td>DEF5611</td><td>30cm x 30cm, Packet of 2</td><td>50</td></tr><tr><td>DEF5613</td><td>30cm x 60cm</td><td>100</td></tr><tr><td>DEF5610</td><td>30cm x 60cm, Packet of 2</td><td>50</td></tr></table> DESSING D TOWEL Made from rayon and polyester.  DEF656 <table><tr><td>DEF656</td><td>32cm x 34cm, Packet of 2</td><td>50</td></tr><tr><td>DEF831</td><td>32cm x 34cm, Packet of 2</td><td>300</td></tr><tr><td>DEF832</td><td>32cm x 34cm, Packet of 10</td><td>100</td></tr><tr><td>DEF894</td><td>34cm x 60cm</td><td>50</td></tr><tr><td>DEF471A</td><td>34cm x 60cm</td><td>300</td></tr><tr><td>DEF910</td><td>34cm x 60cm, Packet of 5</td><td>25</td></tr></table>			DEF5612	30cm x 30cm	100	DEF5614	30cm x 30cm	1000	DEF5611	30cm x 30cm, Packet of 2	50	DEF5613	30cm x 60cm	100	DEF5610	30cm x 60cm, Packet of 2	50	DEF656	32cm x 34cm, Packet of 2	50	DEF831	32cm x 34cm, Packet of 2	300	DEF832	32cm x 34cm, Packet of 10	100	DEF894	34cm x 60cm	50	DEF471A	34cm x 60cm	300	DEF910	34cm x 60cm, Packet of 5	25	ADHESIVE DRAPE TOWEL  DEF5701 <table><tr><td>DEF5701</td><td>48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS</td><td>15</td></tr></table>			DEF5701	48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS	15
DEF5612	30cm x 30cm	100																																							
DEF5614	30cm x 30cm	1000																																							
DEF5611	30cm x 30cm, Packet of 2	50																																							
DEF5613	30cm x 60cm	100																																							
DEF5610	30cm x 60cm, Packet of 2	50																																							
DEF656	32cm x 34cm, Packet of 2	50																																							
DEF831	32cm x 34cm, Packet of 2	300																																							
DEF832	32cm x 34cm, Packet of 10	100																																							
DEF894	34cm x 60cm	50																																							
DEF471A	34cm x 60cm	300																																							
DEF910	34cm x 60cm, Packet of 5	25																																							
DEF5701	48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS	15																																							

CONSUMABLES



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

THEATRE CONSUMABLES

All surgical consumables are pre-packed and sterile, ready to be used in the operating room.

LATEX-FREE BANDS

ID indicates Inner Diameter. BW indicates Band Width.



DEF624A

DEF1590A	#16 – (ID) 60mm x 1.5mm (BW) Salmon Colour, Packet of 2, Double Packed	100
DEF1793A	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 2	25
DEF3542	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 2	100
DEF3560	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 3	100
DEF669A	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 5 Double Packed	25
DEF5797	#33 – (ID) 90mm x 3mm (BW) Blue, Packet of 10	50
DEF619	#34 – (ID) 100mm x 3mm (BW) Clear, Packet of 5	25
DEF600	#34 – (ID) 100mm x 3mm (BW) Clear, Packet of 10	100
DEF624A	#64 – (ID) 90mm x 6mm (BW) Salmon Colour, Packet of 1	100
DEF5528	#64 – (ID) 90mm x 6mm (BW) Clear, Packet of 5	100
DEF626A	#64 – (ID) 90mm x 6mm (BW) Salmon Colour, Packet of 5 Double Packed	100

CLAMPS/CLIPS/HOLDERS

TUBE CLAMPS

DEF1575A	Tube Clamp with Lock Blue Plastic 12cm	50
----------	--	----



DEF1575A

DEF2190	Towel Clamp Green Plastic 11.5cm	25
---------	----------------------------------	----



DEF2190

DEF2488	Towel Clip Stainless Steel 11cm Sharp/Sharp	10
---------	---	----



DEF2488

DEF2298	Tube Clamps 13cm Smooth Jaw	10
---------	-----------------------------	----



DEF2298

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CLAMPS/CLIPS/HOLDERS (CONTINUED)

DRAPE CLIP

General purpose, reusable Drape/Peach Clip, used during surgery for clamping purposes to ensure items remain where placed. A wide easy to use grip surface ensures heavy duty holding power.



DEROYAL-DRAPE-CLIP

DEROYAL-DRAPE-CLIP	Assorted Colours Non Sterile	1
--------------------	------------------------------	---



UTILITY-CLIP-UC-1-P, UTILITY-CLIP-UC-1-R

UTILITY-CLIP-UC-1-P	Utility Clip Pink Non Sterile	1
UTILITY-CLIP-UC-1-R	Utility Clip Red Non Sterile	1
DEF6030	Utility Clip Pink Sterile	50
DEF6031	Utility Clip Red Sterile	50

TUBE HOLDER

Ideal for securing and holding tubes during a procedure. Adhesive backing for securement onto a drape or towel. Adjustable velcro to allow gentle and easy placement.



DEF1889

DEF1889	Tube Holder with Velcro Adhesive 7.5cm x 7.5cm	25
---------	---	----

COMBINE DRESSING

Highly absorbent layer of cellulose blend enclosed in a soft and conformable non woven fabric. Each non woven combine dressing is pre-cut and each end is sealed into a variety of sizes.



DEF711

DEF711	9cm x 20cm	25
DEF719	9cm x 20cm, Packet of 5	100
DEF720	9cm x 20cm, Packet of 10	50
DEF3727	10cm x 15cm	25
DEF697	10cm x 90cm	25
DEF715	20cm x 20cm	25
DEF80	20cm x 60cm	100
DEF1160	20cm x 80cm	50

COMBINE ROLL

DEF714	Non Woven Combine Roll 10cm x 10m, Double Packed	16
DEF1004	Non Woven Combine Roll 20cm x 2.5m	25

COMBINE SPECIALISED



DEF804

DEF804	Non Woven Amputation Dressing 50cm x 70cm	50
--------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DENTAL

BITE BLOCK

100% cotton gauze with a barium sulfate detectable thread.



DEF1828A

DEF1828	Bite Block 7.5cm X-ray with Tail	100
DEF1246	Bite Block 7.5cm X-ray with Tail Non Sterile	100
DEF1232	Bite Block 7.5cm X-ray without Tail Non Sterile	100

TOOTHBRUSH



DEF690

DEF690	Toothbrush	25
--------	------------	----

DENTAL PRESSURE PACK



DEF799

DEF799	X-ray with 30cm Cotton Tail	100
DEF1063	X-ray with 35cm Tape, Packet of 2	100
DEF144	3cm x 4cm X-ray Small with 10cm tail, Packet of 2	50



DEF144

DEF145	4cm x 5cm X-ray Large with 9cm tail, Packet of 2	50
--------	--	----



DEF145

DEF1209	Dental Roll 4mm x 35mm, Packet of 5	100
---------	-------------------------------------	-----



DEF1209

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DEFOGGER

Specially formulated solution that prevents lens fogging. Comes with a non-abrasive, x-ray detectable foam pad. Sterile.



DEROYAL-28-0101


DEROYAL-28-0101	DeFogger Oval Bottle	20
-----------------	----------------------	----

ELECTROSURGICAL

DIATHERMY TIP CLEANER

This coated scratch pad cleaner quickly and easily removes debris from instrument tips and diathermys. The adhesive backing will secure the Tip Cleaner onto draped surfaces and contains a radiopaque strip to ensure X-ray detectability.

DEF1886	Tip Cleaner 5cm x 5cm with Adhesive and Radiopaque Strip	50
---------	--	----



DEF1886

EYE

EYE SPEAR

Highly absorbent.



DEF484

DEF484	Eye Spear, Packet of 2	100
DEF485	Eye Spear, Packet of 5 Double Packed	20
DEF486	Eye Spear, Packet of 6 Double Packed	20
DEF487	Eye Spear, Packet of 10 Double Packed	20

EYE STROLL



DEF112

DEF112	Eye Stroll, Packet of 5	50
DEF081	Eye Stroll, Packet of 10 Double Packed	50

QUIVER

A single use disposable polypropylene Quiver used to hold diathermy during surgical procedures. Quiver is recyclable (metal ring needs to be removed before being placed in recycling container).

DEF2772	5cm x 25cm with Metal Ring	90
---------	----------------------------	----



DEF2772

DEF3068	5cm x 25cm with Cable Tie	90
---------	---------------------------	----



DEF3068

DEF3291	9cm x 35cm with Metal Ring	40
DEF3296	9cm x 35cm with Cable Tie	40
DEF3905	4.5cm x 15.5cm Neuro	100
DEF3906	4.5cm x 20.5cm Neuro	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FOAM PAD DRESSING

Polyurethane foam pads for use in a clinical setting.



DEF792

DEF792	100mm x 100mm x 13mm	100
DEF769	100mm x 100mm x 25mm	50
DEF783	100mm x 180mm x 13mm	50
DEF784	100mm x 180mm x 25mm	40
DEF738	150mm x 150mm x 13mm	50



DEF738

DEF1738	150mm x 315mm x 25mm Non Sterile	50
DEF914	150mm x 200mm x 25mm	25
DEF1549	200mm x 25mm x 25mm	100



DEF1549

DEF913	200mm x 300mm x 25mm	25
DEF739	250mm x 250mm x 13mm	50
DEF786	300mm x 600mm x 13mm Double Packed	20

PREP APPLICATORS

DEF921	Prep Foam 50mm x 60mm x 25mm White No Handle	100
--------	---	-----



DEF921

DEF1846	Prep Foam on a Handle 15cm Blue	100
---------	------------------------------------	-----



DEF1846

DEF1847	Prep Foam on a Handle 15cm Blue, Packet of 2	100
DEF2035	Prep Foam on a Handle 20cm Orange	100



DEF2035

PREP SETS

DEF3106	PREP SET 2 x Gallipot 150ml 3 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow Wrapped in a 91cm x 91cm Field	20
DEF3458A	PREP SET 5 x Gauze Swabs 10cm x 10cm XRD 2 x Sponge Holder 19cm Plastic 2 x Gallipot 150ml 1 x Kidney Dish 700ml Clear 1 x Label Sheet - 1 Set of 8 Wrapped in a 91cm x 115cm Field	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

IDENTILOOPS

These radio-opaque Identiloops are used in surgical procedures for retraction, identification or occlusion of arteries, veins, nerves or tendons. A soft and flexible, inert medical grade silicone which is oval in cross-section, allows for even pressure distribution to minimise damage to tissues. Available in 4 colours – Blue, White, Red and Yellow in mini and maxi sizes. All double packed.

2200220-2	Mini-Red 40cm, Packet of 2	10
2200220-4	Mini-Red 40cm, Packet of 4	10



2200220

2200420-2	Mini-Yellow 40cm, Packet of 2	10
-----------	-------------------------------	----



2200420

2200620-2	Mini-Blue 40cm, Packet of 2	10
2200620-4	Mini-Blue 40cm, Packet of 4	10



2200620

2200920-2	Mini-White 40cm, Packet of 2	10
2200920-4	Mini-White 40cm, Packet of 4	10



2200920

2202220-2	Maxi-Red 40cm, Packet of 2	10
2202220-4	Maxi-Red 40cm, Packet of 4	10



2202220

2202420-2	Maxi-Yellow 40cm, Packet of 2	10
2202420-4	Maxi-Yellow 40cm, Packet of 4	10



2202420

2202620-2	Maxi-Blue 40cm, Packet of 2	10
-----------	-----------------------------	----



2202620

2202920-2	Maxi-White 40cm, Packet of 2	10
-----------	------------------------------	----



2202920

DEF1193	Super Maxi-Blue 40cm, Packet of 2	10
DEF1318	Mixed Maxi 40cm, Packet of 5	10
DEF2182	Super Maxi-Blue 60cm, Packet of 2	10



DEF2182

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
LIGHT HANDLE COVERS Our range of flexible light handle covers are single sterile and come either coloured as blue or green. These adaptable covers are twisted or clicked onto the light handle adapter and are firmly held for duration required. We can supply the compatible light handle adapters for each operating room light. Easy and versatile to use.			SKIN HOOKS		
					
DEF2082			DEF336A		
DEF2082	Light Handle Cover Green Friction Fit	50	DEF336A	5 Small Hooks in a Needle Case Double Packed	20
DEF2919	Light Handle Cover Green Friction Fit, Packet of 2	40	DEF1032A	5 Large Hooks in a Needle Case Double Packed	20
DEF2675	Light Handle Cover Blue, Click On	50			
			DEF2402		
DEF2675			DEF2402	Catspaw Retractor, Packet of 2	10
DEF3170	Light Handle Cover Blue Click On, Packet of 2	40	TOUCH FASTENER		
DEF2557	Light Handle Cover Foam Blue	50	DEF5826	5cm x 16cm Velcro Strip X-ray Detectable, Pack of 2	100
					
DEF2557			DEF5826		
IRRIGATION BULB SYRINGE			YARD STRIPS		
AS011P	Irrigation Bulb Syringe 60ml	50	DEF858	1/2 Yard Strips (53cm) X-ray Packet of 5	50
					
AS011P			DEF858		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

NEEDLE COUNTERS DEROYAL


Excellent selection of Needle Counters from DeRoyal with easy to read numbers and description, facilitating correct counting of blades and needles used in a procedure. Full surface magnetic face provides extra space for sharps, pins, blade or needles.

'Y' codes indicate Needle Counter case is Yellow.

DELUXE NEEDLE COUNTERS

MAGNET

25-0103Y	Single Magnet, Clear Top 10 Count, 20 Capacity	50
----------	---	----



25-0103Y

MAGNET/FOAM COMBINATION

25-0401	Magnet/Foam Block 20 Count, 40 Capacity	20
---------	--	----




25-0401

25-0701	Foam Strips/Magnet 15 Count, 30 Capacity	50
---------	---	----



25-0701

25-0702	Magnet/Foam Strips 30 Count, 60 Capacity	27
---------	---	----



25-0702

FOAM

25-0802	Foam Strips, Clear Top 15 Count	50
---------	------------------------------------	----



25-0802

25-0504	Foam Block/Adhesive 40 Count, 80 Capacity	12
---------	--	----



25-0504

25-0803	Foam Strips 60 Count, 120 Capacity	27
---------	---------------------------------------	----



25-0803

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
DEBLADE NEEDLE COUNTERS					
This range of Needle Counters all feature a blade removal system which provides simple removal of contaminated blades.					
MAGNET					
25-1104Y	Double Magnet with DeBlade 30 Count, 60 Capacity	50			
					
25-1104Y					
25-1108	Double Magnet with DeBlade 60 Count, 120 Capacity	27			
					
25-1108					
MAGNET/FOAM COMBINATION					
25-1401	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20			
25-1401Y	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20			
					
25-1401					
25-1701	Magnet/Foam Strips with DeBlade 15 Count, 30 Capacity	50			
					
25-1701					
25-1702	Magnet/Foam Strips with DeBlade 30 Count, 60 Capacity	27			
					
25-1702					
MARKING PENCIL					
DEF464	Prepacked Marking Pencil – Omnichrom 108-9	24			
					
DEF464					
SAFETY PINS					
DEF071SP	#2 (38mm)	50			
DEF073SP	#4 (57mm)	50			
					
DEF073SP					

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SKIN MARKING PENS

Skin Markers are used during surgical procedures and contain medical grade, non irritating gentian violet ink. They mark a smooth, clear line on a patient's skin, and the ink is non toxic, non tattooing and quick drying.

26-001	Skin Marker Standard Tip with 15cm Ruler	50
--------	--	----



26-001

26-002	Skin Marker Fine Tip with 15cm Ruler	50
--------	--------------------------------------	----



26-002

26-001PR	Skin Marker Standard Tip Prep Resistant Ink with 15cm Ruler	50
----------	---	----

26-003PR	Skin Marker Standard Tip Prep Resistant Ink with 15cm Ruler and 8 labels (55mm x 11mm each)	50
----------	---	----

26-005PR	Skin Marker Standard Tip Prep Resistant Ink	50
----------	---	----



26-001PR

NON-STERILE

26-018NS	Mini Pre-Op Skin Marker Bulk Non Sterile	200
----------	--	-----



26-018NS

UTILITY MARKER

Ideal for writing on instruments, specimens, vials and tubes, basins, bowls, and other objects as needed within the sterile field. This marker is not to be used on skin.

DEF3219	Utility Marker Fine Tip Black Ink	25
---------	-----------------------------------	----



DEF3219

SPONGES & SWABS

SPONGES ABDOMINAL X-RAY

Soft highly absorbent cotton sponges. Each complete with an x-ray radiopaque element for easy detection.



DEF162

DEF162	23cm x 23cm with Tape, XRD	100
DEF160	23cm x 23cm with Tape, XRD Packet of 3	100
DEF168	23cm x 23cm with Tape, XRD Packet of 5	100
DEF1351	30cm x 30cm, XRD, Packet of 5 Double Packed	75
DEF1352	30cm x 48cm, XRD, Packet of 5 Double Packed	60



DEF1351

CODE

DESCRIPTION

BOX QTY.

CODE

DESCRIPTION

BOX QTY.

GAUZE SWABS X-RAY

Excellent absorbency and quality swabs with 8 ply cotton blend. Includes an x-ray radiopaque element for easy detection.



DEF044

DEF202	7.5cm x 7.5cm X-ray, Packet of 5 Double Packed	100
DEF044	10cm x 10cm X-ray, Packet of 5 Double Packed	100
DEF997	10cm x 10cm X-ray, Packet of 5 32ply, Double Packed	50

GAUZE SWABS GREEN



DEF190

DEF190	10cm x 10cm Green, Packet of 2	100
DEF192	10cm x 10cm Green, Packet of 3	100
DEF207	10cm x 10cm Green, Packet of 5 Double Packed	100
DEF225	10cm x 10cm Green, Packet of 20	50
DEF196	10cm x 10cm Green, Bulk Non Sterile	100

GAUZE SWABS FOLDED

DEF347	Gauze Swab 7cm x 21cm X-ray, Folded Packet of 5	50
--------	---	----



DEF347

DEF185	Gauze Swab 10cm x 10cm 32ply X-ray, Packet of 10, Double Packed	100
--------	---	-----



DEF185

DEF1074	Gauze Swab 10cm x 20cm X-ray Folded, Packet of 5	25
---------	--	----

GAUZE HEAD DRESSINGS

DEF611	Head Dressing Kit includes: 8 x Gauze Swabs 10cm x 20cm 20 x Gauze Swabs 10cm x 10cm 2 x 90cm x 90cm Gauze, Folded	35
--------	---	----

GAUZE PREP SPONGES X-RAY

Round sponges with X-ray radiopaque element. Excellent absorbency. Used as a prep sponge prior to a surgical procedure.



DEF251

DEF255	Sponges Prep X-ray, 6cm x 6cm	100
DEF251	Sponges Prep X-ray, 6cm x 6cm Packet of 2	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPONGES & SWABS (CONTINUED)

GAUZE PEANUT SPONGES

Round sponges ideal for delicate surgical sponging or soft tissue dissection. X-ray radiopaque element for easy detection.

30-106	5 in a Foam Holder, XRD	50
--------	-------------------------	----



30-106

SPECIALISED GAUZE



DEF169

DEF986	Donor Site Gauze Dressing, 5cm x 10cm, Double Packed	100
DEF169	Recovery Gauze 5cm x 5cm Plain ,32 Ply	100

NEUROSURGICAL PATTY

Our neurosurgical patties are made from high-quality cotton/rayon fibres to prevent disintegration and linting. They are pliable when wet, allowing patties to conform to irregular contours and provide excellent absorbency. The x-ray detectable string is secured on one side of the patty. A counting card is also included.



UNS1275

UNS55	5mm x 5mm Ultracot Pack of 10	20
UNS1075	10mm x 75mm Ultracot Pack of 10	20
UNS1212	12mm x 12mm Ultracot Pack of 10	20
UNS1225	12mm x 25mm Ultracot Pack of 10	20
UNS1240	12mm x 40mm Ultracot Pack of 10	20
UNS1250	12mm x 50mm Ultracot Pack of 10	20
UNS1275	12mm x 75mm Ultracot Pack of 10	20
UNS2020	20mm x 20mm Ultracot Pack of 10	20
UNS2525	25mm x 25mm Ultracot Pack of 10	20
UNS2575	25mm x 75mm Ultracot Pack of 10	20

SCOPE WARMER

Long, double sided Scope Warmer from that pre-warms endoscopes and laparoscopes to body temperature. Easy to use activate by flexing/ clicking the internal disk. Provides 1 hour of warming.

31SW01	ClearView™ Scope Warmer	10
--------	-------------------------	----



31SW01

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
PVC TUBING WITH XRD Medical grade, clear and transparent tubing. A radiopaque strip along length of tubing ensures X-ray detectability.			ESMARCH BANDAGE Offers excellent compression and elastic properties to achieve exsanguination, whilst delivering protection for latex sensitive patients and health professionals.		
DEF3051	PVC Tubing with XRD 16Fr x 17cm, Packet of 6 Double Packed	25	DEF1006	6.5cm x 2m Green Silicone Double Packed	10
					
DEF3051			DEF1006		
DEF3052	PVC Tubing with XRD 16Fr x 22cm, Packet of 6 Double Packed	25	DEF2764	10cm x 2.7m Blue Neoprene Double Packed	20
					
			DEF2764		
SCALPEL BLADE REMOVER Single handed blade remover that can hold up to three used surgical blades. Simple to use with a quick 'push, click and pull' action to remove and contain used scalpel blades.			35-309	7.6cm x 2.7m Blue Neoprene	20
QSSVCAS-3Y	Quicksmart Blade Cassette	100	35-409	10cm x 2.7m Blue Neoprene	20
			35-609	15cm x 2.7m Blue Neoprene	20
QSSVCAS-3Y					
Please note - blade not included.			35-309		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL RETRACTOR

ONETRAC™

ONETRAC™ is a single-use cordless surgical retractor with a built-in LED light source and smoke evacuation channel.

It provides illumination in the depths of a surgical tissue pocket or cavity. No additional parts or components or fibre optic cables are required. Ideal for breast or reconstructive procedures, as well as other surgical procedures requiring access, illumination and retraction of soft tissue. Reduce cost, increase efficiency and eliminate unnecessary risk of cross-contamination with ONETRAC.



ONETRAC

C100150-5	40mm x 20mm	5
C100170-5	55mm x 8mm	5
C100140-5	90mm x 22mm	5
C100180-5	105mm x 60mm	5
C100165-5	135mm x 20mm	5
C100110-5	135mm x 30mm	5
C100190-5	175mm x 30mm	5



ONETRAC LX

Features integrated multi-LED light sources and dual smoke evacuation channels.

C100270	30mm x 15mm & 50mm x 20mm Dual Blade	1
C100240	90mm x 20mm	1
C100270	90mm x 30mm	1
C100230	90mm x 60mm	1
C100290	135mm x 20mm	1
C100210	135mm x 30mm	1
C100260	135mm x 60mm	1
C100220	175mm x 30mm	1



C100260, C100230

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE TUBING

Our 100% medical grade, clear silicone tubing range is available pre-cut and sterile or as non-sterile rolls. They are double packed unless otherwise specified.

Sizing below is represented as:

ID (Internal Diameter) x OD (Outer Diameter) x Length of tube in cm or m.

DEF2307	0.30mm x 0.64mm x 20cm	75
DEF2522	1mm x 2mm x 50cm	50
DEF477	1.57mm X 2.41mm X 30cm Packet of 2	25
DEF1545	1.57mm x 3.18mm x 1m	15
DEF1428	1.57 mm x 2.41mm x 15cm	100
DEF473	1mm x 3mm x 50cm, Packet of 4	20
DEF3649	1.98mm x 3.18mm x 75cm Packet of 2	20
DEF5712-10	2mm x 4mm x 45cm	10
DEF2314	3.18mm x 6.35mm x 25cm	60
DEF5933	3mm x 5mm x 80cm	20
DEF2038A-10	4mm x 6mm x 45cm	10
DEF1665	5mm x 7mm x 20cm Packet of 2	25
DEF2167B-10	5mm x 7mm x 45cm	10
DEF1253	5mm x 8mm x 1.6m	50
DEF1261A	5mm x 9mm x 1.5m Single Packed	30
DEF5655-10	6mm x 9mm x 45cm	10
DEF5652	6.4mm x 11mm x 3m	8
DEF1818A-10	7mm x 10mm x 45cm	10
DEF2032	7mm x 10mm x 3m	8
DEF1859A-10	8mm x 12mm x 45cm	10
DEF1204	8mm x 14mm x 3m	5
DEF2468	10mm x 16mm x 1.5m	16

SILICONE TUBING WITH CONNECTOR

Flexible pre-cut and sterile silicone tubing with connector attached. All of this range is double packed.




DEF1186

DEF735	5mm x 8mm x 1m	50
DEF1186	5mm x 9mm x 30cm	50
DEF306A	5mm x 9mm x 50cm	50
DEF1712	5mm x 9mm x 60cm	40
DEF1337	5mm x 9mm x 1m	50
DEF1231A	5mm x 9mm x 1.5m	30
DEF305A	7mm x 11mm x 1m	30

SUTURE BOOTS WITH ADHESIVE TRAY

Suture boots offer a secure way to hold sutures without slippage or damage during a surgical procedure. They are fitted onto the distal working end of surgical instruments to provide an improved atraumatic grip and even pressure distribution when clamping an artery or vein. The tray comes with an adhesive strip for securing and easy removal. Radio-opaque element ensures X-ray detection. Available in a Yellow colour.



30-702

DEF1093	2 on an Adhesive Tray Double Packed	25
DEF1094	4 on an Adhesive Tray Double Packed	25
30-702	6 on an Adhesive Tray	50
DEF1096	8 on an Adhesive Tray Double Packed	25
30-701	10 on an Adhesive Tray	50
DEF2284	12 on an Adhesive Tray Double Packed	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

WOUND DRAINAGE

Wound drains and evacuators provide drainage for various surgical procedures, allowing for the removal of blood and excess fluid to improve patient outcome. Our wound drainage systems are made up of evacuators, flat drains, fluted drains, round drains, capillary action drains and easy flow drains.

They are all made from 100% medical grade silicone, which is inert, soft and flexible and designed to reduce the potential of post – operative infections.

WOUND EVACUATOR BULB

Our wound evacuator is a clear bulb with easy to read graduations up to 100ml/cc, and is designed to slowly draw fluid and blood from the wound site. It creates an ideal vacuum to a maximum of 125mmHg pressure and is compatible with our flat drains, fluted drains and round drains.



31300380DS

31300380DS	Evacuator Bulb 100ml/cc	10
------------	-------------------------	----

PERFORATED FLAT WOUND DRAIN

Soft and comfortable silicone flat drains with multiple drainage holes to enhance the fluid removal in a wound. They have an ergonomic flat shape which enables an easier placement. X-ray detectable.



11100200780DS

11100200780DS	7mm x 800mm with a 220mm Channel, White XRD Strip	10
11100100780DS	7mm x 800mm with a 220mm Channel, Blue XRD Strip	10
11100201080DS	10mm x 800mm with a 220mm Channel, White XRD Strip	10
11100101080DS	10mm x 800mm with a 220mm Channel, Blue XRD Strip	10

FLUTED FLAT WOUND DRAIN

Clear soft silicone drain with white x-ray detectable end.

11101400780DS	7mm x 1020m with 200mm Channel and Trocar	10
11101411080DS	10mm x 1020mm with 200mm Channel and Trocar	10

FLUTED ROUND WOUND DRAIN

Fluted drains are made of soft, flexible silicone with or without attached trocar. Designed with a multi channel to allow for great drainage and assisting in minimal tissue in-growth. The inner core of the H design offers drain strength and flexibility. X-ray detectable.



11003021080DS	10Fr x 1100mm with 300mm Channel and Connector, Blue XRD Strip	10
11003121080DS	10Fr x 1100mm with 200mm Channel and Connector, Blue XRD Strip	10
11003121580DS	15Fr x 1200mm with 300mm Channel and Connector, Blue XRD Strip	10
11003121980DS	19Fr x 1200mm with 300mm Channel and Connector, Blue XRD Strip	10
11003122480DS	24Fr x 1100mm with 300mm Channel and Connector, Blue XRD Strip	10

FLUTED ROUND WOUND DRAIN SETS

11003131080DS	10Fr x 1100mm with 200mm Channel and Trocar, Blue XRD Strip	10
11003131580DS	15Fr x 1200mm with 300mm Channel and Trocar, Blue XRD Strip	10
11003131980DS	19Fr x 1200mm with 300mm Channel and Trocar, Blue XRD Strip	10
11003831080DS	10Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10
11003811580DS	15Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10
11003811980DS	19Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10

PERFORATED ROUND WOUND DRAIN

Clear soft silicone drain with a blue x-ray detectable line.

11000301280DS	12Fr 500mm with 5 holes 2x4mm	10
11000301880DS	18Fr 500mm with 5 holes 4x6mm	10
11000302480DS	24Fr 500mm with 5 holes 5x8mm	10
11000302780DS	27Fr 500mm with 5 holes 6x9mm	10
11000303080DS	30Fr 500mm with 5 holes 7x10mm	10
11000303680DS	36Fr 500mm with 5 holes 8x12mm	10
11000600780DS	7Fr 1250mm with 72 Holes	10
11000601080DS	10Fr 1250mm with 72 Holes	10
11000601580DS	15Fr 1250mm with 72 Holes	10
11000601980DS	19Fr 1250mm with 72 Holes	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL CONSUMABLES

ANAESTHETIC FACE MASK

Our anaesthetic face masks are ergonomically designed to optimise the application of anaesthetics. The masks have an inflating air cushion that is adjustable by an inflation valve, for sealing and patient comfort.

The transparent face plate and cushion provides easy patient facial visibility. The masks are for single patient use and are provided non-sterile.



A1701

A1701	Size 0 - Neonate	50
A3011	Size 0 - Neonate, Strawberry Scented	50
A1702	Size 1 - Infant	50
A3012	Size 1 - Infant, Strawberry Scented	50
A1703	Size 2 - Toddler	50
A3013	Size 2 - Toddler, Strawberry Scented	50
A1704	Size 3 - Child / Small Adult	50
A3014	Size 3 - Child / Small Adult Strawberry Scented	50
A1705	Size 4 - Medium Adult	50
A3015	Size 4 - Medium Adult, Strawberry Scented	50
A1706	Size 5 - Large Adult	50
A3016	Size 5 - Large Adult, Strawberry Scented	50
A1707	Size 6 - X-Large Adult	50
A3017	Size 6 - X-Large Adult, Strawberry Scented	50

BURNS

BURNS DRESSING

Specialised products for patients with burns. Non adherent, highly absorbent wound dressings.



DEF199

DEF199	10cm x 20cm, Packet of 20	50
DEF679	50cm x 70cm, Melolin Double Packed	35
DEF3570	50cm x 175cm Melolin Double Packed	25
DEF974	90cm x 2m	50



DEF679

BURNS GAUZE

Soft roll suitable for wounds and burns.

DEF965	10cm x 80cm, Packet of 16	50
DEF2436	25cm x 90cm, Packet of 2	350
DEF1039	60cm x 90cm, Packet of 6	100
DEF342	90cm x 2m	50
DEF344	90cm x 2m, Packet of 4	25
DEF1662	90cm x 5m	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CLIP SEAL BAG



DEF3286

DEF3286	20cm x 25cm, Packet of 20	50
DEF1703	28cm x 38cm	25

FREEZER BAG



DEF639

DEF639	30cm x 45cm	100
DEF66	1 x 30cm x 45cm Bag with 2 x Latex-free Bands #34 Double Packed	100
DEF1336	2 x 30cm x 45cm Bags with 2 x Latex-free Bands #34 Double Packed	100

DRAWSTRING BAG



DEF846

DEF846	28cm x 54.5cm, Double Packed	25
DEF864	44.5cm x 46.5cm	25
DEF861	44.5cm x 46.5cm	100
DEF1687	44.5cm x 46.5cm, Packet of 2	50
DEF848	50cm x 50cm	25
DEF847	50cm x 85cm	25

PLASTIC BAG



DEF1878

DEF1878	14cm x 30cm	100
DEF638	28cm x 32.5cm	100
DEF642	45.5cm x 51cm with Cotton Tape 60cm Long Double Packed	100
DEF808	Heavy Duty 57.5cm x 62cm	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COTTON PRODUCTS

COTTON BALLS

Packaged and sterilised for ease of entry into a sterile field. Ideal for prepping, wound cleaning, perineal care and fluid application.



DEF074

DEF074	3cm, Packet of 5	100
--------	------------------	-----

COTTON WOOL ROLL

Soft and absorbent cotton wool roll used for padding in patient care. Sterilised and ready to be used.



DEF874

DEF874	15cm x 55cm, Double Packed	50
DEF875	30cm x 90cm, Double Packed	25
DEF880	38cm x 165cm, Double Packed	15

SPONGES

SPONGES NASAL BOLSTER

Used after nasal or sinus surgery to absorb any discharge or ooze from the nose. High quality non woven net filled with absorbent cotton padding, and stretch elastic sides to secure around the ears. Extremely soft and comfortable for post surgery patients. Latex free and available as sterile or non sterile.



DEF1381

DEF1381	2cm x 8cm with Ear Loops	20
DEF1380	2cm x 8cm with Ear Loops Bulk Non Sterile	20
DEF1382	2cm x 8cm with Ear Loops	100

NON WOVEN PRODUCTS

NON WOVEN BALLS

DEF1354	2.5cm, Packet of 5	100
---------	--------------------	-----

NON WOVEN GAUZE SWABS - STERILE

DEF3933A	5cm x 5cm, Packet of 2	1000
DEF1131	5cm x 5cm, Packet of 5	100

NON WOVEN GAUZE SWABS - NON STERILE

DEF1869	5cm x 5cm, Bulk	100
DEF1868	7.5cm x 7.5cm, Bulk	100
DEF1867	10cm x 10cm, Bulk	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SWABS

GAUZE SWABS PLAIN - STERILE

Sterile and ready to be used for a wide variety of applications, including as a dressing, for prepping procedures and wound cleaning. Available as a sterile Plain 8 ply gauze.



DEF181

DEF183	7.5cm x 7.5cm, Packet of 3	100
DEF181	7.5cm x 7.5cm, Packet of 5	100
DEF371	7.5cm x 7.5cm, Packet of 20	50
DEF189	10cm x 10cm, Packet of 2	100
DEF186	10cm x 10cm, Packet of 3	100



DEF045

DEF045	10cm x 10cm, Packet of 5	100
DEF208	10cm x 10cm, Packet of 10	50
DEF213	10cm x 10cm, Packet of 10	300
DEF178	10cm x 10cm, Packet of 20	25
DEF177	10cm x 10cm, Packet of 25	25
DEF214	10cm x 10cm, Packet of 25	200
DEF1033	10cm x 20cm, Packet of 4	100
DEF198	10cm x 20cm, Packet of 5	250
DEF195	10cm x 20cm, Packet of 10	50
DEF853	10cm x 80cm, Packet of 4 16 ply Folded	50

GAUZE SWABS PLAIN – NON STERILE

DEF180	5cm x 5cm, Bulk	100
DEF197	7.5 x 7.5cm, Bulk	100
DEF136	10cm x 10cm, Bulk	100

GAUZE DRESSING

DEF571	15cm x 22cm	100
--------	-------------	-----

PURCOTTON SWABS

Made from a non-woven fabric using natural cotton, our purcotton swabs have a high liquid absorbency and are 100% biodegradable. They are also low linting and have an excellent wet-strength.

DEF5753	10cm x 10cm, Pack of 10	100
DEF1062	10cm x 10cm, Non Sterile	100

SANITARY DISPOSAL BAG

DEF709	Packet of 100, Non Sterile	20
--------	----------------------------	----



DEF709

VAGINAL PLUG



DEF850SP

DEF850SP	90mm x 30mm x 30mm	100
----------	--------------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SWAB STICKS/COTTON TIPPED APPLICATORS

Packaged and sterilised, ready to use for taking a swab or specimen from a patient. Paper, plastic or wooden shaft in different length sizes with a single or double ended soft cotton head.

SWAB STICKS WITH WOODEN SHAFT



DEF057

DEF057	7.5cm Single Ended, Packet of 3	100
DEF442	7.5cm Single Ended, Packet of 5	100
DEF439	15cm Single Ended, Packet of 2	100
DEF059	15cm Single Ended, Packet of 3	100
DEF443	15cm Single Ended, Packet of 5	100



DEF439

SWAB STICKS WITH WOODEN SHAFT – NON STERILE

DEF062	15cm Single Ended, Bulk	1000
--------	-------------------------	------

SWAB STICKS WITH PAPER SHAFT



DEF421

DEF058	7.5cm Double Ended, Packet of 3	100
DEF421	7.5cm Double Ended, Packet of 10	100

SWAB STICKS WITH PAPER SHAFT – NON STERILE

DEF061	7.5cm Double Ended, Bulk	1250
--------	--------------------------	------



DEF061

SWAB STICKS JUMBO



DEF401

DEF400	15cm, Packet of 2	150
DEF401	15cm, Packet of 3	150
DEF402	15cm, Packet of 5	150
DEF1593	20cm	200
DEF455	20cm, Packet of 2	150
DEF407	20cm, Packet of 3	150



DEF407

SWAB STICKS JUMBO – NON STERILE

DEF055	15cm Bulk	200
DEF054	20cm Bulk	200

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE TUBING ROLL – NON STERILE

Excellent range of medical grade silicone tubing available in a variety of different metric sizes. Available as non sterile rolls of 15 or 25 metres.



DEF203005009	0.51mm (ID) x 0.94mm (OD) 15m Roll
DEF2110070439	0.7mm (ID) x 1.6mm (OD) 15m Roll
DEF201010020	1mm (ID) x 2mm (OD) 25m Roll
DEF2110100234	1mm (ID) x 3mm (OD) 25m Roll
DEF2110100439	1.02mm (ID) x 2.16mm (OD) 15m Roll
DEF2110150939	1.57mm (ID) x 2.41mm (OD) 15m Roll
DEF2110151039	1.57mm (ID) x 3.18mm (OD) 15m Roll
DEF2110190239	1.98mm (ID) x 3.18mm (OD) 15m Roll
DEF2110200334	2mm (ID) x 4mm (OD) 25m Roll
DEF2110300234	3mm (ID) x 5mm (OD) 25m Roll
DEF2110310239	3.18mm (ID) x 6.35mm (OD) 15m Roll
DEF2110400434	4mm (ID) x 6mm (OD) 25m Roll
DEF2110400734	4mm (ID) x 8mm (OD) 25m Roll
DEF20304707901	4.77mm (ID) x 7.9mm (OD) 15m Roll
DEF2110500334	5mm (ID) x 8mm (OD) 25m Roll
DEF2110500534	5mm (ID) x 9mm (OD) 25m Roll
DEF2110600234	6mm (ID) x 9mm (OD) 25m Roll
DEF2110600334	6mm (ID) x 10mm (OD) 25m Roll
DEF2110630139	6.4mm (ID) x 9.5mm (OD) 15m Roll
DEF2110630239	6.4mm (ID) x 11mm (OD) 15m Roll
DEF2110630339	6.4mm (ID) x 12.7mm (OD) 15m Roll
DEF2110700334	7mm (ID) x 10mm (OD) 25m Roll
DEF2110700534	7mm (ID) x 11mm (OD) 25m Roll

DEF2110790239	7.9mm (ID) x 12.7mm (OD) 15m Roll
DEF2110800334	8mm (ID) x 12mm (OD) 25m Roll
DEF2110800534	8mm (ID) x 14mm (OD) 25m Roll
DEF2110900334	9mm (ID) x 13mm (OD) 25m Roll
DEF2110950139	9.5mm (ID) x 12.7mm (OD) 15m Roll
DEF2110950339	9.5mm (ID) x 15.9mm (OD) 15m Roll



CUSTOM PROCEDURE AND STANDARD PACKS

CUSTOM PROCEDURE PACKS

COST EFFECTIVE AND TIME EFFICIENT: A PACK FOR EVERY PROCEDURE

Our Custom Procedure Pack range covers all areas of the hospital from Cardiac, Orthopaedic, Vascular, ENT and more. They contain all the individual consumables for a surgical procedure, wrapped in a 'field' and sterilised. All products are ready to use when opened, unlike individually wrapped consumables that are touched many times prior to use and require opening separately. Plus there is reduced wastage. Using a Procedure pack saves on setup time for procedures and ensures that all components are readily available for every case. Ordering and inventory control is easier.

CUSTOM TAILOR A PROCEDURE PACK TO MEET YOUR SPECIFIC NEEDS

Select individual consumables such as patient drapes and surgical gowns, to diathermy, suction tubing, abdominal packs, gauze and instruments. Your pack will only include items specifically chosen by you, the surgeon and theatre nurses. Defries assembles these customised packs in cleanrooms according to hospital specifications.

Custom procedure packs can be varied to include cost-effective components as they become available, and allow easy accountability of costs per procedure. As a local Melbourne manufacturer, we are able to respond quickly to any changes requested by a hospital to the configuration of your custom packs.

We'll provide you with ongoing education, listen to feedback and are always available to discuss your specific needs.



READY-FOR-USE STANDARD PROCEDURE PACKS

We have a range of standard procedure packs such as the anaesthetic packs, epidural packs, regional trays, catheter packs, Bronny IV insertion pack, lumbar puncture packs, ICU packs including central line packs, arterial line packs, all designed, made and ready to be ordered.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PACKS

All our anaesthetic packs come sterile, prepacked and ready to use. Indications for use include:

- General and Local Anaesthetic procedures
- Spinal/Epidural/Regional/Lumbar Puncture procedures
- IV Access

NRFIT components are available on request for inclusion into packs.

Select from the packs below or contact us to customise your own anaesthetic pack.

ANAESTHETIC PACKS

DEF1409	1 x Dressing Tray 24cm x 15cm x 4cm 2 x Swabs 10cm x 10cm Green	36
---------	--	----



DEF1409

DEF1162	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 3 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	60
---------	---	----



DEF1162

ANAESTHETIC TRAY

DEF1524	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 2 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	80
DEF406	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Polybacked Towel 40cm x 42cm	60

ANAESTHETIC PACKS

DEF1618	1 x Drape Blue Plastic 75cm x 75cm with 10cm x 15cm Fenestration with Adhesive Tabs on Corners of Fenestration and Corners of Drape 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 5 x Gauze Swabs 10cm x 10cm Green 1 x Forcep Plastic Green 1 x Syringe 5ml Luer Slip 1 x Syringe 3ml Luer Slip 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt Drawing up 1 x Needle 25g x 1" 1 x Dressing Towel Quick Dry 33cm x 33cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
---------	---	----



DEF1618

ANAESTHETIC BASIC TRAY AND SYRINGE KIT

DEF4117E	1 x Anaesthetic Tray with 4 Compartments 22.5cm x 14cm x 2.5cm 1 x Syringe 3ml Luer Slip 1 x Syringe 5ml Luer Slip 1 x Syringe 10ml Luer Slip Eccentric Tip 1 x Syringe 30ml Luer Slip Eccentric Tip 2 x Gauze Swabs 10cm x 10cm Green Wrapped in Crepe Paper 75cm x 75cm	50
----------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PROCEDURE PACK

DEF2255C	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 91cm x 115cm 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 3ml Luer Lock 1 x Syringe 5ml Luer Lock 1 x Needle 25g x 1.5" 2 x Filter Needle 5um 2 x Forceps Pearman Plastic 2 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 1 x Gallipot 60ml 1 x Label Sheet - 1 Set of 16 Wrapped in a Blue Plastic Field 91cm x 115cm	60
----------	--	----

NRFIT REGIONAL ANAESTHETIC PACK

DEF5964	1 x Procedure Drape 120cm x 140cm with Adhesive and 10cm x 15cm Fenestration 1 x Drape 100cm x 100cm with Adhesive 1 x Syringe 1ml Luer Slip 1 x Syringe 5ml Luer Slip 1 x Needle 18g x 1.5" Blunt 1 x NRFIT Needle 19g x 1.5" with Filter 1 x Needle 23g x 1.25" 1 x NRFIT Syringe 1ml Luer Slip 1 x NRFIT Syringe 5ml Luer Slip 5 x Gauze Swabs 10cm x 10cm Green 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow 2 x Gallipot 60ml 1 x Label Sheet - 1 Set of 28 1 x Clipseal bag (for small components) Wrapped in a Back Table Cover 140cm x 190cm	10
---------	--	----

ANAESTHETIC/SPINAL PACK

DEF2232B	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 76cm x 106cm 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 5ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 19g x 1.5" 1 x Needle 25g x 1" 1 x Forceps Plastic Green 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm Green 5 x Prep Swabs Large XRD 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Label Sheet - 1 Set of 7 Wrapped in a Blue Plastic Field 75cm x 75cm	60
----------	--	----



DEF2232B

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
EPIDURAL PACKS			COMPLETE EPIDURAL PACK		
DEF1796	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 2 x Forceps Plastic Green 5 x Gauze Swab 10cm x 10cm Green 1 x Dressing Towel 32cm x 60cm Wrapped in a Blue Plastic Field 91cm x 115cm	50	DEF1480	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Forceps Plastic Green 1 x Syringe 5ml Luer Slip 1 x Syringe 3ml Luer Slip 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1" 5 x Gauze Swab 10cm x 10cm Green 1 x Dressing Towel (Hand Towel) 32cm x 32.5cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1655B	1 x SMS Drape 91cm x 91cm with 12cm x 15cm Fenestration and Adhesive 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Forceps Plastic Green 5 x Gauze Swab 10cm x 10cm Green 1 x Gallipot 60ml 1 x Towel 40cm x 42cm 2ply 1 x Pack Labels - 1 Set of 16 Wrapped in a Blue Plastic Field 91cm x 115cm	60	EPIDURAL/SPINAL TRAY		
DEF3085B	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and ACorner Tabs 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 3ml Luer Slip 1 x Syringe 5ml Luer Slip 2 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1" 1 x Needle 25g x 1.5" 5 x Gauze Swabs 10cm x 10cm Green 1 x Dressing Towel 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	50	DEF2345C	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray with 3 Compartments 20cm x 25cm x 2.5cm 1 x Syringe 10ml Luer Lock 1 x Syringe 5ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1.5" 1 x Filter Needle 5um 1 x Foam Applicator Blue 15cm 5 x Gauze Swabs 10cm x 10cm Green 2 x Dressing Towel (Hand Towel) 32cm x 32.5cm 1 x Gallipot 60ml 1 x Plastic Bag 20cm x 20cm(with small components enclosed) Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF3854	1 x Drape Laminated 120cm x 140cm with 7cm x 10cm Fenestration Blue 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Sponge Holder 19cm Blue Plastic 5 x Gauze Swabs 10cm x 10cm Green 2 x Gallipot 150ml 1 x Gown Surgical T-Cut Unferinforced XL 1 x Huck Towel 40cm x 68cm Blue Wrapped in a Back Table Cover 152cm x 228cm	15			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CATHETER PACKS

All Catheter Procedure Packs have been designed in consultation with the healthcare professionals who use them. They include all the consumables used to insert or change a catheter aseptically. Alternatively please contact us to customise a catheter pack to your own specific requirements.

CATHETER PACKS



DEF1601

DEF1601	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Blue Huck Towel 40cm x 68cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF2074	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Pearman 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 2 x Kidney Dishes 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24

DEF2061A	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Green 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF1262A	1 x Drape SMS 61cm x 61cm with a 10cm x 10cm Fenestration 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Tube Clamp with Lock 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1650B	1 x Drape Azure Minor Procedure 50cm x 60cm with 6.5cm Circular Fenestration and Adhesive around Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1522	3 x Drape Polybacked 40cm x 42cm 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm Plain Wrapped in a Blue Plastic Field 75cm x 75cm	60

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL PACKS

CATHETER PACKS

DEF1036	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Forceps Plastic Green 6 x Cotton Balls Large 1 x Dressing Towel T 31.5cm x 61cm 2 x Kidney Dishes 700ml Clear 2 x Gallipots 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1932A	1 x Drape PolyBacked 50cm x 60cm with 10cm Circular Fenestration 1 x Syringe 30ml Luer Slip Eccentric Tip 1 x Sponge Holder 24cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel 32cm x 32.5cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	30

IV INSERTION PACKS

For the insertion of an IV line.

BRONNY IV PACK

DEF920	1 x Plastic Drape 45cm x 60cm 2 x Gauze Swabs 5cm x 5cm Plain	200
--------	--	-----

BRONNY IV PACK WITH 4 SWABS

DEF1676	1 x Plastic Drape 45cm x 60cm 4 x Gauze Swabs 5cm x 5cm Plain	200
---------	--	-----

MOUTH TOILET PACK

DEF399	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 10 x M & T Swab Sticks 15cm Paper Shaft 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	50
--------	--	----

FINE PROCEDURE PACK

DEF2453A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Scalpel Handle # 3 4 x Gauze Swabs 7.5cm x 7.5cm Plain 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
----------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INTENSIVE CARE PACKS

THYROID PACK

DEF3799B	1 x Forcep Spencer Wells 13cm Straight 1 x Scissors 13cm Sharp/Blunt 1 x Staple Remover Scissor Action with Plastic Handle 1 x Stitch Cutter 10 x Gauze Swabs 10cm x 10cm Plain 1 x Huck Towel 40cm x 68cm Blue Wrapped in a SMS Blue Field 61cm x 61cm	20
----------	---	----

TRACHY PROCEDURE PACK

DEF2445	1 x Tracheal Dilator Trousseau 2 Prong 1 x Retractor Bergmann 16.5cm	1
---------	---	---

LARGE DRESSING PACK

DEF2234A	1 x Dressing Tray 24.5cm x 28cm x 4.5cm 3 x Forceps Plastic Green 10 x Non Woven Balls 10 x Non Woven Swabs 7.5cm x 7.5cm 2 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	10
----------	--	----

CVC INSERTION PACK

DEF2393A	1 x Scissors 13cm Sharp/Blunt Straight 1 x Needle Holder 13cm 1 x Forcep Adson Micro Toothed 2 x Forcep Spencer Wells 13cm Straight 5 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 2 x Towel Cotton 40cm x 68cm Blue 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	25
----------	--	----

LUMBAR PUNCTURE PACK

DEF2746E	1 x Drape SMS 61cm x 61cm with 14cm Circular Fenestration and Adhesive Corner Tabs 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Microlance Needle 25g with 16mm Regular Wall 1 x Syringe 5ml Luer Lock 2 x Forcep Plastic Green 1 x Asguard Transparent Film Dressing 6cm x 7cm High-tack Adhesion 5 x Gauze Swabs 5cm x 5cm Plain 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Labelled Specimen Jar 5ml 4 x Labelled Specimen Tube 10ml with 'V' Bottom Wrapped in a Blue Plastic Field 75cm x 75cm	24
----------	---	----

INTERCOSTAL PACK

DEF1609	1 x Needle Holder 16cm 2 x Forcep Spencer Wells 13cm Straight 2 x Forcep Spencer Wells 13cm Curved 1 x Scissors 12.5cm Sharp/Sharp Straight 1 x Scalpel Handle # 4 1 x Forcep Gillies Toothed 15.5cm 6 x Cotton Balls Large 4 x Gauze Swabs 10cm x 10cm Plain 4 x Dressing Towel (Hand Towel) 31.5cm x 34cm 2 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a SMS Blue Field 91cm x 91cm	24
---------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

STANDARD SURGICAL PACKS

ARTERIAL TRAY

DEF2161A	1 x Drape Polybacked 33cm x 46cm with 8cm Square Fenestration 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Scissors 13cm Sharp/Blunt Straight 1 x Scissors Iris 11.5cm Sharp/Sharp 1 x Needle Holder 13cm 2 x Forceps Artery Mosquito 12.5cm Straight 2 x Forceps Artery Mosquito 12.5cm Curved 1 x Forcep Iris Plain 10cm 1 x Forcep Iris Toothed 10cm 1 x Probe Lachrymal Size 1-2 1 x Probe Lachrymal Size 3-4 6 x Gauze Swabs 7.5cm x 7.5cm Plain 2 x Huck Towel 40cm x 68cm Blue Wrapped in a Blue Plastic Field 75cm x 75cm	5
----------	---	---

UNIVERSAL GENERAL PACK 1

DEF5767	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Unreinforced Gown Level 3 X-Large 1 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 5 x Gauze Swabs 10cm x 10cm X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Foam Applicator 20cm Orange 2 x Light Handle Cover Green 1 x Quiver 1 x Needle Counter Magnet/Foam Strips 30 Count 1 x Labels for Theatre - sheet of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	4
---------	---	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
UNIVERSAL GENERAL PACK 2			UNIVERSAL GENERAL PACK 3		
DEF5768A	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Unreinforced Gown Level 3 X-Large 1 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 1 x Needle 23g x 1.25" 1 x Needle 18g Blunt 1 x Scalpel Blade #15 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Prep Swab Large X-Ray 2 x Light Handle Cover Green 1 x Bowl 1L Blue 2 x Kidney Dish 700ml Clear 3 x Gallipot 150ml 2 x Sharps Tray 500ml Yellow 1 x Sharps Safety Transfer Device 1 x Needle Counter Foam Strips with Clear Lid 15 Count Yellow 1 x Labels for Theatre - sheet of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm Heavy Duty GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	3	DEF5769	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 2 x Fabric Reinforced Gown Level 3 Extra Large 2 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 1 x Needle 23g x 1.25" 1 x Scalpel Blade #15 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Prep Swab Large X-Ray 1 x Needle Counter Foam Strips with Clear Lid 15 Count Yellow 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Fabric Reinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a 91cm x 91cm SMS Wrap	4

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
UNIVERSAL GENERAL PACK 4			INTRAOCULAR SURGICAL PACK 1		
DEF5770B	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 130cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Syringe 20ml Luer Lock 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Gauze Swabs 45cm x 45cm X-Ray 1 x Gauze Prep Ball X-ray 1 x Bowl 500ml Blue 4 x Kidney Dish 700ml Clear 2 x Kidney Dish 700ml Yellow 1 x Kidney Dish 1L Blue 3 x Gallipot 150ml 1 x Needle Counter Foam Strips/Magnetic 30 Count 1 x Quiver 1 x Labels for Theatre - 1 Set of 16 1 x Paper Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm Heavy Duty with Full Reinforcement GOWN PACK 1 x Poly Reinforced Gown Level 3 X-Large 2 x Poly Reinforced Gown Level 3 XX-Large 6 x Hand Towel Blue Paper Wrapped in a SMS Wrap 121cm x 121cm	3	DEF5900A	1 x Ophthalmic Drape 140cm x 160cm with 10cm x 10cm Fenestration 1 x Incise Drape 10cm x 20cm 1 x Plastic Drape 45cm x 60cm with Adhesive Tape 1 x Eye Shield 2 x Steristrips 12mm x 100mm - 0.27mm per 3 1 x Syringe 1ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 26g x 1/2" 1 x Anterior Chamber Cannula Angled Rycroft 30G x 7/8" (0.3 X 22mm) Hub Colour: Yellow 1 x Hydrodissector Curved Kellan 25G x 7/8" (0.5 x 22mm) Hub Colour: Orange 1 x Instrument Wipe 5 x Gauze Swabs 10cm x 10m X-Ray 3 x Gauze Ball XRD 1 x Eye Pad 6cm x 8cm Non Woven 1 x Sponge Holder 19cm Plastic 1 x Bowl 250ml Blue 2 x Bowl 500ml Blue 1 x Kidney Dish 700ml Yellow 1 x Kidney Dish 1L Blue 1 x Gallipot 150ml 2 x Pill Cup/ Measuring Cup 60ml 1 x Label Sheet Ophthalmic - 1 Set of 64 1 x Paper Disposal Bag 30cm x 16.5cm with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Unreinforced Gown Level 3 X-Large 2 x Hand Towel 40cm x 68cm Cotton Blue Wrapped in a SMS Wrap 121cm x 121cm	6

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
CATARACT SURGICAL PACK 1			MINOR LITHOTOMY SURGICAL PACK 1		
DEF4261A	1 x Ophthalmic Drape 140cm x 160cm with 10cm x 10cm Fenestration 1 x Plastic Drape 45cm x 61cm with Adhesive Tape 2 x Steristrips 12mm x 100mm - 0.27mm per 3 1 x Tegaderm Dressing 10cm x 12cm 1 x Syringe 1ml Luer Lock 1 x Syringe 1ml Luer Slip 3 x Syringe 3ml Luer Lock 2 x Needle 18g x 1.5" Blunt 1 x Needle 26g x 1/2" 2 x Anterior Chamber Cannula Angled Rycroft 30G x 7/8" (0.3 X 22mm) Hub Colour: Yellow 5 x Gauze Swabs 10cm x 10m X-Ray 4 x Prep Gauze Swab Small X-Ray 1 x Eye Shield 1 x Eye Pad 6cm x 8cm Non Woven 5 x Eye Spear 17mm x 7mm Pointed Tip 1 x Bowl 250ml Blue 2 x Kidney Dish 700ml Clear 3 x Gallipot 150ml 1 x Needle Counter Foam Strips/Magnetic 30 Count Yellow 1 x Label Sheet Ophthalmic - 1 Set of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 112cm x 229cm GOWN PACK 1 x Surgical Gown Level 3 Large 1 x Surgical Gown Level 3 X-Large 2 x Hand Towel 31.5cm x 61cm Wrapped in a SMS Wrap 121cm x 121cm	6	DEF5774	1 x Lithotomy T Drape 140cm x 195cm x 74cm with 6cm x 14cm Fenestration and Adhesive 2 x Lithotomy Leggings 80cm x 140cm 5 x Gauze Swabs 10cm x 10m X-Ray 1 x Prep Swab Large X-Ray 1 x Bowl 1L Blue 1 x Kidney Dish 700ml Blue 2 x Gallipot 150ml 1 x Labels for Theatre - 1 Set of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Fabri Reinforced Gown Level 3 Large 1 x Fabri Reinforced Gown Level 3 X-Large 2 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	4



MEDICAL DEVICE SYMBOLS

CUSTOM PACK CHECK SHEETS

Most of our Custom Procedure Packs include a Check Sheet listing the pack's componentry. The Check Sheet may also include the medical device symbols below.



MANUFACTURER - indicates the medical device manufacturer and their address.



DATE OF MANUFACTURE - indicates the date when the medical device was manufactured.



USE-BY DATE - indicates the date which the medical device must be used by.



BATCH CODE - indicates the manufacturer's batch code so that the batch or lot can be identified.



CATALOGUE NUMBER - indicates the manufacturer's catalogue number/ product code so that the medical device can be identified.



STERILISED USING ETHYLENE OXIDE - indicates that a medical device has been sterilised using ethylene oxide.



SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE - the outer packaging is a single sterile barrier. Inside, protective packaging is designed to prevent damage to the contents or to help with aseptic presentation.



MEDICAL DEVICE - indicates the item is a medical device.



DO NOT REUSE - indicates a medical device that is intended for one single use only.



DOES NOT CONTAIN NATURAL RUBBER LATEX - indicates that natural rubber latex as a material of construction within the medical device or the packaging of the medical device is not present.



KEEP AWAY FROM SUNLIGHT - indicates a medical device that needs protection from light sources and/or heat.



KEEP DRY - indicates a medical device that needs to be protected from moisture.



PLASTIC REDUCTION - procedure pack contains an item/s that are an alternative to a plastic product. Examples include a bowl made from sugar cane pulp instead of plastic, a product is wrapped with a paper band instead of double packed in plastic or a paper disposal bag has been included in place of a plastic bag.



SUSTAINABLE SOLUTIONS - procedure pack contains an item/s that is more sustainable than a previously used product. For example, our PurCotton swabs are manufactured using a manufacturing process that transforms raw cotton into a finished product in only 8 hours, unlike traditional gauze which can take up to 6 weeks, therefore minimising the bio-burden on the environment.

HOLLOWWARE



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BOWL

A wide selection of bowls that are recyclable and sturdy to use throughout the hospital setting. Graduated marking in ml and cc for easy measurement. Sterile and latex free.

DEF2901	500ml Blue	20
---------	------------	----



DEF2901

DEF2902	1L Blue	20
---------	---------	----



DEF2902

DEF2489	1.5L Blue Double Packed	20
---------	-------------------------	----



DEF2489

DEF2550	1.5L Blue	20
---------	-----------	----



DEF2550

DEF2553	2.5L Blue Double Packed	10
---------	-------------------------	----



NON STERILE

DEF5930NS	500ml Blue	390
DEF5929NS	1L Blue	340

SPLASH BOWL



DEF2650A

DEF2650A	Splash Bowl 6L wrapped in 1 x SMS Field 130cm x 150cm	6
DEF2434A	Splash Bowl 6L wrapped in 2 x SMS Fields 130cm x 150cm	6
DEF5689NS	6L Bowl Blue Non Sterile	70

ENVIRO-SUSTAINABLE SPLASH BOWL

DEF5653	1 x Enviro-Sustainable Splash Bowl 5.5L made from 100% plant-based composite Wrapped in a Blue SMS Field 130cm x 150cm	6
---------	---	---



DEF5653

BOWL DRAPE

This flexible bowl drape is made from a tough plastic material that fits over a re-usable Splash Bowl. Size is 132cm x 132cm, bowl size reinforcement is 36cm x 36cm and wrapped in a laminated field 50cm x 50cm.

80-13506-S	Bowl Drape	30
------------	------------	----



80-13506-S

NON STERILE

DEF5950NS	Bowl Drape	50
-----------	------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BOWL AND SPATULA



DEF1530

DEF1530	Bowl 600ml and Spatula Double Packed	10
DEF1531	Bowl 600ml Double Packed	10
DEF1533NS	Bowl 600ml Bulk Non Sterile	10

GALLIPOT



DEF867

DEF867	60ml, Height 2cm, Diameter 7.5cm	50
DEF860	60ml, Height 2cm, Diameter 7.5cm	250
DEF1147	150ml, Height 5cm, Diameter 7.5cm	25
DEF888	150ml, Height 5cm, Diameter 7.5cm	150
DEF2792	150ml, Height 5cm, Diameter 7.5cm Packet of 3	25



DEF2792

NON STERILE

DEF2362NS	150ml	200
-----------	-------	-----

JUG

These jugs are a clear recyclable plastic for easy viewing of contents. Comfortable handles and are graduated on both sides in ml and cc for easy measurement.

DEF2773	1L Tall Jug Increments in 100ml to 1L	12
---------	---------------------------------------	----



DEF2773

DEF3184A	1.2L Jug Increments in 100ml to 1L	12
DEF2092	1.2L Jug Increments in 100ml to 1L Double Packed	12



DEF2092

NON STERILE

DEF5931NS	1L Jug Increments in 100ml to 1L	100
DEF2092NS	1.2L Jug Increments in 100ml to 1L	70
DEF5928NS	2L Jug Increments in 200ml to 2L	20

MEASURING CUP



DEF459

DEF459	60ml, Height 5cm, Diameter 5cm	50
DEF460	60ml, Height 5cm, Diameter 5cm	250

NON STERILE

DEF1526NS	60ml, Height 5cm, Diameter 5cm	100
-----------	--------------------------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPECIMEN JAR

DEF758A	90ml - supplied with lid screwed on	50
DEF843A	90ml - supplied with lid screwed on	150



DEF758A

DEF842SP	120ml Yellow Lid	150
DEF840SP	250ml Yellow Lid	80
DEF1109SP	250ml Yellow Lid	25
DEF1139SP	500ml White Lid	36



DEF1139SP

SPECIMEN TUBE

DEF5968	10ml Freestanding Vee Bottom Tube with Natural Cap	50
DEF839C	10ml Freestanding Vee Bottom Tube with Natural Cap and Numbered Labels Pack of 3	300



DEF839B

TRAY

DEF2546	Square Fully Deep Tray 28cm x 24.5cm x 4.5cm	20
---------	---	----



DEF2546

DEF321	Dressing Tray 3 compartments 23cm x 9cm x 3cm	80
--------	--	----



DEF321

DEF2690A	Dressing Tray 3 Compartment 20cm x 25cm x 2.5cm	30
----------	--	----



DEF2690A

DEF2942	Sharps Tray 500ml 26.5cm x 9.5cm x 4cm	24
DEF3398	Sharps Tray 500ml 26.5cm x 9.5cm x 4cm Packet of 3	20



DEF2942

NON STERILE

DEF5952NS	Sharps Tray 500ml 26.5cm x 9.5cm x 4cm	200
DEF5953NS	Dressing Tray 24cm x 15cm x 4cm	500
DEF5954NS	Dressing Tray 3 compartments 23cm x 9cm x 3cm	1000

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
KIDNEY DISH					
Our disposable and recyclable polypropylene Kidney Dishes are flexible and easy to hold, with graduated markings up to either 700ml or 1L. Available sterile or non sterile for any area of the hospital.					
DEF2791	Blue 700ml	35			
		DEF2791			
DEF1174	Clear 700ml	35			
		DEF1174			
DEF2109	Yellow 700ml	35			
		DEF2109			
DEF2551	Blue 1L	20			
		DEF2551			
NON STERILE					
DEF5957NS	Blue 700ml	250			
DEF5955NS	Clear 700ml	250			
DEF5956NS	Yellow 700ml	250			
DEF5958NS	Green 700ml	250			

SHARPS SAFETY TRANSFER DEVICE		
Device ensures the safe transfer and passing of sharps between scrub nurses and surgeons by shielding the sharp section of the blade.		
DEF3551	Sharps Safety Transfer Device	48
		DEF3551

SCALPEL HOLDER		
DEF2820	Scalpel Holder	50
		DEF2820

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BOWL AND KIDNEY DISH SETS

PERI-OPERATIVE KIDNEY DISH SET

DEF2541A	2 x Kidney Dish Green 700ml 1 x Kidney Dish Clear 700ml 3 x Gallipot 150ml Wrapped in a Blue SMS Field 75cm x 75cm	20
----------	--	----

KIDNEY DISH SET

DEF2957	1 x Kidney Dish Yellow 700ml 1 x Kidney Dish Blue 700ml 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	20
---------	--	----



DEF2957

KIDNEY DISH AND GALLIPOT PACK

DEF409A	1 x Kidney Dish Clear 700ml 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	35
---------	---	----

KIDNEY DISH SET

DEF2837	1 x Kidney Dish Clear 700ml 1 x Kidney Dish Blue 700ml 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	20
---------	---	----

MAJOR BOWL SET

DEF2549	1 x Bowl 6L 1 x Bowl 2.5L 2 x Kidney Dish Clear 700ml 2 x Gallipot 150ml 1 x Tray with Two Compartments 24cm x 12cm x 4.5cm Wrapped in a Blue SMS Field 130cm x 150cm	5
---------	--	---

STANDARD BOWL SET

DEF3180B	1 x Bowl 6L 2 x Bowl 250ml 1 x Kidney Dish Clear 700ml 1 x Kidney Dish Yellow 700ml 1 x Kidney Dish Blue 1000ml Wrapped in a Blue SMS Field 150cm x 150cm	5
----------	---	---

BOWL SET

DEF2993	1 x Bowl 1.5L 1 x Bowl 1L 2 x Kidney Dish Blue 700ml 1 x Kidney Dish Yellow 700ml 4 x Gallipot 150ml Wrapped in a Blue SMS Field 91cm x 91cm	12
---------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ENVIRO-SUSTAINABLE KIDNEY DISH SET

DEF5671	1 x Enviro-Sustainable Kidney Dish 700ml made from 100% plant-based composite 1 x Gallipot 150ml Wrapped in a Blue SMS Field 75cm x 75cm	35
---------	--	----



ENVIRO-SUSTAINABLE MINOR BOWL SET

DEF5703	2 x Enviro-Sustainable Kidney Dish 700ml made from 100% plant-based composite 2 x Gallipot 150ml Wrapped in a Blue SMS Field 75cm x 75cm	20
---------	--	----



ENVIRO-SUSTAINABLE MINOR BOWL SET

DEF5792	1 x Enviro-Sustainable Bowl 1L made from 100% plant-based composite 2 x Enviro-Sustainable Kidney Dish 700ml made from 100% plant-based composite 2 x Gallipot 150ml 3 x Sharps Tray 500ml Yellow Wrapped in a Blue SMS Field 91cm x 91cm	10
---------	---	----



ENVIRO-SUSTAINABLE MAJOR BOWL SET

DEF5702	1 x Enviro-Sustainable Bowl 5.5L made from 100% plant-based composite 1 x Enviro-Sustainable Bowl 1L made from 100% plant-based composite 2 x Enviro-Sustainable Kidney Dish 700ml made from 100% plant-based composite 2 x Gallipot 150ml Wrapped in a Blue SMS Field 150cm x 200cm	6
---------	--	---



HOLLOWWARE FOR CUSTOM PACKS

ENVIRO-SUSTAINABLE HOLLOWWARE

Made from 100% plant-based composite, this range offers an alternative to plastic bowls used in the operating theatre.

- Designed to hold fluids and medical supplies
- Easy to see graduations of ml and oz
- Coated with biocompatible film
- Easy to hold and stable to use

This range is currently available for inclusion in custom procedure packs or non-sterile.



NGS700EB-250NS	Enviro-Sustainable Kidney Dish 700ml
----------------	--------------------------------------



NGS700EB-250NS

NGS550RB-80NS	Enviro-Sustainable Bowl 5.5L
---------------	------------------------------



NGS550RB-80NS

NGS32SB-400NS	Enviro-Sustainable Bowl 1L
---------------	----------------------------



NGS32SB-400NS

NGS1CTC-250NS	Enviro-Sustainable One Compartment Tray 23.5cm x 12cm x 5cm
---------------	---

NGS2CTC-250NS	Enviro-Sustainable Two Compartment Tray 22cm x 11cm x 4.5cm
---------------	---



NGS2CTC-250NS

The background of the entire page is a close-up, slightly blurred photograph of various surgical instruments, including forceps and scissors, arranged in blue plastic trays. A prominent yellow rectangular box is overlaid on the left side of the image, containing text. The title 'INSTRUMENTS AND INSTRUMENT SETS' is written in large, white, sans-serif capital letters across the top of the image, partially overlapping the yellow box and the background.

INSTRUMENTS AND INSTRUMENT SETS

Our single use sterile instruments will equip you with quality tools for your surgical and medical procedures. The technical design of each instrument is modelled on reusable instruments, whilst providing you the benefit of no repeat sterilisation and a reduced risk of cross infection.

Each instrument has a green dipped handle to visually identify it as single use.

Medical grade stainless steel instruments are made with an exacting mechanical finish.

They can be recycled into meltdown metal or incinerated to become non-toxic file ash.

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CLAMPS



DEF1781

DEF2298	13cm Smooth Jaw	10
DEF1781	18cm Smooth Jaw	1
DEF2775	18cm Smooth Jaw	10
DEF2275	18cm, Smooth Jaw , Packet of 2	5
DEF1575A	Tube Clamp with Lock12cm Blue Plastic	50



DEF1575A

DENTAL INSTRUMENTS KIT

DEF1457	1 x Mirror, 1 x Explorer, 1 x Forcep	100
---------	--------------------------------------	-----



DEF1457

FORCEPS

ADSON FORCECP



DEF2142

DEF2147	Plain 12.5cm	10
DEF2148	Plain 12.5cm	100
DEF2142	Toothed 12.5cm	10
DEF2143	Toothed 12.5cm	100

MICRO ADSON FORCECP



DEF2346A

DEF2135-10	Plain 12.5cm	10
DEF2135	Plain 12.5cm	100
DEF2346A	Toothed Super Fine 12.5cm	10
DEF2134-10	Toothed 12.5cm	10
DEF2134	Toothed 12.5cm	100



DEF2134-10

ROCHESTER PEAN FORCECP



DEF1782B

DEF2920	20cm Curved	10
DEF1782B	20cm Straight	10
DEF3233A	20cm Straight, Packet of 2	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
MOSQUITO ARTERY FORCEP			SPENCER WELLS ARTERY FORCEP		
 DEF2413-10			 DEF2277-10		
DEF2413-10	Straight 12.5cm	10	DEF2277-10	Straight 13cm	10
DEF2413	Straight 12.5cm	100	DEF2277	Straight 13cm	100
DEF5918	Straight 12.5cm Micro	10	DEF2415-10	Curved 13cm	10
DEF2414-10	Curved 12.5cm	10	 DEF2415-10		
 DEF2414-10			DEF2415	Curved 13cm	100
DEF2414	Curved 12.5cm	100	DEF2417-10	Straight 16cm	10
DEF2108	Curved 12.5cm Micro	10	 DEF2417-10		
 DEF2108			DEF2417	Straight 16cm	100
			DEF2416-10	Curved 16cm	10
			 DEF2416-10		
			DEF2416	Curved 16cm	100
			DEF3899	Curved 20cm	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CRILE ARTERY FORCEP



DEF2441

DEF2441	Straight 14cm	10
DEF2246	Straight 14cm	50
DEF1073	Straight 14cm	100
DEF2442	Curved 14cm	10
DEF5844	Curved 14cm	100



DEF2442

DRESSING FORCEP



DEF2152

DEF2152	Blunt 13cm	10
DEF2153	Blunt 13cm	100
DEF1072	Blunt 16cm	100



DEF1304

DEF1304	Pointed with Tip Protector 13cm	10
DEF1264A	Pointed with Tip Protector 13cm	100

GILLIES FORCEP



DEF2418-10

DEF2418-10	Toothed 15.5cm	10
DEF2418	Toothed 15.5cm	100

MCINDOE FORCEP



DEF2677

DEF2677	Non Toothed 15.5cm	10
DEF2948	Non Toothed 15.5cm	100

MAGILL FORCEP



DEF1807

DEF1807	Child 20cm	10
DEF2083	Adult 25cm	10



DEF2083

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GREEN PLASTIC FORCEP



DEF686

DEF686	Green Plastic 11cm	100
--------	--------------------	-----

SPECIALISED FORCEPS

DEF1578	Forcep Wilde 15cm Angled	10
---------	--------------------------	----



DEF1578

DEF2578	Kelman-McPherson Forcep 11cm	10
---------	------------------------------	----



DEF2578

DEF2443	Lister Sinus Forcep 15.5cm	10
---------	----------------------------	----



DEF2443

DEF1579	Bozman Forcep Straight 25cm	5
---------	-----------------------------	---



DEF1579

DEF2238	Packing / Hartmann Croc-Action Forcep 15cm	5
DEF2750	Hartmann Croc-Action Forcep 14cm	10



DEF2750

DEF3141	Fine Croc-Action Forcep 15.5cm Double Packed	10
---------	--	----



DEF3141

DEF3873	Alligator Forcep 21cm	25
---------	-----------------------	----



DEF3873

DEF2695	Iris Forcep 10cm Plain	10
---------	------------------------	----



DEF2695

DEF3205	Iris Forcep Fine 10cm Curved	10
---------	------------------------------	----



DEF3205

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------



DEF2683	Tilley Ear Forcep 15cm	10
DEF2875	Tenaculum Forcep 25cm	10



DEF3596	Roberts Forcep 25cm Curved	10
---------	----------------------------	----



DEF2174	Splinter Forcep Feichenfeldt 11cm	10
---------	-----------------------------------	----



DEF2149	Needle Holder 13cm	10
DEF2136	Needle Holder 13cm	100
DEF2150	Needle Holder 16cm	10



DEF2151	Needle Holder 16cm	100
DEF2707	Nievert Needle Holder 13cm	10



DEF3582B	Weitlaner Retractor 4/3 Blunt with a Enviro-Sustainable Kidney Dish 700ml made from 100% plant-based composite	35
----------	--	----



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SCALPEL

SCALPEL HANDLE



DEF2145

DEF2145	#3	10
DEF2137	#3	50
DEF2146	#3	100

SAFETY SCALPEL AND HANDLE

DEF2701	#11	100
DEF2702	#15	100
DEF3935	#23	100



DEF2702

SCISSORS

SHARP/BLUNT SCISSORS



DEF1870

DEF1870	Sharp/Blunt 13cm	100
---------	------------------	-----



DEF5894

DEF5894	Sharp/Blunt 13cm with Plastic Handles	100
---------	---------------------------------------	-----

SHARP/SHARP SCISSORS



DEF1311

DEF1311	Sharp/Sharp 12.5cm	10
DEF2276	Sharp/Sharp 12.5cm	100

BLUNT/BLUNT SCISSORS

DEF1637	Blunt/Blunt 12.5cm	100
---------	--------------------	-----



DEF1637

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

IRIS SCISSORS



DEF2141

DEF2144	Sharp/Sharp 7.5cm	10
DEF2141	Sharp/Sharp 11.5cm	10
DEF2133	Sharp/Sharp 11.5cm	100
DEF2460	Sharp/Sharp Curved 11.5cm	10

METZENBAUM SCISSORS



DEF2511

DEF2511	Curved 16cm	10
---------	-------------	----

TENOTOMY SCISSORS



DEF2584

DEF2584	Sharp/Sharp 10cm	10
---------	------------------	----

MAYO SCISSORS



DEF2408

DEF2408	Straight 17cm	10
DEF3361	Straight 17cm	50
DEF2117	Curved 17cm	10
DEF3370	Curved 17cm	50
DEF3794	Straight 23cm	50



DEF2117

VANNAS SCISSORS



DEF2973

DEF2972	Straight Fine Tip 9cm	10
DEF2973	Curved Fine Tip 9cm	10

UNIVERSAL SCISSORS



DEF2851

DEF2851	Black Handle 18cm	25
---------	-------------------	----

EPISIOTOMY SCISSORS



DEF5645

DEF5645	14cm	10
---------	------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SPECULUMS

EAR SPECULUM

DEF2240	Ear Speculum	5
---------	--------------	---



DEF2240

VAGINAL SPECULUM



DEF2399

DEF2399	Medium 14cm	20
DEF2400	Large 15.5cm	20

NOSE SPECULUM



DEF1356, DEF1357, DEF1358

DEF1356-10	Thudichum #1 20mm x 70mm	10
DEF1356	Thudichum #1 20mm x 70mm	25
DEF1357-10	Thudichum #2 22mm x 70mm	10
DEF1357	Thudichum #2 22mm x 70mm	25
DEF1358-10	Thudichum #3 25mm x 70mm	10
DEF1358	Thudichum #3 25mm x 70mm	25
DEF2292	Mod Wein 14cm with 22mm Tip	10
DEF2293	Mod Wein 14cm with 30mm Tip	10
DEF2294	Mod Wein 14cm with 34mm Tip	10



DEF2292, DEF2293, DEF2294

VAGINAL SPECULUM WITH LED LIGHT SOURCE

Single use disposable vaginal speculum with built-in LED light source. Each comes ready for use – no assembly required and no extra parts needed. High powered, battery operated LED light source that produces peak light for 30 plus minutes with no heat generated.

STERILE SPECULUM WITH LIGHT SOURCE

C020300	21mm x 116mm Small White	12
C020310	29mm x 136mm Medium Green	12
C020320	33mm x 148mm Large Blue	12

NON STERILE SPECULUM WITH LIGHT SOURCE

ER-SPEC-SMALL	21mm x 116mm Small White	18
ER-SPEC-MEDIUM	29mm x 136mm Medium Green	18
ER-SPEC-LARGE	33mm x 148mm Large Blue	16



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAL SPECULUM NON STERILE

ANOSPEC WITH LED LIGHT SOURCE

Single use disposable Anoscope with built-in LED light source from Cooper Surgical. Designed to aid in the examination of the anus to screen, diagnose and evaluate problems of the anal canal. This ready-to-use device eliminates the time and expense of reprocessing. Each Anospec comes complete with its own obturator.

ANO-SPEC-SMALL	90mm x 18mm – Slotted	20
ANO-SPEC-MEDIUM	96mm x 23mm – Slotted	16
ANO-SPEC-LARGE	103mm x 18mm – Bevelled	20

SPONGE HOLDER

DEF2609	Ramsey 24cm Stainless Steel	10
---------	-----------------------------	----



DEF2609

DEF1923	Plastic Blue 19cm	10
---------	-------------------	----



DEF1923

DEF2071	Plastic Blue 24cm	10
---------	-------------------	----

DEF3447	Plastic Blue 24cm, Packet of 2	50
---------	--------------------------------	----



DEF2071

STAPLE REMOVER



DEF1826

DEF1826	Staple Remover Plastic Handle	50
---------	-------------------------------	----

SUCTION INSTRUMENTS

YANKAUER SUCTION INSTRUMENT

Popular suction instrument with a large lumen providing rapid aspiration whilst resisting clogging. Removes fluid and debris from a surgical site or body orifice.



DEF2329, DEF2925 DEF3030

DEF2329	With Tip 30cm	10
DEF2925	No Tip Medium Bore 30cm	10
DEF3030	No Tip Straight 30cm	10
DEF3190	Plastic Vented	50
DEF3192	Plastic Non Vented	50



DEF3190, DEF3192

DEF3594	Plastic with On/Off Switch	50
---------	----------------------------	----



DEF3594

POOLE SUCTION INSTRUMENT

Designed to evacuate pooled blood, fluid and debris during deep abdominal surgeries and caesarean sections.

DEF2554	25cm	10
---------	------	----



DEF2554

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FRAZIER SUCTION INSTRUMENT

Used in open surgical procedures where fine, delicate suctioning is necessary to remove debris from the surgical site. Comes with an obturator to clear the instrument of obstructions during use.

FRAZIER SUCTION INSTRUMENT STANDARD



DEF2325, DEF2326, DEF2327, DEF2328

DEF2325	6 Fr Gauge 17cm	10
DEF2326	8 Fr Gauge 17cm	10
DEF2327	10 Fr Gauge 17cm	10
DEF2328	12 Fr Gauge 17cm	10

FRAZIER SUCTION INSTRUMENT LONG



DEF3482, DEF3483, DEF3484

DEF3562	6 Fr Gauge 20.5cm	50
DEF3482	8 Fr Gauge 20.5cm	50
DEF3483	10 Fr Gauge 20.5cm	50
DEF3484	12 Fr Gauge 20.5cm	50
DEF3485	14 Fr Gauge 20.5cm	50

FRAZIER SUCTION INSTRUMENT EXTRA LONG

The suction tube connector opening allows for a secure fit of a 20ml Luer Slip syringe for flushing of the suction tip.

DEF3080A	8 Fr Gauge 23cm	50
DEF3081A	10 Fr Gauge 23cm	50
DEF3082A	12 Fr Gauge 23cm	50
DEF3486	14 Fr Gauge 23cm	50

FUJITA SUCTION INSTRUMENT

Used in neurosurgical procedures to clear surgical sites of fluids. This instrument has an integrated tear-drop tapered finger control.

FUJITA SUCTION INSTRUMENT SHORT



DEF3500, DEF3502, DFE3501

DEF3500	4 Fr Gauge 18cm	50
DEF3501	5 Fr Gauge 18cm	50
DEF3502	7 Fr Gauge 18cm	50
DEF5915	9 Fr Gauge 18cm	50

FUJITA SUCTION INSTRUMENT STANDARD

DEF3503	4 Fr Gauge 23cm	50
DEF3504	7 Fr Gauge 23cm	50
DEF5916	9 Fr Gauge 23cm	50

PAEDIATRIC SUCTION INSTRUMENT

DEF2497	Small 27cm	10
DEF2498	Medium 29cm	10
DEF5798	Plastic	25



DEF2498, DEF2497

ADSON FRAZIER SUCTION INSTRUMENT

DEF3857	6 Fr Gauge 21cm	20
DEF2817	8 Fr Gauge 21cm	20
DEF2818	10 Fr Gauge 21cm	20
DEF2819	12 Fr Gauge 21cm	20



DEF3857, DEF2817, DEF2818, DEF2819

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CARDIAC SUCTION INSTRUMENT

DEF3721	34cm	25
		
DEF3721		
DEF5502	34cm	25
		
DEF5502		
DEF5557	Fine Tip Narrow	25
DEF5558	Fine Tip Standard	25
		
DEF5557		

TRACHEAL DILATOR

DEF2373	3 Prong	1
		
DEF2373		
DEF1467	Trousseau 2 Prong	1
DEF2585	Trousseau 2 Prong	5

RHOTON SUCTION TUBE

		
DEF2807		
DEF3787	5FG x 10cm	10
DEF3786	7FG x 10cm	10
DEF3871	5FG x 13cm	10
DEF3858	7FG x 13cm	10
DEF2807A	10FG x 13cm	10

WIRE CUTTER

		
DEF1800		
DEF1800	Universal 12cm	50
DEF2003	Universal 12cm	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INSTRUMENT PACKS

PLASTIC SET

DEF1328A	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep McIndoe 15.5cm	10
----------	--	----

SUTURE KIT

DEF1811	1 x Safety Scalpel # 11 with Handle 1 x Needle Holder 16cm	25
---------	---	----

FINE SUTURE SET

DEF1656	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain	25
---------	--	----

SUTURE SET SMALL

DEF1560	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 13cm	10
---------	--	----



DEF1560

SUTURE SET LARGE

DEF1561	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Needle Holder 16cm	10
---------	--	----

SUTURE KIT

DEF1494	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Blunt 13cm Plain 1 x Needle Holder 16cm	25
---------	--	----

FINE FACIAL MICRO SUTURE SET

DEF1636-10	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm	10
DEF1636	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm	100

FINE PLASTIC SUTURE TRAY

DEF2051	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Scalpel Handle # 3 2 x Non Woven Swabs 7.5cm x 7.5cm 5 x Gauze Swabs Plain 10cm x 10cm Wrapped in a Blue Plastic Field 50cm x 50cm	10
---------	---	----

PAEDIATRIC SUTURE SET

DEF2955	1 x Scissors 13cm Sharp/Blunt 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm	25
---------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
FINE PLASTIC SUTURE TRAY			SUTURE TRAY		
DEF2507B	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm 1 x Scalpel Handle # 3 5 x Gauze Swabs 10cm x 10cm Plain 2 x Non Woven Swabs 7.5cm x 7.5cm Wrapped in a Blue Plastic Field 75cm x 75cm	10	DEF2050A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Towel Blue Cotton 40cm x 68cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Scissor 13cm Sharp/Blunt 1 x Forcep Adson 12.5cm Plain 1 x Needle Holder 13cm 5 x Gauze Swab 10cm x 10cm Plain 2 x Non Woven Swabs 7.5cm x 7.5cm Wrapped in a Blue Plastic Field 50cm x 50cm	10
SUTURE KIT			FINE PLASTIC SUTURE TRAY WITH DRAPE		
DEF2835A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Scissor 13cm Sharp/Blunt 1 x Forcep Adson 12.5cm Micro Toothed 1 x Forcep Green Plastic 1 x Needle Holder 13cm 1 x Gauze Swab 10cm x 10cm Plain Wrapped in a Blue Plastic Field 75cm x 75cm	10	DEF3468	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Super Fine Toothed 1 x Needle Holder Nievert 13cm 5 x Gauze Swabs 10cm x 10cm Plain Placed into a plastic bag 20cm x 20cm	40
SUTURE SET			FINE SUTURE SET		
DEF2422	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Forcep Green Plastic 1 x Needle Holder 16cm 5 x Non Woven Swabs 7.5cm x 7.5cm 5 x Non Woven Balls Wrapped in a Blue Plastic Field 50cm x 50cm	60	DEF2226	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel T Small 31.5cm x 34cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Forcep Green Plastic 1 x Needle Holder 13cm 10 x Gauze Swabs 10cm x 10cm Plain 1 x Syringe 5ml Luer Slip 1 x Needle 25g x 1" Wrapped in a Blue Plastic Field 75cm x 75cm	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SUTURE KIT

DEF1671D	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Forcep Blunt 13cm 1 x Forcep Green Plastic 1 x Needle Holder 16cm 1 x Scalpel Blade #15 and Handle 6 x Gauze Swabs 10cm x 10cm Plain Wrapped in a Blue Plastic Field 50cm x 50cm	10
----------	---	----

SUTURE PACK

DEF5926A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Scalpel Blade #15 and Handle 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Forcep 13cm Pointed 1 x Forcep Green Plastic 1 x Needle Holder 13cm Nievert 3 x Gauze Swabs 10cm x 10cm Plain 5 Non Woven Balls Wrapped in a Blue Plastic Field 50cm x 50cm	40
----------	--	----

SUTURE KIT

DEF2270C	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Towel White 2 Ply 40cm x 42cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Plain 1 x Forcep Adson 12.5cm Micro Toothed 1 x Forcep Green Plastic 1 x Needle Holder 13cm 10 x Gauze Swabs 10cm x 10cm X-ray Wrapped in a Blue Plastic Field 50cm x 50cm	10
----------	---	----

INSTRUMENT DRESSING PACK

DEF1677A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Scissors 13cm Sharp/Blunt 1 x Forcep 13cm Blunt 1 x Forcep 13cm Pointed 1 x Forcep Green Plastic 5 x Gauze Swabs 7.5cm x 7.5cm Plain Wrapped in a Blue Plastic Field 75cm x 75cm	60
----------	---	----

INSTRUMENT PACK

DEF1913B	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Toothed 1 x Forcep Adson 12.5cm Plain 1 x Needle Holder 13cm	10
----------	---	----

MATERNITY PRODUCTS



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BABY CUP WITH LID

DEF2129SP	Baby Cup with Lid	100
-----------	-------------------	-----



DEF2129SP

MATERNITY PAD



DEF721

DEF721	10cm x 34cm	50
--------	-------------	----

NON STERILE

DEF238NS	10cm x 34cm, Packet of 10	25
----------	---------------------------	----

BABY RECEIVING BLANKET

Soft, cotton flannel blanket used for wrapping a baby after birth.



DEF3227

DEF3227	Flannel 75cm x 100cm	40
---------	----------------------	----

DEF3621	Flannel 75cm x 100cm Packet of 2	40
---------	-------------------------------------	----

PLACENTA BASIN

DEF2822	26cm x 18.5cm x 7cm	12
---------	---------------------	----



DEF2822

NON STERILE

DEF5949NS	26cm x 18.5cm x 7cm	100
-----------	---------------------	-----

BABY CAP

DEF5969NS	One Size Non Sterile (suitable for newborn to 6 months)	160
-----------	--	-----



DEF5969NS

UMBILICAL CORD CUTTING DEVICE

DEF2535	Umbilical Cord Cutting Scissor	10
---------	--------------------------------	----



DEF2535

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MATERNITY PACKS

MATERNITY SUTURE SET

DEF2131A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Needle Holder 18cm 1 x Forcep Gillies 15.5cm Toothed 1 x Forcep Mosquito Artery 12.5cm Straight 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	10
----------	--	----

DELIVERY SUTURE SET

DEF2118	1 x Scissors Mayo 17cm Straight 1 x Needle Holder 16cm 1 x Forcep Bonney 18cm ½ tooth 2 x Forcep Crile 14cm Curved	25
---------	---	----

EPISIOTOMY SUTURE SET

DEF2671A	1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Needle Holder 16cm 1 x Forcep Gillies 15.5cm Toothed Wrapped in a Blue SMS field 91cm x 91cm	10
----------	---	----

LABOUR PACK

DEF2526B	2 x Forcep Rochester Pean 20cm Straight 1 x Bowl 600ml Green 1 x Bowl 2.5L 1 x Gallipot 150ml Wrapped in a Blue SMS Field 91cm x 91cm	10
----------	---	----

IUD INSERTION KIT

This procedure Pack is for an IUD insertion procedure.

The kit is sterile and disposable, ensuring all components are available for every insertion.

Reduced Packaging waste, plus the benefit of no repeat sterilisation costs.

DEF3256	2 x Swab Sticks 20cm Paper Shaft 2 x Gallipot 150ml 1 x Kidney Dish 1000ml Blue 1 x Vaginal Speculum Large 1 x Scissor Mayo 23cm Straight 1 x Forcep Tenaculum 25cm 1 x Sponge Holder 19cm Plastic 1 x Uterine Sound (Hysterometer CH14) 4.67mm with Nozzle Wrapped in a Blue Plastic Field 75cm x 75cm	10
---------	---	----



DEF3256

DELIVERY BOWL SET

DEF3471	1 x Placenta Basin 26cm x 18.5cm x 7cm 1 x Gallipot 150ml 1 x Kidney Dish 1000ml Blue Wrapped in a Blue SMS Drape 91cm x 91cm	10
---------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DELIVERY SET

DEF3931C	1 x Placenta Basin 26cm x 18.5cm x 7cm 3 x Forcep Rochester Pean 20cm Curved 1 x Scissor Mayo 17cm Curved 1 x Adhesive Absorbent Pad 55cm x 70cm Wrapped in a Blue SMS Drape 91cm x 91cm	12
----------	--	----

DELIVERY BUNDLE

DEF3410A	1 x Kidney Dish 700ml Blue 1 x Kidney Dish 1000ml Blue 2 x Forcep Rochester Pean 22cm Straight 1 x Scissor Mayo 17cm Straight 1 x Scissor Mayo 17cm Curved 3 x Umbilical Cord Clamp 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue SMS Drape 91cm x 91cm	20
----------	--	----

DELIVERY PACK

DEF4101G	1 x Gallipot 150ml 1 x Kidney Dish 700ml Clear 2 x Kidney Dish 1000ml Blue 1 x Bowl 1 Litre Blue 10 x Abdominal Sponges 23cm x 23cm X-ray Detectable with Tape 1 x Square Drape 98cm x 98cm Wrapped in a Back Table Cover 112cm x 228cm GOWN PACK 1 x Unreinforced Gown Level 3 Extra Large 1 x Hand Towel Blue Cotton 40cm x 68cm Wrapped in a SMS Wrap 91cm x 91cm	5
----------	--	---

MEDICAL IMAGING AND CATH LAB PRODUCTS




CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANGIO PROCEDURAL DRAPES

FEMORAL ANGIO RADIAL DRAPE

We also have a wide range of designs of Angiography Drapes available as a component in our Custom Procedure Packs.

DEF3801	240cm x 370cm with Airlaid Fabric PE Panels and Reinforcement Dual 7cm Round Fenestration with Built-in Cut-Out Incise Dual 7.5cm x 6cm Oval Fenestration with Adhesive Wrapped in a SMS 91cm x 91cm Field	12
---------	---	----



PROBE COVER



DEF301

DEF301	Probe Cover 10cm x 100cm with a Sealed End	100
DEF648	Drawstring Probe Cover 23cm x 100cm	25
DEF1482	Probe Cover 23cm x 100cm with a Sealed End	100
DEF2630	Drawstring Probe Cover 23cm x 150cm	25
DEF1088	Ultrasound Probe Cover 23cm x 80cm with a Sealed End	100
DEF643	Ultrasound Probe Cover 23cm x 150cm with a Sealed End	25
DEF2409	Ultrasound Probe Cover 25cm x 150cm with a Sealed End and a Latex Free Band #64	100

EQUIPMENT COVERS

IMAGE INTENSIFIER COVER



DEF764

DEF764	Cap 61cm x 91cm, Double Packed	40
DEF770	Cover 76cm x 106cm, Double Packed	50
DEF3247A	Cover 90cm x 135cm	25
DEF767	Cover 91cm x 115cm	50
DEF774	Cover 120cm x 140cm	40
DEF1348	Cover 102cm x 150cm with 2 Securing Tapes, Cuff and Labels, Double Packed	15
DEF1045	Cover 105cm x 150cm with 1 Securing Tape, Cuff and Labels, Double Packed	15

LEAD SCREEN COVER



DEF1245

DEF1245	80cm x 80cm with Elastic Opening	50
DEF3281	75cm x 95cm	50
DEF2073	Fluoroscope/Lead Screen Cover 76cm x 91cm	50
DEF1161	80cm x 110cm	50
DEF805	90cm x 107cm	50
DEF3669	91cm x 112cm	50
DEF1762	104cm x 105cm	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GUIDEWIRE BOWL

Designed with a built up lip and recessed bottom to hold guidewire in place and allows for easy removal.

DEF2579	Guidewire Bowl 2500ml	10
---------	-----------------------	----



DEF2579

GUIDEWIRE DRAPES

GUIDEWIRE BOOK

For storage of multiple wires for quick and easy access. Made from low linting absorbent material. Guidewire Book has open outer edges. Guidewire Wallet has sewn edges.



DEF3577

DEF3577	6 Page Book with Adhesive Tab	50
DEF3576	6 Page Wallet	50
DEF5962	9 Page Wallet	60

SPLASH STOP WASTE COLLECTION

An innovative alternative to messy, open basins.

A closed bowl with a secure, non-splatter lid which minimises staff and patient exposure to blood-borne pathogens.

Contains a secure, non-splatter lid with a tapered insertion site and a one-step biohazard warning label.

77-301370	Waste Collection Bowl	20
77-301371	Waste Collection Bowl with Solidifier	20



71-301371

MEDICAL IMAGING & CATH LAB PACKS

We have an extensive range of packs to cover Angio, Pacemaker, TAVI and hybrid-theater requirements. Contact Defries to discuss your custom pack requirements.

X-RAY INSTRUMENT PACK

DEF2485	1 x Scissors 13cm Sharp/Blunt 1 x Forcep Gillies 15.5cm Toothed 1 x Forcep Crile 14cm Straight 1 x Needle Holder 16cm 1 x Forcep Artery Mosquito 12.5cm Curved 1 x Forcep McIndoe 15.5cm Plain	10
---------	---	----

ANGIO INSTRUMENT SET

DEF2768	1 x Scissors 12.5cm Sharp/Sharp 1 x Forcep Artery Mosquito 12.5cm Straight 1 x Forcep Artery Mosquito 12.5cm Curved 1 x Sponge Holder 24cm Plastic Blue	10
---------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANGIO PACK

DEF2903M	1 x Angio Drape 343cm x 221cm with Femoral Fenestration 3 x Adhesive Drape 55cm x 70cm with Adhesive 2 x Absorbent Drape 100cm x 150cm 1 x Plastic Drape 120cm x 140cm 3 x Absorbent Foam Applicator 15cm Blue 1 x Needle Arterial 18G x 7cm Pink Hub 1 x Needle 18g x 1.5" Blunt 2 x Needle 23g x 1.25" 1 x Syringe 3ml Luer Lock 1 x Syringe 5ml Luer Lock 2 x Syringe 10ml Luer Lock 2 x Syringe 20ml Luer Slip Eccentric Tip 1 x Forcep Artery Mosquito 12.5cm Straight 1 x Scissor 13cm Sharp/Blunt 1 x Safety Scalpel #11 25 x Non-Woven Swabs 10cm x 10cm 1 x Fluoroscope/Lead Screen Cover 76cm x 91cm 1 x Image Intensifier Cover 61cm x 91cm 1 x Pure Cotton Hand Towel 32cm x 60cm 2 x Bowl 500ml Blue 1 x Bowl 1.5 Litre Blue 2 x Gallipot 150ml 1 x Kidney Dish 700ml Clear 1 x Label Sheet - 1 Set of 92 Wrapped in a Back Table Cover 112cm x 228cm GOWN PACK 1 x Comfort Gown Reinforced Fabric Large 1 x Comfort Gown Reinforced Fabric XL 4 x Pure Cotton Hand Towel 32cm x 60cm Wrapped in a SMS Drape 91cm x 91cm	4
----------	---	---

PACEMAKER PACK

DEF5732	1 x Pacemaker Drape 208cm x 292cm with Reinforced Dual Fenestration, Built-in Incise and Pouch 5 x Absorbent Foam Applicator 15cm Blue 1 x Needle 18g x 1.5" Blunt 2 x Needle 23g x 1.25" 1 x Syringe 10ml Luer Lock 3 x Syringe 20ml Luer Lock 1 x Scissor 13cm Sharp/Blunt 1 x Scalpel Blade #15 and Handle 15 x Swabs 10cm x 10cm X-ray 5 x Abdominal Composite Sponges 30cm x 30cm X-ray 1 x Lead Screen Cover 90cm x 107cm 1 x Image Intensifier Cover 61cm x 91cm 1 x Splash Stop Plus with Solidifier Waste Collection System 3 x Bowl 1.5 Litre Blue 3 x Gallipot 150ml 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow 1 x Label Sheet - 1 Set of 75 1 x Clipseal Bag (with small components enclosed) Wrapped in a Back Table Cover 152cm x 228cm	4
---------	--	---



MEDICAL LABELS

Please contact Customer Service on 1300 550 278 to discuss your label requirements. Our durable pre-printed Sterile and Non Sterile labels are used throughout the hospital in different wards such as the Emergency Department, General Wards, Sterilising Dept, Operating Theatre and Intensive Care.

Our labels are designed to clearly and easily identify items in a hospital on a professional, easy to read label.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

GENERAL USE NON STERILE LABELS

WHITE BACKGROUND

These labels are available with a white background and a coloured print in either blue, black, red, yellow, green and magenta. Rolls of 100 or 1000 labels.

STANDARD SIZE	ANAESTHETIC SIZE
---------------	------------------

21mm x 64mm	12mm x 40mm
-------------	-------------

32mm x 50mm	15mm x 30mm
-------------	-------------

32mm x 70mm	
-------------	--

50mm x 70mm	
-------------	--

40mm x 102mm	
--------------	--

40mm x 170mm	
--------------	--

80 x 102mm – provided as 2 rolls of 500	
---	--

FULL COLOUR BACKGROUND

Any colour combination can be made for the background or print for this type of label. Rolls of 100 or 1000 labels.

STANDARD SIZE	ANAESTHETIC SIZE
---------------	------------------

21mm x 64mm	12mm x 40mm
-------------	-------------

32mm x 50mm	15mm x 30mm
-------------	-------------

32mm x 70mm	
-------------	--

50mm x 70mm	
-------------	--

40mm x 102mm	
--------------	--

80 x 102mm – provided as 2 rolls of 500	
---	--

PIGGY BACK LABELS

These labels are made up of two labels where one labels sits on top of (or Piggy backs) the other label. The top label can be removed and adhered. Available in rolls or 100 or 1000 labels.

WHITE BACKGROUND

Available with a coloured print in either blue, black, red, yellow, green and magenta.

PBLABELS – BLANK	PBLABELS – PRINTED
------------------	--------------------

25mm x 50mm	25mm x 50mm
-------------	-------------

32mm x 50mm	32mm x 50mm
-------------	-------------

36mm x 105mm	36mm x 105mm
--------------	--------------

50mm x 70mm	50mm x 70mm
-------------	-------------

COLOURED BACKGROUND

Black print. Background colour as indicated.

25mm x 50mm - CYAN	
--------------------	--

25mm x 50mm - YELLOW	
----------------------	--

32mm x 50mm - YELLOW	
----------------------	--

50mm x 70mm - YELLOW	
----------------------	--

FLUORESCENT LABELS

Black print on all backgrounds. Red print also available on Yellow background. Available in Orange, Yellow, Red, Green, Blue or Pink.

32mm x 50mm – ONLY Fluorescent Red	
------------------------------------	--

30mm x 70mm	
-------------	--

50mm x 70mm – ONLY Fluorescent Green	
--------------------------------------	--

80mm x 102mm – PLUS Fluorescent Purple	
--	--

STERILE ANAESTHESIA DRUG LABELS

Used in the operating theatre to label syringes during procedures to help prevent medication mix ups. These easy to use labels come pre-printed and pre-packed in a peel down sterile packet. Available in any colour print on a white background. Supplied in boxes of 100.

15mm x 30mm	12mm x 40mm
-------------	-------------

Prepacked in 1's	Prepacked in 1's
------------------	------------------

Prepacked in 2's	Prepacked in 2's
------------------	------------------

Prepacked in 3's	Prepacked in 3's
------------------	------------------

Prepacked in 6's	Prepacked in 7's
------------------	------------------

Prepacked in 12's	Prepacked in 14's
-------------------	-------------------

Prepacked in 18's	Prepacked in 21's
-------------------	-------------------

Prepacked in 24's	Prepacked in 28's
-------------------	-------------------

Prepacked in 30's	
-------------------	--

These labels can have any colour combination for the background or print. Available in boxes of 100.

15mm x 30mm	12mm x 40mm
-------------	-------------

Prepacked in 1's	Prepacked in 1's
------------------	------------------

Prepacked in 2's	Prepacked in 2's
------------------	------------------

Prepacked in 3's	Prepacked in 3's
------------------	------------------

Prepacked in 6's	Prepacked in 7's
------------------	------------------

Prepacked in 12's	Prepacked in 14's
-------------------	-------------------

Prepacked in 18's	Prepacked in 21's
-------------------	-------------------

Prepacked in 24's	Prepacked in 28's
-------------------	-------------------

Prepacked in 30's	
-------------------	--

PIGGY BACK LABEL ROLL

36mm x 105mm labels on a blank cylinder enabling health services to print off their own labels. 3200 labels on the roll.

LAB3612-ROLL

A healthcare worker, likely a nurse or doctor, is shown from the chest up, wearing full personal protective equipment (PPE). She is wearing a blue surgical cap, a light blue surgical mask covering her nose and mouth, and a blue surgical gown. She is looking directly at the camera with a slight smile. The background is a blurred clinical setting with white walls and a door.

PERSONAL PROTECTIVE APPAREL

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ISOLATION GOWN

Full coverage level 2 isolation gown. Made from a two-layered polypropylene/polyethylene laminate material. With hook and loop neck closure, elastic wrist and thermal bonded sleeve seams.



60-3109

60-3109	Large, Length 117cm, Sleeve 58cm Chest 74cm, Yellow, Knitted Cuff Non Sterile	100
60-3110	Extra Large, Length 122cm Sleeve 61cm, Chest 79cm, Yellow Knitted Cuff, Non Sterile	100

CHEMOTHERAPY GOWNS

Chemoguard Chemotherapy Gowns are a lightweight gown made from a poly-coated spunbound fabric with a closed-back design for full coverage. Adjustable tie neck closure for ease of donning and doffing. A generous fit with long sleeves and fitted knit cuffs.

60-3200	Large Blue Non Sterile	60
60-3201	Extra Large Blue Non Sterile	60

SPEyes EYEWEAR



Our popular DeRoyal SPEyes eyewear is designed to minimise exposure to pathogens from waste fluid. They combine innovation, quality, value, and comfort for clinicians where personal protective equipment is a vital part of self-care.

- Latex free

- Wrap-around disposable lens for complete protection
- UltraClear™ technology prevents glare, fog, and static
- Frame fits completely against forehead for splash protection from above
- Sealed lens mounting for added fluid protection
- Clip-On/Off design for quick and easy changing
- Coloured top strip allows finger-print-free handling



23-505

23-505	SPEyes Eye Shields Pre-assembled in Assorted Colours, Non Sterile	100
23-506	Dispenser box with pre-assembled SPEyes Eye Shields, Non Sterile (20 per box 6 boxes per case)	120



23-506

23-501	SPEyes Eye Shield Replacement Lens, Non Sterile	100
23-502	SPEyes Eye Shield Frame Non Sterile	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FACE SHIELD

Our Face Shields feature an anti-fog visor with vented foam. Each mask has a peel away protective cover with tabs to prevent scratches. Each shield is 23cm long and 30cm wide.



DEF5526		
DEF5526	Face Shield Clear Non Sterile	10
FACE-SHIELD-FULL	Face Shield Clear Non Sterile	300

OVERSHOES



DEF1145		
DEF1145	Non Skid Universal Overshoes Polypropylene Non Sterile	100



SURGICAL GOWNS

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL GOWNS

Our single-use surgical gowns have specifically been developed for the Australian and New Zealand environment. Our range encompasses gowns for short and dry procedures through to high fluid and lengthy procedures.

OUR COMFORT RANGE – OUR SOFTEST AND MOST COMFORTABLE GOWNS

- Comfort Surgical Gown – Level 2.
- Comfort Unreinforced with Glove Assist – Level 3 gown with a reinforced sleeve to allow extra strength during gloving, but still cool, soft and ideal for short and dry procedures.
- Comfort Plus Fabric Reinforced Surgical Gown – Level 3 gown includes fluid resistant fabric reinforced panels yet remains lightweight and durable.
- Comfort Plus Poly Reinforced Surgical Gown – Level 3 gown which includes poly reinforced impervious panels for added protection.
- Comfort Zoned Protection – Fully impervious zoned panel on front and sleeves. Made from high tech BVB (breathable viral barrier) fabric, but still soft, comfortable and cool which meets AAMI Level 4 fabric standards.

OUR CLASSIC RANGE:

- Classic Standard Surgical Gown – suitable for most surgical procedures.
- Classic Plus Reinforced Surgical Gown – Includes a reinforced panel on the front and sleeves, for surgical cases requiring extra protection.
- Coolmax Surgical Gown – breathable and cool, these gowns are notable for their superb comfort and being fully impervious. With AAMI Level 4 fabric.

Designed in close co-operation with hospital staff, our Surgical Gowns are cool and breathable with an easy fit, while offering excellent protection against strike through from microbial and liquid penetration.




CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COMFORT SURGICAL GOWN

LEVEL 2 COMFORT UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 35gsm T-cut Surgical Gown ideal for dry and short procedures.




Comfort Unreinforced

DEF5316	Medium 119cm	25
DEF5317	Large 127cm	25
DEF5318	X-Large 137cm	25
DEF5319	XL-XLong 149cm	25

LEVEL 3 COMFORT UNREINFORCED WITH GLOVE ASSIST

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41gsm T-cut Surgical Gown ideal for dry and short procedures. Has innovative Glove Assist fabric reinforcement to sleeve area to ensure extra strength. Comfortable and cool to wear.




Comfort Unreinforced With Glove Assist

DEF5301	Medium 119cm	24
DEF5302	Large 127cm	24
DEF5303	X-Large 137cm	24
DEF5339	X-Large Twin Pack	24
DEF5304	XL-XLong 150cm	24
DEF5305	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS FABRIC REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41gsm T-cut Surgical Gown with fabric reinforcement on front and sleeves. Fluid resistant on panels. Comfortable to wear.




Comfort Plus Fabric Reinforced

DEF5306	Medium 119cm	24
DEF5307	Large 127cm	24
DEF5308	X-Large 137cm	24
DEF5328	X-Large Twin Pack	24
DEF5309	XL-XLong 150cm	24
DEF5310	XX-Large 155cm	24
DEF5329	XX-Large Twin Pack	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 3 COMFORT PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41 gsm T-cut Surgical Gown which includes poly reinforced panels. Durable and impervious to fluid on panels. Comfortable to wear.



Comfort Plus Poly Reinforced

DEF5311	Medium 119cm	24
DEF5312	Large 127cm	24
DEF5327	Large & X-Large Twin Pack	24
DEF5313	X-Large 137cm	24
DEF5326	X-Large Twin Pack	24
DEF5314	XL-XLong 150cm	24
DEF5315	XX-Large 155cm	24
DEF5330	XX-Large Twin Pack	24

LEVEL 4 COMFORT ZONED PROTECTION

Surgical Gown T-cut wrapped in a 91cm x 91cm SMS field with two hand towels. Made from a 41gsm SMMS material with a poly-ethylene breathable reinforcement on the sleeves and the front of the gown providing protection against fluids, bacteria and viruses.



Comfort Zoned Protection

DEF5334	Medium 109cm	24
DEF5331	Large 117cm	24
DEF5338	Large & X-Large Twin Pack	24
DEF5332	X-Large 127cm	24
DEF5336	X-Large Twin Pack	24
DEF5333	XX-Large 133cm	24
DEF5759	XL-XLong 135cm	24
DEF5337	XX-Large Twin Pack	24
DEF5335	XXX-Large 135cm	24



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 4 COMFORT MAX

Cool and comfortable T-cut surgical gown made from a 60gsm breathable impervious trilaminate. Provides protection from fluid, bacteria and viruses. Wrapped in a SMS field with two hand towels.




Comfort Max

DEF6020	Medium 119cm	24
DEF6021	Large 127cm	24
DEF6022	X-Large 137cm	24
DEF6023	X-Large Long 150cm	24
DEF5024	XX-Large 155cm	24

LEVEL 3 CLASSIC PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown with impervious poly reinforcement. Colour coded neck band – purple for level 3 Gown.




Classic Plus Poly Reinforced

DEF5005A	Medium 119cm	25
DEF5006A	Large 130cm	25
DEF5007A	X-Large 140cm	25
DEF5009A	XL-XLong 157cm	25
DEF5054	XX-Large with Long Arms 157cm	25

CLASSIC SURGICAL GOWN

LEVEL 3 CLASSIC UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown ideal for short and dry procedures.




Classic Unreinforced Gown

DEF5000A	Medium 119cm	25
DEF5001A	Large 130cm	25
DEF5002A	X-Large 140cm	25
DEF5004A	XL-XLong 157cm	25

LEVEL 4 COOLMAX

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile Level 4, trilaminate material which is a BVB fabric – offers a Breathable Viral Barrier, providing excellent protection against fluids, bacteria and viruses.

High quality T-cut Surgical Gown ideal for high fluid procedures or lengthy procedures. Colour Coded neck band – blue for level 4 Gown.



Coolmax

DEF5010B	Medium 119cm	25
DEF5011B	Large 130cm	25
DEF5012B	X-Large 140cm	25
DEF5013B	XL-XLong 157cm	25
DEF5014B	XXL-Large 157cm	25

CODE

DESCRIPTION

BOX QTY.

CODE


DESCRIPTION

BOX QTY.

EXACT MEDICAL SURGICAL GOWN

LEVEL 3 STANDARD

40gsm SMMS Surgical Gown wrapped in a 91cm x 91cm SMS field with one hand towel.




Standard

DEF5940	Medium 104cm	44
DEF5941	Large 110.5cm	44
DEF5942	X-Large 122cm	44
DEF5943	XX-Large 130.5cm	36

LEVEL 3 POLY REINFORCED

35gsm SMMS Surgical Gown with poly reinforced panels on the front and sleeves, wrapped in a 91cm x 91cm SMS field with one hand towel.



Poly Reinforced

DEF5944	Medium 104cm	40
DEF5945	Large 110.5cm	36
DEF5946	X-Large 122cm	32
DEF5947	XX-Large 130.5cm	28



WOUND CARE

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAEMOSTATIC BANDAGES

First aid and haemostatic products developed in Japan to simplify haemostasis treatment. These adhesive bandages protect and secure punctured skin after an IV injection or catheter removal. Used for various situations such as covering a punctured site following an injection, to stop the flow of blood from an intravenous line using pressure immobilisation, and to treat bleeding of a peripheral vein.

IV PRESSURE PAD

These IV Pressure pads are suitable to use following the removal of an IV, after giving an intravenous injection or after drawing blood or immunisation injections. Comprises of a 36mm round adhesive with a 16mm haemostatic pad.

The 3mm thick round pad is made of highly absorbent non-woven cellulose and creates a compression effect, resulting in bleeding stopping quickly. After bleeding has stopped, remove the pressure pad slowly and carefully.

DEF929	IV Pressure Pad 36mm White	100
--------	----------------------------	-----



DEF929

STEPTY HAEMOSTATIC BANDAGE

These pressure pads comprise of a rectangular 4 way stretch adhesive with a thick compression haemostatic pad. Designed as a protective covering of peripheral veins after blood sampling, a blood transfusion or an infusion.

Two products available with different pad thicknesses.

DEF1251	Stepty Pressure Pad 39mm x 80mm, 6mm Thick, White	50
---------	--	----



DEF1251

DEF1252	Stepty P Pressure Pad 40mm x 120mm, 9mm Thick, White	30
---------	---	----



DEF1252

CHUSHAVAN – FIRST AID ADHESIVE BANDAGE

These adhesive bandages are ideal as a protective covering after an injection. They are 27mm x 27mm adhesive size with a pad size of 13mm x 13mm with 1mm thickness.

DEF1250	Chushavan Pressure Pad Medium, Packet of 8	Box of 60 packets
---------	--	-------------------



DEF1250

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BANDAGES

Excellent range of sterile bandages that can be used in the operating room.

HOSPITAL WEIGHT CREPE

Made from hospital quality, porous cotton with non-fraying edges. Ideally suited for light compression and holding a dressing in place.



DEF822

DEF3782	2.5cm x 1.5m 108gsm Single Packed	48
DEF819	2.5cm x 1.5m 108gsm Double Packed	48
DEF5871	5cm x 1.5m 75gsm Single Packed	80
DEF3783	5cm x 1.5m 108gsm Single Packed	48
DEF820	5cm x 1.5m 108gsm Double Packed	48
DEF5872	7.5cm x 1.5m 75gsm Single Packed	80
DEF3784	7.5cm x 1.5m 108gsm Single Packed	48
DEF821	7.5cm x 1.5m 108gsm Double Packed	48
DEF5873	10cm x 1.5m 75gsm Single Packed	80
DEF3780	10cm x 1.5m 108gsm Single Packed	48
DEF822	10cm x 1.5m 108gsm Double Packed	48
DEF5874	15cm x 1.5m 75gsm Single Packed	80
DEF3781	15cm x 1.5m 108gsm Single Packed	36
DEF823	15cm x 1.5m 108gsm Double Packed	36

MEDIUM WEIGHT CREPE

Medium weight elastic crepe bandage provides support for soft tissue injuries and will secure a dressing in place. Non fraying edges and a soft, comfortable cotton.



DEF268

DEF268	10cm x 1.6m Single Packed	48
DEF269	15cm x 1.6m Single Packed	48

HEAVY DUTY CREPE

Sterile, heavy weight supporting crepe bandage that provides a firm support to joint injuries. Has a centre guide for assisting application and is comfortable to wear.



DEF272

DEF270	5cm x 2.3m Double Packed	48
DEF271	7.5cm x 2.3m Double Packed	48
DEF272	10cm x 2.3m Double Packed	48
DEF273	15cm x 2.3m Double Packed	24

MONTEX SUPER SPONGE BANDAGE

Super absorbent bandage made of prewashed fluff-dried woven cotton gauze with a crinkle-weave pattern for loft and bulk to cushion and protect wound areas. Finished edges are designed to reduce loose ends and li



DEF5544

DEF5544	11.4cm x 3.7m 6ply	120
---------	--------------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ORTHOPAEDIC PADDING

CONFORMING BANDAGE

Soft non adhesive bandage ideal as a secondary dressing for high body movement areas such as joints. Provides light two-way stretch and conforms easily.



DEF260

DEF1140	2.5cm x 1.5m, Double Packed	48
DEF258	5cm x 1.5m Double Packed	48
DEF259	7.5cm x 1.5m Double Packed	48
DEF260	10cm x 1.5m Double Packed	48
DEF261	15cm x 1.5m Double Packed	48

WATER RESISTANT UNDERCAST PADDING

Comfortable and smooth orthopaedic padding featuring high cushioning and protection for a fracture. Easily conformable and is resistant to moisture absorbency.



DEF824SP

DEF829	5cm x 2.7m Double Packed	48
DEF830SP	7.5cm x 2.7m Single Packed	48
DEF824SP	10cm x 2.7m Single Packed	48
DEF825SP	15cm x 2.7m Single Packed	24

COHESIVE BANDAGE

Sterile, lightweight, highly conformable bandage that is non slipping. This bandage is intended for general support and compression of soft tissue injuries.

DEF5888	10cm x 2m Single Packed	48
DEF1222A	10cm x 2m Double Packed	18
DEF3152	15cm x 2m Double Packed	12



DEF1222A

COTTON UNDERCAST PADDING

100% cotton undercast padding to provide soft cushioning and excellent absorbency.



DEF5708

DEF5710	5cm x 2.75m Double Packed	200
DEF5707	7.5cm x 2.75m Double Packed	150
DEF5708	10cm x 2.75m Double Packed	150
DEF5709	15cm x 2.75m Double Packed	80

CODE

DESCRIPTION

BOX QTY.

CODE

DESCRIPTION

BOX QTY.

ABSORBENT UNDERCAST PADDING

Absorbent padding which is multilayered to keep the skin dry and comfortable. Provides excellent comfort and cushioning.



DEF275

DEF828SP	7.5cm x 2.75m Single Packed	48
DEF275	10cm x 2.75m Single Packed	48
DEF276	15cm x 2.75m Single Packed	24

SURGICAL TAPE

Used to secure a dressing onto a wound, surgical tapes have a low allergenic adhesive which is designed to hold firmly onto the skin or dressing materials but removes easily without damaging the skin.

ADHESIVE SURGICAL FIXATION TAPE – F

Sterile surgical tape that is used for easy fixation in the operating theatre. Coated with a hypoallergenic adhesive whilst the backing is a stretch non woven polyester. Suitable for fixation of wound dressings or fixation of tubes, catheters or instruments.



DEF1504

DEF1504	10cm x 12cm	100
DEF1505	10cm x 20cm	100
DEF1506	10cm x 45cm	100
DEF1614	15cm x 70cm	50
DEF1507	20cm x 20cm	100
DEF3943	30cm x 45cm	100

ADHESIVE SURGICAL FIXATION TAPE – H

Excellent fixation tape which offers stretch and conformability and a low allergenic adhesive. Non woven fabric dressing sheet designed for rapid and quick fixation. Pre-cut and sterilised to a variety of sizes for immediate use in the operating theatre.



DEF1044

DEF1630	2.5cm x 15cm	100
DEF1623	2.5cm x 30cm	50
DEF1647	5cm x 8cm	50
DEF1995	5cm x 15cm	66
DEF1461	5cm x 24cm	50
DEF1599	5cm x 25cm	50
DEF1624	5cm x 30cm	50
DEF1797	5cm x 50cm	60
DEF1044	10cm x 10cm	100
DEF943	10cm x 20cm	50
DEF944	10cm x 25cm	50
DEF1057	10cm x 36.5cm	100
DEF1798	10cm x 50cm	60
DEF2195	15cm x 7cm	140
DEF1600	15cm x 30cm	50
DEF946	15cm x 40cm	25
DEF1026	15cm x 60cm	100
DEF1025	15cm x 100cm	100
DEF2160	20cm x 8cm	125
DEF937	20cm x 25cm	25
DEF678	20cm x 30cm	100
DEF475	20cm x 75cm Double Packed	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
TRACHEOSTOMY SPONGE Soft, absorbent gauze ideal to use around any drainage area to prevent chaffing.			COTTON TAPE FOR RETENTION 100% cotton tape pre-cut and sterilised non-medical device. Ready for use where a sterile securing tape is required.		
DEF053	10cm x 10cm Plain	100			
			DEF167		
			DEF173	1.5cm x 60cm Double Packed	100
			DEF2106	1.5cm x 1.5m	100
			DEF918	1.2cm x 500m, Non Sterile	1
			DEF919	2.5cm x 250m, Non Sterile	1
			POVIDONE-IODINE DRESSING MEGHdin® is a non-adherent povidone-iodine dressing designed to provide antimicrobial protection while helping to maintain a moist wound environment. Its low-adherent construction helps minimise trauma during dressing changes, making it suitable for a range of acute and chronic wounds.		
			MEGH DIN DRESSING MEGHdin® combines broad-spectrum antimicrobial protection with a low-adherent wound contact layer for the management of acute and chronic wounds.		
DEF2165	10cm x 10cm Y Cut Non woven, Packet of 2	100			
			P05-MEGHDIN		
			P05-MEGHDIN	5cm x 5cm Non Adherent Povidone Iodine Dressing	25
CR3900-ADVAZORB-T	8cm x 8cm Standard Hydrophilic Foam Tracheostomy Dressing	20	P95-MEGHDIN	9.5cm x 9.5cm Non Adherent Povidone Iodine Dressing	10
					
			CR3900-ADVAZORB-T		

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MANUKA HONEY DRESSINGS

The Advancis Medical Manuka Honey Dressing Range includes a selection of wound care dressings containing medical-grade Manuka honey. Suitable for use across a variety of wound types, the range is designed to support wound management by helping to maintain an appropriate wound healing environment and assisting clinicians in meeting individual patient needs.

ACTIVON TUBE®

Activon Tube® is a medical-grade Manuka honey dressing presented in an easy-to-apply tube format. It is suitable for use on acute and chronic wounds and can be applied directly to the wound bed or used in conjunction with an appropriate secondary dressing.



CR4493-ACTIVON-TUBE

CR4493-ACTIVON-TUBE	20g Manuka Honey	1
CR4493-12-ACTIVON-TUBE	20g Manuka Honey	12

ACTIVON TULLE® DRESSING

Activon Tulle® is a non-adherent viscose mesh dressing impregnated with 100% medical-grade Manuka honey. Designed for direct contact with the wound bed, it helps maintain a moist wound healing environment and can be used in the management of a wide range of acute and chronic wounds. Its flexible structure allows it to conform closely to the wound surface while requiring a secondary dressing.



CR3761-ACTIVON-TULLE

CR3761-ACTIVON-TULLE	5cm x 5cm Honey Coated Tulle	5
CR3658-ACTIVON-TULLE	10cm x 10cm Honey Coated Tulle	5

ALGIVON® ALGINATE DRESSING

Algivon® is a soft alginate dressing containing 100% medical-grade Manuka honey. Designed for use on moderate to heavily exuding wounds, it combines the absorbent properties of alginate with the benefits of Manuka honey to help maintain a moist wound healing environment. Its conformable structure makes it suitable for cavity wounds and areas requiring gentle packing.



CR3831-ALGIVON

CR3831-ALGIVON	5cm x 5cm Alginate with Honey	5
CR3659-ALGIVON	10cm x 10cm Alginate with Honey	5

ALGIVON PLUS® ALGINATE DRESSING

Algivon® Plus is a non-adherent alginate dressing containing 100% medical-grade Manuka honey. Designed for use on wounds with light to moderate exudate, it helps maintain a moist wound healing environment while allowing the dressing to conform closely to the wound bed. The dressing can be cut to size and is suitable for a variety of acute and chronic wounds.



CR4230-ALGIVON-PLUS

CR4230-ALGIVON-PLUS	5cm x 5cm Reinforced Alginate with Honey	5
CR4225-ALGIVON-PLUS	10cm x 10cm Reinforced Alginate with Honey	5

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE DRESSINGS

Advancis Medical's silicone dressing range provide gentle and versatile solutions for wound protection and fixation. Designed to minimise trauma during dressing changes, these silicone-based products support patient comfort while addressing a variety of wound care and skin protection needs.

SILTAPE® SILICONE TAPE

A soft silicone perforated fixation tape designed to secure dressings while minimising skin trauma upon removal. It is particularly suitable for patients with fragile or sensitive skin and can be repositioned during application.

CR3938-SILTAPE	2cm x 3m, Non Sterile	1
CR3938-SILTAPE-12	2cm x 3m, Non Sterile	12



CR3938-SILTAPE

CR3939-SILTAPE	4cm x 1.5m, Non Sterile	1
CR3939-SILTAPE-12	4cm x 1.5m, Non Sterile	12

SILFLEX® SILICONE DRESSING

A soft silicone wound contact dressing designed to protect the wound bed while allowing exudate to pass through to a secondary absorbent dressing. Its gentle adherence helps minimise trauma and pain during dressing changes.



CR3922-SILFLEX

CR3922-SILFLEX	5cm x 7cm Silicone Wound Contact Layer	10
CR3923-SILFLEX	8cm x 10cm Silicone Wound Contact Layer	10
CR3924-SILFLEX	12cm x 15cm Silicone Wound Contact Layer	10
CR3925-SILFLEX	20cm x 30cm Silicone Wound Contact Layer	10
CR4006-SILFLEX	35cm x 60cm Silicone Wound Contact Layer Sterile	5

VELLAFILM® FILM DRESSING

Vellafilm® is a transparent, atraumatic, polyurethane film dressing with Silfix® soft silicone adherence, designed to protect fragile skin and superficial wounds. Its waterproof, breathable construction allows monitoring of the wound site without removing the dressing.



CR3919-VELLAFILM

CR3919-VELLAFILM	12cm x 12cm Adherent Silicone Film Dressing	10
CR3921-VELLAFILM	15cm x 20cm Adherent Silicone Film Dressing	10
CR3917-VELLAFILM	12cm x 35cm Adherent Silicone Film Dressing	10

FOAM DRESSINGS

Advancis Medical's foam dressing range offers absorbent solutions for the management of a wide variety of acute and chronic wounds. Designed to manage exudate and maintain a moist wound healing environment, the range provides options to meet different anatomical and clinical needs.

ADVAZORB SACRAL® FOAM DRESSING

Advazorb Sacral® is a hydrophilic foam dressing designed to contour to the sacral area, for wounds with low to medium levels of exudate. With a soft silicone wound contact layer for gentle adherence, Advazorb Sacral® will only adhere to dry skin surrounding the wound and can be easily lifted for adjustment without losing its adherent properties.



CR4401-ADVAZORB-SACRAL

CR4401-ADVAZORB-SACRAL	17cm x 19cm Hydrophilic Foam Sacral Dressing with Adhesive and Border	10
CR4402-ADVAZORB-SACRAL	22cm x 23cm Hydrophilic Foam Sacral Dressing with Adhesive and Border	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ADVAZORB® FOAM DRESSING

A polyurethane foam dressing with a breathable film backing designed for the management of moderate to heavily exuding wounds. It provides absorbency while helping to maintain a moist wound healing environment.



CR4165-ADVAZORB

CR4165-ADVAZORB	7.5cm x 7.5cm Hydrophilic Foam Dressing with Film Backing	10
CR4166-ADVAZORB	10cm x 10cm Hydrophilic Foam Dressing with Film Backing	10
CR4167-ADVAZORB	12.5cm x 12.5cm Hydrophilic Foam Dressing with Film Backing	10
CR4168-ADVAZORB	15cm x 15cm Hydrophilic Foam Dressing with Film Backing	10
CR4169-ADVAZORB	10cm x 20cm Hydrophilic Foam Dressing with Film Backing	10
CR4170-ADVAZORB	20cm x 20cm Hydrophilic Foam Dressing with Film Backing	10

ADVAZORB SILFIX® DRESSING

Advazorb® Silfix is a soft silicone foam dressing that combines the absorbency with a gentle, atraumatic silicone adhesive layer. Designed for wounds with low to moderate levels of exudate, it helps maintain a moist wound healing environment while minimising trauma during dressing changes, making it suitable for patients with fragile or sensitive skin.



CR4177-ADVAZORB-SILFIX

CR4177-ADVAZORB-SILFIX	7.5cm x 7.5cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4178-ADVAZORB-SILFIX	10cm x 10cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4179-ADVAZORB-SILFIX	12.5cm x 12.5cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4180-ADVAZORB-SILFIX	15cm x 15cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4181-ADVAZORB-SILFIX	10cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4182-ADVAZORB-SILFIX	20cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive	10

ADVAZORB HEEL® FOAM DRESSING

Advazorb Heel® is a shaped foam dressing designed specifically for the heel. Suitable for wounds with low to moderate exudate, it combines absorbent foam with an anatomical design to provide a secure fit, helping to manage exudate and maintain a moist wound healing environment.



CR4228-ADVAZORB-HEEL

CR4228-ADVAZORB-HEEL	17cm x 21cm Hydrophilic Foam Dressing with Film Backing	5
----------------------	---	---

CODE

DESCRIPTION

BOX QTY.

CODE

DESCRIPTION

BOX QTY.

ADVAZORB BORDER® DRESSING

Advazorb Border® is a self-adhesive foam dressing with a soft silicone wound contact layer, designed for wounds with low to moderate exudate. The absorbent foam and perforated silicone layer allow exudate to pass into the dressing while protecting the wound bed, helping to maintain a moist wound environment and minimise trauma during dressing changes.



CR4190-ADVAZORB-BORDER

CR4190-ADVAZORB-BORDER	7.5cm x 7.5cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4191-ADVAZORB-BORDER	10cm x 10cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4192-ADVAZORB-BORDER	12.5cm x 12.5cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4193-ADVAZORB-BORDER	15cm x 15cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4194-ADVAZORB-BORDER	10cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4195-ADVAZORB-BORDER	20cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4196-ADVAZORB-BORDER	10cm x 30cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10

SUPER ABSORBENT DRESSINGS

The Advancis Medical Eclipse® range offers highly absorbent dressings for the management of moderate to heavily exuding wounds. Designed to provide fluid handling and retention, the range helps protect peri-wound skin and maintain an optimal wound healing environment.

ECLYPSE® DRESSING

Eclipse® is a super absorbent secondary dressing designed for the management of moderate to heavily exuding wounds. Its high-capacity absorbent core rapidly wicks, absorbs and locks away exudate, helping to reduce the risk of leakage and maceration while supporting an optimal wound healing environment. A water-resistant backing helps prevent strike-through and prolong wear time.



CR3818-ECLYPSE

CR3818-ECLYPSE	10cm x 10cm Super Absorbent Dressing	20
CR3769-ECLYPSE	15cm x 15cm Super Absorbent Dressing	20
CR3743-ECLYPSE	20cm x 30cm Super Absorbent Dressing	20
CR3808-ECLYPSE	60cm x 40cm Super Absorbent Dressing	10

ECLYPSE BOOT® DRESSING

Eclipse Boot® is a super absorbent secondary dressing specifically designed for use on the foot and ankle. Its anatomically shaped design provides coverage while the high-capacity absorbent core rapidly wicks, absorbs and locks away exudate, helping to reduce the risk of leakage and maceration. A water-resistant backing helps prevent strike-through and prolong wear time.



CR4237-ECLYPSE-BOOT

CR4237-ECLYPSE-BOOT	55cm x 47cm Super Absorbent Leg Wrap	5
CR4014-ECLYPSE-BOOT	60cm x 70cm Super Absorbent Leg Wrap	5

INDEX



A

ABDOMINAL SPONGES.....	23
ANAESTHETIC FACE MASK	30
ANAESTHETIC PACKS	
Epidural pack.....	40
Procedure pack.....	38-39
Spinal pack/tray.....	39-40
ANAESTHETIC TRAY.....	38
ANAESTHESIA DRUG LABELS	81
ANAL SPECULUM.....	66
ANGIO INSTRUMENT SET.....	78
ANGIO PACK.....	79
ARTERIAL TRAY.....	44

B

BABY	
Blanket	73
Cap	73
Cup.....	73
BACK TABLE COVER	5
BAGS	
Clip seal bag.....	31
Drawstring bag.....	31
Freezer bag.....	31
Plastic bag	31
Sanitary disposal bag.....	33
BANDAGES	
Crepe	93
Cohesive.....	94
Conforming	94
Esmarch	26
Haemostatic.....	92
Heavy duty.....	93
Hospital weight.....	93
Medium weight.....	93
Super sponge, Montex	93
Undercast	94-95
BANDS	
Latex free bands	14

BITE BLOCK.....	16
BLANKET – BABY RECEIVING	73
BOWL	
Bowl and spatula.....	52
Drape.....	51
Guidewire.....	78
Kidney dish.....	54
Plant-based	51,56
Plastic	51
Set - Bowl / Kidney Dish / Gallipot	55
Splash	51
Splash Stop	78
BOWL DRAPE	51
BOWL & KIDNEY DISH SET.....	55-56, 74
BRONNY IV PACKS	42
BULB IRRIGATION SYRINGE	20
BURNS	
Dressing.....	30
Gauze	30

C

C-ARM DRAPE.....	5
CAMERA CORD COVER.....	5
CAPS BERET	84
CARDIAC SUCTION INSTRUMENT	68
CATH LAB PACKS	78-79
CATHETER PACKS	41-42
CATS PAW RETRACTOR.....	20
CHAIR COVER.....	5
CHEMOTHERAPY GOWN.....	83
CHUSHAVEN PRESSURE PAD	92
CLAMPS	
Tube clamp blue plastic	14, 58
Tube clamp smooth jaw	14, 58
CLIP SEAL BAG	31

CLIPS		CUSTOM PROCEDURE PACKS	36
Drape clip	15	CVC INSERTION PACK	43
Towel clip green plastic	14		
Towel clip stainless steel	14	D	
COHESIVE BANDAGE	94	DELIVERY SUTURE SET	74
COMBINE		DELIVERY SET	74,75
Dressing	15	DEFOGGER	17
Roll	15	DENTAL	
Specialised	15	Bite block	16
CONFORMING BANDAGE	94	Instruments kit	58
CORD CLAMP CUTTING DEVICE	73	Pressure pack	16
COTTON TIPPED APPLICATORS/SWAB STICKS		Roll	16
Mouth and throat jumbo	34	Toothbrush	16
Paper shaft	34	DIATHERMY	
Wooden shaft	34	Quiver	17
COTTON		Tip cleaner	17
Balls	32	DRAINS	
Roll	32	Evacuator bulb	29
Tape	96	Fluted flat wound drain	29
Wool	32	Fluted round wound drain	29
COVERS		Fluted round wound drain set	29
Back table	5	Perforated flat wound drain	29
C-Arm	5	Perforated round wound drain	29
Camera cord	5	DRAPES	
Chair	5	Adhesive towel	12
Drill sleeve	5	Aqua drape sheets	9
Handle sleeve	6	Azure minor procedure	9
Image intensifier	77	Bowl	51
Lead screen	77	C-Arm	5
Light handle cover	20	Epidural	8
Mayo stand	6	Femoral radial angiography	77
Plastic sleeve	6	Guidewire book and wallet	78
Probe	77	Leggings	10
Protective	6	Minor procedure	9
Trolley	6	O.R. table /C-Arm protector	5
Ultrasound	77	Plastic	8
CREPE BANDAGE	93	Polybacked	10
CUPS		PurCotton	11
Baby cup with lid	73	SMS	10
Measuring cup	52	Stockinette	10
		Utility	10

DRAPE PACKS

Aqua universal drape pack heavy duty long sides	9
Aqua universal drape pack heavy duty medium sides	9
Aqua universal drape pack long sides.....	9
Aqua universal drape pack long sides without mayo cover.....	9

DRAPE CLIP	15
------------------	----

DRAWSTRING BAG	31
----------------------	----

DRESSING PACK LARGE	43
---------------------------	----

DRESSINGS

Burns.....	30
Combine	15
Combine roll.....	15
Foam	98-99
Manuka honey.....	97
Povidone-Iodine	96
Silicone	98
Super absorbent	100

DRESSING TRAY	53
---------------------	----

DRILL SLEEVE COVER	5
--------------------------	---

E

EAR SPECULUM	65
--------------------	----

ELECTROSURGICAL

Diathermy tip cleaner.....	17
Quiver	17

ENVIRO-SUSTAINABLE

Bowl and bowl set	51,56
Kidney dish and kidney dish set.....	56

EPIDURAL DRAPES	8
-----------------------	---

EPIDURAL/SPINAL/ANAESTHETIC PACKS	38-40
---	-------

EPISIOTOMY SUTURE SET	74
-----------------------------	----

ESMARCH BANDAGE

Green silicone.....	26
Blue neoprene.....	26

EYE

Eye spear	17
Eye stroll.....	17
Face shield.....	84
SPEye protective eyewear	83

EVACUATOR BULB	29
----------------------	----

F

FACE MASK ANAESTHETIC	30
-----------------------------	----

FACE SHIELD FULL.....	84
-----------------------	----

FEMORAL RADIAL ANGIOGRAPHY DRAPE	77
--	----

FOAM DRESSING	98-99
---------------------	-------

FOAM PAD	18
----------------	----

FORCEPS

Adson.....	58
Alligator	61
Bozman.....	61
Crile.....	60
Croc action.....	61
Dressing	60
Gillies	60
Hartmann croc action	61
Iris forcep.....	61
Kelman-Mcpherson.....	61
Lister sinus	61
Magill.....	60
McIndoe.....	60
Micro Adson	58
Mosquito.....	59
Plastic green	61
Roberts.....	62
Rochester Pean.....	58
Spencer Wells.....	59
Splinter Feichenfeldt.....	62
Tenaculum.....	62
Tilley ear	62
Wilde	61

FRAZIER SUCTION INSTRUMENT	67
----------------------------------	----

FREEZER BAG	31
-------------------	----

FUJITA SUCTION INSTRUMENT	67
---------------------------------	----

G

GALLIPOT.....	52
---------------	----

GAUZE SWABS/SPONGES		INSTRUMENTS.....	57
Gauze burns.....	30	INSTRUMENT PACKS.....	69-71, 78
Gauze folded.....	24	INSTRUMENT DRESSING PACK.....	71
Gauze green.....	24	INTERCOSTAL PACK.....	43
Gauze head dressing.....	24	IRRIGATION BULB SYRINGE.....	20
Gauze non woven.....	32	ISOLATION GOWN.....	83
Gauze plain.....	33	IUD INSERTION KIT.....	74
Gauze x-ray.....	24	IV INSERTION PACK.....	42
PurCotton.....	33	IV PRESSURE PAD.....	92
Sponges abdominal.....	23		
Sponges peanut.....	25	J	
Yard Strips.....	20	JUG.....	52
GOWNS		K	
Chemotherapy.....	83	KIDNEY DISH	
Classic surgical.....	89	Plant-based.....	56
Comfort surgical.....	87-89	Plastic.....	54
Coolmax surgical.....	89	KIDNEY DISH SETS.....	55
Exact Medical surgical.....	90		
Isolation.....	83	L	
Surgical.....	87-90	LABELS	
GUIDEWIRE BOWL.....	78	Anaesthesia.....	81
GUIDEWIRE BOOK.....	78	General medical.....	81
GUIDEWIRE WALLET.....	78	Piggybacked.....	81
H		LABOUR PACK.....	74
HANDLE SLEEVE COVER.....	6	LEAD SCREEN COVER.....	77
HAND DRYING		LEGGINGS.....	10
Dressing D towel.....	12	LIGHT HANDLE COVER.....	20
Dressing detowel.....	12	LUMBER PUNCTURE PACK.....	43
Huck towel blue.....	11		
Huck towel green.....	11	M	
Huck towel white.....	11	MANUKA HONEY DRESSING.....	97
PurCotton towel.....	11	MARKING PEN.....	23
HOLLOWWARE.....	50	MARKING PENCIL.....	22
HOOKS		MAYO STAND COVER.....	6
Skin.....	20	MATERNITY PAD.....	73
I			
IDENTILOOPS.....	19		
IMAGE INTENSIFIER COVER.....	77		

MATERNITY SUTURE SET	74	PREP APPLICATORS	
MEASURING CUP	52	Gauze prep sponges X-ray	24
MEDICAL DEVICE SYMBOLS	49	Prep foam	18
MEDICAL IMAGING PACKS	78-79	Prep foam on a handle	18
MEDICAL LABELS	80	PREP SET	18
MOUTH & THROAT JUMBO SWAB STICKS	34	PROBE COVER	77
MOUTH TOILET PACK	42	PROTECTIVE COVER	6
N		PRESSURE PAD	
NASAL BOLSTER	32	Chushavan	92
NASAL SPECULUM		I.V.	92
Mod Wein	65	Stepty	92
Thudichum	65	Stepty P	92
NEEDLE COUNTER	21-22	PURCOTTON	
NEEDLE HOLDER	62	Drape	11
NEUROSURGICAL PATTY	25	Towel	11
NON-WOVEN		Swabs	33
Balls	32	PVC TUBING WITH XRD	26
Combine	15	Q	
Gauze	32	QUIVER	17
O		R	
ONETRAC RETRACTOR	27	RETRACTOR	
OVERSHOES	84	Cats paw	20
P		Onetrac	27
PACEMAKER PACK	79	Weitlaner with kidney dish	62
PADDING UNDERCAST	94	S	
PEANUT SPONGE	25	SAFETY PINS	22
PENCIL	22	SANITARY DISPOSABLE BAG	33
PERSONAL PROTECTIVE APPAREL	82	SCALPEL	
PLACENTA BASIN	73	Blade Remover	26
PLASTIC BAG	31	Handle	63
POVIDONE-IODINE DRESSING	96	Holder	54
		Scalpel blade & handle	63

SCISSORS		SPONGES	
Blunt/blunt.....	63	Abdominal.....	23
Dressing.....	63	Nasal bolster.....	32
Episiotomy.....	64	Peanut.....	25
Iris.....	64	Tracheostomy.....	96
Mayo.....	64	STAPLE REMOVER.....	66
Metzenbaum.....	64	STEPTY & STEPTY P PRESSURE PAD.....	92
Sharp/blunt.....	63	STOCKINETTE.....	10
Sharp/sharp.....	63	SUCTION INSTRUMENT	
Tenotomy.....	64	Adson Frazier.....	67
Umbilical cord cutting.....	73	Cardiac.....	68
Universal.....	64	Frazier.....	67
Vannas.....	64	Fujita.....	67
SCOPE WARMER.....	25	ONETRAC.....	67
SHARPS TRANSFER DEVICE.....	54	Paediatric.....	67
SHARPS TRAY.....	53	Poole.....	66
SILICONE DRESSING.....	98	Rhoton.....	68
SILICONE TUBING		Yankauer.....	66
Pre-cut sterile.....	28	SURGICAL GOWN	
Roll non sterile.....	35	Classic surgical.....	89
With connector.....	28	Comfort surgical.....	87-89
SILTAPE.....	98	Coolmax surgical.....	89
SKIN HOOK.....	20	Exact Medical surgical.....	90
SKIN MARKING PEN.....	23	SURGICAL PACKS - STANDARD	
SPEYES EYEWEAR PROTECTION.....	83	Cataract.....	47
SPECIMEN JAR.....	53	Intraocular.....	46
SPECIMEN TUBE.....	53	Minor Lithotomy.....	47
SPECULUMS		Universal General.....	44-46
Anal.....	66	SURGICAL RETRACTOR	
Ear.....	65	Onetrac.....	27
Nasal.....	65	SUTURE BOOTS.....	28
Vaginal.....	65	SUTURE SETS.....	69-71, 74
SPLASH BOWL.....	51	SYMBOLS	
SPLASH STOP.....	78	Medical device.....	49
SPONGE HOLDER.....	66		

SWABS/GAUZE

Gauze head dressing.....	24
Gauze non-woven.....	32
Gauze peanut.....	25
Gauze prep.....	24
Swabs green.....	24
Swabs folded.....	24
Swabs plain.....	33
Swabs PurCotton.....	33
Swabs x-ray.....	24

SWAB STICKS/COTTON TIPPED APPLICATOR

Mouth and throat jumbo.....	34
Paper shaft.....	34
Wooden shaft.....	34

T

TAPE

Cotton.....	96
Siltape silicone.....	98
Surgical tape F.....	95
Surgical tape H.....	95

THYROID PACK.....43

TIP CLEANER.....17

TOOTHBRUSH.....16

TOUCH FASTENER.....20

TOWEL

Clip.....	14
Dressing D towel.....	12
Dressing towel detowel.....	12
Huck towel blue.....	11
Huck towel green.....	11
Huck towel white.....	11
PurCotton towel.....	11

TRACHEAL DILATOR.....68

TRACHY PROCEDURE PACK.....43

TRACHEOSTOMY SPONGE.....96

TRAY

Dressing tray.....	53
Sharps.....	53
Square fully deep.....	53

TROLLEY COVER.....6

TUBE CLAMP.....14, 58

TUBE HOLDER.....15

TUBING.....26, 28, 35

U

UNDERCAST PADDING.....94-95

ULTRASOUND PROBE COVER.....77

UMBILICAL CORD CUTTING SCISSORS.....73

UNIVERSAL DRAPE PACKS.....9

UNIVERSAL GENERAL PACKS.....44-46

UTILITY MARKER.....23

V

VAGINAL PLUG.....33

VAGINAL SPECULUM.....65

W

WOUND CARE.....91

WOUND DRAIN

Evacuator bulb.....	29
Fluted flat.....	29
Fluted round.....	29
Fluted round set.....	29
Perforated flat.....	29
Perforated round.....	29

WIRE CUTTER.....68

X

X-RAY INSTRUMENT PACK.....78

Y

YANKAUER SUCTION INSTRUMENT.....66



Head Office

1-7 Vision Street
Dandenong South VIC 3175
t: 61 3 9706 3600
f: 61 3 9706 3619

freecall: 1300 550 278
e: sales@defries.com.au
www.defries.com.au

Manufacturing Site

18-24 Villas Road
Dandenong South VIC 3175