

CLINICAL DRAPES AND EQUIPMENT COVERS



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

EQUIPMENT COVERS

BACK TABLE COVER

DEF5016AQ	152cm x 228cm (60" x 90") Aqua	40
DEF3544	150cm x 190cm Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3545	150cm x 230cm Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3678	150cm x 300cm Sapphire, Heavy Duty, Impervious and Absorbent	12
DEF5802	152cm x 228cm Aqua, Zoned, Heavy Duty	20

CHAIR COVER

		
DEF772		
DEF772	2 x 57.5cm x 62cm with 2 x Latex-free Bands #33	100

CAMERA CORD COVER

		
DEF1201		
DEF1201	Camera Cord Cover 15cm x 230cm with Cuff Sealed End & Tapes, Double Packed	50

DRILL SLEEVE COVER

		
DEF889		
DEF1224	5cm x 80cm with 3 Labels	100
DEF1368	5cm x 110cm with Tape	100
DEF1516	7.5cm x 100cm with Adhesive both Ends	100
DEF1154	7.5cm x 150cm with 2 x Latex-free Bands #31	100
DEF889	10cm x 150cm	100
DEF618	10cm x 152cm with Tape Double Packed	25

C-ARM DRAPE

Large, one-piece clear poly cover for complete coverage of C-Arm with white poly adhesive straps to secure the drape in place. Packaged with a telescoping fold. Adhesive straps included to secure excess drape material.

30-0201	104cm x 188cm (41" x 74")	20
---------	---------------------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HANDLE SLEEVE COVER



DEF1053

DEF1053	14cm x 30cm with Adhesive	100
DEF1466	14cm x 30cm with Adhesive Packet of 2	100

PROTECTIVE COVER

DEF613	15cm x 100cm, Plastic	50
DEF647	20cm x 20cm, Plastic	100
DEF650	30.5cm x 39cm, Plastic	100
DEF763	104cm x 150cm, Clear Plastic	40
DEF1573	150cm x 208cm, Clear Plastic	50

TROLLEY COVER



DEF641

DEF641	Trolley Bag Cover 61cm x 102cm	50
DEF1239	Shower Trolley Drape 2m x 3m	25

MAYO STAND COVER

DEF5019AQ	61cm x 135cm, Aqua	30
DEF761	62cm x 91cm Clear Plastic with Sealed End	100
DEF775	70cm x 137cm, Blue Plastic with Cuff	50
DEF6100	80cm x 145cm, Sapphire, Heavy Duty, Blue Plastic with Reinforcement	25

PLASTIC SLEEVE COVER



DEF454

DEF454	15cm x 100cm, with 2 x Latex-free Bands #31, Double Packed	100
DEF1744	15cm x 200cm, with 1 x Latex-free Band #31 and Sealed End, Double Packed	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------



CUSTOMISED PROCEDURE DRAPES

THE FLEXIBILITY OF CHOICE

At Defries we have a comprehensive range of patient drapes which cover all procedures in the operating room. We provide our customers with superior products that match their clinical needs.

Our drapes are made to specific requirements such as size, fabric, fenestrations and extra details like velcro taping and clear plastic anaesthetic screens.

CHOICE OF FABRIC IS IMPORTANT

At Defries we use a variety of drape fabrics:

Azure drapes use a viscose impervious fabric with a polyethylene laminate, suitable for general surgery with minimal fluid or short procedures. This range includes the Azure minor procedure drapes.

Sapphire drapes use a bi-laminate or tri-laminate fabric with extra reinforcement providing excellent barrier protection. This range includes back table covers, mayo covers, some procedure drapes and bespoke designed drapes to meet specific customer requirements.

Aqua drapes are our newest drape range available and use an innovative drape fabric that is a tri-laminate – a three layer drape which offers an absorbent top level, an impervious middle coating and a soft and comfortable patient layer. Soft and light, this material feels the closest to linen when draping the patient.

This high performance fabric is soft and drapable with additional features such as reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. We have an excellent standard range of designs covering all procedures from general surgery and universal packs, to orthopaedics and specialised surgery such as neuro.

Supreme drapes are made from a tri-laminate fabric and are specifically designed for high level performance whilst being soft, drapable and easy to use. These drapes provide great absorbency from fabric reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. The Supreme range includes drape sheets and specialty procedure drapes.

Our drape design includes the instructions for correct aseptic application with markings and directions on the drape to assist ease of use.

All drape fabrics are conformable to patients' shapes. They are low linting and have excellent barrier properties to prevent strikethrough of airborne and blood borne bacteria, meeting regulatory standards.

Please contact your Defries representative for further information.

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

DRAPES FOR EVERY PROCEDURE

Our wide range of drapes include:

- **ANCILLARY**

Absorbent Pad, Adhesive Drape Towel, Drape Tape

- **ANGIOGRAPHY & CATH LAB**

Angio, Femoral Radial Angio, Biliary, Hickman, Pacemaker, Radial

- **CARDIOVASCULAR**

Cardiac & Cardiac Splits, CV Groin, Cardiac Saddle

- **ENT**

Minor Procedure, Oral, Turban, Split, Paediatric ENT Split

- **EQUIPMENT COVERS**

Back Table Cover, Mayo Stand Covers, Utility

- **GENERAL DRAPES**

Drape Sheets, Fenestrated Minor Procedure, Laparoscopy, Laparotomy, PICC/Arterial, Central Line, Universal Drape Packs

- **GYNAECOLOGY & OBSTETRICS**

Caesarean, Endoscopy, Laparoscopy, Lithotomy, Leggings, Under Buttocks

- **NEPHROLOGY**

PERC

- **NEURO SPINAL**

Neuro-craniotomy, Isolation, Transsphenoidal

- **PAEDIATRICS**

Neonatal Split, Neonatal, Paediatric Drape Sheets

- **OPHTHALMIC**

Eye

- **ORTHOPAEDICS**

Extremity Foot, Hand, Hip, Knee, Shoulder Arthroscopy, Shoulder Extremity, Shoulder Beach Chair, Split & U, Stockinette, Vertical Isolation

- **UROLOGY**

TUR

GENERAL DRAPES AND PACKS

CLEAR PLASTIC DRAPE



DEF742

DEF742	45cm x 61cm	200
DEF079	61cm x 91cm	200
DEF749	70cm x 120cm	100
DEF745	76cm x 106cm	100
DEF3910	76cm x 152cm, Double Packed	100
DEF151	91cm x 115cm	100
DEF2430A	Procedure Adhesive Drape 91cm x 115cm with 2.5cm x 115cm adhesive down length of drape	50
DEF748	120cm x 140cm	100
DEF752	120cm x 140cm, Double Packed	100
DEF5713	200cm x 200cm	100

EPIDURAL DRAPE



DEF1481

DEF1481	76cm x 106cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100
DEF748FA-S	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	10
DEF3394	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

UNIVERSAL DRAPE PACKS

AQUA UNIVERSAL DRAPE PACK – LONG SIDES

DEF5400	1 x Top Drape, 150cm x 276cm with Reinforcement and Adhesive 2 x Side Drapes, 99cm x 145cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 193cm x 193cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 80cm x 145cm with Reinforcement Wrapped in a Back Table Cover, 152cm x 248cm with Reinforcement	5
---------	---	---

AQUA UNIVERSAL DRAPE PACK – HEAVY DUTY - MEDIUM SIDES

DEF5800A	1 x Top Drape, 150cm x 276cm with Reinforcement and Adhesive 2 x Side Drapes, 99cm x 145cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 193cm x 193cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 80cm x 145cm with Reinforcement Wrapped in a Back Table Cover, 152cm x 248cm with Reinforcement	6
----------	---	---

AQUA UNIVERSAL DRAPE PACK – HEAVY DUTY - LONG SIDES

DEF5801-4	1 x Top Drape, 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drapes, 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement and Adhesive 1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement	4
-----------	--	---

AQUA UNIVERSAL DRAPE PACK – LONG SIDES WITHOUT MAYO STAND COVER

DEF5243AQ	1 x Top Drape, 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drapes, 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement and Adhesive Wrapped in a Back Table Cover, 127cm x 230cm with Reinforcement	5
-----------	---	---

AZURE MINOR PROCEDURE DRAPE

Soft viscose fabric combined with a polyethylene laminate to create a drape ideal for minor procedures.

DEF5094	38cm x 45cm Azure Circumcision Drape with a 3cm Circular Fenestration	50
DEF5090	50cm x 60cm Azure Minor Procedure Drape with a 6.5cm Circular Fenestration and Full Adhesive	25
DEF5201	90cm x 90cm Fenestrated Drape with a 15cm x 20cm Rectangular Fenestration	25
DEF5091	91cm x 91cm Azure Minor Procedure Drape with a 10cm x 15cm Oval Fenestration and Full Adhesive	40
DEF2923	112cm x 153cm ICU Procedure Drape with a 5.5cm x 8cm Rectangular Fenestration and Full Adhesive	10
DEF5092	200cm x 300cm Azure Minor Procedure Drape with a 15cm x 25cm Fenestration and Full Adhesive	12


AQUA DRAPE SHEETS

Innovative tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.

DEF5141AQ	Square Drape, 100cm x 100cm	40
DEF5025AQ	Medium Drape, 100cm x 180cm	20
DEF5024AQ	Half Drape, 110cm x 150cm	20
DEF5026AQ	3/4 Drape, 150cm x 200cm	20
DEF5028AQ	Large Drape, 180cm x 215cm	20
DEF5030AQ	Top/Head Drape 135cm x 250cm with Reinforcement and Adhesive	15
DEF5828	Top/Head Drape 150cm x 276cm with Reinforcement Adhesive and Tube Holders	15
DEF5031AQ	Side Drape 100cm x 130cm with Reinforcement and Adhesive	40
DEF5023AQ	Side Drape 100cm x 195cm with Reinforcement and Adhesive	20
DEF5032AQ	Bottom/Foot Drape 190cm x 185cm with Reinforcement and Adhesive	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

UTILITY DRAPE



DEF5130

DEF5130	65cm x 65cm with Adhesive wrapped in a 91cm x 91cm SMS	50
---------	--	----

FABRIC SMS DRAPE

SMS fabric is a 3 layer polypropylene material bonded together to form a drape that is fluid repellent. Ideal for minor procedures.



DEF1594

DEF1122	30.5cm x 61cm	100
DEF1594	61cm x 61cm	50
DEF3806	61cm x 61cm, Non Sterile	100
DEF2300A	91cm x 91cm	100
DEF3807	91cm x 91cm, Non Sterile	100
DEF3805	150cm x 200cm, Non Sterile	50

POLYBACKED TOWEL / DRAPE

An absorbent drape with a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor procedures.



DEF833

DEF833	Polybacked Towel, 33cm x 46cm White	750
DEF1287	Polybacked Towel, 40cm x 60cm Blue	100
DEF1288A	Polybacked Towel, 40cm x 60cm with a 10cm Circular Fenestration Blue	100
DEF2072	Polybacked Procedure Drape 75cm x 105cm with a 10cm x 15cm Fenestration and a 5cm Adhesive length along the 10cm sides of Fenestration, Blue	25

STOCKINETTE



DEF5060AQ

9503-36	7.6cm x 91cm (3" x 36"), Single Ply	20
DEF5060AQ	20cm x 96cm, Medium, Aqua	25
DEF5061AQ	23cm x 122cm, Large, Aqua	25
DEF5062AQ	30cm x 122cm, XLarge, Aqua	25

LEGGINGS

DEF2481AQ	Aqua Lithotomy Leggings 80cm x 140cm, Packet of 2	26
DEF5681	Clear Lithotomy Leggings 82cm x 125cm, Packet of 2	26

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SUPPLEMENTARY

HAND DRYING AND STERILE FIELDS

BLUE HUCK TOWEL

High absorbency cotton towels.



DEF1300

DEF1976	40cm x 68cm, Blue	30
DEF1300	40cm x 68cm, Blue	100
DEF3274	40cm x 68cm, Blue, Packet of 2	50
DEF2018	40cm x 68cm, Blue, Packet of 5, Double Packed	40

GREEN HUCK TOWEL



DEF2619

DEF2619	40cm x 68cm, Green, Packet of 2	100
DEF2140	40cm x 68cm, Green, Packet of 5	40
DEF5556	60cm x 90cm, Green	100

WHITE HUCK TOWEL

DEF2753	40cm x 68cm, White	30
---------	--------------------	----



DEF2753

PURCOTTON TOWEL

This PurCotton towel is a superior alternative to traditional Operating Room hand towels. They are a non-woven fabric made with a natural cotton, have a soft and comfortable feel, high liquid absorbency and are 100% biodegradable. They are also low linting and have an excellent wet-strength.



DEF3536A

DEF3536A	42cm x 65cm, Blue	100
DEF2521A	61cm x 61cm, White, Packet of 4	100
DEF3674	90cm x 100cm, White	100



DEF2917A

DEF2905A	40cm x 60cm White Non Sterile	100
DEF2917A	42cm x 35cm, White Non Sterile	200
DEF5824	90cm x 100cm, White Non Sterile	100

PURCOTTON DRAPE



DEF2404A

DEF2404A	76cm x 106cm wrapped in a 91cm x 91cm SMS	30
DEF2412	76cm x 106cm, Packet of 2 wrapped in a 91cm x 91cm SMS	50
DEF2663D	120cm x 140cm with a 7cm circular fenestration	20

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.																																							
<div><div>DESSING TOWEL (DETOWEL)</div><div>Made from airlaid paper with cotton yarn.</div><div></div><div>DEF5612</div><table><tr><td>DEF5612</td><td>30cm x 30cm</td><td>100</td></tr><tr><td>DEF5614</td><td>30cm x 30cm</td><td>1000</td></tr><tr><td>DEF5611</td><td>30cm x 30cm, Packet of 2</td><td>50</td></tr><tr><td>DEF5613</td><td>30cm x 60cm</td><td>100</td></tr><tr><td>DEF5610</td><td>30cm x 60cm, Packet of 2</td><td>50</td></tr></table><div>DESSING D TOWEL</div><div>Made from rayon and polyester.</div><div></div><div>DEF656</div><table><tr><td>DEF656</td><td>32cm x 34cm, Packet of 2</td><td>50</td></tr><tr><td>DEF831</td><td>32cm x 34cm, Packet of 2</td><td>300</td></tr><tr><td>DEF832</td><td>32cm x 34cm, Packet of 10</td><td>100</td></tr><tr><td>DEF894</td><td>34cm x 60cm</td><td>50</td></tr><tr><td>DEF471A</td><td>34cm x 60cm</td><td>300</td></tr><tr><td>DEF910</td><td>34cm x 60cm, Packet of 5</td><td>25</td></tr></table></div>						DEF5612	30cm x 30cm	100	DEF5614	30cm x 30cm	1000	DEF5611	30cm x 30cm, Packet of 2	50	DEF5613	30cm x 60cm	100	DEF5610	30cm x 60cm, Packet of 2	50	DEF656	32cm x 34cm, Packet of 2	50	DEF831	32cm x 34cm, Packet of 2	300	DEF832	32cm x 34cm, Packet of 10	100	DEF894	34cm x 60cm	50	DEF471A	34cm x 60cm	300	DEF910	34cm x 60cm, Packet of 5	25	<div><div>ADHESIVE DRAPE TOWEL</div><div></div><div>DEF5701</div><table><tr><td>DEF5701</td><td>48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS</td><td>15</td></tr></table></div>			DEF5701	48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS	15
DEF5612	30cm x 30cm	100																																										
DEF5614	30cm x 30cm	1000																																										
DEF5611	30cm x 30cm, Packet of 2	50																																										
DEF5613	30cm x 60cm	100																																										
DEF5610	30cm x 60cm, Packet of 2	50																																										
DEF656	32cm x 34cm, Packet of 2	50																																										
DEF831	32cm x 34cm, Packet of 2	300																																										
DEF832	32cm x 34cm, Packet of 10	100																																										
DEF894	34cm x 60cm	50																																										
DEF471A	34cm x 60cm	300																																										
DEF910	34cm x 60cm, Packet of 5	25																																										
DEF5701	48cm x 75cm, Pack of 4 wrapped in a 91cm x 91cm SMS	15																																										