



CUSTOM PROCEDURE AND STANDARD PACKS

CUSTOM PROCEDURE PACKS

COST EFFECTIVE AND TIME EFFICIENT: A PACK FOR EVERY PROCEDURE

Our Custom Procedure Pack range covers all areas of the hospital from Cardiac, Orthopaedic, Vascular, ENT and more. They contain all the individual consumables for a surgical procedure, wrapped in a 'field' and sterilised. All products are ready to use when opened, unlike individually wrapped consumables that are touched many times prior to use and require opening separately. Plus there is reduced wastage. Using a Procedure pack saves on setup time for procedures and ensures that all components are readily available for every case. Ordering and inventory control is easier.

CUSTOM TAILOR A PROCEDURE PACK TO MEET YOUR SPECIFIC NEEDS

Select individual consumables such as patient drapes and surgical gowns, to diathermy, suction tubing, abdominal packs, gauze and instruments. Your pack will only include items specifically chosen by you, the surgeon and theatre nurses. Defries assembles these customised packs in cleanrooms according to hospital specifications.

Custom procedure packs can be varied to include cost-effective components as they become available, and allow easy accountability of costs per procedure. As a local Melbourne manufacturer, we are able to respond quickly to any changes requested by a hospital to the configuration of your custom packs.

We'll provide you with ongoing education, listen to feedback and are always available to discuss your specific needs.



READY-FOR-USE STANDARD PROCEDURE PACKS

We have a range of standard procedure packs such as the anaesthetic packs, epidural packs, regional trays, catheter packs, Bronny IV insertion pack, lumbar puncture packs, ICU packs including central line packs, arterial line packs, all designed, made and ready to be ordered.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PACKS

All our anaesthetic packs come sterile, prepacked and ready to use. Indications for use include:

- General and Local Anaesthetic procedures
- Spinal/Epidural/Regional/Lumbar Puncture procedures
- IV Access

NRFIT components are available on request for inclusion into packs.

Select from the packs below or contact us to customise your own anaesthetic pack.

ANAESTHETIC PACKS

DEF1409	1 x Dressing Tray 24cm x 15cm x 4cm 2 x Swabs 10cm x 10cm Green	36
---------	--	----



DEF1409

DEF1162	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 3 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	60
---------	---	----



DEF1162

ANAESTHETIC PACKS

DEF1618	1 x Drape Blue Plastic 75cm x 75cm with 10cm x 15cm Fenestration with Adhesive Tabs on Corners of Fenestration and Corners of Drape 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 5 x Gauze Swabs 10cm x 10cm Green 1 x Forcep Plastic Green 1 x Syringe 5ml Luer Slip 1 x Syringe 3ml Luer Slip 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt Drawing up 1 x Needle 25g x 1" 1 x Dressing Towel Quick Dry 33cm x 33cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
---------	---	----



DEF1618

ANAESTHETIC TRAY

DEF1524	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 2 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	80
DEF406	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Polybacked Towel 40cm x 42cm	60

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PROCEDURE PACK

DEF2255B	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 91cm x 115cm 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 3ml Luer Lock 1 x Syringe 5ml Luer Lock 1 x Needle 25g x 1.5" Hypodermic 2 x Filter Needle 5um 2 x Forceps Pearman Plastic 2 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 91cm x 115cm	60
----------	---	----

ANAESTHETIC/SPINAL PACK

DEF2232A	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 76cm x 106cm 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 5ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt Hypodermic 1 x Needle 19g x 1.5" Hypodermic 1 x Needle 25g x 1" Hypodermic 1 x Forceps Plastic Green 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm Green 5 x Prep Swabs Large XRD 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
----------	--	----

ANAESTHETIC BASIC TRAY AND SYRINGE KIT

DEF4117E	1 x Anaesthetic Tray with 4 Compartments 22.5cm x 14cm x 2.5cm 1 x Syringe 3ml Luer Slip 1 x Syringe 5ml Luer Slip 1 x Syringe 10ml Luer Slip Eccentric Tip 1 x Syringe 30ml Luer Slip Eccentric Tip 2 x Gauze Swabs 10cm x 10cm Green Wrapped in Crepe Paper 75cm x 75cm	50
----------	---	----



DEF2232A

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
EPIDURAL PACKS					
DEF1796	<p>1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs</p> <p>2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm</p> <p>2 x Forceps Plastic Green</p> <p>5 x Gauze Swab 10cm x 10cm Green</p> <p>1 x Dressing Towel (Hand Towel) Large 32cm x 60cm</p> <p>Wrapped in a Blue Plastic Field 91cm x 115cm</p>	50			
DEF1655B	<p>1 x SMS Drape 91cm x 91cm with 12cm x 15cm Fenestration and Adhesive</p> <p>1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm</p> <p>1 x Forceps Plastic Green</p> <p>5 x Gauze Swab 10cm x 10cm Green</p> <p>1 x Gallipot 60ml</p> <p>1 x Towel 40cm x 42cm 2ply</p> <p>Wrapped in a Blue Plastic Field 91cm x 115cm</p>	60			
DEF3085A	<p>1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and ACorner Tabs</p> <p>2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm</p> <p>1 x Syringe 5ml Luer Slip</p> <p>1 x Syringe 2ml Luer Slip</p> <p>2 x Needle 18g x 1.5" Blunt Hypodermic</p> <p>1 x Needle 25g x 1" Hypodermic</p> <p>1 x Needle 25g x 1.5" Hypodermic</p> <p>5 x Gauze Swabs 10cm x 10cm Green</p> <p>1 x Dressing Towel (Hand Towel) 31.5cm x 34cm</p> <p>Wrapped in a Blue Plastic Field 75cm x 75cm</p>	50			
DEF3854	<p>1 x Drape Laminated 120cm x 140cm with 7cm x 10cm Fenestration Blue</p> <p>1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm</p> <p>1 x Clamp 19cm Blue Plastic</p> <p>5 x Gauze Swabs 10cm x 10cm Green</p> <p>2 x Measuring Cup 150ml Plastic</p> <p>1 x Gown Surgical T-Cut Unferinforced XL</p> <p>1 x Huck Towel 43cm x 69cm Blue</p> <p>Wrapped in a Complex Tablecloth 115cm x 150cm</p>	15			
COMPLETE EPIDURAL PACK					
DEF1480	<p>1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs</p> <p>1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm</p> <p>1 x Forceps Plastic Green</p> <p>1 x Syringe 5ml Luer Slip</p> <p>1 x Syringe 3ml Luer Slip</p> <p>1 x Needle 19g x 1.5" Hypodermic</p> <p>1 x Needle 18g x 1.5" Blunt Hypodermic</p> <p>1 x Needle 25g x 1" Hypodermic</p> <p>5 x Gauze Swab 10cm x 10cm Green</p> <p>1 x Dressing Towel (Hand Towel) 31.5cm x 34cm</p> <p>Wrapped in a Blue Plastic Field 75cm x 75cm</p>	60			
EPIDURAL/SPINAL TRAY					
DEF2345C	<p>1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs</p> <p>1 x Dressing Tray with 3 Compartments 20cm x 25cm x 2.5cm</p> <p>1 x Syringe 10ml Luer Lock</p> <p>1 x Syringe 5ml Luer Lock</p> <p>1 x Syringe 3ml Luer Lock</p> <p>1 x Needle 18g x 1.5" Blunt Hypodermic</p> <p>1 x Needle 25g x 1.5" Hypodermic</p> <p>1 x Filter Needle 5um</p> <p>1 x Foam Applicator Blue 15cm</p> <p>5 x Gauze Swabs 10cm x 10cm Green</p> <p>2 x Dressing Towel (Hand Towel) 31.5cm x 34cm</p> <p>1 x Gallipot 60ml</p> <p>Wrapped in a Blue Plastic Field 91cm x 115cm</p>	24			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CATHETER PACKS

All Catheter Procedure Packs have been designed in consultation with the healthcare professionals who use them. They include all the consumables used to insert or change a catheter aseptically. Alternatively please contact us to customise a catheter pack to your own specific requirements.

CATHETER PACKS



DEF1601

DEF1601	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Blue Huck Towel 40cm x 68cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF2074	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Pearman 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 2 x Kidney Dishes 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24

DEF2061A	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Green 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF1262	1 x Drape SMS 61cm x 61cm with a 10cm x 10cm Fenestration 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Tube Clamp with Lock 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1333TA	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip Concentric Tip 2 x Forceps Plastic Green 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1650B	1 x Drape Azure Minor Procedure 50cm x 60cm with 6.5cm Circular Fenestration and Adhesive around Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 2 x Forceps Plastic Green 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
DEF2253E	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 2 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Cotton Balls Large 1 x Dressing Towel T 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60			
DEF1522	3 x Drape Polybacked 33cm x 46cm 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm Plain Wrapped in a Blue Plastic Field 75cm x 75cm	60			
DEF1036	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Forceps Plastic Green 6 x Cotton Balls Large 1 x Dressing Towel T 31.5cm x 61cm 2 x Kidney Dishes 700ml Clear 2 x Gallipots 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	30			
DEF1038A	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 1 x Tube Clamp with Lock 1 x Kidney Dish 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	30			
DEF1932A	1 x Drape PolyBacked 50cm x 60cm with 10cm Circular Fenestration 1 x Syringe 30ml Luer Slip Eccentric Tip 1 x Sponge Holder 24cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	30			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL PACKS

IV INSERTION PACKS

For the insertion of an IV line.

BRONNY IV PACK

DEF920	1 x Plastic Drape 45cm x 60cm 2 x Gauze Swabs 5cm x 5cm Plain	200
--------	--	-----

BRONNY IV PACK WITH 4 SWABS

DEF1676	1 x Plastic Drape 45cm x 60cm 4 x Gauze Swabs 5cm x 5cm Plain	200
---------	--	-----

THYROID PACK

DEF3799B	1 x Forcep Spencer Wells 13cm Straight 1 x Scissors 13cm Sharp/Blunt 1 x Staple Remover Scissor Action with Plastic Handle 1 x Stitch Cutter 10 x Gauze Swabs 10cm x 10cm Plain 1 x Huck Towel 40cm x 68cm Blue Wrapped in a SMS Blue Field 61cm x 61cm	20
----------	---	----

MOUTH TOILET PACK

DEF399	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 10 x M & T Swab Sticks 15cm Paper Shaft 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	50
--------	--	----

PLASTICS PACK

DEF219A	1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 12 x Gauze Swabs 10cm x 10cm Plain 2 x Dressing Towel (Hand Towel) 31.5cm x 34cm 2 x Swab Sticks 15cm with Wooden Shaft 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 50cm x 50cm	60
---------	---	----

FINE PROCEDURE PACK

DEF2453A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Scalpel Handle # 3 4 x Gauze Swabs 7.5cm x 7.5cm Plain 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
----------	--	----

LARGE DRESSING PACK

DEF2234A	1 x Dressing Tray 24.5cm x 28cm x 4.5cm 3 x Forceps Plastic Green 10 x Non Woven Balls 10 x Non Woven Swabs 7.5cm x 7.5cm 2 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	10
----------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LUMBAR PUNCTURE PACK

DEF2746D	1 x Drape SMS 61cm x 61cm with 14cm Circular Fenestration and Adhesive Corner Tabs 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Microlance Needle 25g with 16mm Regular Wall 1 x Syringe 5ml Luer Lock 2 x Forcep Plastic Green 1 x IV 3000 Dressing 6cm x 8.5cm 5 x Gauze Swabs 5cm x 5cm Plain 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Labelled Specimen Jar 5ml 4 x Labelled Specimen Tube 10ml with 'V' Bottom Wrapped in a Blue Plastic Field 75cm x 75cm	24
----------	---	----

INTERCOSTAL PACK

DEF1609	1 x Needle Holder 16cm 2 x Forcep Spencer Wells 13cm Straight 2 x Forcep Spencer Wells 13cm Curved 1 x Scissors 12.5cm Sharp/Sharp Straight 1 x Scalpel Handle # 4 1 x Forcep Gillies Toothed 15.5cm 6 x Cotton Balls Large 4 x Gauze Swabs 10cm x 10cm Plain 4 x Dressing Towel (Hand Towel) 31.5cm x 34cm 2 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a SMS Blue Field 91cm x 91cm	24
---------	---	----

INTENSIVE CARE PACKS

TRACHY PROCEDURE PACK

DEF2445	1 x Tracheal Dilator Trousseau 2 Prong 1 x Retractor Bergmann 16.5cm	1
---------	---	---

CVC INSERTION PACK

DEF2393A	1 x Scissors 13cm Sharp/Blunt Straight 1 x Needle Holder 13cm 1 x Forcep Adson Micro Toothed 2 x Forcep Spencer Wells 13cm Straight 5 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 2 x Huck Towel 48cm x 60cm Blue 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	25
----------	--	----

ARTERIAL TRAY

DEF2161A	1 x Drape Polybacked 33cm x 46cm with 8cm Square Fenestration 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Scissors 13cm Sharp/Blunt Straight 1 x Scissors Iris 11.5cm Sharp/Sharp 1 x Needle Holder 13cm 2 x Forceps Artery Mosquito 12.5cm Straight 2 x Forceps Artery Mosquito 12.5cm Curved 1 x Forcep Iris Plain 10cm 1 x Forcep Iris Toothed 10cm 1 x Probe Lachrymal Size 1-2 1 x Probe Lachrymal Size 3-4 6 x Gauze Swabs 7.5cm x 7.5cm Plain 2 x Huck Towel 48cm x 60cm Blue Wrapped in a Clear Plastic Field 75cm x 75cm	5
----------	--	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

STANDARD SURGICAL PACKS

C-SECTION STANDARD SURGICAL PACK

DEF5620	<p>1 x Caesarean Drape 249cm x 312cm x 196cm with Pouch and Circular Fenestration and Built In Incise</p> <p>1 x Baby Blanket 75cm x 100cm</p> <p>1 x Yankauer Suction Handle Bulbous Tip</p> <p>1 x Surgical Blade #23</p> <p>1 x Umbilical Cord Clamp</p> <p>10 x Gauze Swabs 10cm x 10m X-Ray</p> <p>5 x Sponges 30cm x 50cm with X-Ray</p> <p>1 x Maternity Pad 10cm x 34cm</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Mayo Cover 128cm x 60cm</p> <p>1 x Placenta Basin 3L</p> <p>1 x Lid for Placenta Basin</p> <p>1 x Quiver</p> <p>1 x Bowl 1L Blue</p> <p>3 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>3 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	---	---

ENT STANDARD SURGICAL PACK

DEF5621	<p>1 x ENT Split Drape 200cm x 260cm with 7.5cm x 65cm Split</p> <p>1 x Head Turban Drape 75cm x 100cm</p> <p>1 x Surgical Blade #15</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Bowl 1L Blue</p> <p>1 x Kidney Dish 700ml Clear</p> <p>2 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	--	---

HAND STANDARD SURGICAL PACK

DEF5622	<p>1 x Hand Drape with Armboards 370cm x 290cm x 194cm with Reinforcement and Elastic Aperture with 5cm Round Fenestration</p> <p>1 x Syringe 10ml Luer Lock</p> <p>1 x Needle 18g Blunt</p> <p>1 x Needle 23g</p> <p>1 x Surgical Blade #15</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>2 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	--	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FOOT STANDARD SURGICAL PACK

DEF5623	<p>1 x Foot Extremity Drape 305cm x 330cm x 230cm with Reinforcement and Elastic Aperture with 5cm Round Fenestration</p> <p>1 x Poly U Drape 152cm x 183cm</p> <p>1 x Syringe 20ml Luer Lock</p> <p>1 x Needle 18g Blunt</p> <p>1 x Needle 23g</p> <p>1 x Surgical Blade #10</p> <p>1 x Surgical Blade #15</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Bowl 1L Blue</p> <p>1 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>2 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	---	---

HIP STANDARD SURGICAL PACK

DEF5625	<p>1 x Hip Drape 290cm x 350cm x 217cm with Bilateral Pouches, Reinforcement and Elastic Aperture with Tear Drop Fenestration</p> <p>1 x Poly U Drape 152cm x 213cm</p> <p>1 x Stockinette 30cm x 122cm</p> <p>1 x Surgical Blade #10</p> <p>1 x Yankauer Suction Handle Bulbous Tip</p> <p>1 x Suction Tubing 35mm (ID) x 9.91mm (OD) x 3.6m</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>5 x Sponges 50cm x 30cm with X-Ray</p> <p>5 x Sponges 30cm x 30cm with X-Ray</p> <p>1 x Crepe Bandage 15cm x 4.5m</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Mayo Stand Cover Heavy Duty 76cm x 145cm</p> <p>1 x Quiver</p> <p>1 x Bowl 1L Blue</p> <p>3 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>3 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	--	---

STANDARD SURGICAL SET UP KIT

DEF5624	<p>2 x Prep Foam Sticks 16cm Blue</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Mayo Cover 128cm x 60cm</p> <p>2 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>2 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p>	9
---------	--	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

KNEE STANDARD SURGICAL PACK

DEF5626	<p>1 x Extremity Drape with Armboards 315cm x 365cm x 230cm with Reinforcement and Elastic Aperture with 6cm Round Fenestration</p> <p>1 x Poly U Drape 152cm x 183cm</p> <p>1 x Stockinette 30cm x 122cm</p> <p>1 x Yankauer Suction Handle Bulbous Tip</p> <p>1 x Suction Tubing 35mm (ID) x 9.91mm (OD) x 3.6m</p> <p>1 x Surgical Blade #10</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>5 x Sponges 50cm x 30cm with X-Ray</p> <p>5 x Sponges 30cm x 30cm with X-Ray</p> <p>2 x Crepe Bandage 15cm x 4.5m</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Quiver</p> <p>1 x Bowl 1L Blue</p> <p>3 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>3 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	--	---

LITHOTOMY STANDARD DAY STAY PACK

DEF5627	<p>1 x Lithotomy T-Drape 240cm x 195cm x 75cm with Reinforcement and 7.5cm x 15cm Fenestration</p> <p>1 x Under Buttocks Drape 90cm x 120cm with Cuff</p> <p>2 x Lithotomy Legging 75cm x 150cm</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Bowl 1L Blue</p> <p>1 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>2 x Gallipot 150ml</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p>	4
---------	---	---

MINOR STANDARD SURGICAL PACK

DEF5628	<p>1 x Head Sheet 135cm x 250cm</p> <p>1 x Bottom Sheet 185cm x 190cm</p> <p>2 x Side Sheet 100cm x 110cm</p> <p>1 x Surgical Blade #10</p> <p>1 x Surgical Blade #15</p> <p>5 x Gauze Swabs 10cm x 10m X-Ray</p> <p>1 x Blue Huck Towel 40cm x 68cm</p> <p>1 x Bowl 1L Blue</p> <p>1 x Kidney Dish 700ml Clear</p> <p>1 x Kidney Dish 700ml Yellow</p> <p>2 x Gallipot 150ml</p> <p>1 x Needle Counter 20 Count</p> <p>1 x Paper Rubbish Bag 30cm x 16cm with Adhesive</p> <p>Wrapped in a Back Table Cover 150cm x 230cm</p> <p>GOWN PACK</p> <p>1 x Unreinforced Gown Level 3 Large</p> <p>2 x Paper Hand Towel 30cm x 60cm</p> <p>Wrapped in a SMS Wrap</p>	4
---------	---	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MAJOR STANDARD SURGICAL PACK		
DEF5629	1 x Head Sheet 135cm x 250cm 1 x Bottom Sheet 185cm x 190cm 2 x Side Sheet 110cm x 183cm 2 x Drape Tape 50cm x 10cm 1 x Yankauer Suction Handle Bulbous Tip 1 x Suction Tubing 35mm (ID) x 9.91mm (OD) x 3.6m 1 x Surgical Blade #15 1 x Surgical Blade #20 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Blue Huck Towel 40cm x 68cm 1 x Mayo Cover 128cm x 60cm 1 x Quiver 1 x Bowl 1L Blue 2 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow 3 x Gallipot 150ml 1 x Needle Counter 20 Count 1 x Paper Rubbish Bag 30cm x 16cm with Adhesive Wrapped in a Back Table Cover 150cm x 230cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 2 x Paper Hand Towel 30cm x 60cm Wrapped in a SMS Wrap	4

TURP STANDARD SURGICAL PACK		
DEF5630	1 x TURP Drape All-In-One with Pouch 160cm x 166cm x 160cm 1 x Suction Tubing 35mm (ID) x 9.91mm (OD) x 3.6m 5 x Gauze Swabs 10cm x 10m X-Ray 2 x Blue Huck Towel 40cm x 68cm 2 x Kidney Dish 700ml Clear 2 x Gallipot 150ml 1 x Paper Rubbish Bag 30cm x 16cm with Adhesive Wrapped in a Back Table Cover 150cm x 230cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 2 x Paper Hand Towel 30cm x 60cm Wrapped in a SMS Wrap	6

MEDICAL DEVICE SYMBOLS

CUSTOM PACK CHECK SHEETS

Most of our Custom Procedure Packs include a Check Sheet listing the pack’s componentry. The Check Sheet may also include the medical device symbols below.



MANUFACTURER - indicates the medical device manufacturer and their address.



DATE OF MANUFACTURE - indicates the date when the medical device was manufactured.



USE-BY DATE - indicates the date which the medical device must be used by.



BATCH CODE - indicates the manufacturer’s batch code so that the batch or lot can be identified.



CATALOGUE NUMBER - indicates the manufacturer’s catalogue number/ product code so that the medical device can be identified.



STERILISED USING ETHYLENE OXIDE - indicates that a medical device has been sterilised using ethylene oxide.



SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE - the outer packaging is a single sterile barrier. Inside, protective packaging is designed to prevent damage to the contents or to help with aseptic presentation.



MEDICAL DEVICE - indicates the item is a medical device.



DO NOT REUSE - indicates a medical device that is intended for one single use only.



DOES NOT CONTAIN NATURAL RUBBER LATEX - indicates that natural rubber latex as a material of construction within the medical device or the packaging of the medical device is not present.



KEEP AWAY FROM SUNLIGHT - indicates a medical device that needs protection from light sources and/or heat.



KEEP DRY - indicates a medical device that needs to be protected from moisture.



PLASTIC REDUCTION - procedure pack contains an item/s that are an alternative to a plastic product. Examples include a bowl made from sugar cane pulp instead of plastic, a product is wrapped with a paper band instead of double packed in plastic or a paper disposal bag has been included in place of a plastic bag.



SUSTAINABLE SOLUTIONS - procedure pack contains an item/s that is more sustainable than a previously used product. For example, our PurCotton swabs are manufactured using a manufacturing process that transforms raw cotton into a finished product in only 8 hours, unlike traditional gauze which can take up to 6 weeks, therefore minimising the bio-burden on the environment.