



CUSTOM PROCEDURE AND STANDARD PACKS

CUSTOM PROCEDURE PACKS

COST EFFECTIVE AND TIME EFFICIENT: A PACK FOR EVERY PROCEDURE

Our Custom Procedure Pack range covers all areas of the hospital from Cardiac, Orthopaedic, Vascular, ENT and more. They contain all the individual consumables for a surgical procedure, wrapped in a 'field' and sterilised. All products are ready to use when opened, unlike individually wrapped consumables that are touched many times prior to use and require opening separately. Plus there is reduced wastage. Using a Procedure pack saves on setup time for procedures and ensures that all components are readily available for every case. Ordering and inventory control is easier.

CUSTOM TAILOR A PROCEDURE PACK TO MEET YOUR SPECIFIC NEEDS

Select individual consumables such as patient drapes and surgical gowns, to diathermy, suction tubing, abdominal packs, gauze and instruments. Your pack will only include items specifically chosen by you, the surgeon and theatre nurses. Defries assembles these customised packs in cleanrooms according to hospital specifications.

Custom procedure packs can be varied to include cost-effective components as they become available, and allow easy accountability of costs per procedure. As a local Melbourne manufacturer, we are able to respond quickly to any changes requested by a hospital to the configuration of your custom packs.

We'll provide you with ongoing education, listen to feedback and are always available to discuss your specific needs.



READY-FOR-USE STANDARD PROCEDURE PACKS

We have a range of standard procedure packs such as the anaesthetic packs, epidural packs, regional trays, catheter packs, Bronny IV insertion pack, lumbar puncture packs, ICU packs including central line packs, arterial line packs, all designed, made and ready to be ordered.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PACKS

All our anaesthetic packs come sterile, prepacked and ready to use. Indications for use include:

- General and Local Anaesthetic procedures
- Spinal/Epidural/Regional/Lumbar Puncture procedures
- IV Access

NRFIT components are available on request for inclusion into packs.

Select from the packs below or contact us to customise your own anaesthetic pack.

ANAESTHETIC PACKS

DEF1409	1 x Dressing Tray 24cm x 15cm x 4cm 2 x Swabs 10cm x 10cm Green	36
---------	--	----



DEF1409

DEF1162	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 3 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	60
---------	---	----



DEF1162

ANAESTHETIC TRAY

DEF1524	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 2 x Gauze Swabs 10cm x 10cm Green 1 x Towel 40cm x 42cm 2ply	80
DEF406	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 5 x Gauze Swabs 10cm x 10cm Green 1 x Polybacked Towel 40cm x 42cm	60

ANAESTHETIC PACKS

DEF1618	1 x Drape Blue Plastic 75cm x 75cm with 10cm x 15cm Fenestration with Adhesive Tabs on Corners of Fenestration and Corners of Drape 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 5 x Gauze Swabs 10cm x 10cm Green 1 x Forcep Plastic Green 1 x Syringe 5ml Luer Slip 1 x Syringe 3ml Luer Slip 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt Drawing up 1 x Needle 25g x 1" 1 x Dressing Towel Quick Dry 33cm x 33cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
---------	---	----



DEF1618

ANAESTHETIC BASIC TRAY AND SYRINGE KIT

DEF4117E	1 x Anaesthetic Tray with 4 Compartments 22.5cm x 14cm x 2.5cm 1 x Syringe 3ml Luer Slip 1 x Syringe 5ml Luer Slip 1 x Syringe 10ml Luer Slip Eccentric Tip 1 x Syringe 30ml Luer Slip Eccentric Tip 2 x Gauze Swabs 10cm x 10cm Green Wrapped in Crepe Paper 75cm x 75cm	50
----------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ANAESTHETIC PROCEDURE PACK

DEF2255C	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 91cm x 115cm 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 3ml Luer Lock 1 x Syringe 5ml Luer Lock 1 x Needle 25g x 1.5" 2 x Filter Needle 5um 2 x Forceps Pearman Plastic 2 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 1 x Gallipot 60ml 1 x Label Sheet - 1 Set of 16 Wrapped in a Blue Plastic Field 91cm x 115cm	60
----------	--	----

NRFIT REGIONAL ANAESTHETIC PACK

DEF5964	1 x Procedure Drape 120cm x 140cm with Adhesive and 10cm x 15cm Fenestration 1 x Drape 100cm x 100cm with Adhesive 1 x Syringe 1ml Luer Slip 1 x Syringe 5ml Luer Slip 1 x Needle 18g x 1.5" Blunt 1 x NRFIT Needle 19g x 1.5" with Filter 1 x Needle 23g x 1.25" 1 x NRFIT Syringe 1ml Luer Slip 1 x NRFIT Syringe 5ml Luer Slip 5 x Gauze Swabs 10cm x 10cm Green 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow 2 x Gallipot 60ml 1 x Label Sheet - 1 Set of 28 1 x Clipseal bag (for small components) Wrapped in a Back Table Cover 140cm x 190cm	10
---------	--	----

ANAESTHETIC/SPINAL PACK

DEF2232B	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Drape Plastic 76cm x 106cm 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 5ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 19g x 1.5" 1 x Needle 25g x 1" 1 x Forceps Plastic Green 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm Green 5 x Prep Swabs Large XRD 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Label Sheet - 1 Set of 7 Wrapped in a Blue Plastic Field 75cm x 75cm	60
----------	--	----



DEF2232B

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
EPIDURAL PACKS			COMPLETE EPIDURAL PACK		
DEF1796	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 2 x Forceps Plastic Green 5 x Gauze Swab 10cm x 10cm Green 1 x Dressing Towel 32cm x 60cm Wrapped in a Blue Plastic Field 91cm x 115cm	50	DEF1480	1 x Drape Epidural 76cm x 106cm Clear Plastic with Full 2.5cm Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Forceps Plastic Green 1 x Syringe 5ml Luer Slip 1 x Syringe 3ml Luer Slip 1 x Needle 19g x 1.5" 1 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1" 5 x Gauze Swab 10cm x 10cm Green 1 x Dressing Towel (Hand Towel) 32cm x 32.5cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1655B	1 x SMS Drape 91cm x 91cm with 12cm x 15cm Fenestration and Adhesive 1 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Forceps Plastic Green 5 x Gauze Swab 10cm x 10cm Green 1 x Gallipot 60ml 1 x Towel 40cm x 42cm 2ply 1 x Pack Labels - 1 Set of 16 Wrapped in a Blue Plastic Field 91cm x 115cm	60	EPIDURAL/SPINAL TRAY		
DEF3085B	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and ACorner Tabs 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 3ml Luer Slip 1 x Syringe 5ml Luer Slip 2 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1" 1 x Needle 25g x 1.5" 5 x Gauze Swabs 10cm x 10cm Green 1 x Dressing Towel 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	50	DEF2345C	1 x Drape Epidural 120cm x 140cm with 2.5cm Full Adhesive around the 10cm x 15cm Fenestration and Adhesive on Corner Tabs 1 x Dressing Tray with 3 Compartments 20cm x 25cm x 2.5cm 1 x Syringe 10ml Luer Lock 1 x Syringe 5ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 25g x 1.5" 1 x Filter Needle 5um 1 x Foam Applicator Blue 15cm 5 x Gauze Swabs 10cm x 10cm Green 2 x Dressing Towel (Hand Towel) 32cm x 32.5cm 1 x Gallipot 60ml 1 x Plastic Bag 20cm x 20cm(with small components enclosed) Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF3854	1 x Drape Laminated 120cm x 140cm with 7cm x 10cm Fenestration Blue 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Sponge Holder 19cm Blue Plastic 5 x Gauze Swabs 10cm x 10cm Green 2 x Gallipot 150ml 1 x Gown Surgical T-Cut Unferinforced XL 1 x Huck Towel 40cm x 68cm Blue Wrapped in a Back Table Cover 152cm x 228cm	15			

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CATHETER PACKS

All Catheter Procedure Packs have been designed in consultation with the healthcare professionals who use them. They include all the consumables used to insert or change a catheter aseptically. Alternatively please contact us to customise a catheter pack to your own specific requirements.

CATHETER PACKS



DEF1601

DEF1601	1 x Drape Polybacked 50cm x 60cm with a 10cm Circular Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Blue Huck Towel 40cm x 68cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF2074	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Pearman 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 2 x Kidney Dishes 700ml Clear 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24

DEF2061A	1 x Drape Polybacked 75cm x 105cm with 5cm x 15cm Fenestration 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Sponge Holder 19cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Mediwipe Towel 35cm x 60cm 1 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Green 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	24
DEF1262A	1 x Drape SMS 61cm x 61cm with a 10cm x 10cm Fenestration 1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 1 x Syringe 10ml Luer Slip Concentric Tip 1 x Forcep Plastic Pearman 1 x Tube Clamp with Lock 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1650B	1 x Drape Azure Minor Procedure 50cm x 60cm with 6.5cm Circular Fenestration and Adhesive around Fenestration 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm Wrapped in a Blue Plastic Field 75cm x 75cm	60
DEF1522	3 x Drape Polybacked 40cm x 42cm 2 x Dressing Tray with 3 Compartments 23cm x 9cm x 3cm 1 x Syringe 10ml Luer Slip 2 x Forceps Plastic Green 5 x Gauze Swabs 10cm x 10cm Plain Wrapped in a Blue Plastic Field 75cm x 75cm	60

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MEDICAL PACKS

CATHETER PACKS

DEF1036	1 x Drape SMS 91cm x 91cm with a 10cm Circular Fenestration 2 x Forceps Plastic Green 6 x Cotton Balls Large 1 x Dressing Towel T 31.5cm x 61cm 2 x Kidney Dishes 700ml Clear 2 x Gallipots 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
DEF1932A	1 x Drape PolyBacked 50cm x 60cm with 10cm Circular Fenestration 1 x Syringe 30ml Luer Slip Eccentric Tip 1 x Sponge Holder 24cm Plastic 5 x Gauze Swabs 10cm x 10cm X-ray 1 x Dressing Towel 32cm x 32.5cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	30

IV INSERTION PACKS

For the insertion of an IV line.

BRONNY IV PACK

DEF920	1 x Plastic Drape 45cm x 60cm 2 x Gauze Swabs 5cm x 5cm Plain	200
--------	--	-----

BRONNY IV PACK WITH 4 SWABS

DEF1676	1 x Plastic Drape 45cm x 60cm 4 x Gauze Swabs 5cm x 5cm Plain	200
---------	--	-----

MOUTH TOILET PACK

DEF399	1 x Dressing Tray 23cm x 9cm x 3cm with 3 Compartments 10 x M & T Swab Sticks 15cm Paper Shaft 1 x Gallipot 150ml Wrapped in a Blue Plastic Field 75cm x 75cm	50
--------	--	----

FINE PROCEDURE PACK

DEF2453A	1 x Drape Poly Backed 50cm x 60cm with a 10cm Circular Fenestration and Adhesive Tabs in Drape Corners 1 x Scissor Iris 11.5cm Sharp/Sharp 1 x Forcep Adson 12.5cm Micro Toothed 1 x Needle Holder 13cm 1 x Forcep Mosquito Artery 12.5cm Curved 1 x Scalpel Handle # 3 4 x Gauze Swabs 7.5cm x 7.5cm Plain 1 x Kidney Dish 700ml Clear 1 x Gallipot 60ml Wrapped in a Blue Plastic Field 75cm x 75cm	30
----------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

INTENSIVE CARE PACKS

THYROID PACK

DEF3799B	1 x Forcep Spencer Wells 13cm Straight 1 x Scissors 13cm Sharp/Blunt 1 x Staple Remover Scissor Action with Plastic Handle 1 x Stitch Cutter 10 x Gauze Swabs 10cm x 10cm Plain 1 x Huck Towel 40cm x 68cm Blue Wrapped in a SMS Blue Field 61cm x 61cm	20
----------	---	----

TRACHY PROCEDURE PACK

DEF2445	1 x Tracheal Dilator Trousseau 2 Prong 1 x Retractor Bergmann 16.5cm	1
---------	---	---

LARGE DRESSING PACK

DEF2234A	1 x Dressing Tray 24.5cm x 28cm x 4.5cm 3 x Forceps Plastic Green 10 x Non Woven Balls 10 x Non Woven Swabs 7.5cm x 7.5cm 2 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	10
----------	--	----

CVC INSERTION PACK

DEF2393A	1 x Scissors 13cm Sharp/Blunt Straight 1 x Needle Holder 13cm 1 x Forcep Adson Micro Toothed 2 x Forcep Spencer Wells 13cm Straight 5 x Non Woven Balls 5 x Gauze Swabs 10cm x 10cm Plain 2 x Towel Cotton 40cm x 68cm Blue 1 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a Blue Plastic Field 91cm x 115cm	25
----------	--	----

LUMBAR PUNCTURE PACK

DEF2746E	1 x Drape SMS 61cm x 61cm with 14cm Circular Fenestration and Adhesive Corner Tabs 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Microlance Needle 25g with 16mm Regular Wall 1 x Syringe 5ml Luer Lock 2 x Forcep Plastic Green 1 x Asguard Transparent Film Dressing 6cm x 7cm High-tack Adhesion 5 x Gauze Swabs 5cm x 5cm Plain 1 x Dressing Towel (Hand Towel) 31.5cm x 34cm 1 x Labelled Specimen Jar 5ml 4 x Labelled Specimen Tube 10ml with 'V' Bottom Wrapped in a Blue Plastic Field 75cm x 75cm	24
----------	---	----

INTERCOSTAL PACK

DEF1609	1 x Needle Holder 16cm 2 x Forcep Spencer Wells 13cm Straight 2 x Forcep Spencer Wells 13cm Curved 1 x Scissors 12.5cm Sharp/Sharp Straight 1 x Scalpel Handle # 4 1 x Forcep Gillies Toothed 15.5cm 6 x Cotton Balls Large 4 x Gauze Swabs 10cm x 10cm Plain 4 x Dressing Towel (Hand Towel) 31.5cm x 34cm 2 x Kidney Dish 700ml Clear 2 x Gallipot 150ml Wrapped in a SMS Blue Field 91cm x 91cm	24
---------	---	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

STANDARD SURGICAL PACKS

ARTERIAL TRAY

DEF2161A	1 x Drape Polybacked 33cm x 46cm with 8cm Square Fenestration 1 x Dressing Tray with 3 Compartments 25.5cm x 20.5cm x 2.5cm 1 x Scissors 13cm Sharp/Blunt Straight 1 x Scissors Iris 11.5cm Sharp/Sharp 1 x Needle Holder 13cm 2 x Forceps Artery Mosquito 12.5cm Straight 2 x Forceps Artery Mosquito 12.5cm Curved 1 x Forcep Iris Plain 10cm 1 x Forcep Iris Toothed 10cm 1 x Probe Lachrymal Size 1-2 1 x Probe Lachrymal Size 3-4 6 x Gauze Swabs 7.5cm x 7.5cm Plain 2 x Huck Towel 40cm x 68cm Blue Wrapped in a Blue Plastic Field 75cm x 75cm	5
----------	---	---

UNIVERSAL GENERAL PACK 1

DEF5767	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Unreinforced Gown Level 3 X-Large 1 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 5 x Gauze Swabs 10cm x 10cm X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Foam Applicator 20cm Orange 2 x Light Handle Cover Green 1 x Quiver 1 x Needle Counter Magnet/Foam Strips 30 Count 1 x Labels for Theatre - sheet of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	4
---------	---	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
UNIVERSAL GENERAL PACK 2			UNIVERSAL GENERAL PACK 3		
DEF5768A	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Unreinforced Gown Level 3 X-Large 1 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 1 x Needle 23g x 1.25" 1 x Needle 18g Blunt 1 x Scalpel Blade #15 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Prep Swab Large X-Ray 2 x Light Handle Cover Green 1 x Bowl 1L Blue 2 x Kidney Dish 700ml Clear 3 x Gallipot 150ml 2 x Sharps Tray 500ml Yellow 1 x Sharps Safety Transfer Device 1 x Needle Counter Foam Strips with Clear Lid 15 Count Yellow 1 x Labels for Theatre - sheet of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm Heavy Duty GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	3	DEF5769	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 195cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 2 x Fabric Reinforced Gown Level 3 Extra Large 2 x Hand Towel Blue Paper 1 x Syringe 20ml Luer Lock 1 x Needle 23g x 1.25" 1 x Scalpel Blade #15 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Sponges 30cm x 30cm with X-Ray 1 x Prep Swab Large X-Ray 1 x Needle Counter Foam Strips with Clear Lid 15 Count Yellow 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Fabric Reinforced Gown Level 3 Large 1 x Hand Towel Blue Paper Wrapped in a 91cm x 91cm SMS Wrap	4

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
UNIVERSAL GENERAL PACK 4			INTRAOCULAR SURGICAL PACK 1		
DEF5770B	1 x Top Drape 135cm x 250cm with Reinforcement and Adhesive 2 x Side Drape 100cm x 130cm with Reinforcement and Adhesive 1 x Bottom Drape 185cm x 190cm with Reinforcement and Adhesive 1 x Syringe 20ml Luer Lock 5 x Gauze Swabs 10cm x 10m X-Ray 5 x Gauze Swabs 45cm x 45cm X-Ray 1 x Gauze Prep Ball X-ray 1 x Bowl 500ml Blue 4 x Kidney Dish 700ml Clear 2 x Kidney Dish 700ml Yellow 1 x Kidney Dish 1L Blue 3 x Gallipot 150ml 1 x Needle Counter Foam Strips/Magnetic 30 Count 1 x Quiver 1 x Labels for Theatre - 1 Set of 16 1 x Paper Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm Heavy Duty with Full Reinforcement GOWN PACK 1 x Poly Reinforced Gown Level 3 X-Large 2 x Poly Reinforced Gown Level 3 XX-Large 6 x Hand Towel Blue Paper Wrapped in a SMS Wrap 121cm x 121cm	3	DEF5900A	1 x Ophthalmic Drape 140cm x 160cm with 10cm x 10cm Fenestration 1 x Incise Drape 10cm x 20cm 1 x Plastic Drape 45cm x 60cm with Adhesive Tape 1 x Eye Shield 2 x Steristrips 12mm x 100mm - 0.27mm per 3 1 x Syringe 1ml Luer Lock 1 x Syringe 3ml Luer Lock 1 x Needle 18g x 1.5" Blunt 1 x Needle 26g x 1/2" 1 x Anterior Chamber Cannula Angled Rycroft 30G x 7/8" (0.3 X 22mm) Hub Colour: Yellow 1 x Hydrodissector Curved Kellan 25G x 7/8" (0.5 x 22mm) Hub Colour: Orange 1 x Instrument Wipe 5 x Gauze Swabs 10cm x 10m X-Ray 3 x Gauze Ball XRD 1 x Eye Pad 6cm x 8cm Non Woven 1 x Sponge Holder 19cm Plastic 1 x Bowl 250ml Blue 2 x Bowl 500ml Blue 1 x Kidney Dish 700ml Yellow 1 x Kidney Dish 1L Blue 1 x Gallipot 150ml 2 x Pill Cup/ Measuring Cup 60ml 1 x Label Sheet Ophthalmic - 1 Set of 64 1 x Paper Disposal Bag 30cm x 16.5cm with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Unreinforced Gown Level 3 Large 1 x Unreinforced Gown Level 3 X-Large 2 x Hand Towel 40cm x 68cm Cotton Blue Wrapped in a SMS Wrap 121cm x 121cm	6

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
CATARACT SURGICAL PACK 1			MINOR LITHOTOMY SURGICAL PACK 1		
DEF4261A	1 x Ophthalmic Drape 140cm x 160cm with 10cm x 10cm Fenestration 1 x Plastic Drape 45cm x 61cm with Adhesive Tape 2 x Steristrips 12mm x 100mm - 0.27mm per 3 1 x Tegaderm Dressing 10cm x 12cm 1 x Syringe 1ml Luer Lock 1 x Syringe 1ml Luer Slip 3 x Syringe 3ml Luer Lock 2 x Needle 18g x 1.5" Blunt 1 x Needle 26g x 1/2" 2 x Anterior Chamber Cannula Angled Rycroft 30G x 7/8" (0.3 X 22mm) Hub Colour: Yellow 5 x Gauze Swabs 10cm x 10m X-Ray 4 x Prep Gauze Swab Small X-Ray 1 x Eye Shield 1 x Eye Pad 6cm x 8cm Non Woven 5 x Eye Spear 17mm x 7mm Pointed Tip 1 x Bowl 250ml Blue 2 x Kidney Dish 700ml Clear 3 x Gallipot 150ml 1 x Needle Counter Foam Strips/Magnetic 30 Count Yellow 1 x Label Sheet Ophthalmic - 1 Set of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 112cm x 229cm GOWN PACK 1 x Surgical Gown Level 3 Large 1 x Surgical Gown Level 3 X-Large 2 x Hand Towel 31.5cm x 61cm Wrapped in a SMS Wrap 121cm x 121cm	6	DEF5774	1 x Lithotomy T Drape 140cm x 195cm x 74cm with 6cm x 14cm Fenestration and Adhesive 2 x Lithotomy Leggings 80cm x 140cm 5 x Gauze Swabs 10cm x 10m X-Ray 1 x Prep Swab Large X-Ray 1 x Bowl 1L Blue 1 x Kidney Dish 700ml Blue 2 x Gallipot 150ml 1 x Labels for Theatre - 1 Set of 64 1 x Paper Disposal Bag with small components enclosed Wrapped in a Back Table Cover 152cm x 228cm GOWN PACK 1 x Fabri Reinforced Gown Level 3 Large 1 x Fabri Reinforced Gown Level 3 X-Large 2 x Hand Towel Blue Paper Wrapped in a SMS Wrap 91cm x 91cm	4