



**PERSONAL
PROTECTIVE
APPAREL**

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ISOLATION GOWN

Full coverage level 2 isolation gown. Made from a two-layered polypropylene/polyethylene laminate material. With hook and loop neck closure, elastic wrist and thermal bonded sleeve seams.



60-3109

60-3109	Large, Length 117cm, Sleeve 58cm Chest 74cm, Yellow, Knitted Cuff Non Sterile	100
60-3110	Extra Large, Length 122cm Sleeve 61cm, Chest 79cm, Yellow Knitted Cuff, Non Sterile	100

CHEMOTHERAPY GOWNS

Chemoguard Chemotherapy Gowns are a lightweight gown made from a poly-coated spunbound fabric with a closed-back design for full coverage. Adjustable tie neck closure for ease of donning and doffing. A generous fit with long sleeves and fitted knit cuffs.

60-3200	Large Blue Non Sterile	60
60-3201	Extra Large Blue Non Sterile	60

SPEyes EYEWEAR



Our popular DeRoyal SPEyes eyewear is designed to minimise exposure to pathogens from waste fluid. They combine innovation, quality, value, and comfort for clinicians where personal protective equipment is a vital part of self-care.

- Latex free

- Wrap-around disposable lens for complete protection
- UltraClear™ technology prevents glare, fog, and static
- Frame fits completely against forehead for splash protection from above
- Sealed lens mounting for added fluid protection
- Clip-On/Off design for quick and easy changing
- Coloured top strip allows finger-print-free handling



23-505

23-505	SPEyes Eye Shields Pre-assembled in Assorted Colours, Non Sterile	100
23-506	Dispenser box with pre-assembled SPEyes Eye Shields, Non Sterile (20 per box 6 boxes per case)	120



23-506

23-501	SPEyes Eye Shield Replacement Lens, Non Sterile	100
23-502	SPEyes Eye Shield Frame Non Sterile	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

FACE SHIELD

Our Face Shields feature an anti-fog visor with vented foam. Each mask has a peel away protective cover with tabs to prevent scratches. Each shield is 23cm long and 30cm wide.



DEF5526

DEF5526	Face Shield Clear Non Sterile	10
FACE-SHIELD-FULL	Face Shield Clear Non Sterile	300

BERET CAPS



DEF1142

DEF1142	53cm (21") Blue Non Sterile	250
---------	-----------------------------	-----

OVERSHOES



DEF1145

DEF1145	Non Skid Universal Overshoes Polypropylene Non Sterile	100
---------	--	-----

PROTECTIVE GLASSES

Anti-fog glasses with high impact polycarbonate frames.



DEF1755

DEF1755	Disposable Glasses Clear Polycarbonate Non Sterile	1
DEF1755-250	Disposable Glasses Clear Polycarbonate Non Sterile	250



SURGICAL GOWNS

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL GOWNS

Our single-use surgical gowns have specifically been developed for the Australian and New Zealand environment. Our range encompasses gowns for short and dry procedures through to high fluid and lengthy procedures.

OUR COMFORT RANGE – OUR SOFTEST AND MOST COMFORTABLE GOWNS

- Comfort Surgical Gown – Level 2.
- Comfort Unreinforced with Glove Assist – Level 3 gown with a reinforced sleeve to allow extra strength during gloving, but still cool, soft and ideal for short and dry procedures.
- Comfort Plus Fabric Reinforced Surgical Gown – Level 3 gown includes fluid resistant fabric reinforced panels yet remains lightweight and durable.
- Comfort Plus Poly Reinforced Surgical Gown – Level 3 gown which includes poly reinforced impervious panels for added protection.
- Comfort Zoned Protection – Fully impervious zoned panel on front and sleeves. Made from high tech BVB (breathable viral barrier) fabric, but still soft, comfortable and cool which meets AAMI Level 4 fabric standards.

OUR CLASSIC RANGE:

- Classic Standard Surgical Gown – suitable for most surgical procedures.
- Classic Plus Reinforced Surgical Gown – Includes a reinforced panel on the front and sleeves, for surgical cases requiring extra protection.
- Coolmax Surgical Gown – breathable and cool, these gowns are notable for their superb comfort and being fully impervious. With AAMI Level 4 fabric.

Designed in close co-operation with hospital staff, our Surgical Gowns are cool and breathable with an easy fit, while offering excellent protection against strikethrough from microbial and liquid penetration.




CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COMFORT SURGICAL GOWN

LEVEL 2 COMFORT UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 35gsm T-cut Surgical Gown ideal for dry and short procedures.




Comfort Unreinforced

DEF5316	Medium 119cm	25
DEF5317	Large 127cm	25
DEF5318	X-Large 137cm	25
DEF5319	XL-XLong 149cm	25

LEVEL 3 COMFORT UNREINFORCED WITH GLOVE ASSIST

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41gsm T-cut Surgical Gown ideal for dry and short procedures. Has innovative Glove Assist fabric reinforcement to sleeve area to ensure extra strength. Comfortable and cool to wear.




Comfort Unreinforced With Glove Assist

DEF5301	Medium 119cm	24
DEF5302	Large 127cm	24
DEF5303	X-Large 137cm	24
DEF5339	X-Large Twin Pack	24
DEF5304	XL-XLong 150cm	24
DEF5305	XX-Large 155cm	24

LEVEL 3 COMFORT PLUS FABRIC REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41gsm T-cut Surgical Gown with fabric reinforcement on front and sleeves. Fluid resistant on panels. Comfortable to wear.




Comfort Plus Fabric Reinforced

DEF5306	Medium 119cm	24
DEF5307	Large 127cm	24
DEF5308	X-Large 137cm	24
DEF5328	X-Large Twin Pack	24
DEF5309	XL-XLong 150cm	24
DEF5310	XX-Large 155cm	24
DEF5329	XX-Large Twin Pack	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 3 COMFORT PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with two hand towels. Sterile 41 gsm T-cut Surgical Gown which includes poly reinforced panels. Durable and impervious to fluid on panels. Comfortable to wear.




Comfort Plus Poly Reinforced

DEF5311	Medium 119cm	24
DEF5312	Large 127cm	24
DEF5327	Large & X-Large Twin Pack	24
DEF5313	X-Large 137cm	24
DEF5326	X-Large Twin Pack	24
DEF5314	XL-XLong 150cm	24
DEF5315	XX-Large 155cm	24
DEF5330	XX-Large Twin Pack	24

LEVEL 4 COMFORT ZONED PROTECTION

Surgical Gown T-cut wrapped in a 91cm x 91cm SMS field with two hand towels. Made from a 41gsm SMMS material with a polyethylene breathable reinforcement on the sleeves and the front of the gown providing protection against fluids, bacteria and viruses.



Comfort Zoned Protection


DEF5334	Medium 109cm	24
DEF5331	Large 117cm	24
DEF5338	Large & X-Large Twin Pack	24
DEF5332	X-Large 127cm	24
DEF5336	X-Large Twin Pack	24
DEF5333	XX-Large 133cm	24
DEF5759	XL-XLong 135cm	24
DEF5337	XX-Large Twin Pack	24
DEF5335	XXX-Large 135cm	24



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

LEVEL 4 COMFORT MAX

Cool and comfortable T-cut surgical gown made from a 60gsm breathable impervious trilaminate. Provides protection from fluid, bacteria and viruses. Wrapped in a SMS field with two hand towels.




Comfort Max

DEF6020	Medium 119cm	24
DEF6021	Large 127cm	24
DEF6022	X-Large 137cm	24
DEF6023	X-Large Long 150cm	24
DEF5024	XX-Large 155cm	24

LEVEL 3 CLASSIC PLUS POLY REINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown with impervious poly reinforcement. Colour coded neck band – purple for level 3 Gown.




Classic Plus Poly Reinforced

DEF5005A	Medium 119cm	25
DEF5006A	Large 130cm	25
DEF5007A	X-Large 140cm	25
DEF5009A	XL-XLong 157cm	25
DEF5054	XX-Large with Long Arms 157cm	25

CLASSIC SURGICAL GOWN

LEVEL 3 CLASSIC UNREINFORCED

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile 50gsm T-cut Surgical Gown ideal for short and dry procedures.




Classic Unreinforced Gown

DEF5000A	Medium 119cm	25
DEF5001A	Large 130cm	25
DEF5002A	X-Large 140cm	25
DEF5004A	XL-XLong 157cm	25

LEVEL 4 COOLMAX

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile Level 4, trilaminate material which is a BVB fabric – offers a Breathable Viral Barrier, providing excellent protection against fluids, bacteria and viruses.

High quality T-cut Surgical Gown ideal for high fluid procedures or lengthy procedures. Colour Coded neck band – blue for level 4 Gown.



Coolmax

DEF5010B	Medium 119cm	25
DEF5011B	Large 130cm	25
DEF5012B	X-Large 140cm	25
DEF5013B	XL-XLong 157cm	25
DEF5014B	XXL-Large 157cm	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SCRUB SUIT

Scrub Suits are made from a comfortable and soft non-woven fabric, with ultrasonic seam stitching. Waist tie to adjust pants. The top also features two pockets. Each top and pants set is non-sterile, individually wrapped. Ideal for nurses in the operating room, ER Departments, GP Clinics, Dentists, Vets and Students.



SCRUB-SUIT-LGE-GREEN

SCRUB-SUIT-SML-GREEN	Small Top 71cm Long x 54cm Wide Pants - 105cm Long x 51cm Waist	50
SCRUB-SUIT-MED-GREEN	Medium Top 74cm Long x 57cm Wide Pants 106cm Long x 53m Waist	50
SCRUB-SUIT-LGE-GREEN	Large Top 76cm Long x 61cm Wide Pants 108cm Long x 57m Waist	50
SCRUB-SUIT-XL-GREEN	Extra Large Top 78cm Long x 63cm Wide Pants 112cm Long x 61m Waist	50
SCRUB-SUIT-XXXL-GREEN	Extra Extra Extra Large Top 82cm Long x 77cm Wide Pants 121cm Long x 72cm Waist	50