



WOUND CARE

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAEMOSTATIC BANDAGES

First aid and haemostatic products developed in Japan to simplify haemostasis treatment. These adhesive bandages protect and secure punctured skin after an IV injection or catheter removal. Used for various situations such as covering a punctured site following an injection, to stop the flow of blood from an intravenous line using pressure immobilisation, and to treat bleeding of a peripheral vein.

IV PRESSURE PAD

These IV Pressure pads are suitable to use following the removal of an IV, after giving an intravenous injection or after drawing blood or immunisation injections. Comprises of a 36mm round adhesive with a 16mm haemostatic pad.

The 3mm thick round pad is made of highly absorbent non-woven cellulose and creates a compression effect, resulting in bleeding stopping quickly. After bleeding has stopped, remove the pressure pad slowly and carefully.

DEF929	IV Pressure Pad 36mm White	100
--------	----------------------------	-----



DEF929

STEPTY HAEMOSTATIC BANDAGE

These pressure pads comprise of a rectangular 4 way stretch adhesive with a thick compression haemostatic pad. Designed as a protective covering of peripheral veins after blood sampling, a blood transfusion or an infusion.

Two products available with different pad thicknesses.

DEF1251	Stepty Pressure Pad 39mm x 80mm, 6mm Thick, White	50
---------	--	----



DEF1251

DEF1252	Stepty P Pressure Pad 40mm x 120mm, 9mm Thick, White	30
---------	---	----



DEF1252

CHUSHAVAN – FIRST AID ADHESIVE BANDAGE

These adhesive bandages are ideal as a protective covering after an injection. They are 27mm x 27mm adhesive size with a pad size of 13mm x 13mm with 1mm thickness.

DEF1250	Chushavan Pressure Pad Medium, Packet of 8	Box of 60 packets
---------	---	-------------------



DEF1250

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BANDAGES

Excellent range of sterile bandages that can be used in the operating room.

MEDIUM WEIGHT CREPE

Medium weight elastic crepe bandage provides support for soft tissue injuries and will secure a dressing in place. Non fraying edges and a soft, comfortable cotton.



DEF268

DEF268	10cm x 1.6m Single Packed	48
DEF269	15cm x 1.6m Single Packed	48

HEAVY DUTY CREPE

Sterile, heavy weight supporting crepe bandage that provides a firm support to joint injuries. Has a centre guide for assisting application and is comfortable to wear.



DEF272

DEF270	5cm x 2.3m Double Packed	48
DEF271	7.5cm x 2.3m Double Packed	48
DEF272	10cm x 2.3m Double Packed	48
DEF273	15cm x 2.3m Double Packed	24

HOSPITAL WEIGHT CREPE

Made from hospital quality, porous cotton with non-fraying edges. Ideally suited for light compression and holding a dressing in place.



DEF822

DEF3782	2.5cm x 1.5m Single Packed	48
DEF819	2.5cm x 1.5m Double Packed	48
DEF3783	5cm x 1.5m Single Packed	48
DEF820	5cm x 1.5m Double Packed	48
DEF3784	7.5cm x 1.5m Single Packed	48
DEF821	7.5cm x 1.5m Double Packed	48
DEF3780	10cm x 1.5m Single Packed	48
DEF822	10cm x 1.5m Double Packed	48
DEF3781	15cm x 1.5m Single Packed	36
DEF823	15cm x 1.5m Double Packed	36

COHESIVE BANDAGE

Sterile, lightweight, highly conformable bandage that is non slipping. This bandage is intended for general support and compression of soft tissue injuries.



DEF1222A

DEF1222A	10cm x 2m Double Packed	18
DEF3152	15cm x 2m Double Packed	12

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ORTHOPAEDIC PADDING

CONFORMING BANDAGE

Soft non adhesive bandage ideal as a secondary dressing for high body movement areas such as joints. Provides light two-way stretch and conforms easily.



DEF260

DEF1140	2.5cm x 1.5m, Double Packed	48
DEF258	5cm x 1.5m Double Packed	48
DEF260	10cm x 1.5m Double Packed	48
DEF261	15cm x 1.5m Double Packed	48

WATER RESISTANT UNDERCAST PADDING

Comfortable and smooth orthopaedic padding featuring high cushioning and protection for a fracture. Easily conformable and is resistant to moisture absorbency.



DEF824SP

DEF829	5cm x 2.7m Double Packed	48
DEF830SP	7.5cm x 2.7m Single Packed	48
DEF824SP	10cm x 2.7m Single Packed	48
DEF825SP	15cm x 2.7m Single Packed	24

MONTEX SUPER SPONGE BANDAGE

Super absorbent bandage made of prewashed fluff-dried woven cotton gauze with a crinkle-weave pattern for loft and bulk to cushion and protect wound areas. Finished edges are designed to reduce loose ends and lint.



DEF5544

DEF5544	11.4cm x 3.7m 6ply	120
---------	--------------------	-----

COTTON UNDERCAST PADDING

100% cotton undercast padding to provide soft cushioning and excellent absorbency.



DEF5708

DEF5710	5cm x 2.75m Double Packed	200
DEF5707	7.5cm x 2.75m Double Packed	150
DEF5708	10cm x 2.75m Double Packed	150
DEF5709	15cm x 2.75m Double Packed	80

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ABSORBENT UNDERCAST PADDING

Absorbent padding which is multilayered to keep the skin dry and comfortable. Provides excellent comfort and cushioning.



DEF275


DEF828	7.5cm x 2.75m Double Packed	48
DEF275	10cm x 2.75m Single Packed	48
DEF276	15cm x 2.75m Single Packed	24

SURGICAL TAPE

Used to secure a dressing onto a wound, surgical tapes have a low allergenic adhesive which is designed to hold firmly onto the skin or dressing materials but removes easily without damaging the skin.

ADHESIVE SURGICAL FIXATION TAPE – F

Sterile surgical tape that is used for easy fixation in the operating theatre. Coated with a hypoallergenic adhesive whilst the backing is a stretch non woven polyester. Suitable for fixation of wound dressings or fixation of tubes, catheters or instruments.




DEF1504

DEF1504	10cm x 12cm	100
DEF1505	10cm x 20cm	100
DEF1506	10cm x 45cm	100
DEF1614	15cm x 70cm	50
DEF1507	20cm x 20cm	100
DEF3943	30cm x 45cm	100

ADHESIVE SURGICAL FIXATION TAPE – H

Excellent fixation tape which offers stretch and conformability and a low allergenic adhesive. Non woven fabric dressing sheet designed for rapid and quick fixation. Pre-cut and sterilised to a variety of sizes for immediate use in the operating theatre.



DEF1044

DEF1630	2.5cm x 15cm	100
DEF1623	2.5cm x 30cm	50
DEF1647	5cm x 8cm	50
DEF1995	5cm x 15cm	66
DEF1461	5cm x 24cm	50
DEF1599	5cm x 25cm	50
DEF1624	5cm x 30cm	50
DEF1797	5cm x 50cm	60
DEF1044	10cm x 10cm	100
DEF943	10cm x 20cm	50
DEF944	10cm x 25cm	50
DEF1057	10cm x 36.5cm	100
DEF1798	10cm x 50cm	60
DEF2195	15cm x 7cm	140
DEF1600	15cm x 30cm	50
DEF946	15cm x 40cm	25
DEF1026	15cm x 60cm	100
DEF1025	15cm x 100cm	100
DEF2160	20cm x 8cm	125
DEF937	20cm x 25cm	25
DEF678	20cm x 30cm	100
DEF475	20cm x 75cm Double Packed	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

COTTON TAPE FOR RETENTION

100% pure cotton tape pre-cut and sterilised ready for use in the operating room. Ideal for retention of tubes or retractors during a procedure.



DEF167

DEF173	1.5cm x 60cm Double Packed	100
DEF2106	1.5cm x 1.5m	100
DEF918	1.2cm x 500m, Non Sterile	1
DEF919	2.5cm x 250m, Non Sterile	1

TRACHEOSTOMY SPONGE

Soft, absorbent gauze ideal to use around any drainage area to prevent chaffing.

DEF053	10cm x 10cm Plain	100
--------	-------------------	-----



DEF053

DEF2165	10cm x 10cm Y Cut Non woven, Packet of 2	100
---------	---	-----



DEF2165

CR3900-ADVAZORB-T	8cm x 8cm Standard Hydrophilic Foam Tracheostomy Dressing	20
-------------------	---	----



CR3900-ADVAZORB-T

SILICONE TAPE – SILTAPE ROLL

Soft silicone perforated tape roll by Advancis Medical.

Gentle on the skin and easy removal for inspection of wound or dressings. Ideal for use on fragile or sensitive skin.

CR3938-SILTAPE	2cm x 3m, Non Sterile	1
CR3938-SILTAPE-12	2cm x 3m, Non Sterile	12



CR3938-SILTAPE

CR3939-SILTAPE	4cm x 1.5m, Non Sterile	1
CR3939-SILTAPE-12	4cm x 1.5m, Non Sterile	12