



WOUND CARE

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAEMOSTATIC BANDAGES

First aid and haemostatic products developed in Japan to simplify haemostasis treatment. These adhesive bandages protect and secure punctured skin after an IV injection or catheter removal. Used for various situations such as covering a punctured site following an injection, to stop the flow of blood from an intravenous line using pressure immobilisation, and to treat bleeding of a peripheral vein.

IV PRESSURE PAD

These IV Pressure pads are suitable to use following the removal of an IV, after giving an intravenous injection or after drawing blood or immunisation injections. Comprises of a 36mm round adhesive with a 16mm haemostatic pad.

The 3mm thick round pad is made of highly absorbent non-woven cellulose and creates a compression effect, resulting in bleeding stopping quickly. After bleeding has stopped, remove the pressure pad slowly and carefully.

DEF929	IV Pressure Pad 36mm White	100
--------	----------------------------	-----



DEF929

STEPTY HAEMOSTATIC BANDAGE

These pressure pads comprise of a rectangular 4 way stretch adhesive with a thick compression haemostatic pad. Designed as a protective covering of peripheral veins after blood sampling, a blood transfusion or an infusion.

Two products available with different pad thicknesses.

DEF1251	Stepty Pressure Pad 39mm x 80mm, 6mm Thick, White	50
---------	--	----



DEF1251

DEF1252	Stepty P Pressure Pad 40mm x 120mm, 9mm Thick, White	30
---------	---	----



DEF1252

CHUSHAVAN – FIRST AID ADHESIVE BANDAGE

These adhesive bandages are ideal as a protective covering after an injection. They are 27mm x 27mm adhesive size with a pad size of 13mm x 13mm with 1mm thickness.

DEF1250	Chushavan Pressure Pad Medium, Packet of 8	Box of 60 packets
---------	--	-------------------



DEF1250

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BANDAGES

Excellent range of sterile bandages that can be used in the operating room.

HOSPITAL WEIGHT CREPE

Made from hospital quality, porous cotton with non-fraying edges. Ideally suited for light compression and holding a dressing in place.



DEF822

DEF3782	2.5cm x 1.5m 108gsm Single Packed	48
DEF819	2.5cm x 1.5m 108gsm Double Packed	48
DEF5871	5cm x 1.5m 75gsm Single Packed	80
DEF3783	5cm x 1.5m 108gsm Single Packed	48
DEF820	5cm x 1.5m 108gsm Double Packed	48
DEF5872	7.5cm x 1.5m 75gsm Single Packed	80
DEF3784	7.5cm x 1.5m 108gsm Single Packed	48
DEF821	7.5cm x 1.5m 108gsm Double Packed	48
DEF5873	10cm x 1.5m 75gsm Single Packed	80
DEF3780	10cm x 1.5m 108gsm Single Packed	48
DEF822	10cm x 1.5m 108gsm Double Packed	48
DEF5874	15cm x 1.5m 75gsm Single Packed	80
DEF3781	15cm x 1.5m 108gsm Single Packed	36
DEF823	15cm x 1.5m 108gsm Double Packed	36

MEDIUM WEIGHT CREPE

Medium weight elastic crepe bandage provides support for soft tissue injuries and will secure a dressing in place. Non fraying edges and a soft, comfortable cotton.



DEF268

DEF268	10cm x 1.6m Single Packed	48
DEF269	15cm x 1.6m Single Packed	48

HEAVY DUTY CREPE

Sterile, heavy weight supporting crepe bandage that provides a firm support to joint injuries. Has a centre guide for assisting application and is comfortable to wear.



DEF272

DEF270	5cm x 2.3m Double Packed	48
DEF271	7.5cm x 2.3m Double Packed	48
DEF272	10cm x 2.3m Double Packed	48
DEF273	15cm x 2.3m Double Packed	24

MONTEX SUPER SPONGE BANDAGE

Super absorbent bandage made of prewashed fluff-dried woven cotton gauze with a crinkle-weave pattern for loft and bulk to cushion and protect wound areas. Finished edges are designed to reduce loose ends and li



DEF5544

DEF5544	11.4cm x 3.7m 6ply	120
---------	--------------------	-----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ORTHOPAEDIC PADDING

CONFORMING BANDAGE

Soft non adhesive bandage ideal as a secondary dressing for high body movement areas such as joints. Provides light two-way stretch and conforms easily.



DEF260

DEF1140	2.5cm x 1.5m, Double Packed	48
DEF258	5cm x 1.5m Double Packed	48
DEF259	7.5cm x 1.5m Double Packed	48
DEF260	10cm x 1.5m Double Packed	48
DEF261	15cm x 1.5m Double Packed	48

WATER RESISTANT UNDERCAST PADDING

Comfortable and smooth orthopaedic padding featuring high cushioning and protection for a fracture. Easily conformable and is resistant to moisture absorbency.



DEF824SP

DEF829	5cm x 2.7m Double Packed	48
DEF830SP	7.5cm x 2.7m Single Packed	48
DEF824SP	10cm x 2.7m Single Packed	48
DEF825SP	15cm x 2.7m Single Packed	24

COHESIVE BANDAGE

Sterile, lightweight, highly conformable bandage that is non slipping. This bandage is intended for general support and compression of soft tissue injuries.

DEF5888	10cm x 2m Single Packed	48
DEF1222A	10cm x 2m Double Packed	18
DEF3152	15cm x 2m Double Packed	12



DEF1222A

COTTON UNDERCAST PADDING

100% cotton undercast padding to provide soft cushioning and excellent absorbency.



DEF5708

DEF5710	5cm x 2.75m Double Packed	200
DEF5707	7.5cm x 2.75m Double Packed	150
DEF5708	10cm x 2.75m Double Packed	150
DEF5709	15cm x 2.75m Double Packed	80

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ABSORBENT UNDERCAST PADDING

Absorbent padding which is multilayered to keep the skin dry and comfortable. Provides excellent comfort and cushioning.



DEF275


DEF828SP	7.5cm x 2.75m Single Packed	48
DEF275	10cm x 2.75m Single Packed	48
DEF276	15cm x 2.75m Single Packed	24

SURGICAL TAPE

Used to secure a dressing onto a wound, surgical tapes have a low allergenic adhesive which is designed to hold firmly onto the skin or dressing materials but removes easily without damaging the skin.

ADHESIVE SURGICAL FIXATION TAPE – F

Sterile surgical tape that is used for easy fixation in the operating theatre. Coated with a hypoallergenic adhesive whilst the backing is a stretch non woven polyester. Suitable for fixation of wound dressings or fixation of tubes, catheters or instruments.




DEF1504

DEF1504	10cm x 12cm	100
DEF1505	10cm x 20cm	100
DEF1506	10cm x 45cm	100
DEF1614	15cm x 70cm	50
DEF1507	20cm x 20cm	100
DEF3943	30cm x 45cm	100

ADHESIVE SURGICAL FIXATION TAPE – H

Excellent fixation tape which offers stretch and conformability and a low allergenic adhesive. Non woven fabric dressing sheet designed for rapid and quick fixation. Pre-cut and sterilised to a variety of sizes for immediate use in the operating theatre.



DEF1044

DEF1630	2.5cm x 15cm	100
DEF1623	2.5cm x 30cm	50
DEF1647	5cm x 8cm	50
DEF1995	5cm x 15cm	66
DEF1461	5cm x 24cm	50
DEF1599	5cm x 25cm	50
DEF1624	5cm x 30cm	50
DEF1797	5cm x 50cm	60
DEF1044	10cm x 10cm	100
DEF943	10cm x 20cm	50
DEF944	10cm x 25cm	50
DEF1057	10cm x 36.5cm	100
DEF1798	10cm x 50cm	60
DEF2195	15cm x 7cm	140
DEF1600	15cm x 30cm	50
DEF946	15cm x 40cm	25
DEF1026	15cm x 60cm	100
DEF1025	15cm x 100cm	100
DEF2160	20cm x 8cm	125
DEF937	20cm x 25cm	25
DEF678	20cm x 30cm	100
DEF475	20cm x 75cm Double Packed	50

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

TRACHEOSTOMY SPONGE

Soft, absorbent gauze ideal to use around any drainage area to prevent chaffing.

DEF053	10cm x 10cm Plain	100
--------	-------------------	-----



DEF053

DEF2165	10cm x 10cm Y Cut Non woven, Packet of 2	100
---------	---	-----



DEF2165

CR3900-ADVAZORB-T	8cm x 8cm Standard Hydrophilic Foam Tracheostomy Dressing	20
-------------------	---	----



CR3900-ADVAZORB-T

COTTON TAPE FOR RETENTION

100% cotton tape pre-cut and sterilised non-medical device. Ready for use where a sterile securing tape is required.



DEF167


DEF173	1.5cm x 60cm Double Packed	100
DEF2106	1.5cm x 1.5m	100
DEF918	1.2cm x 500m, Non Sterile	1
DEF919	2.5cm x 250m, Non Sterile	1

POVIDONE-IODINE DRESSING

MEGHdin® is a non-adherent povidone-iodine dressing designed to provide antimicrobial protection while helping to maintain a moist wound environment. Its low-adherent construction helps minimise trauma during dressing changes, making it suitable for a range of acute and chronic wounds.

MEGHDIN DRESSING

MEGHdin® combines broad-spectrum antimicrobial protection with a low-adherent wound contact layer for the management of acute and chronic wounds.



P05-MEGHDIN

P05-MEGHDIN	5cm x 5cm Non Adherent Povidone Iodine Dressing	25
P95-MEGHDIN	9.5cm x 9.5cm Non Adherent Povidone Iodine Dressing	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

MANUKA HONEY DRESSINGS

The Advancis Medical Manuka Honey Dressing Range includes a selection of wound care dressings containing medical-grade Manuka honey. Suitable for use across a variety of wound types, the range is designed to support wound management by helping to maintain an appropriate wound healing environment and assisting clinicians in meeting individual patient needs.

ACTIVON TUBE®

Activon Tube® is a medical-grade Manuka honey dressing presented in an easy-to-apply tube format. It is suitable for use on acute and chronic wounds and can be applied directly to the wound bed or used in conjunction with an appropriate secondary dressing.




CR4493-ACTIVON-TUBE

CR4493-ACTIVON-TUBE	20g Manuka Honey	1
CR4493-12-ACTIVON-TUBE	20g Manuka Honey	12

ACTIVON TULLE® DRESSING

Activon Tulle® is a non-adherent viscose mesh dressing impregnated with 100% medical-grade Manuka honey. Designed for direct contact with the wound bed, it helps maintain a moist wound healing environment and can be used in the management of a wide range of acute and chronic wounds. Its flexible structure allows it to conform closely to the wound surface while requiring a secondary dressing.



CR3761-ACTIVON-TULLE

CR3761-ACTIVON-TULLE	5cm x 5cm Honey Coated Tulle	5
CR3658-ACTIVON-TULLE	10cm x 10cm Honey Coated Tulle	5

ALGIVON® ALGINATE DRESSING

Algivon® is a soft alginate dressing containing 100% medical-grade Manuka honey. Designed for use on moderate to heavily exuding wounds, it combines the absorbent properties of alginate with the benefits of Manuka honey to help maintain a moist wound healing environment. Its conformable structure makes it suitable for cavity wounds and areas requiring gentle packing.



CR3831-ALGIVON

CR3831-ALGIVON	5cm x 5cm Alginate with Honey	5
CR3659-ALGIVON	10cm x 10cm Alginate with Honey	5

ALGIVON PLUS® ALGINATE DRESSING

Algivon® Plus is a non-adherent alginate dressing containing 100% medical-grade Manuka honey. Designed for use on wounds with light to moderate exudate, it helps maintain a moist wound healing environment while allowing the dressing to conform closely to the wound bed. The dressing can be cut to size and is suitable for a variety of acute and chronic wounds.



CR4230-ALGIVON-PLUS

CR4230-ALGIVON-PLUS	5cm x 5cm Reinforced Alginate with Honey	5
CR4225-ALGIVON-PLUS	10cm x 10cm Reinforced Alginate with Honey	5

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE DRESSINGS

Advancis Medical's silicone dressing range provide gentle and versatile solutions for wound protection and fixation. Designed to minimise trauma during dressing changes, these silicone-based products support patient comfort while addressing a variety of wound care and skin protection needs.

SILTAPE® SILICONE TAPE

A soft silicone perforated fixation tape designed to secure dressings while minimising skin trauma upon removal. It is particularly suitable for patients with fragile or sensitive skin and can be repositioned during application.

CR3938-SILTAPE	2cm x 3m, Non Sterile	1
CR3938-SILTAPE-12	2cm x 3m, Non Sterile	12



CR3938-SILTAPE

CR3939-SILTAPE	4cm x 1.5m, Non Sterile	1
CR3939-SILTAPE-12	4cm x 1.5m, Non Sterile	12

SILFLEX® SILICONE DRESSING

A soft silicone wound contact dressing designed to protect the wound bed while allowing exudate to pass through to a secondary absorbent dressing. Its gentle adherence helps minimise trauma and pain during dressing changes.



CR3922-SILFLEX

CR3922-SILFLEX	5cm x 7cm Silicone Wound Contact Layer	10
CR3923-SILFLEX	8cm x 10cm Silicone Wound Contact Layer	10
CR3924-SILFLEX	12cm x 15cm Silicone Wound Contact Layer	10
CR3925-SILFLEX	20cm x 30cm Silicone Wound Contact Layer	10
CR4006-SILFLEX	35cm x 60cm Silicone Wound Contact Layer Sterile	5

VELLAFILM® FILM DRESSING

Vellafilm® is a transparent, atraumatic, polyurethane film dressing with Silfix® soft silicone adherence, designed to protect fragile skin and superficial wounds. Its waterproof, breathable construction allows monitoring of the wound site without removing the dressing.



CR3919-VELLAFILM

CR3919-VELLAFILM	12cm x 12cm Adherent Silicone Film Dressing	10
CR3921-VELLAFILM	15cm x 20cm Adherent Silicone Film Dressing	10
CR3917-VELLAFILM	12cm x 35cm Adherent Silicone Film Dressing	10

FOAM DRESSINGS

Advancis Medical's foam dressing range offers absorbent solutions for the management of a wide variety of acute and chronic wounds. Designed to manage exudate and maintain a moist wound healing environment, the range provides options to meet different anatomical and clinical needs.

ADVAZORB SACRAL® FOAM DRESSING

Advazorb Sacral® is a hydrophilic foam dressing designed to contour to the sacral area, for wounds with low to medium levels of exudate. With a soft silicone wound contact layer for gentle adherence, Advazorb Sacral® will only adhere to dry skin surrounding the wound and can be easily lifted for adjustment without losing its adherent properties.



CR4401-ADVAZORB-SACRAL

CR4401-ADVAZORB-SACRAL	17cm x 19cm Hydrophilic Foam Sacral Dressing with Adhesive and Border	10
CR4402-ADVAZORB-SACRAL	22cm x 23cm Hydrophilic Foam Sacral Dressing with Adhesive and Border	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ADVAZORB® FOAM DRESSING

A polyurethane foam dressing with a breathable film backing designed for the management of moderate to heavily exuding wounds. It provides absorbency while helping to maintain a moist wound healing environment.



CR4165-ADVAZORB

CR4165-ADVAZORB	7.5cm x 7.5cm Hydrophilic Foam Dressing with Film Backing	10
CR4166-ADVAZORB	10cm x 10cm Hydrophilic Foam Dressing with Film Backing	10
CR4167-ADVAZORB	12.5cm x 12.5cm Hydrophilic Foam Dressing with Film Backing	10
CR4168-ADVAZORB	15cm x 15cm Hydrophilic Foam Dressing with Film Backing	10
CR4169-ADVAZORB	10cm x 20cm Hydrophilic Foam Dressing with Film Backing	10
CR4170-ADVAZORB	20cm x 20cm Hydrophilic Foam Dressing with Film Backing	10

ADVAZORB SILFIX® DRESSING

Advazorb® Silfix is a soft silicone foam dressing that combines the absorbency with a gentle, atraumatic silicone adhesive layer. Designed for wounds with low to moderate levels of exudate, it helps maintain a moist wound healing environment while minimising trauma during dressing changes, making it suitable for patients with fragile or sensitive skin.



CR4177-ADVAZORB-SILFIX

CR4177-ADVAZORB-SILFIX	7.5cm x 7.5cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4178-ADVAZORB-SILFIX	10cm x 10cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4179-ADVAZORB-SILFIX	12.5cm x 12.5cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4180-ADVAZORB-SILFIX	15cm x 15cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4181-ADVAZORB-SILFIX	10cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive	10
CR4182-ADVAZORB-SILFIX	20cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive	10

ADVAZORB HEEL® FOAM DRESSING

Advazorb Heel® is a shaped foam dressing designed specifically for the heel. Suitable for wounds with low to moderate exudate, it combines absorbent foam with an anatomical design to provide a secure fit, helping to manage exudate and maintain a moist wound healing environment.



CR4228-ADVAZORB-HEEL

CR4228-ADVAZORB-HEEL	17cm x 21cm Hydrophilic Foam Dressing with Film Backing	5
----------------------	---	---

CODE

DESCRIPTION

BOX QTY.

CODE

DESCRIPTION

BOX QTY.

ADVAZORB BORDER® DRESSING

Advazorb Border® is a self-adhesive foam dressing with a soft silicone wound contact layer, designed for wounds with low to moderate exudate. The absorbent foam and perforated silicone layer allow exudate to pass into the dressing while protecting the wound bed, helping to maintain a moist wound environment and minimise trauma during dressing changes.



CR4190-ADVAZORB-BORDER

CR4190-ADVAZORB-BORDER	7.5cm x 7.5cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4191-ADVAZORB-BORDER	10cm x 10cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4192-ADVAZORB-BORDER	12.5cm x 12.5cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4193-ADVAZORB-BORDER	15cm x 15cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4194-ADVAZORB-BORDER	10cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4195-ADVAZORB-BORDER	20cm x 20cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10
CR4196-ADVAZORB-BORDER	10cm x 30cm Hydrophilic Foam Dressing with Silicone Adhesive and Border	10

SUPER ABSORBENT DRESSINGS

The Advancis Medical Eclipse® range offers highly absorbent dressings for the management of moderate to heavily exuding wounds. Designed to provide fluid handling and retention, the range helps protect peri-wound skin and maintain an optimal wound healing environment.

ECLYPSE® DRESSING

Eclipse® is a super absorbent secondary dressing designed for the management of moderate to heavily exuding wounds. Its high-capacity absorbent core rapidly wicks, absorbs and locks away exudate, helping to reduce the risk of leakage and maceration while supporting an optimal wound healing environment. A water-resistant backing helps prevent strike-through and prolong wear time.



CR3818-ECLYPSE

CR3818-ECLYPSE	10cm x 10cm Super Absorbent Dressing	20
CR3769-ECLYPSE	15cm x 15cm Super Absorbent Dressing	20
CR3743-ECLYPSE	20cm x 30cm Super Absorbent Dressing	20
CR3808-ECLYPSE	60cm x 40cm Super Absorbent Dressing	10

ECLYPSE BOOT® DRESSING

Eclipse Boot® is a super absorbent secondary dressing specifically designed for use on the foot and ankle. Its anatomically shaped design provides coverage while the high-capacity absorbent core rapidly wicks, absorbs and locks away exudate, helping to reduce the risk of leakage and maceration. A water-resistant backing helps prevent strike-through and prolong wear time.



CR4237-ECLYPSE-BOOT

CR4237-ECLYPSE-BOOT	55cm x 47cm Super Absorbent Leg Wrap	5
CR4014-ECLYPSE-BOOT	60cm x 70cm Super Absorbent Leg Wrap	5