

Face Shield

CLINICAL PROTECTIVE APPAREL



Specifications

Length: 23cm

Width: 30cm

Vented foam with strap

Peel away protective cover with tabs to prevent scratches

Anti-fog

Non-Sterile, Single use

Latex Free

ARTG 334977

Meets AS/ NZS 1337

Product Number	UOM
DEF5526	Box of 10
FACE-SHIELD-FULL	Box of 300



Please direct any questions to your **Sales Account Manager** or contact our **Customer Service Team** on **1300 550 278**.

FShield_v2 Jul 21