

SURGICAL GOWNS

May 2025



Designed for our local Australian conditions, our surgical gowns have been developed in collaboration with clinicians to protect staff in critical surgical environments.

We have a wide range of gowns suitable for all surgical cases in the operating room environment, from quick and short procedures, through to high fluid complex procedures.


Gowns are available single sterile or for inclusion in custom procedure packs.

SURGICAL GOWNS

FOR EVERY PROCEDURE

SOFT POLYESTER CUFFS | T-CUT DESIGN | VELCRO NECKLINE FOR EASY ADJUSTMENT
INCLUDES PAPER/HUCK TOWELS WITHIN GOWN PACK | WRAPPED IN STERILE FIELD
LEVEL CONFORMS TO AMMI PB70.2013 FOR LIQUID BARRIER PERFORMANCE

ESSENTIAL PROTECTION



ULTIMATE PROTECTION

LEVEL 2 - ideal for low-risk settings and dry and short procedures


Comfort Unreinforced Surgical Gown



Comfortable T-Cut Surgical Gown.

35gsm alcohol-repellent SMMS non-woven fabric.

Wrapped in a SMS field with two purcotton hand towels.



Level 2

Order Code	Size	Length	Width	UOM Qty
DEF5316	Medium	119cm	65cm	25
DEF5317	Large	127cm	67cm	25
DEF5318	X-Large	137cm	75cm	25
DEF5319	XX-Large	149cm	78cm	25

LEVEL 3 - suitable for low-risk settings and dry and short procedures

Comfort Unreinforced Surgical Gown with Glove Assist




Comfortable and light T-Cut Surgical Gown.

41gsm anti-static and alcohol repellent treated SMS non-woven fabric with additional 35gsm alcohol repellent Spunmelt fabric on the sleeves.

Soft and cool to wear.

Has innovative Glove Assist fabric reinforcement on sleeve area for extra strength.

Wrapped in a SMS field with two blue paper hand towels.



Level 3

Order Code	Size	Length	Width	UOM Qty
DEF5301	Medium	119cm	62 cm	24
DEF5302	Large	127cm	67cm	24
DEF5303	X-Large	137cm	72cm	24
DEF5304	X-Large X-Long	150cm	72cm	24
DEF5305	XX-Large	155cm	72cm	24

LEVEL 3 - suitable for low-risk settings and dry and short procedures


Classic Unreinforced Surgical Gown



Comfortable T-Cut Surgical Gown.

50 gsm non-woven SMMS fabric.

Wrapped in a SMS field with one blue huck towel.



Order Code	Size	Length	Width	UOM Qty
DEF5000A	Medium	119cm	71cm	25
DEF5001A	Large	130cm	74cm	25
DEF5002A	X-Large	140cm	82cm	25
DEF5003A	XX-Large	157cm	82cm	25
DEF5004A	X-Large X-Long *	155cm	72cm	25

* Sleeves are 10cm longer than XXL gown

LEVEL 3 - suitable for medium or moderate fluid procedures

Comfort Fabric Reinforced Surgical Gown




Comfortable T-Cut Surgical Gown.

41gsm anti-static and alcohol repellant treated SMS non-woven fabric with an additional 41gsm anti-static and alcohol repellant treated SMS on front panel and sleeves.

Fluid resistant on panels.

Wrapped in a SMS field with two blue paper hand towels.



Order Code	Size	Length	Width	UOM Qty
DEF5306	Medium	119cm	62cm	24
DEF5307	Large	127cm	67cm	24
DEF5308	X-Large	137cm	72cm	24
DEF5328	X-Large Twin Pack	137cm	72cm	24
DEF5309	X-Large X-Long	150cm	72cm	24
DEF5310	XX-Large	155cm	72cm	24
DEF5329	XX-Large Twin Pack	155cm	72cm	24

LEVEL 3 - suitable for medium or moderate fluid procedures

Comfort Poly Reinforced Surgical Gown




Comfortable T-Cut Surgical Gown.

41gsm anti-static and alcohol repellent treated SMS non-woven fabric with 40gsm poly film reinforcement on front and sleeves.

Durable and impervious to fluid on panels.

Wrapped in a SMS field with two blue paper hand towels.



Level 3

Order Code	Size	Length	Width	UOM Qty
DEF5311	Medium	119cm	62cm	24
DEF5312	Large	127cm	67cm	24
DEF5327	Large & X Large Twin Pack	127cm 137cm	67cm 72cm	24
DEF5313	X-Large	137cm	72cm	24
DEF5326	X-Large Twin Pack	137cm	72cm	24
DEF5314	X-Large X-Long	150cm	72cm	24
DEF5315	XX-Large	155cm	72cm	24
DEF5330	XX-Large Twin Pack	155cm	72cm	24

Classic Poly Reinforced Surgical Gown




Comfortable T-Cut Surgical Gown.

43gsm alcohol-repellent SMMS fabric with 40gsm breathable poly reinforcement on front and sleeves.

Durable and impervious to fluid on panels.

Wrapped in a SMS field with one huck towel.



Level 3

Order Code	Size	Length	Width	UOM Qty
DEF5005A	Medium	119cm	71cm	25
DEF5006A	Large	130cm	74cm	25
DEF5007A	X-Large	140cm	82cm	25
DEF5009A	X-Large X-Long	157cm	82cm	25
DEF5054	XX-Large *	157cm	82cm	25

* Sleeves are 10cm longer than XLXL gown

LEVEL 4 - long, complicated or high fluid procedures

Comfort Zoned Protection Surgical Gown




Cool and comfortable T-Cut Surgical Gown.

35gsm SMS with anti-static and alcohol-repellent treatment. With 67gsm breathable sleeves and 34gsm poly reinforcement on the front panel.

Durable and impervious to fluid on panels.

Provides protection from fluid, bacteria and viruses.

Wrapped in a SMS field with two white paper hand towels.



Order Code	Size	Length	Width	UOM Qty
DEF5334	Medium	109cm	63cm	24
DEF5331	Large	117cm	72cm	24
DEF5338	Large & X Large Twin Pack	117cm 127cm	72cm 75cm	24
DEF5332	X-Large	127cm	75cm	24
DEF5336	X-Large Twin Pack	127cm	75cm	24
DEF5333	XX-Large	133cm	80cm	24
DEF5337	XX-Large Twin Pack	133cm	80cm	24
DEF5759	X-Large X-Long	135cm	75cm	24
DEF5335	XXX-Large	135cm	80cm	24

Comfort Max Surgical Gown




Cool and comfortable T-Cut Surgical Gown.

60gsm breathable impervious trilaminate.

Provides protection from fluid, bacteria and viruses.

Wrapped in a SMS field with two blue paper hand towels.



Order Code	Size	Length	Width	UOM Qty
DEF6020	Medium	119cm	62cm	24
DEF6021	Large	127cm	67cm	24
DEF6022	X-Large	137cm	72cm	24
DEF6023	X-Large Long	150cm	72cm	24
DEF6024	XX-Large	155cm	72cm	24

LEVEL 4 - long, complicated or high fluid procedures

Coolmax Surgical Gown




Comfortable T-Cut Surgical Gown.

64gsm BVB trilaminate fabric offering a breathable viral barrier.

Provides protection from fluid, bacteria and viruses.


Wrapped in a SMS field with one blue huck towel.



Level 4

Order Code	Size	Length	Width	UOM Qty
DEF5010B	Medium	119cm	71cm	25
DEF5011B	Large	130cm	74cm	25
DEF5012B	X-Large	140cm	82cm	25
DEF5013B	X-Large X-Long	157cm	82cm	25
DEF5014B	XX-Large *	157cm	82cm	25

* Sleeves are 10cm longer than XLXL gown



Model is 174cm tall and wears a size Large gown.

Gown_May25

Head Office

1-7 Vision Street
Dandenong South VIC 3175
t: 61 3 9706 3600
f: 61 3 9706 3619

freecall: 1300 550 278
e: sales@defries.com.au
www.defries.com.au

Manufacturing Site

18-24 Villas Road
Dandenong South VIC 3175