

SURGICAL GOWNS

November 2024



Designed for our local Australian conditions, our surgical gowns have been developed in collaboration with clinicians to protect staff in critical surgical environments.

We have a wide range of gowns suitable for all surgical cases in the operating room environment, from quick and short procedures, through to high fluid complex procedures.


Gowns are available single sterile or for inclusion in custom procedure packs.

SURGICAL GOWNS

FOR EVERY PROCEDURE

SOFT POLYESTER CUFFS | T-CUT DESIGN | VELCRO NECKLINE FOR EASY ADJUSTMENT
INCLUDES PAPER/HUCK TOWELS WITHIN GOWN PACK | WRAPPED IN STERILE FIELD
LEVEL CONFORMS TO AMMI PB70.2013 FOR LIQUID BARRIER PERFORMANCE

ESSENTIAL PROTECTION



ULTIMATE PROTECTION

LEVEL 2 - ideal for low-risk settings and dry and short procedures


Comfort Unreinforced Surgical Gown



Comfortable T-Cut Surgical Gown.

35gsm alcohol-repellent SMMS non-woven fabric.

Wrapped in a SMS field with two purcotton hand towels.



Level 2

| Order Code | Size | Length | Width | UOM Qty |
|------------|----------|--------|-------|---------|
| DEF5316 | Medium | 119cm | 65cm | 25 |
| DEF5317 | Large | 127cm | 67cm | 25 |
| DEF5318 | X-Large | 137cm | 75cm | 25 |
| DEF5319 | XX-Large | 149cm | 78cm | 25 |

LEVEL 3 - suitable for low-risk settings and dry and short procedures

Comfort Unreinforced Surgical Gown with Glove Assist



Comfortable and light T-Cut Surgical Gown.

41gsm anti-static and alcohol repellent treated SMS non-woven fabric with additional 35gsm alcohol repellent Spunmelt fabric on the sleeves.

Soft and cool to wear.

Has innovative Glove Assist fabric reinforcement on sleeve area for extra strength.

Wrapped in a SMS field with two blue paper hand towels.



Level 3


| Order Code | Size | Length | Width | UOM Qty |
|------------|----------------|--------|-------|---------|
| DEF5301 | Medium | 119cm | 62 cm | 24 |
| DEF5302 | Large | 127cm | 67cm | 24 |
| DEF5303 | X-Large | 137cm | 72cm | 24 |
| DEF5304 | X-Large X-Long | 150cm | 72cm | 24 |
| DEF5305 | XX-Large | 155cm | 72cm | 24 |

LEVEL 3 - suitable for low-risk settings and dry and short procedures

Classic Unreinforced Surgical Gown



Comfortable T-Cut Surgical Gown.
50 gsm non-woven SMMS fabric.
Wrapped in a SMS field with one blue huck towel.



Level 3

| Order Code | Size | Length | Width | UOM Qty |
|------------|------------------|--------|-------|---------|
| DEF5000A | Medium | 119cm | 71cm | 25 |
| DEF5001A | Large | 130cm | 74cm | 25 |
| DEF5002A | X-Large | 140cm | 82cm | 25 |
| DEF5003A | XX-Large | 157cm | 82cm | 25 |
| DEF5004A | X-Large X-Long * | 155cm | 72cm | 25 |


* Sleeves are 10cm longer than XXL gown

LEVEL 3 - suitable for medium or moderate fluid procedures

Comfort Fabric Reinforced Surgical Gown



Comfortable T-Cut Surgical Gown.
41 gsm anti-static and alcohol repellent treated SMS non-woven fabric with an additional 41 gsm anti-static and alcohol repellent treated SMS on front panel and sleeves.
Fluid resistant on panels.
Wrapped in a SMS field with two blue paper hand towels.



Level 3

| Order Code | Size | Length | Width | UOM Qty |
|------------|--------------------|--------|-------|---------|
| DEF5306 | Medium | 119cm | 62cm | 24 |
| DEF5307 | Large | 127cm | 67cm | 24 |
| DEF5308 | X-Large | 137cm | 72cm | 24 |
| DEF5328 | X-Large Twin Pack | 137cm | 72cm | 24 |
| DEF5309 | X-Large X-Long | 150cm | 72cm | 24 |
| DEF5310 | XX-Large | 155cm | 72cm | 24 |
| DEF5329 | XX-Large Twin Pack | 155cm | 72cm | 24 |

LEVEL 3 - suitable for medium or moderate fluid procedures

Comfort Poly Reinforced Surgical Gown




Comfortable T-Cut Surgical Gown.

41gsm anti-static and alcohol repellent treated SMS non-woven fabric with 40gsm poly film reinforcement on front and sleeves.

Durable and impervious to fluid on panels.

Wrapped in a SMS field with two blue paper hand towels.



Level 3

| Order Code | Size | Length | Width | UOM Qty |
|------------|---------------------------|----------------|--------------|---------|
| DEF5311 | Medium | 119cm | 62cm | 24 |
| DEF5312 | Large | 127cm | 67cm | 24 |
| DEF5327 | Large & X Large Twin Pack | 127cm 137cm | 67cm 72cm | 24 |
| DEF5313 | X-Large | 137cm | 72cm | 24 |
| DEF5326 | X-Large Twin Pack | 137cm | 72cm | 24 |
| DEF5314 | X-Large X-Long | 150cm | 72cm | 24 |
| DEF5315 | XX-Large | 155cm | 72cm | 24 |
| DEF5330 | XX-Large Twin Pack | 155cm | 72cm | 24 |

Classic Reinforced Surgical Gown



Comfortable T-Cut Surgical Gown.

43gsm alcohol-repellent SMMS fabric with 40gsm breathable poly reinforcement on front and sleeves.

Durable and impervious to fluid on panels.

Wrapped in a SMS field with one huck towel.



Level 3

| Order Code | Size | Length | Width | UOM Qty |
|------------|----------------|--------|-------|---------|
| DEF5005A | Medium | 119cm | 71cm | 25 |
| DEF5006A | Large | 130cm | 74cm | 25 |
| DEF5007A | X-Large | 140cm | 82cm | 25 |
| DEF5009A | X-Large X-Long | 157cm | 82cm | 25 |
| DEF5054 | XX-Large * | 157cm | 82cm | 25 |

* Sleeves are 10cm longer than XLXL gown

LEVEL 4 - long, complicated or high fluid procedures

Comfort Zoned Protection Surgical Gown




Cool and comfortable T-Cut Surgical Gown.

35gsm SMS with anti-static and alcohol-repellent treatment. With 67gsm breathable sleeves and 34gsm poly reinforcement on the front panel.

Durable and impervious to fluid on panels.

Provides protection from fluid, bacteria and viruses.

Wrapped in a SMS field with two white paper hand towels.



| Order Code | Size | Length | Width | UOM Qty |
|------------|---------------------------|----------------|--------------|---------|
| DEF5334 | Medium | 109cm | 63cm | 24 |
| DEF5331 | Large | 117cm | 72cm | 24 |
| DEF5338 | Large & X Large Twin Pack | 117cm 127cm | 72cm 75cm | 24 |
| DEF5332 | X-Large | 127cm | 75cm | 24 |
| DEF5336 | X-Large Twin Pack | 127cm | 75cm | 24 |
| DEF5333 | XX-Large | 133cm | 80cm | 24 |
| DEF5337 | XX-Large Twin Pack | 133cm | 80cm | 24 |
| DEF5759 | X-Large X-Long | 135cm | 75cm | 24 |
| DEF5335 | XXX-Large | 135cm | 80cm | 24 |

Comfort Max Surgical Gown




Cool and comfortable T-Cut Surgical Gown.

60gsm breathable impervious trilaminate.

Provides protection from fluid, bacteria and viruses.

Wrapped in a SMS field with two blue paper hand towels.



| Order Code | Size | Length | Width | UOM Qty |
|------------|--------------|--------|-------|---------|
| DEF6020 | Medium | 119cm | 62cm | 24 |
| DEF6021 | Large | 127cm | 67cm | 24 |
| DEF6022 | X-Large | 137cm | 72cm | 24 |
| DEF6023 | X-Large Long | 150cm | 72cm | 24 |
| DEF6024 | XX-Large | 155cm | 72cm | 24 |

LEVEL 4 - long, complicated or high fluid procedures

Coolmax Surgical Gown




Comfortable T-Cut Surgical Gown.

64gsm BVB trilaminate fabric offering a breathable viral barrier.

Provides protection from fluid, bacteria and viruses.


Wrapped in a SMS field with one blue huck towel.



Level 4

| Order Code | Size | Length | Width | UOM Qty |
|------------|----------------|--------|-------|---------|
| DEF5010B | Medium | 119cm | 71cm | 25 |
| DEF5011B | Large | 130cm | 74cm | 25 |
| DEF5012B | X-Large | 140cm | 82cm | 25 |
| DEF5013B | X-Large X-Long | 157cm | 82cm | 25 |
| DEF5014B | XX-Large * | 157cm | 82cm | 25 |

* Sleeves are 10cm longer than XLXL gown



Model is 174cm tall and wears a size Large gown.

Gown_Nov24

Head Office

1-7 Vision Street
Dandenong South VIC 3175
t: 61 3 9706 3600
f: 61 3 9706 3619

freecall: 1300 550 278
e: sales@defries.com.au
www.defries.com.au

Manufacturing Site

18-24 Villas Road
Dandenong South VIC 3175