

# Plant-Based Products

SUSTAINABLY DESIGNED MEDICAL CONSUMABLES



Replace your single-use plastic medical products with our plant-based alternatives, clinically developed for use in the operating room and custom procedure packs.

- Made with 100% plant-based materials\*
- An environmentally friendly and biodegradable range
- Coated with a biocompatible film
- Designed for clinical performance
- Sustainable manufacturing process

\*(Needle counters are 95% plastic free – contains EPP foam and magnet)

## SMART SUSTAINABLE DESIGN

By taking an agricultural sugarcane by-product and upcycling it into a sustainable healthcare solution, these products are environmentally friendly and biodegradable.



100% Australian & family owned since 1980

Distributed throughout Australia and New Zealand by:



# Plant-Based Products

SUSTAINABLY DESIGNED MEDICAL CONSUMABLES



Defries Industries is proud to bring this innovative, environmentally friendly product range exclusively to the Australian and New Zealand healthcare market.

## NEWGEN PLANT-BASED HOLLOWWARE

Made from 100% plant-based composite, this range offers an alternative to plastic bowls used in the operating theatre.

- Designed to hold fluids and medical supplies
- Easy to see graduations of ml and oz
- Coated with biocompatible film
- Easy to hold and stable to use
- Latex free

### PLANT-BASED KIDNEY DISH FOR CUSTOM PACKS

| CODE           | DESCRIPTION              |
|----------------|--------------------------|
| NGS700EB-250NS | NEWGEN KIDNEY DISH 700ML |



### PLANT-BASED 6 L BOWL FOR CUSTOM PACKS

| CODE           | DESCRIPTION    |
|----------------|----------------|
| NGS5500RB-80NS | NEWGEN BOWL 6L |



### PLANT-BASED 1 L BOWL FOR CUSTOM PACKS

| CODE          | DESCRIPTION    |
|---------------|----------------|
| NGS32SB-400NS | NEWGEN BOWL 1L |



## NEWGEN PLANT-BASED NEEDLE COUNTERS

The outer box is made with bagasse, a post-agricultural by-product from sugarcane production, that is upcycled to create the fibre pulp. The hinges, latch, foam and magnet are similar to those found in existing needle counters and made from medical grade materials.

- Easy-to-read numbers for securing of needles
- Secure latch with audible closure
- Adhesive tabs secure the needle counter to any workspace
- Full-surface magnet provides space for larger needles and blades
- ASTM F2132 compliant for puncture resistance.

### PLANT-BASED NEEDLE COUNTERS FOR CUSTOM PACKS

| CODE             | DESCRIPTION   |
|------------------|---|
| NGSNC40FMMA-90NS | NEWGEN NEEDLE COUNTER, 40 COUNT, FOAM BLOCK, DOUBLE MAGNET, ADHESIVE TABS |



| CODE              | DESCRIPTION   |
|-------------------|---|
| NGSNC80FFMMA-90NS | NEWGEN NEEDLE COUNTER, 40-80 COUNT, DOUBLE FOAM, DOUBLE MAGNET, ADHESIVE TABS |



NewGen\_Prod\_v3 Sep 21

Distributed throughout Australia and New Zealand by:

T: 03 9706 3600 | Freecall: 1300 550 278  
Email: sales@defries.com.au

