

# Siltape®

Soft silicone perforated tape roll



Siltape® is a soft silicone fixation tape delivered on a handy roll with perforations for ease of use, no scissors are required. Siltape® can be used in any situation where a regular adhesive tape would be used. It is particularly useful where friable or sensitive skin is an issue. It is ideal for use on paediatrics and neonates.



## Features

- ✔ Atraumatic soft silicone, reduced pain on change, adjustment or removal
- ✔ Soft & conformable
- ✔ Adheres only to dry skin not to a moist wound
- ✔ Easily adjustable with no loss of adherence

## Indications

- ✔ Hold dressings in place
- ✔ Taping down eyelids (in theatre)
- ✔ Securing IV lines
- ✔ Over small incision sites
- ✔ Under oxygen masks on bridges of noses
- ✔ Securing fistula needles
- ✔ Protection from pressure damage or sores from tubing

Size	Stock code	Pack size
2cm x 3m	CR3938	1
4cm x 1.5m	CR3939	1

Distributed throughout Australia by:



t: 03 9706 3600 freecall: 1300 550 278  
e: sales@defries.com.au www.defries.com.au

