

> SURGICAL DRAPES





Designed by Australians for Australian hospital conditions.

> SURGICAL DRAPES



basic drapes

Categories for drapes include:

general drapes

equipment covers

- paediatrics
- cardiac
- vascular
- angiography
- EENT
- neuro
- plastics
- ophthalmic
- orthopaedics including general ortho, extremity, shoulder, hip, arthroscopy
- gynaecological
- obstetrics
- urology.

Designed by Australians for Australian hospital conditions.



Head Office 18-24 Villas Road Dandenong South VIC 3175 t: 61 3 9706 3600

f: 61 3 9706 3619 freecall: 1300 550 278 e: sales@defries.com.au

www.defries.com.au

Manufacturing Site 17-20 & 22 Kempson Court Keysborough VIC 3173





The Supreme Drape fabric is unique

This purpose-designed trilaminate fabric offers:

- an absorbent layer
- an impermeable layer
- a patient comfort layer.

The fabric has excellent barrier properties to prevent strikethrough of airborne and blood-borne bacteria. It is both drapeable and low linting.

Additional features for high level performance include:

- Velcro line holders
- armboard cover adhesive
- integrated pouches of various sizes
- custom made to your specifications.

These drapes meet or exceed EN13795, and they meet fire retardant standard CPSC 16 CFR Part 1610.

The Azure Minor Procedure Drapes

- A comfortable fabric combined with a PE laminate
- Appropriate for most types of minor procedures.



