

SURGICAL GOWNS




CLASSIC AND COOLMAX RANGES

Surgical Gowns with varying degrees of protection.

ALL OF OUR CLASSIC & COOLMAX GOWNS:

- are T-cut with velcro neckline tabs for easy adjustment
- have high quality cuffs made with a soft cotton fabric
- meet AAMI according to Level 3 or 4 surgical gown standards

- come wrapped in a 91cm x 91cm SMS field with 40cm x 68cm huck towel
- are sterile and latex free
- come in cartons of 25.



AAMI LEVEL 3

AAMI LEVEL 4

	CLASSIC UNREINFORCED				CLASSIC POLY REINFORCED				COOLMAX SURGICAL GOWN			
SIZE	M	L	XL	XLXL	M	L	XL	XLXL	M	L	XL	XLXL
CODE	DEF5000A	DEF5001A	DEF5002A	DEF5004A	DEF5005A	DEF5006A	DEF5007A	DEF5009A	DEF5010A	DEF5011A	DEF5012A	DEF5013A
FEATURES	<ul style="list-style-type: none"> • Comfortable and light • Soft material and cool to wear 				<ul style="list-style-type: none"> • Poly reinforcement on front panel and sleeves • Durable and impervious to fluid on panels • Comfortable to wear 				<ul style="list-style-type: none"> • Trilaminare Breathable Viral Barrier (BVB) material • High-tech breathable fabric which is fully impervious • Comfortable and cool to wear 			
THICKNESS	50 gsm				50 gsm + PE Zoned				68 gsm			
PROCEDURE	Dry and short				Medium or moderate fluid				High-fluid or lengthy			

MATERIALS KEY



Unreinforced



Poly Reinforced



Trilaminare BVB

Surg_Gown_v2May21