





# SURGICAL GOWNS

## COMFORT RANGE

Surgical Gowns with varying degrees of protection.

ALL OF OUR GOWNS:

- are T-cut with velcro neckline tabs for easy adjustment
- have high quality cuffs made with a soft cotton fabric
- meet AAMI standards according to Level 2, 3 or 4 and exceed EN 13795
- come wrapped in a 91cm x 91cm SMS field with two 40cm x 58cm hand towels
- are sterile and latex free
- come in cartons of 24



### AAMI LEVEL 3

### AAMI LEVEL 4

	COMFORT UNREINFORCED WITH GLOVE ASSIST					COMFORT PLUS FABRIC REINFORCED					COMFORT PLUS POLY REINFORCED					COMFORT ZONED PROTECTION			
SIZE	M	L	XL	XLXL	XXL	M	L	XL	XLXL	XXL	M	L	XL	XLXL	XXL	M	L	XL	XXL
CODE	DEF5301	DEF5302	DEF5303	DEF5304	DEF5305	DEF5306	DEF5307	DEF5308	DEF5309	DEF5310	DEF5311	DEF5312	DEF5313	DEF5314	DEF5315	DEF5334	DEF5331	DEF5332	DEF5333
FEATURES	<ul style="list-style-type: none"><li>• Glove Assist fabric-reinforcement on wrist and sleeve area for extra strength</li><li>• Comfortable and light</li><li>• Soft material and cool to wear</li></ul>					<ul style="list-style-type: none"><li>• Fabric-reinforcement on front panel and sleeves</li><li>• Fluid resistant on panels</li><li>• Comfortable to wear</li></ul>					<ul style="list-style-type: none"><li>• Poly-reinforcement on front panel and sleeves</li><li>• Durable and impervious to fluid on panels</li><li>• Comfortable to wear</li></ul>					<ul style="list-style-type: none"><li>• Critical area zoned gown</li><li>• Main gown is a SMMS material with a PE poly-ethylene breathable reinforcement on sleeves and front of gown</li><li>• Excellent barrier and viral zoned protection</li></ul>			
THICKNESS	41 gsm					41 gsm + Fabric reinforcement Zoned					41 gsm + PE Zoned					41 gsm + PE Zoned			
PROCEDURE	Dry and short					Medium or moderate-fluid					High-fluid					High-fluid or lengthy			

MATERIALS     Unreinforced     Glove Assist     Fabric Reinforced     Poly Reinforced     Zoned