



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

HAEMOSTATIC BANDAGES

First aid and Haemostatic products developed in Japan to simplify haemostasis treatment. These adhesive bandages protect and secure punctured skin after an IV injection or catheter removal. Used for various situations such as covering a punctured site following an injection, to stop the flow of blood from an intravenous line using pressure immobilisation, and to treat bleeding of a peripheral vein.

CHUSHAVAN MEDIUM SIZE – FIRST AID ADHESIVE BANDAGE

These adhesive bandages are ideal as a protective covering after an injection. They are 27mm x 27mm adhesive size with a pad size of 13mm x 13mm with 1mm thickness.

DEF1250	Chushavan Pressure Pad Medium	Box of 60 packets – each packet contains 8 pressure pads
---------	-------------------------------	----------------------------------------------------------



DEF1250

IV PRESSURE PAD

These IV Pressure pads are suitable to use following the removal of an IV, after giving an intravenous injection or after drawing blood or immunisation injections. Comprises of a 36mm round adhesive with a 16mm haemostatic pad. The 3mm thick round pad is made of highly absorbent non-woven cellulose and creates a compression effect, resulting in bleeding stopping quickly. After bleeding has stopped, remove the pressure pad slowly and carefully.

DEF929	IV Pressure Pad White	100
--------	-----------------------	-----



DEF929

STEPTY HAEMOSTATIC BANDAGE

These pressure pads comprise of a rectangular 4 way stretch adhesive with a thick compression haemostatic pad. Designed as a protective covering of peripheral veins after blood sampling, a blood transfusion or an infusion. Two products available with different pad thicknesses.

DEF1251	Stepty Pressure Pad 6mm thick, length 39mm x 80mm, white	50
---------	----------------------------------------------------------	----



DEF1251

DEF1252	Stepty P Pressure Pad 9mm thick, length 40mm x 120mm, white	30
---------	-------------------------------------------------------------	----



DEF1252

ISLAND AND PAD DRESSING

Quality island and pad dressing by Nichiban, Japan. These dressings are gentle, secure and have a long lasting urethane adhesive. The pad is a high quality absorbent cotton, whilst the clear adhesive is breathable, vapour permeable and skin friendly. Ideal for use as a post surgical dressing or for cuts and abrasions. Sterile and latex free.

CPSPO508	Cathereplus Film & Pad Island Dressing 5x8cm (Pad 2.5x5cm)	20
CPSPO915	Cathereplus Film & Pad Island Dressing 9x15cm (Pad 5x10cm)	20
CPSPO921	Cathereplus Film & Pad Island Dressing 9x21cm (Pad 5x16cm)	20

CATHEREEPLUS ROLL

A new urethane gel adhesive which offers high permeability and breathability. Comfortable and gentle on the skin, and easy to remove. Pull out necessary length and cut using scissors.

CPSR1010	100mm x 10m	1
----------	-------------	---

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

BANDAGES

Excellent range of sterile bandages that can be used in the operating room.

MEDIUM WEIGHT CREPE

Medium weight elastic crepe bandage provides support for soft tissue injuries and will secure a dressing in place. Non fraying edges and a soft, comfortable cotton.



DEF268

DEF266	5cm x 1.6m Single Packed	48
DEF267	7.5cm x 1.6m Single Packed	48
DEF268	10cm x 1.6m Single Packed	48
DEF269	15cm x 1.6m Single Packed	48

HOSPITAL WEIGHT CREPE

Made from hospital quality, porous cotton with non-fraying edges. Ideally suited for light compression and holding a dressing in place.



DEF822

DEF819	2.5cm x 1.5m Double Packed	48
DEF820	5cm x 1.5m Double Packed	48
DEF821	7.5cm x 1.5m Double Packed	48
DEF822	10cm x 1.5m Double Packed	48
DEF823	15cm x 1.5m Double Packed	36

COTTON/WOOL CREPE BANDAGE (VIC)

Blended wool & cotton elastic crepe bandage with light support properties. Ideal for providing warmth and absorbency.



DEF517

DEF516	7.5cm x 2.2m Double Packed	48
DEF517	10cm x 2.2m Double Packed	48
DEF518	15cm x 2.2m Double Packed	24

HEAVY DUTY

Sterile, heavy weight supporting crepe bandage that provides a firm support to joint injuries. Has a centre guide for assisting application and is comfortable to wear.



DEF272

DEF270	5cm x 2.3m Double Packed	48
DEF271	7.5cm x 2.3m Double Packed	48
DEF272	10cm x 2.3m Double Packed	48
DEF273	15cm x 2.3m Double Packed	48

COHESIVE BANDAGE

Sterile, lightweight, highly conformable bandage that is non slipping. It will grip against the skin, without any adhesive, so ideal for sensitive skin.

DEF1056	6cm x 1.5m Double Packed	25
---------	--------------------------	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

CONFORMING BANDAGE

Soft, cotton 2ply woven bandage ideal as a secondary dressing for difficult to cover areas of the body. It provides a light, two way stretch and conforms easily.

DEF1140	2.5cm x 1.5m Double Packed	48
DEF1100	5cm x 3.5m Double Packed	48

H CONFORMING GAUZE BANDAGE

Conforming, non adhesive gauze bandage and is a light weight elastic bandage secondary dressing. Ideal for high body movement areas such as joints.



DEF260

DEF258	5cm x 1.5m Double Packed	48
DEF259	7.5cm x 1.5m Double Packed	48
DEF260	10cm x 1.5m Double Packed	48
DEF261	15cm x 1.5m Double Packed	48

ORTHOPAEDIC PADDING

NATURAL UNDERCAST PADDING

Soft, undercast orthopaedic padding made from rayon fleece and used under casts to provide comfort and cushioning for a fracture. Pulls apart easily for a quick application.



DEF825N

DEF825N	15cm x 2.7m Double Packed	24
---------	---------------------------	----

WATER RESISTANT UNDERCAST PADDING

Comfortable and smooth orthopaedic padding featuring high cushioning and protection for a fracture. Easily conformable and is resistant to moisture absorbency.



DEF824

DEF829	5cm x 2.7m Double Packed	48
DEF830	7.5cm x 2.7m Double Packed	48
DEF824	10cm x 2.7m Double Packed	48
DEF825	15cm x 2.7m Double Packed	24

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL TAPE

Surgical or medical tapes are used to secure a dressing onto a wound. They have a low allergenic adhesive which is designed to hold firmly onto the skin or dressing materials but removes easily without damaging the skin.

ABSORBENT UNDERCAST PADDING

Absorbent padding which is multilayered to keep the skin dry and comfortable. Provides excellent comfort and cushioning.



DEF826

DEF828	7.5cm x 2.75m Double Packed	48
DEF826	10cm x 2.75m Double Packed	48
DEF827	15cm x 2.75m Double Packed	24

COTTON UNDERCAST PADDING

100% cotton undercast padding to provide soft cushioning and excellent absorbency. Has a unique crinkle weave pattern for loft and bulk. Can be used with either synthetic or plaster casting material.



DEF675

DEF673	5cm x 3.6m Double Packed	48
DEF674	7.5cm x 3.6m Double Packed	48
DEF675	10cm x 3.6m Double Packed	48
DEF676	15cm x 3.6m Double Packed	24

ADHESIVE SURGICAL FIXATION TAPE – F

Sterile surgical tape that is used for easy fixation in the operating theatre. Coated with a hypoallergenic adhesive whilst the backing is a stretch non woven polyester. Suitable for fixation of wound dressings or fixation of tubes, catheters or instruments.




DEF1504

DEF1504	10cm x 12cm	100
DEF1505	10cm x 20cm	100
DEF1506	10cm x 45cm	100
DEF1614	15cm x 70cm	50
DEF1507	20cm x 20cm	100
DEF3943	30cm x 45cm	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

ADHESIVE SURGICAL FIXATION TAPE – H

Excellent fixation tape which offers stretch and conformability and a low allergenic adhesive. Non woven fabric dressing sheet designed for rapid and quick fixation. Pre-cut and sterilised to a variety of sizes for immediate use in the operating theatre.



DEF1044

DEF1630	2.5cm x 15cm	100
DEF1623	2.5cm x 30cm	50
DEF3174	2.5cm x 100cm	100
DEF1647	5cm x 8cm	50
DEF1995	5cm x 15cm	66
DEF1461	5cm x 24cm	50
DEF1599	5cm x 25cm	50
DEF1624	5cm x 30cm	50
DEF1797	5cm x 50cm	60
DEF1044	10cm x 10cm	100
DEF943	10cm x 20cm	50
DEF944	10cm x 25cm	50
DEF1057	10cm x 36.5cm	100
DEF1798	10cm x 50cm	60
DEF2195	15cm x 7cm	140
DEF1600	15cm x 30cm	50
DEF946	15cm x 40cm	25
DEF1026	15cm x 60cm	100
DEF1025	15cm x 100cm	100
DEF2160	20cm x 8cm	125
DEF937	20cm x 25cm	25
DEF678	20cm x 30cm	100
DEF475	20cm x 75cm Double Packed	50
DEF3942	30cm x 45cm	100

SURGICAL TAPE NICHIMESH

Meshpore Tape from Japan, it is made of an elastic non-woven mesh with a gentle adhesive. Properties include excellent flexibility, softness, breathability and moisture permeability. Ideal for protective covering and securing.



DEF947

DEF947	5cm x 30cm, packet of 10, Double Packed	50
--------	-----------------------------------------	----

COTTON TAPE FOR RETENTION

100% pure cotton tape pre-cut and sterilised ready for use in the operating room. Ideal for retention of tubes or retractors during a procedure.



DEF167

DEF173	1.5cm x 60cm Double Packed	100
DEF2106	1.25cm x 1.5m	100
DEF167	2.5cm x 2m Rondo	100

COTTON TAPE FOR RETENTION ROLL NON STERILE



DEF919

DEF918	1.2cm x 500m	1
DEF919	2.5cm x 250m	1
DEF1591	6mm x 1000m	1

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SILICONE TAPE – SILTAPE ROLL

Soft silicone perforated tape roll by Advancis Medical.
Gentle on the skin and easy removal and inspection of wound or dressings. Ideal use on fragile or sensitive skin.

* CR3938-SILTAPE	Siltape 2cmx 3m	1
------------------	-----------------	---



CR3938

* CR3939	Siltape 4cm x 1.5m	1
----------	--------------------	---



CR3939

* Not available in New Zealand

TRACHEOSTOMY SPONGE

Soft, absorbent gauze ideal to use around any drainage area to prevent chaffing.

DEF053	10cm x 10cm Plain	100
--------	-------------------	-----



DEF053

DEF2165	10cm x 10cm Y Cut non woven, packet of 2	100
---------	------------------------------------------	-----



DEF2165

CR3900	Ultra Thin Hydrophilic Foam Tracheostomy Dressing 8cm x 8cm	20
--------	-------------------------------------------------------------	----

CR4212	Standard Hydrophilic Foam Tracheostomy Dressing 8cm x 8cm	20
--------	-----------------------------------------------------------	----

